Aquadien™ DNA Extraction and Purification Kit
iQ-Check™ Screen *Legionella* spp. Kit
iQ-Check Quanti *Legionella* spp. Kit
iQ-Check Screen *L. pneumophila* Kit
iQ-Check Quanti *L. pneumophila* Kit

User Guide

Tests for the DNA extraction, purification, and real-time PCR detection and quantification of *Legionella* spp. and *Legionella pneumophila* in water samples

Catalog #3578121, Aquadien DNA Extraction and Purification Kit
Catalog #3578104, iQ-Check Screen *Legionella* spp. Kit
Catalog #3578102, iQ-Check Quanti *Legionella* spp. Kit
Catalog #3578105, iQ-Check Screen *L. pneumophila* Kit
Catalog #3578103, iQ-Check Quanti *L. pneumophila* Kit

Table of Contents

Section 1	Introduction	1
Section 2	The iQ-Check <i>Legionella</i> Technology	1
Section 3	Kit Components.....	2
Section 4	Shelf Life and Storage	3
Section 5	Materials Required but Not Supplied.....	3
	Equipment.....	3
	Supplies	3
Section 6	Safety Precautions and Recommendations for Best Results	4
Section 7	Protocol.....	5
	A. Sampling and Transportation	5
	B. Water Filtration and DNA Extraction	5
	C. Real-Time PCR.....	10
	D. Data Analysis.....	11
Section 8	Confirmation of Positive Results	17
Section 9	Confirmation of Single Colonies Using the iQ-Check Kit.....	18
Section 10	Test Performance and Validations	18
Section 11	References.....	19
Section 12	Revision History	19
	Appendix 1 — Membrane Folding Guide	20
	Appendix 2 — PCR Mix Calculation Guide	21

Section 1

Introduction

Legionella are gram-negative bacteria present in all aquatic environments. Infection can cause acute pneumonia, Legionnaires' disease, and a milder form of pulmonary infection, Pontiac fever. In Europe, the number of declared cases of Legionnaires' disease increases each year. In the United States, the Centers for Disease Control and Prevention estimates the number of cases of Legionellosis at between 10,000 and 20,000 per year. *Legionella pneumophila* species is responsible for approximately 90% of all clinical cases.

Regular control of the presence of *Legionella* in water supply systems such as air conditioning cooling towers, spa pools, fountains, hot and cold-water systems, is the only way of preventing the disease. Detection of *Legionella* is compulsory or strongly recommended in most developed countries. Conventional culture methods for detection of *Legionella* spp. and *L. pneumophila* in water samples present several disadvantages, including low sensitivity and long incubation periods (10–14 days for positive samples).

The Aquadien Kit allows an optimal DNA extraction and purification from bacteria present in water samples for subsequent real-time polymerase chain reaction (PCR) detection. The principle of the extraction is based on alkaline lysis of bacteria in the presence of thermal shock and DNA purification using ultrafiltration.

The iQ-Check *Legionella* spp. and iQ-Check *L. pneumophila* Kits are formulated for rapid detection or quantification of *Legionella* (spp. or *pneumophila*) in water samples. Quantification is possible with the use of the calibrated iQ-Check *Legionella* DNA Quantification Standards in the amplification step. The results are obtained in less than 2 hr following the filtration and DNA extraction step.

Section 2

iQ-Check *Legionella* Technology

The iQ-Check *Legionella* Kits are tests based on gene amplification and detection by real-time PCR. Ready-to-use PCR reagents contain DNA primers and a DNA probe specific for *Legionella* spp. or *L. pneumophila*, as well as DNA polymerase and nucleotides. Detection and data analysis are optimized for use with a Bio-Rad real-time PCR instrument, such as the CFX96 Touch Deep Well and CFX Opus Real-Time PCR Systems.


PCR is a powerful technique used to generate many copies of target DNA. During the PCR reaction, several cycles of heating and cooling denature the DNA. Primers then anneal to the target region, where the DNA polymerase uses the primers and deoxynucleotide triphosphates (dNTPs) to extend the DNA, creating copies of the target DNA. These copies are called amplicons.

In real-time PCR, specific probes detect the DNA during amplification by hybridizing to the amplicons. These probes are linked to a fluorophore that fluoresces only when hybridized to the target sequence. FAM is the fluorophore linked to the probe in this kit that hybridizes to the *Legionella* spp.– or *L. pneumophila*–specific DNA sequence. In the absence of target DNA, no fluorescence will be detected. As the number of amplicons increases with each round of amplification, fluorescence intensity also increases. At the annealing step of each PCR cycle, the optical module or detector measures this fluorescence. The software associated with the instrument plots fluorescence intensity against number of cycles. This method allows a simple determination of the presence or absence of *Legionella* spp. or *L. pneumophila* in a sample.

Section 3 Kit Components

A synthetic DNA internal control is included in the reaction mix. This control is amplified with a specific probe at the same time as the *Legionella* spp. or *L. pneumophila* target DNA sequence and is detected by a second fluorophore (HEX). This allows for the validation of any negative result.

The iQ-Check kits can detect or quantify *Legionella* spp. or *L. pneumophila* in water samples in four main steps:

* Please refer to the user guide of the iQ-Check Free DNA Removal Solution (#10000058391) for conditions of use.

Section 3 Kit Components

The Aquadien DNA Extraction and Purification Kit contains sufficient reagents for 96 tests (94 samples).

Reagent ID	Reagent	Quantity Provided, ml
R1	Lysis solution	1 bottle, 100
R2	Elution buffer	1 tube, 25
	Purification column	96 columns
	Collector vials	192 vials

The iQ-Check Screen *Legionella* spp. and *L. pneumophila* Kits contain sufficient reagents for 96 tests (94 samples). The Quanti *Legionella* spp. and *L. pneumophila* Kits additionally contain quantification standards sufficient to quantify up to 43 samples.

Reagent ID	Reagent	Quantity Provided, ml
B	Fluorescent probes	1 tube, 0.55
C	Amplification mix	1 tube, 4.4
D	PCR negative control	1 tube, 0.5
E	PCR positive control	1 tube, 0.25
Qs1	PCR standard 1	1 tube, 0.08 (white cap)
Qs2	PCR standard 2	1 tube, 0.08 (yellow cap)
Qs3	PCR standard 3	1 tube, 0.08 (orange cap)
Qs4	PCR standard 4	1 tube, 0.08 (red cap)
	Reference material ¹	1 tube, 0.08 (blue cap)
	Calibration card	1 card

¹ The reference material contains calibrated DNA connected to the standard reference material, independent from the quantification standards.

Section 4

Shelf Life and Storage

Once received, the kit must be stored at 2–8°C. Reagents stored at this temperature can be used until the expiration date indicated on the tubes. Do not freeze the reagents.

Section 5

Materials Required but Not Supplied

Equipment

- Bunsen burner
- Filtration apparatus mounted either on an air pump or vacuum flask
- Biosafety cabinet — Class II
- Water bath, preferably with lid, capable of maintaining $95 \pm 5^\circ\text{C}$
- Vortexer
- Magnetic stir plate
- Heating thermoshaker* capable of maintaining 95°C , with a mixing speed of 1,300 rpm
- Centrifuge with fixed angle rotor, preferably refrigerated
 - Capable of holding 1.5–2.0 ml tubes with a rotation capacity of 6,000 x g
 - Capable of holding 1.5–2.0 ml tubes with a rotation capacity of 12,000 x g for dirty/clogging water samples
- 20, 200, and 1,000 µl micropipets
- Bio-Rad real-time PCR system;* for example, the CFX96 Touch Deep Well (catalog #3600037) or CFX Opus 96 (catalog # 17007992) or CFX Opus Deepwell (catalog # 17007991) Real-Time PCR Systems

Note: We recommend using an uninterrupted power supply (UPS) with the thermal cycler.

* Contact Bio-Rad Technical Support for information on recommended instruments.

Supplies

- Aquadien Reagent W2, 5 ml, in case of purification column clogging (catalog #3578119)
- Aquadien Lysis Solution R1, 100 ml (catalog #3578125)
- CryoTubes, sterile, 4.5 ml, internal thread, for Aquadien and Aquadien W2 protocols

Section 6 Safety Precautions and Recommendations for Best Results

- *Legionella* DNA Free Water, 1 L (catalog #12006823)
- iQ-Check Free DNA Removal Solution (catalog #3594970)
- Polycarbonate membrane filter, porosity of 0.45 µm
- Disposable sterile funnel, 250 ml
- Sterile metallic inox tweezers with rounded ends
- Sterile filter tips adaptable to 20, 200, and 1,000 µl micropipets
- Sterile tubes, 2 ml
- PCR plates, tubes, sealing tape, and caps
- Powder-free gloves
- Single-use respiratory mask
- Bleach, 5%
- Alcohol, 70%
- Cleaning agent such as DNA away or RNase AWAY

Section 6

Safety Precautions and Recommendations for Best Results

- This test must be performed by trained personnel
- Water samples must be handled as potentially infectious material and discarded according to local rules and regulations
- All potentially infectious material should be autoclaved before disposal
- The quality of results depends on strict compliance with Good Laboratory Practices (for example, the ISO 8199 and ISO/TS 12869 standards), especially concerning PCR:
 - Never circulate laboratory equipment (pipets, tubes, etc.) from one workstation to another
 - Decontaminate small laboratory material used for the filtration (for example, the metal tweezers) after each extraction with alcohol and sterilize it by burning on the Bunsen burner
 - Check that tweezers are not hot before use
 - Always use a positive control and a negative control for each series of amplification reactions
 - A negative control should be used during DNA extraction
 - Do not use reagents after their expiration date
 - Vortex reagents from the kit before using them to ensure homogeneity

- Periodically verify the accuracy and precision of pipets, as well as correct functioning of the instruments
- Change gloves often, especially if you suspect they are contaminated
- Clean work spaces periodically with 5% bleach then rinse with sterile *Legionella* DNA Free Water and 70% alcohol or with other decontaminating agents, such as DNA AWAY
- Use powder-free gloves and avoid fingerprints and writing on tube caps. Both will interfere with data acquisition
- For quantification, it is necessary to test each sample, control, and all quantification standards in duplicate. Reference material should be tested according to the NF T90-471 and ISO/TS 12869 standards. Do not mix quantification standards from different batches

Section 7

Protocol

It is strongly recommended to read the entire protocol before starting the test.

A. Sampling and Transportation

Water samples should be collected according to the general standards for bacteria detection and enumeration (NF T90-471 and ISO/TS 12869). The quality of results depends on strict compliance with good practices of water sampling for microbiological analysis (for example, the ISO 19458 standard), especially concerning PCR.

B. Water Filtration and DNA Extraction

General recommendations

- Switch on the water bath or the heating thermoshaker and set to $95 \pm 5^{\circ}\text{C}$. Check the level of the water bath for adequate submersion of the 4.5 ml cryotubes
- Prepare the number of tubes corresponding to the number of samples for DNA extraction. Pipet x ml of R1 in each tube, depending on the chosen protocol

Note: Pipet R1 while stirring at medium speed on a stir plate in order to keep the resin in suspension. Use a pipet tip with a wide opening (for example, use a 1000 μl pipet with the corresponding tip).

- Prepare the number of purification columns required by placing a column in each collector vial

Water filtration

1. Rinse the filtration ramp with 100 ml of *Legionella* DNA Free Water and then decontaminate the ramp by burning with alcohol. Make sure the ramp is dry and cool before positioning the membrane filter. This operation should be repeated after each filtered sample in order to avoid DNA or bacterial contamination between sample filtrations.
2. Place the membrane filter on the filtration apparatus. Place the sterile funnels.
3. Filter 100 ml – 1 L of water sample over the membrane filter.

Aquadien DNA extraction and purification

1. Carefully fold the membrane in half three times to obtain a cone (see Appendix 1).
2. Using tweezers, place the membrane in a cryotube containing 2 ml of R1. The point of the cone formed by the membrane should be positioned toward the top of the cryotube.
3. Vortex for 20 sec.

Note: Check that the membrane is still totally immersed in R1 solution.

4. Incubate in a covered water bath at $95 \pm 5^{\circ}\text{C}$ for 15 min.
5. Vortex for 20 sec.
6. Using a 1 ml sterile filtered pipet tip, carefully remove the membrane by pressing it to the wall of the cryotube in order to recover all the lysate. Throw away the membrane.
7. Allow the cryotubes to sit at room temperature for 20 min. The resin of the R1 solution will pellet at the bottom of the cryotube. The supernatant contains 1.6 ml of extracted DNA. It is also possible to centrifuge the cryotubes at $900 \times g$ for 3 min in a centrifuge adapted for 4.5 ml tubes. Verify that the tubes are at room temperature before starting the purification step.
8. Transfer 500 μl of the supernatant to a purification column. Do not vortex the lysate before collection.
9. Seal each column with the collector vial cap.

Note: Do not pipet the resin pellet. If this occurs, wait 5 min for the resin to pellet again or centrifuge at $900 \times g$ for 3 min.

10. Centrifuge at $6,000 \times g$ for 10 min. If the centrifuge temperature is programmable, set to 20°C .

Note: If clogging of the purification column occurs, increase the centrifugation time without increasing the centrifugation speed.

11. Throw away the liquid contained in the collector vial.
12. Repeat step 8 by pipetting an additional 500 μl of the supernatant in the same purification column and seal with the collector cap.
13. Centrifuge at $6,000 \times g$ for 10 min. If the centrifuge temperature is programmable, set to 20°C .

Note: If clogging of the purification column occurs, increase the centrifugation time without increasing the centrifugation speed. Check that all supernatant is filtered through the column. If not, centrifuge again. If clogging still occurs after a second centrifugation, follow the protocol for purification column clogging below.

14. Add 100 μl of R2 solution in the purification column.
15. Throw away the collector vial. Cover the purification column with a new clean collector vial and turn the whole unit upside down.
16. Centrifuge at $1,000 \times g$ for 3 min. If the centrifuge temperature is programmable, set to 20°C . Throw away the purification column.

Note: In this step, the cap cannot be closed.

17. Store the collector vial containing the 100 µl of purified DNA solution. Use 5 µl of the DNA solution for each real-time PCR reaction. The DNA solution can be stored for several months at –20°C.
The calculation of the Z factor ($Z = 32$) is detailed in Section 7 Protocol, D. Data Analysis, Calculation of the Analyzed Fraction of Sample.

Aquadien DNA extraction and purification —short protocol

1. Carefully fold the membrane in half three times to obtain a cone (see Appendix 1).
2. Using tweezers, place the membrane in a tube containing 1 ml of R1. The point of the cone formed by the membrane should be positioned toward the top of the tube.
3. Incubate at 95°C for 15 min at 1,300 rpm in a heating thermoshaker.
4. Carefully take out the membrane by pressing it to the wall of the tube in order to recover all the lysate. Throw away the membrane.
5. Centrifuge at 900 x g for 3 min.
6. Transfer 500 µl of the supernatant to a purification column. Do not vortex the lysate before collection.
7. Seal each column with the collector vial cap.

Note: Do not pipet the resin pellet. If this occurs, wait 5 min for the resin to pellet again or centrifuge at 900 x g for 3 min.

8. Centrifuge at 6,000 x g for 10 min. If the centrifuge temperature is programmable, set to 20°C.

Note: If clogging of the purification column occurs, increase the centrifugation time without increasing the centrifugation speed.

9. Add 100 µl of R2 solution in the purification column.
10. Throw away the collector vial. Cover the purification column with a new clean collector vial and turn the whole unit upside down.
11. Centrifuge at 1,000 x g for 3 min. If the centrifuge temperature is programmable, set to 20°C. Throw away the purification column.

Note: In this step, the cap cannot be closed.

12. Store the collector vial containing the 100 µl of purified DNA solution. Use 5 µl of the DNA solution for each real-time PCR reaction. The DNA solution can be stored for several months at –20°C.
The calculation of the Z factor ($Z = 32$) is detailed in Section 7 Protocol, D. Data Analysis, Calculation of the Analyzed Fraction of Sample.

Aquadien DNA extraction and purification — in case of purification column clogging – W2 protocol

1. Carefully fold the membrane in half three times to obtain a cone (see Appendix 1).
2. Using tweezers, place the membrane in the cryotube containing 2 ml of R1. The point of the cone formed by the membrane should be positioned toward the top of the cryotube.
3. Vortex 20 sec.

Note: Check that the membrane is still totally submerged in R1 solution.

Section 7 Protocol

4. Incubate in a covered water bath at $95 \pm 5^{\circ}\text{C}$ for 15 min.
5. Vortex 20 sec.
6. Using a 1 ml sterile filtered pipet tip, carefully remove the membrane by pressing it to the wall of the cryotube in order to recover all the lysate. Throw away the membrane.
7. Transfer the sample, including the resin, to a new 2 ml tube. Throw away the cryotube. The lysate can be stored at $2-8^{\circ}\text{C}$ for 24–72 hr.
8. Add 200 μl of cool (4°C) W2 buffer. Vortex 5 sec.
9. Allow the tube to sit in a refrigerator at $4 \pm 2^{\circ}\text{C}$ for 15 min.
10. Centrifuge at $12,000 \times g$ at 4°C for 15 min.
11. Transfer 500 μl of the supernatant to a purification column. Do not vortex the lysate before collection. Seal each column with the collector vial cap.

Note: Do not pipet the resin pellet. If this occurs, wait 5 min for the resin to pellet again or centrifuge at $900 \times g$ for 3 min.

12. Centrifuge at $6,000 \times g$ at 20°C for 10 min.

Note: If clogging of the purification column occurs, increase the centrifugation time without increasing the centrifugation speed.

13. Throw away the liquid contained in the collector vial.
14. Repeat step 11 by pipetting an additional 500 μl of the supernatant in the same purification column and seal with the collector cap.
15. Centrifuge at $6,000 \times g$ at 20°C for 10 min.

Note: If clogging of the purification column occurs, increase the centrifugation time without increasing the centrifugation speed. Check that all supernatant is filtered through the column. If not, centrifuge again.

16. Throw away the liquid contained in the collector vial.
17. Add 250 μl of R2 solution in the purification column.
18. Centrifuge at $6,000 \times g$ at 20°C for 5 min.
19. Add 100 μl of R2 solution in the purification column.
20. Throw away the collector vial and cover the purification column with a new clean collector vial and turn the whole unit upside down.
21. Centrifuge at $1,000 \times g$ at 20°C for 3 min. Throw away the purification column.

Note: In this step, the cap cannot be closed.

22. Store the collector vial containing the 100 µl of purified DNA solution. Use 5 µl of the DNA solution for each real-time PCR reaction. The DNA solution can be stored for several months at –20°C. The calculation of the Z factor ($Z = 36$) is detailed in Section 7 Protocol, D. Data Analysis, Calculation of the Analyzed Fraction of Sample.

Aquadien DNA extraction and purification —in case of purification column clogging – W2 short protocol

1. Carefully fold the membrane in half three times to obtain a cone (see Appendix 1).
2. Using tweezers, place the membrane in a tube containing 1 ml of R1. The point of the cone formed by the membrane should be positioned toward the top of the tube.
3. Incubate at 95°C for 15 min at 1,300 rpm in a thermoshaker.
4. Carefully take out the membrane by pressing it to the wall of the tube in order to recover all the lysate. Throw away the membrane.
5. Add 100 µl of cool (4°C) W2 buffer. Vortex 5 sec.
6. Allow the tube to sit in a refrigerator at $4 \pm 2^\circ\text{C}$ for 15 min.
7. Centrifuge at 12,000 x g at 4°C for 15 min.
8. Transfer 500 µl of the supernatant to a purification column. Do not vortex the lysate before collection. Seal each column with the collector vial cap.

Note: Do not pipet the resin pellet. If this occurs, wait 5 min for the resin to pellet again or centrifuge at 900 x g for 3 min.

9. Centrifuge at 6,000 x g at 20°C for 10 min.

Note: If clogging of the purification column occurs, increase the centrifugation time without increasing the centrifugation speed.

10. Add 125 µl of R2 solution in the purification column.
11. Centrifuge at 6,000 x g at 20°C for 5 min.
12. Add 100 µl of R2 solution in the purification column.
13. Throw away the collector vial and cover the purification column with a new clean collector vial and turn the whole unit upside down.
14. Centrifuge at 1,000 x g at 20°C for 3 min. Throw away the purification column.

Note: In this step, the cap cannot be closed.

15. Store the collector vial containing the 100 µl of purified DNA solution. Use 5 µl of the DNA solution for each real-time PCR reaction. The DNA solution can be stored for several months at –20°C. The calculation of the Z factor ($Z = 36$) is detailed in Section 7 Protocol, D. Data Analysis, Calculation of the Analyzed Fraction of Sample.

Aquadien DNA extraction and purification — Free DNA Removal Solution (FDRS) short protocol

1. Carefully fold the membrane in half three times to obtain a cone (see Appendix 1).
2. Using tweezers, place the membrane in a tube containing 460 µl of *Legionella* DNA Free Water and 40µl of activated FDRS. The point of the cone formed by the membrane should be positioned toward the top of the tube.
3. Invert the tube several times to homogenize. Do not vortex.
4. Incubate at 37°C for 30 min.
5. Add 500 µl R1 to inactivate FDRS for DNA extraction.
6. Vortex for 10 sec.
7. Incubate at 95°C for 15 min at 1,300 rpm in a thermoshaker.
8. Carefully take out the membrane by pressing it to the wall of the tube in order to recover all the lysate. Throw away the membrane.
9. Centrifuge at 900 x g for 3 min.
10. Transfer 500 µl of the supernatant to a purification column. Do not vortex the lysate before collection.
11. Centrifuge at 6,000 x g for 10 min.
12. Add 100 µl of R2 solution and throw away the collector vial.
13. Throw away the collector vial. Cover the purification column with a new clean collector vial and turn the whole unit upside down.
14. Centrifuge at 1,000 x g for 3 min. If the centrifuge temperature is programmable, set to 20°C. Throw away the purification column.

Note: In this step, the cap cannot be closed.

15. Store the collector vial containing the 100 µl of purified DNA solution. Use 5 µl of the DNA solution for each real-time PCR reaction. The DNA solution can be stored for several months at –20°C. The calculation of the Z factor ($Z = 36$) is detailed in Section 7 Protocol, D. Data Analysis, Calculation of the Analyzed Fraction of Sample.

C. Real-Time PCR

Instrument and Software Setup

For instrument and software setup, follow instructions in the real-time PCR system user guide for iQ-Check Kits.

PCR Mix Preparation

1. Prepare PCR mix containing the amplification solution (reagent C) and the fluorescent probes (reagent B). Volume of PCR mix needed depends on the number of samples and controls to be analyzed. At least one positive and one negative control must be included in each PCR run. Use the pipetting table in the Appendix 2 to find the correct volumes to use for each reagent.

Note: For quantification, it is necessary to test each sample, control, and all quantification standards in duplicate. Do not mix quantification standards from different batches.

Note: When running both iQ-Check *Legionella* spp. and iQ-Check *L. pneumophila* on the same plate, because the amplification mixes are specific to each method, it is necessary to run two independent series of quantification standards. In addition, the results of each should be analyzed separately.

2. Use the PCR mix (reagents B + C) immediately after preparation. It is stable for 1 hr maximum at 2–8°C.
3. Pipet 45 µl of PCR mix into each well according to your plate setup.
4. Add 5 µl of DNA extract, reagent D (negative control), and reagent E (positive control) or quantification standards. Do not vortex the sample before pipetting. Hermetically seal the wells of the plate or strips. It is important to avoid bubbles at the bottom of the wells by pipetting carefully. As an optional step to eliminate any bubbles, centrifuge the sealed PCR plate or the PCR strips (quick spin).
5. Place the plate or strips in the thermal cycler. Be sure to place the plate with the A1 well at the upper left corner. Close the reaction module.

Run PCR

To start the PCR run, follow instructions in the real-time PCR system user guide for iQ-Check Kits.

1. Define the plate setup in the CFX Manager Industrial Diagnostic Edition (IDE) Software.
2. For iQ-Check *Legionella* Screen Kits, samples and controls are tested singly.
3. For iQ-Check *Legionella* Quanti kits, samples, controls, and quantification standards are tested in duplicate. Indicate the quantity value for each standard and for the reference material. This value is found on the calibration card in each box of the iQ-Check *Legionella* Quantification Standards and on the certificate of analysis for each product batch.

Note: For quantification assays, the calibration curve generated by one batch can be saved and reused until the expiration date of that batch. The option 'Use reduced set of QS if possible' needs to be activated before the PCR run. At the end of the run, select “Tools” then click on “Save Standard Curve” to save the curve. Subsequently, a single point of the calibration range needs to be run, in duplicate, at each analysis to verify conformity.

4. Select the appropriate thermal protocol, “Legio. spp.” or “Legio. pneumo.,” for qualitative (Screen) kits or quantitative (Quanti) kits and click Run.

D. Data Analysis

Data can be analyzed directly at the end of the PCR run or at a later time by opening the stored data file. Follow instructions in the corresponding CFX Manager IDE Software user manual for opening data files and setting the data analysis parameters.

Calculation of the Analyzed Fraction of Sample (Z Factor)

Any analysis method that includes an extraction step followed by a detection step must be accompanied by a calculation of the fraction of the processed sample that was actually analyzed during the final detection. This value is taken into account in the calculation of the detection limit and the quantification limit of the overall method and is used to give a final result in a quantitative test. To do this, calculate the remaining fraction with respect to the initial sample at each step in the protocol (concentration, elimination, etc.). The Z factor corresponds to the denominator of the fraction analyzed and is specific to each extraction protocol. Z represents the value F/V expressed in the result expression table in ISO/TS 12869 and AFNOR T90-471 standards.

Section 7 Protocol

The CFX Manager Software, IDE automatically calculates the Z factor (Z_2 and Z_3). The Z_1 needs to be modified by the user in the software if the filtered volume is different than 1 L.

1. When 1 L of water is sampled and 5 μ l of DNA is analyzed using PCR, the Z value of the Aquadien protocol and Aquadien short protocol is 32.

It is calculated as follows:

- Sample filtration step: 1,000 ml of water is filtered out of the 1,000 ml sampled. The filtered fraction is 1/1 so $Z_1 = 1$
- DNA extraction and purification step:
 - Aquadien protocol = 1 ml of R1 supernatant is processed through the column out of the 1.6 ml generated. The purified fraction is 1 from 1.6 so $Z_2 = 1.6$
 - Aquadien short protocol = 500 μ l of R1 supernatant is processed through the column out of the 800 μ l generated. The purified fraction is 500 μ l from 800 μ l so $Z_2 = 1.6$
- PCR analysis step: 5 μ l is analyzed by PCR out of 100 μ l extracted DNA. The analyzed fraction is $5/100 = 1/20$ so $Z_3 = 20$

The overall Z factor for the Aquadien protocols is $Z = Z_1 \times Z_2 \times Z_3 = 1 \times 1.6 \times 20 = 32$.

The raw PCR result should be multiplied by 32 to obtain the final quantity of bacteria contained in the initial water sample, expressed in genomic units (GU) per liter of water sample. If the filtered water volume is different than 1 L, it must be captured in these calculations.

2. When 1 L of water is sampled and 5 μ l of DNA is analyzed using PCR, the Z value of the Aquadien W2 protocol and Aquadien W2 short protocol in case of purification column clogging is 36.

It is calculated as follows:

- Sample filtration step: 1,000 ml of the water is filtered out of the 1,000 ml sampled. The filtered fraction is 1/1 so $Z_1 = 1$
- DNA extraction and purification step:
 - Aquadien protocol = 1 ml of R1 supernatant is processed through the column out of the 1.8 ml generated (1.6 ml of R1 + 0.2 ml of W2). The purified fraction is 1 from 1.8 so $Z_2 = 1.8$
 - Aquadien short protocol = 500 μ l of R1 supernatant is processed through the column out of the 900 μ l generated (800 μ l of R1 + 100 μ l of W2). The purified fraction is 500 μ l from 900 μ l so $Z_2 = 1.8$
- PCR analysis step: 5 μ l is analyzed by PCR out of 100 μ l extracted DNA. The analyzed fraction is $5/100 = 1/20$ so $Z_3 = 20$.

The overall Z factor from the Aquadien W2 protocols in case of purification column clogging is:

$$Z = Z_1 \times Z_2 \times Z_3 = 1 \times 1.8 \times 20 = 36.$$

The raw PCR result should be multiplied by 36 to obtain the final quantity of bacteria contained in the initial water sample, expressed in GU per liter of water sample. If the filtered water volume is different than 1 L, it must be captured in these calculations.

3. When 1 L of water is sampled and 5 μ l of DNA is analyzed using PCR, the Z value of the Aquadien FDRS short protocol is 36.

It is calculated as follows:

- Sample filtration step: 1,000 ml of the water is filtered out of the 1,000 ml sampled. The filtered fraction is 1/1 so $Z_1 = 1$
- DNA extraction and purification step: 500 µl supernatant is processed through the column out of the 900 µl generated (460 µl of water and 40 µl activated FDRS with 400 µl of R1). The purified fraction is 900/500 so $Z_2 = 1.8$.
- PCR analysis step: 5 µl is analyzed by PCR out of 100 µl extracted DNA. The analyzed fraction is $5/100 = 1/20$ so $Z_3 = 20$.

The overall Z factor from the Aquadien FDRS short protocol is: $Z = Z_1 \times Z_2 \times Z_3 = 1 \times 1.8 \times 20 = 36$.

The raw PCR result should be multiplied by 36 to obtain the final quantity of bacteria contained in the initial water sample, expressed in GU per liter of water sample. If the filtered water volume is different than 1 L, it must be captured in these calculations.

Interpreting Results – iQ-Check Screen *Legionella*

Once the data analysis parameters have been set, results are interpreted by analyzing the Cq values of each sample (the cycle at which the amplification curve crosses the threshold).

CFX Manager IDE Software allows complete automated analysis for Bio-Rad real-time PCR detection systems.

Controls

Verify the positive and negative controls before interpreting sample results.

For the experiment to be valid, the controls must have the following results, as summarized in the table below. Otherwise, the PCR reaction must be repeated.

	<i>Legionella</i> Detection (FAM Channel)	Internal Control Detection (HEX Channel)
Negative control	Cq = N/A*	$28 \leq Cq \leq 40$
Positive control	$30 \leq Cq \leq 40$	N/A

* The software indicates a Cq value (the cycle at which the amplification curve crosses the threshold) of N/A (not applicable) when the fluorescence of a sample does not rise significantly above the background noise and hence does not cross the threshold.

If results of negative and positive controls differ from those in the table above (invalid control), repeat the run and analysis starting from Section 7 Protocol, C. Real-Time PCR, Instrument and Software Setup.

Samples

Interpretation of sample results is summarized in the following table:

<i>Legionella</i> Detection (FAM Channel)	Internal Control Detection (HEX Channel)	Interpretation
$Cq \geq 10$	N/A	Positive
$Cq > \text{Intercept value}^1$	$28 \leq Cq \leq 40$	Negative

Section 7 Protocol

Cq = N/A	Cq = N/A	Inhibition ²
Cq = N/A	Cq > 40	Inhibition
Cq = N/A	Cq < 28	Analytical issue ³

¹ The intercept value to be considered in the context of a qualitative analysis is Cq = 43 (this value is used to give an indication of the intercept value chosen during the validation). An intercept value can be determined by the user by averaging intercept values obtained and adding two standard deviations.

² When both sample and internal control detection give a Cq value = N/A, the sample must be diluted 1:10 in DNA extraction elution buffer and tested again.

³ Verify that the threshold was correctly placed, or that the curve as raw data is a regular amplification curve. If the curve does not have a characteristic shape, repeat the PCR test.

An invalid interpretation can be given when validation criteria are not met. Check the raw data and proceed as if the sample was inhibited.

Detection limit

The detection limit of the PCR step (Ld PCR) for a *Legionella* PCR detection and quantification method is defined as the lowest number of GU generating a positive result at a 90% confidence limit (NF T90-471). The detection limit of the PCR step (Ld PCR) is 5 GU per 5 µl of extracted DNA. To calculate the theoretical detection limit of the total method (Ld) for a 1 L sample, take into consideration the Z factor. The theoretical Ld for the total method takes into consideration a filtration and extraction yield of 100%.

$$Ld = \frac{5 \times Z \times D}{\text{Volume of filtered sample (L)}}$$

Z = Fraction of the sample deposited in each PCR well

D = Dilution factor if the DNA has been diluted before the PCR run

Interpreting Results — iQ-Check Quanti *Legionella*

Verify that the raw data curve has typical amplification (flat base line followed by an exponential increase of fluorescence and then a flattening out). Check the validity of the quantification standards series and of the controls prior to reading the mean quantity of GU calculated for each sample.

Controls and standards

Verify the controls and standards before interpreting sample results.

For the experiment to be valid, the controls and standards must have the following results, as summarized in the table below. Otherwise, the PCR reaction must be repeated.

	<i>Legionella</i> Detection (FAM Channel)	Internal Control Detection (HEX Channel)
Negative control	Cq = N/A*	$28 \leq Cq \leq 40$
PCR efficiency	$75\% \leq \text{PCR efficiency} \leq 125\%$	N/A
Correlation coefficient (r^2)	≥ 0.99	N/A
Qs1–Qs4	$20 \leq Cq \leq 40$	$28 \leq Cq \leq 40$

* The software indicates a Cq value (the cycle at which the amplification curve crosses the threshold) of N/A (not applicable) when the fluorescence of a sample does not rise significantly above the background noise and hence does not cross the threshold.

If one of the quantification standards is significantly outside of the standard curve, it can be eliminated to optimize the results and reach the above parameters. Only one standard replicate can be eliminated. The reference material must follow NF T90-471 and ISO/TS 12869 standards requirements. If the results of the negative controls and standards differ from those given above, repeat the PCR test.

Samples

Interpretation of sample results is summarized in the following table:

<i>Legionella</i> Detection (FAM Channel)	Internal Control Detection (HEX Channel)	Interpretation
Cq ≥ 10	N/A	Positive
Cq ≥ 10	$Cq \leq \text{Mean Cq Qs} + 3\sigma^1$	Positive
Cq ≥ 10	$Cq > \text{Mean Cq Qs} + 3\sigma$	Inhibition ²
Cq = N/A	$\text{Mean Cq Qs} - 3\sigma \leq Cq \leq \text{Mean Cq Qs} + 3\sigma$	Negative
Cq = N/A	N/A	Inhibition
Cq = N/A	$Cq > \text{Mean Cq Qs} + 3\sigma$	Inhibition
Cq = N/A	$Cq < \text{Mean Cq Qs} - 3\sigma$	Analytical issue ³

¹Mean Cq Qs is the average of the HEX Cq values for all Qs; σ represents standard deviation.

² When both sample and internal control detection give a Cq value = N/A, the sample must be diluted 1:10 in DNA extraction elution buffer and tested again.

³ Verify that the threshold was correctly placed, or that the curve as raw data is a regular amplification curve. If the curve does not have a characteristic shape, repeat the PCR test.

If a sample has a very high concentration of *Legionella* ($Cq_{\text{FAM}} < Cq_{\text{FAM}} \text{ Qs4}$), the Cq_{HEX} might be higher than the Cq_{HEX} of Qs4. This sample is outside the dynamic range of quantification. If precise quantification is required, dilute the sample and repeat PCR.

A well with a Cq value > b intercept value and with no inhibition is not considered positive.

If one positive and one negative PCR result is obtained for a sample, the sample is considered positive.

Calculation of *Legionella* concentration

The values given for each sample in the report correspond to the initial quantity of *Legionella* GU present in 5 µl of DNA extract. To obtain the concentration of *Legionella* in GU/1 L of water sample, take into account the quantification mean value ("mean GU/well") calculated by the software for each sample.

$$X = \frac{[\text{Mean quantity in 5 } \mu\text{l}] \times Z \times D}{\text{Volume of filtered sample (L)}}$$

X = *Legionella* GU contained in 1 L of water sample

Z = Fraction of the sample deposited in each PCR well

D = Dilution factor if the DNA has been diluted before the PCR run

Detection limit

The detection limit of the PCR step (Ld PCR) for a *Legionella* PCR detection and quantification method is defined as the lowest number of GU generating a positive result at a 90% confidence limit (NF T90-471). The detection limit of the PCR step (Ld PCR) is 5 GU per 5 µl of extracted DNA. To calculate the theoretical detection limit of the total method (Ld) for a 1 L sample, take into consideration the Z factor divided by 2, since samples are tested in duplicate. The theoretical Ld for the total method takes into consideration a filtration and extraction yield of 100%.

$$Ld = \frac{5 \times Z \times D}{\text{Volume of filtered sample (L) } \times 2}$$

Z = Fraction of the sample deposited in each PCR well

D = Dilution factor if the DNA has been diluted before the PCR run

Quantification limit

The quantification limit of the PCR step (Lq PCR) of a *Legionella* PCR quantification method is defined as the lowest number of copies allowing quantification repeatability and accuracy as described in the NF T90-471 standard. Qs1, the first point on the quantification standard, corresponds to the quantification limit for each batch of iQ-Check Quanti *Legionella*. The exact value of Qs1 for each batch of quantification standards is indicated on the calibration card and on the certificate of analysis.

$$Lq = Qs1 \times Z$$

If a volume of water different from 1 L is filtered or if the DNA extract has to be diluted, the Lq for the method changes.

$$Lq = \frac{Qs1 \times Z \times D}{\text{Volume of filtered sample (L)}}$$

Upper quantification limit

The upper quantification limit (UQL) of the method corresponds to the value (for 1 L) given by the highest point of the quantification standards, the Qs4 standard. The exact Qs4 value for each batch of quantification standards is indicated on the calibration card and the certificate of analysis. To obtain the UQL of the method, the Qs4 value must be multiplied by the fraction of analyzed sample Z and by the dilution factor D.

$$UQL = \frac{Qs4 \times Z \times D}{\text{Volume of filtered sample (L)}}$$

Expressing Results

The final result is expressed in GU/L with 2 significant figures (for example, 1,300 GU/L for a calculated concentration of 1,256 GU/L).

Genomic Units/PCR Well	Interpretation	Conclusion
<5	<Ld	Target not detected at the detection limit
[5, Qs1]	<Lq	Target present but not quantifiable
[Qs1–Qs4]	X	Target present and quantifiable
> Qs4	UQL	High amount of target*

* Outside the dynamic range of quantification. If precise quantification is required, dilute the sample and repeat PCR.

Section 8

Confirmation of Positive Results

Refer to local regulations.

Section 9

Confirmation of Single Colonies Using iQ-Check Kit

iQ-Check *Legionella* kits may also be used for confirming single isolated *Legionella* colonies on GVPC or BCYE agar plates.

1. Pick an isolated colony with a toothpick, sterile loop, or other adapted consumable (for example, a pipet tip).
2. Resuspend the colony in 200 µl *Legionella* DNA Free Water in a microcentrifuge tube. Vortex for 20 sec.
3. Divide the suspension into two 100 µl aliquots (tube 1 and tube 2).
4. Incubate tube 1 at 95°C for 10 min. Use 5 µl of the suspension with 45 µl of PCR mix (see Section 7 Protocol, C. Real-Time PCR) to run PCR.
5. If $C_q > \text{intercept value}$, the result is negative and the colony cannot be confirmed or identified.
6. If $C_q < 30$, the result is positive and the colony is confirmed and identified. In this case, tube 1 with extracted DNA should be stored at –20°C for up to 3 months, if PCR must be repeated. In addition, tube 2 must be stored at 4°C for possible culturing.
7. If PCR is inhibited, dilute the tube 1 DNA extract in sterile *Legionella* DNA Free Water (1:100) and repeat PCR.
8. If $30 < C_q < \text{intercept value}$, confirm the positive result with another method or repeat the protocol with another isolated colony.

Section 10

Test Performance and Validations

BRD 07/15 – 12/07
BRD 07/16 – 12/07

METHOD FOR WATER
ANALYSIS

<http://nf-validation.afnor.org>

NF Validation

Aquadien DNA Extraction and Purification Kit is used in the scope of the iQ-Check *Legionella* spp. and *L. pneumophila* methods certified NF Validation according to the validation protocol based on the NF T90-471 Standard (June 2015) for the detection and quantification of *Legionella* and/or *Legionella pneumophila* with iQ-Check Real-Time PCR Kits. The validation scope has been extended to the ISO/TS 12869 (April 2019): Water quality – Detection and quantification of *Legionella* spp. and/or *L. pneumophila* by concentration and genic amplification by quantitative polymerase chain reaction (qPCR). Certificate number: iQ-Check *Legionella* spp. BRD 07/15 – 12/07 and iQ-Check *L. pneumophila* BRD 07/16 – 12/07. Valid until: Refer to the certificate available on the AFNOR Certification website.

Section 11

References

AFNOR NF T90-471. Detection and quantification of *Legionella* and/or *Legionella pneumophila* by concentration and PCR amplification.

Engleberg NC et al. (1989). DNA sequence of mip, a *Legionella pneumophila* gene associated with macrophage infectivity. Infect Immun 57,1263–1270.

ISO 11731. Water quality — Enumeration of *Legionella*.

ISO/TS 12869. Detection and quantification of *Legionella* spp. and/or *Legionella pneumophila* by concentration and genic amplification by quantitative polymerase chain reaction (qPCR)

MacDonell MT and Colwell RR. (1987). The nucleotide sequence of the 5S rRNA from *Legionella pneumophila*. Nucleic Acids Res 15, 1335.

NF T 90-431. Water quality — Detection and enumeration of *Legionella* spp. and *Legionella pneumophila* — Method by direct inoculation and after concentration by membrane filtration or centrifugation.


Specht T et al. (1986). Compilation of 5S rRNA and 5S rRNA gene sequences. Nucleic Acids Res 18 Suppl, 2215–2230.

Section 12


Revision History

Release date	Document number	Change
October 2020	10000134675 Ver A	- New document design and update of references and content - Document number change (previous versions were 881116 and 881117)
June 2023	10000134675 Ver B	- Renewal and extension of NF Validation: - CFX Opus 96 Real-Time PCR System - Calibration curve generated by one batch can be saved and reused until the end of that batch - Update of references and content
October 2023	10000134675 Ver C	- NF Validation extension on CFX Opus Deepwell Real-Time PCR System - Update Aquadien kit components - General content updates

Appendix 1 — Membrane Folding Guide

Carefully fold the membrane in half three times to obtain a cone.

Using tweezers, place the membrane in a tube containing 1 ml of R1. The point of the cone formed by the membrane should be positioned toward the top of the tube.

Appendix 2 — PCR Mix Calculation Guide

To find the correct volumes to use when preparing the PCR mix, add the total number of samples and controls to be analyzed and find the corresponding volumes of reagent B and reagent C in the table.

Total Number of Samples and Controls	Probes Reagent B, μ l	Amplification Mix Reagent C, μ l	Total Number of Samples and Controls	Probes Reagent B, μ l	Amplification Mix Reagent C, μ l	Total Number of Samples and Controls	Probes Reagent B, μ l	Amplification Mix Reagent C, μ l
1	5	40	33	178	1,400	65	351	2,800
2	11	86	34	184	1,500	66	356	2,900
3	16	130	35	189	1,500	67	362	2,900
4	22	173	36	194	1,600	68	367	2,900
5	27	216	37	200	1,600	69	373	3,000
6	32	259	38	205	1,600	70	378	3,000
7	38	302	39	211	1,700	71	383	3,100
8	43	346	40	216	1,700	72	389	3,100
9	49	389	41	221	1,800	73	394	3,200
10	54	432	42	227	1,800	74	400	3,200
11	59	475	43	232	1,900	75	405	3,200
12	65	518	44	238	1,900	76	410	3,300
13	70	562	45	243	1,900	77	416	3,300
14	76	605	46	248	2,000	78	421	3,400
15	81	648	47	254	2,000	79	427	3,400
16	86	691	48	259	2,100	80	432	3,500
17	92	734	49	265	2,100	81	437	3,500
18	97	778	50	270	2,200	82	443	3,500
19	103	821	51	275	2,200	83	448	3,600
20	108	864	52	281	2,200	84	454	3,600
21	113	907	53	286	2,300	85	459	3,700
22	119	950	54	292	2,300	86	464	3,700
23	124	994	55	297	2,400	87	470	3,800
24	130	1,000	56	302	2,400	88	475	3,800
25	135	1,100	57	308	2,500	89	481	3,800
26	140	1,100	58	313	2,500	90	486	3,900
27	146	1,200	59	319	2,500	91	491	3,900
28	151	1,200	60	324	2,600	92	497	4,000
29	157	1,300	61	329	2,600	93	502	4,000
30	162	1,300	62	335	2,700	94	508	4,100
31	167	1,300	63	340	2,700	95	513	4,100
32	173	1,400	64	346	2,800	96	518	4,100

Visit bio-rad.com/legionella for more information.

BIO-RAD is a trademark of Bio-Rad Laboratories, Inc.

AQUADIEN and IQ-CHECK are trademarks of Bio-Rad Europe GMBH in certain jurisdictions.

All trademarks used herein are the property of their respective owner.

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BIO-RAD es una marca registrada de Bio-Rad Laboratories, Inc.
AQUADIEN e IQ-CHECK son marcas registradas de Bio-Rad Europe GmbH en diversos países.