

SE

Bio-Prep SE 100/17 and 1000/17 Size Exclusion Chromatography Columns

- Prepacked Biocompatible Cross-Linked Agarose Column
- High Mechanical and Chemical Stability
- Fractionation Ranges 5,000–100,000 Da, 10,000–1,000,000 Da
- High Resolution Separations

Summary

Size exclusion chromatography remains one of the most widely used and accepted chromatographic techniques available. This technique is based on the separation of bio-molecules according to their molecular weight, and has been extensively used for analytical, preparative, and process scale separations.

Product Description

The Bio-Prep SE 100/17 and SE 1000/17 glass columns are based on biocompatible, 17 μ m spherical, crosslinked agarose beads. The beads are produced from highly purified agarose utilizing a proprietary crosslinking process which allows the production of a mechanically and chemically stable material.

Two fractionation ranges are available, 5,000–100,000 Da using the Bio-Prep SE 100/17 column, and 10,000–1,000,000 Da using the Bio-Prep SE 1000/17 column.

Sample volumes of 0.6–2% of the bed volume and flow-rates from 0.1 to 2 ml/min are recommended for normal operation. Larger sample volumes of up to 5% of the bed volume are possible where lower resolution is acceptable *e.g.* desalting.

Typically, lower flow rates and smaller sample volumes yield superior separations. However, the robust nature of these materials makes superior separations at elevated flow rates possible if process throughput is a concern, or during operations such as sanitization, equilibration, etc.

Figure 1. Kd Bio-Prep SE 100/17 gel.




Figure 2. Separation of a gel filtration standard using the SE 100/17 column.




Figure 3. Kd Bio-Prep SE 1000/17 column.




Figure 4. Separation of a gel filtration standard using the Bio-Prep SE 1000/17 column.





Figure 5. Pressure/Flow of Bio-Prep SE100/17 and SE1000/17, 8 x 300 mm Columns.



Ordering Information

	Bio-Prep SE100/17	Bio-Prep SE1000/17
Catalog Number	732-1501	732-1502
Column Dimension	8 x 300 mm	8 x 300 mm
Column Housing	Glass	Glass
Bed Matrix	Cross-linked agarose	Cross-linked agarose
Gel Volume	15.0 ml	15.0 ml
Particle Size	15 -20 µm	15 -20 µm
pH Range	1 to 14	1 to 14
Maximum Pressure	580 psi, 40 bar	430 psi, 30 bar
Recommended Flow-rate	0.1 -1.0 ml/min	0.1 -1.0 ml/min
Maximum Flow-rate	2.0 ml/min	2.0 ml/min
Linear Separation Range	10,000 - 100,000 Da	10,000 - 1,000,000 Da
Exclusion Limit	~ 200,000 Da	~ 1,500,000 Da
Cleaning Agent	0.5 m NaOH	0.5 m NaOH

Applications

Applications range from the separation of protein aggregates (dimers, trimers, etc.) from purified antibody solutions to the separation of complex mixtures of bio-molecules directly from cell culture feedstreams. Size exclusion chromatography is commonly used as a final polishing step in the preparation of many bio-pharmaceuticals and has numerous applications as a QC tool for assessing product purity, concentration, etc.

BIO-RAD

**Bio-Rad
Laboratories**

Life Science
Group

Website www.bio-rad.com **Bio-Rad Laboratories Main Office** 2000 Alfred Nobel Drive, Hercules, CA 94547, Ph. (510) 741-1000, Fx. (510) 741-5800
Also in: **Australia** Ph. 02-9914-2800, Fx. 02-9914-2889 **Austria** Ph. (1)-877 89 01, Fx. (1) 876 56 29 **Belgium** Ph. 09-385 55 11, Fx. 09-385 65 54
Canada Ph. (905) 712-2771, Fx. (905) 712-2990 **China** Ph. (86-10) 2046622, Fx. (86-10) 2051876 **Denmark** Ph. 39 17 9947, Fx. 39 27 1698
Finland Ph. 90 804 2200, Fx. 90 804 1100 **France** Ph. (1) 43 90 46 90, Fx. (1) 46 71 24 67 **Germany** Ph. 089 31884-0, Fx. 089 31884-100
Hong Kong Ph. 7893300, Fx. 7891257 **India** Ph. 91-11-461-0103, Fx. 91-11-461-0765 **Israel** Ph. 03 951 4127, Fx. 03 951 4129
Italy Ph. 02-21609.1, Fx. 02-21609.399 **Japan** Ph. 03-5811-6270, Fx. 03-5811-6272 **The Netherlands** Ph. 0313 18-540666, Fx. 0313 18-542216
New Zealand Ph. 09-443 3099, Fx. 09-443 3097 **Singapore** Ph. (65) 272-9877, Fx. (65) 273-4838 **Spain** Ph. (91) 661 70 85, Fx. (91) 661 96 98
Sweden Ph. 46 (0) 8 627 50 00, Fx. 46 (0) 8 627 54 00 **Switzerland** Ph. 01-809 55 55, Fx. 01-809 55 00
United Kingdom Ph. 0800 181134, Fx. 01442 259118