

* We recommend adapting this protocol to use the Gene Pulser electroporation buffer (catalog #165-2676, 165-2677), which increases cell viability and transfection efficiency in mammalian cell lines.

Cell Type	Mammalian, suspension	Molecules Electroporated	DNA: pRSVneo, linearized, 5.6 kB
Species Used	Human, K562, chronic myeloid leukemia		

Before the Pulse

Cell Growth Medium	RPMI + 20% Fetal Calf Serum (GIBCO/BRL,Sigma)	Growth Phase at Harvest	Log
Wash Solution	RPMI	Pre-pulse Incubation	Not given

The Pulse

Instruments Used Gene Pulser® apparatus & Capacitance

Electroporation Temperature	Room temperature	Cuvette Gap	0.4cm
Electroporation Medium*	RPMI	Voltage	up to 2 kV
Cell Density	5 x10 (6) cells / ml	Field Strength	up to 5 kV/cm
Volume of Cells	0.4 to 0.8 ml	Capacitor	up to 960 µF
DNA Concentration	5 µg / 800 µl	Resistor	(Pulse Controller) Ω none
DNA Resuspension Buffer	Not given	Time Constant	0.4 msec
Volume of DNA	up to 50 µl		
After the Pulse			
Outgrowth Medium	RPMI + 20% Fetal Calf Serum		

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.

Outgrowth Temperature	37 °C
Length of Incubation	48 hr.
Selection Method or Assay Used	G418
Electroporation Efficiency	0.3% clonogenic cells / µg
Per Cent Survival	Not given

Name of Submitter Ms. Genevieve M. Croaker

Institution Address Royal Prince Alfred Hospital
Kanematsu Laboratories
Missenden Road
Camperdown, NSW 2050
Australia

Survey Number
103