

Gene Pulser® Electroprotocols

* We recommend adapting this protocol to use the Gene Pulser electroporation buffer (catalog #165-2676, 165-2677), which increases cell viability and transfection efficiency in mammalian cell lines.

Cell Type Mammalian

Molecules Electroporated DNA: pZ189

Species Used Human, K562, chronic myeloid leukemia

Before the Pulse

Cell Growth Medium RPMI 1640 (GIBCO/BRL,Sigma)

Growth Phase at Harvest 10 (6) cells / ml

Pre-pulse Incubation Not given

Wash Solution Phosphate Buffered Saline

The Pulse

Instruments Used Not given

Electroporation Temperature 20 to 25 °C (but pre-cooled on ice)

Electroporation Medium* Phosphate Buffered Saline

Cuvette Gap 0.4 cm

Cell Density 1 x 10 (7) cells / ml

Voltage 0.6 kV

Volume of Cells 0.4 ml

Field Strength 1.50 V/cm

DNA Concentration 400 µg / ml

DNA Resuspension Buffer TE

Capacitor 25 µF

Volume of DNA 40 µg

Resistor (Pulse Controller) Ω none

After the Pulse

Time Constant 1.5 to 1.9 msec

Outgrowth Medium RPMI 1640

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.

PBS: 1x = 8g NaCl, 0.2g KCl, 0.2g KH₂PO₄, 1.15g Na₂HPO₄

Outgrowth Temperature 37 °C, 5% CO₂

Length of Incubation 1 day

Selection Method or Assay Used Not given

Electroporation Efficiency Not given

Per Cent Survival 50%

Name of Submitter Fumio Yatagai, Scientist

Institution Address RIKEN (The Inst. of Phy. & Chem. Res.)
Radiation Biology
2-1 Hirosawa
Wako-shi, Saitama 351-01
JAPAN

Survey Number

104