

Gene Pulser® Electroprotocols

* We recommend adapting this protocol to use the Gene Pulser electroporation buffer (catalog #165-2676, 165-2677), which increases cell viability and transfection efficiency in mammalian cell lines.

Cell Type	Mammalian, suspension	Molecules Electroporated	DNA: linearized (5kB) n-rasCAT plasmids
Species Used	Human, K562, chronic myeloid leukemia.		

Before the Pulse

Cell Growth Medium	RPMI + 20% Fetal Calf Serum (GIBCO/BRL, Sigma)	Growth Phase at Harvest	log phase
		Pre-pulse Incubation	None

Wash Solution RPMI

The Pulse

Instruments Used Gene Pulser® apparatus & Capacitance

Electroporation Temperature	Room temperature		
Electroporation Medium*	RPMI	Cuvette Gap	0.4 cm
Cell Density	10 (7) cells / cuvette	Voltage	2.0 kV
Volume of Cells	850 µl	Field Strength	5 kV/cm
DNA Concentration	50 µg DNA / cuvette	Capacitor	25 µF
DNA Resuspension Buffer	Not given	Resistor	(Pulse Controller) Ω none. NOT
Volume of DNA	50 to 300 µl	Time Constant	0.4 msec
After the Pulse			
Outgrowth Medium	RPMI + 20% Fetal Calf Serum		

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.

It is NOT RECOMMENDED to use high voltage without the Pulse Controller.

Outgrowth Temperature	37 °C
Length of Incubation	2 days
Selection Method or Assay Used	G418
Electroporation Efficiency	Not given
Per Cent Survival	Not given

Name of Submitter Ms. Jacqueline Thorn

Institution Address Royal Prince Alfred Hospital
Missenden Road,
Camperdown, NSW 2050
AUSTRALIA

Survey Number
105