

Gene Pulser® Electroprotocols

* We recommend adapting this protocol to use the Gene Pulser electroporation buffer (catalog #165-2676, 165-2677), which increases cell viability and transfection efficiency in mammalian cell lines.

Cell Type	Mammalian, suspension	Molecules Electroporated	Electrofusion
Species Used	B-cell line, unspecified, and heteromyeloma		

Before the Pulse

Cell Growth Medium	RPMI+ 10% Fetal Calf Serum+ 2 mM L-Glutamine (GIBCO/ BRL, Sigma)	Growth Phase at Harvest	Exponential growth
		Pre-pulse Incubation	No

Wash Solution RPMI

The Pulse

Instruments Used Gene Pulser® apparatus & Capacitance

Electroporation Temperature	Room temperature		
Electroporation Medium*	RPMI	Cuvette Gap	0.2 cm
Cell Density	10 (8) cells / μ l	Voltage	0.160 kV
Volume of Cells	200 μ l	Field Strength	0.800 kV/cm
DNA Concentration	Not given	Capacitor	960 μ F
DNA Resuspension Buffer	Not given	Resistor	(Pulse Controller) none Ω
Volume of DNA	10 (8) cells / ml	Time Constant	10 msec
After the Pulse			
Outgrowth Medium	HAT medium and azaserine & ouabain		

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.

Outgrowth Temperature	37 °C
Length of Incubation	Not given
Selection Method or Assay Used	Not given
Electroporation Efficiency	Not given
Per Cent Survival	No fusion

Name of Submitter Berta Sanchez, Ph.D.

Institution Address Hospital Virgen Del Rocio
Immunologia
Avda Manuez Siurot s/n
41013 Sevilla
Spain

Survey Number

121