

Gene Pulser® Electroprotocols

* We recommend adapting this protocol to use the Gene Pulser electroporation buffer (catalog #165-2676, 165-2677), which increases cell viability and transfection efficiency in mammalian cell lines.

Cell Type	Mammalian, adherent	Molecules Electroporated	DNA: CMV b-gal , 6 kB, supercoiled
Species Used	Monkey, COS-7, kidney, SV-40 transformed		

Before the Pulse

Cell Growth Medium	DMEM + 10% Fetal Bovine Serum (GIBCO/BRL, Sigma)	Growth Phase at Harvest	Log phase, 70 to 80% confluent
		Pre-pulse Incubation	Room temperature
Wash Solution	Trypsinize		

The Pulse

Instruments Used Gene Pulser® apparatus & Capacitance

Electroporation Temperature	Room temperature		
Electroporation Medium*	DMEM + 10% Fetal Bovine Serum + 5mM BES, pH 7.2	Cuvette Gap	0.4 cm
Cell Density	2 x 10 ⁷ (7) cells / ml	Voltage	0.170 kV
Volume of Cells	10 µg / 250 µl DNA	Field Strength	0.425 kV/cm
DNA Concentration	10 mM Tris, 1 mM EDTA, pH 8.0)	Capacitor	960 µF
DNA Resuspension Buffer	Not given	Resistor	(Pulse Controller) Ω none
Volume of DNA	Not given	Time Constant	45 msec
After the Pulse			
Outgrowth Medium	Not given		

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.

Outgrowth Temperature	Room temperature	
Length of Incubation	10 min., then pellet& resuspend in DMEM	See reference: Ustav, M., and Stenlund, A. 1991. <i>EMBO J.</i> 10 (2):449-457.
Selection Method or Assay Used	Not given	
Electroporation Efficiency	25%	
Per Cent Survival	about 100%	

Name of Submitter Cheng-Ming Chiang

Institution Address University of Rochester Medical Center
Biochemistry, Box 607
601 Elmwood Ave
Rochester, NY 14642

Survey Number

127