

Gene Pulser® Electroprotocols

* We recommend adapting this protocol to use the Gene Pulser electroporation buffer (catalog #165-2676, 165-2677), which increases cell viability and transfection efficiency in mammalian cell lines.

Cell Type	Mammalian, suspension	Molecules Electroporated	DNA: pSV2 neo, 5 kB
Species Used	Mouse, erythroleukemia cells		

Before the Pulse

Cell Growth Medium	RPMI ,10% Fetal Calf Serum (GIBCO/BRL, Sigma)	Growth Phase at Harvest	Log phase
		Pre-pulse Incubation	10 to 20 min.
Wash Solution	Not given		

The Pulse

Instruments Used Gene Pulser® apparatus

Electroporation Temperature	4 °C	Cuvette Gap	0.4 cm
Electroporation Medium*	Not given	Voltage	0.4 kV
Cell Density	10 (7) cells / 50 µl	Field Strength	1 kV/cm
Volume of Cells	50 µl	Capacitor	25 µF
DNA Concentration	Not given	Resistor	(Pulse Controller) none Ω
DNA Resuspension Buffer	Not given	Time Constant	10 msec
Volume of DNA	Not given		
After the Pulse			
Outgrowth Medium	RPMI,10% Fetal Calf Serum		

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.

**Note: Transfection efficiemnnces ofthese cells may be increased with the use of the Capacitance Extender (providing longer time constants).

Outgrowth Temperature	37 °C
Length of Incubation	10 to 30 minutes
Selection Method or Assay Used	
Electroporation Efficiency	Low (**see notes)
Per Cent Survival	20 to 70%

Name of Submitter Iliana Coccia Rese

Institution Address Istituto Superiore Sanita
Dept of Virology
Visle Regina Elena 299
Rome , ITALY 00164

Survey Number
131