

* We recommend adapting this protocol to use the Gene Pulser electroporation buffer (catalog #165-2676, 165-2677), which increases cell viability and transfection efficiency in mammalian cell lines.

Cell Type	Mammalian, suspension	Molecules Electroporated	DNA: linear and circular plasmids, 4-12 kB in size
Species Used	Mouse, 32d, J558L, myeloma		

Before the Pulse

Cell Growth Medium	RPMI, 10% Fetal Calf Serum GIBCO/BRL, Sigma)	Growth Phase at Harvest	log
		Pre-pulse Incubation	10 minutes ice
Wash Solution	log		

The Pulse

Instruments Used Gene Pulser® apparatus & Capacitance

Electroporation Temperature	(ice) 0 °C		
Electroporation Medium*	RPMI, 10% Nu-serum	Cuvette Gap	0.4 cm
Cell Density	10 (7) cells / µl	Voltage	0.20 to 0.40 kV
Volume of Cells	800 µl	Field Strength	0.5 to 1.0 kV/cm
DNA Concentration	1 mg / ml		
DNA Resuspension Buffer	TE (10 mM Tris, 1 mM EDTA, pH 8.0)	Capacitor	960 µF
Volume of DNA	20 to 40 µl	Resistor	(Pulse Controller) none Ω
		Time Constant	17 msec

After the Pulse

Outgrowth Medium Not given

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.

Outgrowth Temperature	Not given
Length of Incubation	Not given
Selection Method or Assay Used	G418; Hygromycin B
Electroporation Efficiency	Not given
Per Cent Survival	Not given

Name of Submitter Alan Matsumoto

Institution Address Johns Hopkins University
725 N. Wolfe St.
Baltimore, MD 21205

Survey Number
134