

Gene Pulser® Electroprotocols

* We recommend adapting this protocol to use the Gene Pulser electroporation buffer (catalog #165-2676, 165-2677), which increases cell viability and transfection efficiency in mammalian cell lines.

Cell Type	Mammalian, suspension	Molecules Electroporated	DNA: immunoglobulin genes in pSV vector
Species Used	Mouse, NSO, myeloma cells		

Before the Pulse

Cell Growth Medium	RPMI + 10% Fetal Calf Serum (GIBCO/BRL, Sigma)	Growth Phase at Harvest	Not given
		Pre-pulse Incubation	10 min. on ice
Wash Solution	Not given		

The Pulse

Instruments Used Not given

Electroporation Temperature	0 °C	Cuvette Gap	0.2 cm
Electroporation Medium*	Not given	Voltage	0.250 kV
Cell Density	12 x 10 (6) / ml	Field Strength	1.25 kV/cm
Volume of Cells	400 µl	Capacitor	960 µF
DNA Concentration	Not given	Resistor	(Pulse Controller) Ω none
DNA Resuspension Buffer	Not given	Time Constant	0.4 to 0.8 msec
Volume of DNA	20 µl		
After the Pulse			
Outgrowth Medium	Not given		

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.

Outgrowth Temperature	0 °C for 10 min. (ice)
Length of Incubation	Not given
Selection Method or Assay Used	mycophenolic acid
Electroporation Efficiency	None
Per Cent Survival	40 to 50 %

Name of Submitter Rosaria Orlandi

Institution Address Istituto Nazionale Tumori
Oncologia Sperimentale
Via Venezian 1
Milano, 20100
ITALY

Survey Number
140