

Gene Pulser® Electroprotocols

* We recommend adapting this protocol to use the Gene Pulser electroporation buffer (catalog #165-2676, 165-2677), which increases cell viability and transfection efficiency in mammalian cell lines.

Cell Type	Mammalian, adherent, suspension	Molecules Electroporated	DNA: ccc <i>Bam</i> Z (5 kB), <i>Eco</i> A cosmid > 40 kB, <i>Sal</i> A cosmid > 40 kB
Species Used	Human, B lymphomas: BJAB, P3HR-1, B95-8; Rat-1		

Before the Pulse

Cell Growth Medium	RPMI-10 (GIBCO/BRL, Sigma)	Growth Phase at Harvest	Log
Wash Solution	None	Pre-pulse Incubation	10 min.

The Pulse

Instruments Used Gene Pulser® apparatus, Capacitance

Electroporation Temperature	4°C	Cuvette Gap	0.4 cm
Electroporation Medium*	RPMI-10 + 15% serum	Voltage	0.20 to 0.25 kV
Cell Density	5 x 10 (6) or 10 (7) per 350 µl	Field Strength	0.08 to 0.1 kV/cm
Volume of Cells	350 µl	Capacitor	960 µF
DNA Concentration	10 to 50 µg / cuvette	Resistor	(Pulse Controller) Ω none
DNA Resuspension Buffer	RPMI-10 or Phosphate Buffered Saline (PBS) or water	Time Constant	28 to 40 msec
Volume of DNA	10 µl		
After the Pulse			
Outgrowth Medium	Not given		

Outgrowth Temperature	37°C	Relevant Publications and/or Comments Note: exponential values designated in parentheses. <i>PNAS</i> 88:1546-1550 (1991). PBS: 1x = 8g NaCl, 0.2g KCl, 0.2g KH ₂ PO ₄ , 1.15g Na ₂ HPO ₄
Length of Incubation	Variable	
Selection Method or Assay Used	Protein expression by immunoblot, immunofluorescence, generation of recombinants.	
Electroporation Efficiency	> 20%	
Per Cent Survival	Variable %	

Name of Submitter Sankar Swaminathan, M.D.

Institution Address Harvard Medical School
Brigham & Women's Hospital
Dept of Infectious Diseases
20 Shattuck St.
Boston, MA 02115

Survey Number
157