

Gene Pulser® Electroprotocols

* We recommend adapting this protocol to use the Gene Pulser electroporation buffer (catalog #165-2676, 165-2677), which increases cell viability and transfection efficiency in mammalian cell lines.

Cell Type	Mammalian, adherent, suspension	Molecules Electroporated	DNA: pTS1-envIF, 11 kB, supercoiled
Species Used	Mouse, NIH/3T3, embryo; Human T-cell line, (PEER)		

Before the Pulse

Cell Growth Medium	MDEM (GIBCO/BRL, Sigma)	Growth Phase at Harvest	Not given
		Pre-pulse Incubation	MDEM
Wash Solution	No		

The Pulse

Instruments Used Gene Pulser® apparatus & Capacitance

Electroporation Temperature	25 °C	Cuvette Gap	0.4 cm
Electroporation Medium*	MDEM	Voltage	0.270 kV
Cell Density	75%	Field Strength	0.675 kV/cm
Volume of Cells	0.8 ml	Capacitor	960 µF
DNA Concentration	1 µg / µl	Resistor	(Pulse Controller) Ω none
DNA Resuspension Buffer	TE (10 mM Tris, 1 mM EDTA, pH 8.0)	Time Constant	12.0 msec
Volume of DNA	35 µl		
After the Pulse			
Outgrowth Medium	MDEM		

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.

Outgrowth Temperature	25 °C
Length of Incubation	10 min.
Selection Method or Assay Used	G418
Electroporation Efficiency	Very good
Per Cent Survival	80%

Name of Submitter Dr. Yuejin Yu, Postdoctoral Fellow

Institution Address University of Texas
M.D. Anderson Cancer Center
Science Park Research Division
Park Road 1C, P.O. Box 389
Smithville, Tx 78957

Survey Number
159