

Gene Pulser® Electroprotocols

* We recommend adapting this protocol to use the Gene Pulser electroporation buffer (catalog #165-2676, 165-2677), which increases cell viability and transfection efficiency in mammalian cell lines.

Cell Type	Mammalian, adherent	Molecules Electroporated	DNA: varies, fibronectin β -gal, genomic DNA, CMUB, etc.
Species Used	Human, fibroblast; Human Hep3b2, hepatocyte; Mouse, L-cells.		

Before the Pulse

Cell Growth Medium	MEM + 10% Fetal Calf Serum (GIBCO/BRL, Sigma)	Growth Phase at Harvest	Log
		Pre-pulse Incubation	5 minutes

Wash Solution Phosphate Buffered Saline and Trypsin

The Pulse **Instruments Used** Gene Pulser® apparatus & Capacitance

Electroporation Temperature	25 °C	Cuvette Gap	0.4 cm
Electroporation Medium*	MEM	Voltage	0.320 kV
Cell Density	1 to 10 x 10 (6) cells / pulse	Field Strength	0.8 kV/cm
Volume of Cells	500 μ l	Capacitor	500 μ F
DNA Concentration	40 to 100 μ g DNA per pulse	Resistor	(Pulse Controller) Ω none
DNA Resuspension Buffer	TE (10 mM Tris, 1 mM EDTA, pH 8.0)	Time Constant	Not given
Volume of DNA	20 to 40 μ l		

After the Pulse
Outgrowth Medium MEM

Relevant Publications and/or Comments
Note: exponential values designated in parentheses.

Outgrowth Temperature	37 °C
Length of Incubation	48 hours to 1 month
Selection Method or Assay Used	GPT and Neomycin
Electroporation Efficiency	60%
Per Cent Survival	50%

Name of Submitter Rane Taylor - Student

Institution Address Baylor College of Medicine
Molecular Genetics
1 Baylor Plaza
Houston, Texas 77030

Survey Number
163