

Gene Pulser® Electroprotocols

* We recommend adapting this protocol to use the Gene Pulser electroporation buffer (catalog #165-2676, 165-2677), which increases cell viability and transfection efficiency in mammalian cell lines.

Cell Type	Mammalian, adherent	Molecules Electroporated	DNA: linearized 8.5 kB plasmid
Species Used	Hybrid, mouse/human, A9 fibroblast, hybrid containing human chromosomes		

Before the Pulse

Cell Growth Medium	DMEM + 10% Fetal Calf Serum + mycophenolic acid + xanthine	Growth Phase at Harvest	10 (7) cells / ml
		Pre-pulse Incubation	10 min. on ice, (0°C)
Wash Solution	HEPES Buffered Saline (HBS)		

The Pulse

Instruments Used Gene Pulser® apparatus & Capacitance

Electroporation Temperature	Room temperature	Cuvette Gap	0.4 cm
Electroporation Medium*	HEPES Buffered Saline (HBS)	Voltage	0.250, 0.450 kV
Cell Density	Approximately 10 (7) cells / ml	Field Strength	0.625, 1.125 kV/cm
Volume of Cells	0.8 ml	Capacitor	125, 250, 500, 960 µF
DNA Concentration	10 µg	Resistor	(Pulse Controller) Ω none
DNA Resuspension Buffer	Not given	Time Constant	Not given
Volume of DNA	Not given		
After the Pulse			
Outgrowth Medium	Not given		

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.

Outgrowth Temperature	37 °C	HBS:	10mM HEPES, pH 7.2, 150 mM NaCl, 5 mM CaCl ₂
Length of Incubation	Not given		
Selection Method or Assay Used	Hygromycin B		
Electroporation Efficiency	Not given		
Per Cent Survival	Not given		

Name of Submitter Keith A. Rosenbach

Institution Address UMDNJ New Jersey Medical School
Microbiology & Molecular Genetics
185 South Orange Ave.
Newark, NJ 07103

Survey Number
165