

Gene Pulser® Electroporation

Cell Type Fungal / Yeast

Molecules Electroporated DNA: supercoiled on linear plasmids containing URA5 gene

Species Used *Cryptococcus neoformans*, ma5 mutants

Before the Pulse

Cell Growth Medium YEPD

Growth Phase at Harvest Logarithmic, O.D.(650) = ~ 1.

Pre-pulse Incubation None

Wash Solution 270 mM Sucrose, 1 mM MgCl₂, 10 mM Tris-HCl, pH 7.5, 4 mM DTT

The Pulse

Instruments Used Gene Pulser® apparatus

Electroporation Temperature Room temperature

Electroporation Medium 270 mM Sucrose, 10 mM Tris HCl, pH 7.5 (no DTT)

Cuvette Gap 0.2 cm

Cell Density Cells concentrated, 100 fold

Voltage 0.470 kV

Volume of Cells 450 µl

Field Strength 2.35 kV/cm

DNA Concentration 0.1 to 1.0 µg, in 1 to 10 µl TE

DNA Resuspension Buffer TE Buffer (10 mM Tris, 1 mM EDTA, pH 8.0)

Capacitor 25 µF

Volume of DNA 1 to 10 µl

Resistor (Pulse Controller) none Ω

After the Pulse

Time Constant 18 to 22 msec

Outgrowth Medium None

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.

Outgrowth Temperature Not given

Length of Incubation Not given

Selection Method or Assay Used SD media (lacking uracil): 6.7 g yeast nitrogen base per liter without amino acids and 20 g glucose / liter

Note: the protocol described in the paper by Edman and Kwon-Chung, *Mol. & Cell. Biol.*, **10**:4538-4544 (1990) is not the one described above. The one above gives 10-100x greater efficiency.

Electroporation Efficiency 2000 transfectants / µl

Per Cent Survival Not tested

Name of Submitter Jeffrey C. Edman, Asst. Prof.

Institution Address Univ. of California
Hormone Research Institute
Box 0534
San Francisco, CA 94143-0534

Survey Number

173