

Cell Type Fungal / Yeast

Molecules Electroporated DNA: supercoiled and linear

Species Used *Dictyostelium discoideum*

Before the Pulse

Cell Growth Medium HL5 (peptone, yeast extract, glucose)
(ATCC Media #671)

Growth Phase at Harvest 1×10^6 to 1×10^7 cells / ml

Pre-pulse Incubation 5 min.

Wash Solution Hepes Buffered Saline

The Pulse

Instruments Used Gene Pulser® apparatus & Pulse Controller

Electroporation Temperature 4 °C

Electroporation Medium Hepes Buffered Saline

Cuvette Gap 0.4 cm

Cell Density 5×10^6 / ml

Voltage 1.25 kV

Volume of Cells 1 ml

Field Strength 3.125 kV/cm

DNA Concentration 20 ng to 20 µg

DNA Resuspension Buffer TE (10 mM Tris, 1 mM EDTA, pH 8.0)

Capacitor 25 µF

Volume of DNA Any volume

Resistor (Pulse Controller) Ω none. NOT

After the Pulse

Time Constant 0.5 to 0.7 msec

Outgrowth Medium HL5

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.

**It is NOT RECOMMENDED to use high voltage with out the Pulse Controller.

HBS: 10mM HEPES, pH 7.2, 150 mM NaCl, 5 mM CaCl₂

Outgrowth Temperature 22 °C

Length of Incubation overnight

Selection Method or Assay Used G418

Electroporation Efficiency 10^{-3} transformants/ cell;
 6×10^3 transformants / µg DNA.

Per Cent Survival >90%

Name of Submitter David Knecht

Institution Address Univ of Connecticut
Dept of Molecular & Cell Biology
U-125
Storrs, CT 06269

Survey Number

175