

Gene Pulser® Electroprotocol

Cell Type Fungal / Yeast

Molecules Electroporated DNA: supercoiled plasmids, linear fragments for integration.

Species Used *Schizosaccharomyces pombe*

Before the Pulse

Cell Growth Medium YE or dropout media (recipes in paper)

Growth Phase at Harvest 1 x 10⁷ (7) cells / ml

Pre-pulse Incubation none

Wash Solution 1.2 M sorbitol (ice cold, filter sterilized)

The Pulse

Instruments Used Gene Pulser® apparatus & Pulse Controller

Electroporation Temperature (Ice cold) 0 °C

Electroporation Medium 1.2 M sorbitol

Cuvette Gap 0.2 cm

Cell Density 1 x 10⁴ (4) cells / ml

Voltage 2.25 kV

Volume of Cells 200 µl

Field Strength 1.125 kV/cm

DNA Concentration 1 ng to 1 µg DNA per pulse

DNA Resuspension Buffer 1.2 M sorbitol

Capacitor 25 µF

Volume of DNA <10 µl

Resistor (Pulse Controller) 200 Ω

After the Pulse

Time Constant 5 msec

Outgrowth Medium SD + necessary nutrients

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.

Reference: Prentice, H. *Nucleic Acid Res.* **20** (3):621.

Outgrowth Temperature 30 °C

Length of Incubation 4 to 6 days

Selection Method or Assay Used auxotrophy

Electroporation Efficiency 1 x 10⁵ (5) to 1 x 10⁶ (6) for autonomous plasmids

Per Cent Survival ~50% (seeTable 1 in reference)

Name of Submitter Holly Prentice

Institution Address Mass Gen Hosp./ Harvard Med School
Molecular Biology / Genetics Dept.
Wellman 10
Boston, MA 02114

Survey Number

176