

Gene Pulser® Electroprotocol

Cell Type Fungal / Yeast

Molecules Electroporated DNA: pUC19 and Bluescript™ -based plasmids.

Species Used *Saccharomyces cerevisiae* - lines derived from S288C

Before the Pulse

Cell Growth Medium YEPD (ATCC#1202/1245)

Growth Phase at Harvest 0.5 to 1.0 at O.D.600; ~100 mls log phase culture concentrated 200x, in water.

Pre-pulse Incubation None required

Wash Solution Water

The Pulse

Instruments Used Gene Pulser® apparatus & Pulse Controller

Electroporation Temperature Room temperature

Electroporation Medium Water

Cuvette Gap 0.2 cm

Voltage 0.6 kV

Cell Density 6 x 10 (10) cells / ml

Volume of Cells 40 to 50 µl

Field Strength 3.0 kV/cm

DNA Concentration Total DNA = 0.2 to 1.0 µg

DNA Resuspension Buffer Not given

Capacitor 25 µF

Volume of DNA ≤5µl

Resistor (Pulse Controller) 200 Ω

After the Pulse

Time Constant Not given

Outgrowth Medium Add water to final vol. 150-200 µl. Plate on selective synthetic media (SD)

Relevant Publications and/or Comments

Note: exponential values designated in parentheses. Washed and concentrated cells can be stored by adding 80% glycerol to a final concentration of 20% and freezing at -80°C. To electroporate thawed cells: pellet cells, remove glycerol containing supernatant, wash cells in 0.5-1.0 ml water, and resuspend cells in water to original volume. Efficiency is not affected significantly by freezing but you must remove glycerol. **Comment:** removing glycerol after thawing cells may help by simply removing lysed cell contents that would alter media conductivity - this could impact efficiency (alters the time constant, t).

Outgrowth Temperature 30 °C

Length of Incubation 2 days

Selection Method or Assay Used Amino acid or nucleotide dropout

Electroporation Efficiency 2 to 10,000 transfectants / µg DNA

Per Cent Survival Not given

Name of Submitter Joachim Li

Institution Address Univ of California /San Francisco
Dept. of Biochem & Biophysics, S-960
513 Parnassus
San Francisco, CA 94143-0448

Survey Number

179