

Cell Type	Plant	Molecules Electroporated	DNA: ds DNA mRNA
Species Used	Maize, Black mexican sweet		

Before the Pulse

Cell Growth Medium	BMS culture medium (see reference in notes)	Growth Phase at Harvest	Log
		Pre-pulse Incubation	5 min. at 42° C

Wash Solution Not given

The Pulse

Instruments Used Promega X-Cell

Electroporation Temperature	0 °C		
Electroporation Medium	10 mM HEPES 1mM CaCl ₂	Cuvette Gap	0.4 cm
Cell Density	2 x 10 ⁶ (6) / ml	Voltage	0. 200 kV
Volume of Cells	0.8 ml	Field Strength	0.5 kV/cm
DNA Concentration	10 µg / ml	Capacitor	1250 µF
DNA Resuspension Buffer	Not given	Resistor	(Pulse Controller) Do not calculate
Volume of DNA	Not given	Time Constant	Do not calculate
After the Pulse			
Outgrowth Medium	Plant growth media		

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.
 All of our recent work uses electroporation.
Ref: Fromm, *et. al.*, *Meth. Enzymol.* **153**:351-366 (1987).

Outgrowth Temperature	10 min. at 4 °C, then 25° C
Length of Incubation	up to 3 days
Selection Method or Assay Used	b-glucuronidase (GUS), luciferase
Electroporation Efficiency	50%
Per Cent Survival	90%

Name of Submitter Dr. Virginia Walbot

Institution Address Stanford University
 Dept. of Biological Sciences
 Stanford, CA 94305-5020

Survey Number
 181