

Cell Type Bacterial, gram negative

Molecules Electroporated DNA: M13 recombinant RF (double stranded) and single-stranded DNA

Species Used *E. coli*, MC1061, modified to contain F' (host for M13)

Before the Pulse

Cell Growth Medium LB

Growth Phase at Harvest O.D. (600) = 0.4 - 0.6

Pre-pulse Incubation none; done as fast as possible

Wash Solution 10% glycerol, 4 washes at 4°C

The Pulse

Instruments Used Gene Pulser® apparatus

Electroporation Temperature 0 to 4 °C

Electroporation Medium 10% glycerol

Cuvette Gap 0.2 cm

Cell Density As high as possible: just vortexed drained

Voltage 2.450 kV

Volume of Cells 40 µl

Field Strength 12.2 kV/cm

DNA Concentration 7.5 µg/µl

DNA Resuspension Buffer 10% glycerol

Capacitor 25 µF

Volume of DNA 5 µl + 40µl cells

Resistor 400 Ω (Pulse Controller)

After the Pulse

Time Constant 9.1 to 9.6 msec

Outgrowth Medium SOC

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.

Transformation efficiency for M13 is about 3x as high as pUC plasmid.

SOC: 2% Bacto tryptone, 0.5% Bacto yeast extract, 10mM NaCl, 2.5mM KCl, 10 mM MgCl₂, 10 mM MgSO₄, 20 mM glucose.

LB: 1% Bacto tryptone, 0.5% Bacto yeast extract, 0.5% NaCl.

Outgrowth Temperature 25 °C

Length of Incubation 1 to 10 min., plated immediately

Selection Method or Assay Used M13 plaques, no drugs

Electroporation Efficiency 1-5x10¹⁰ (RF); 1-50x10⁸ (single-stranded) transformants / µg DNA

Per Cent Survival Not measured

Name of Submitter Ted Gurney, Assoc. Professor

Institution Address University of Utah
Biology Department
Salt Lake City, UT 84112

Survey Number

036