

Cell Type Bacterial, gram negative

Molecules Electroporated DNA: plasmids, not supercoiled, pUC 13, pT3T7, pPL-lambda

Species Used *E. coli*, MC1061, NM522

Before the Pulse

Cell Growth Medium 10 g/l Tryptone, 10 g/l Yeast Extract, 10 g/l NaCl

Growth Phase at Harvest O.D. (600) = 0.9

Pre-pulse Incubation 1 min., 4°C

Wash Solution 10% glycerol

The Pulse

Instruments Used Gene Pulser® apparatus

Electroporation Temperature 4 °C

Electroporation Medium 10% glycerol

Cuvette Gap 0.1 cm

Voltage 2.0 kV

Cell Density 1/500 volume of initial culture

Volume of Cells 50 µl

Field Strength 20 kV/cm

DNA Concentration usually < or = 1 ng / µl

DNA Resuspension Buffer water

Capacitor 25 µF

Volume of DNA 10 µl

Resistor 200 Ω (Pulse Controller)

After the Pulse

Time Constant 4.82 msec

Outgrowth Medium 10g/l Tryptone, 10 g/l NaCl, 5 g/l Yeast Extract

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.

Outgrowth Temperature 37 °C

Length of Incubation 30 min

Selection Method or Assay Used ampicillin

Electroporation Efficiency 10 (8) to 10 (9) transformants / µg DNA

Per Cent Survival unknown

Name of Submitter Dr. Bruce Ross

Institution Address Fairfield Hospital, Clinical Pathology,
Yara Bend Road
Fairfield 3078, Victoria
Australia

Survey Number

037