

Cell Type Bacterial, gram negative
Species Used *E. coli*, MV1184, JM 109

Molecules Electroporated DNA: pUC-119, *Mollusum contagiosum* fragment

Before the Pulse

Cell Growth Medium SOC without glucose, LB

Growth Phase at Harvest O.D. (600) =0.55

Pre-pulse Incubation 0°C

Wash Solution water

The Pulse

Instruments Used Gene Pulser® apparatus

Electroporation Temperature 0°C

Electroporation Medium water

Cell Density 2 x 10 (11) / µl

Volume of Cells 1x10 (10) cells / ml

DNA Concentration Not given

DNA Resuspension Buffer Not given

Volume of DNA 5 µl

Cuvette Gap 0.2 cm

Voltage 2.5 kV

Field Strength 12.5 kV/cm

Capacitor 25 µF

Resistor 400 Ω (Pulse Controller)

Time Constant 7.1 msec

After the Pulse

Outgrowth Medium SOC without glucose
LB

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.
SOC: 2% Bacto tryptone, 0.5% Bacto yeast extract, 10mM NaCl, 2.5mM KCl, 10 mM MgCl₂, 10 mM MgSO₄, 20 mM glucose.
LB: 1% Bacto tryptone, 0.5% Bacto yeast extract, 0.5% NaCl.

Outgrowth Temperature 37 °C

Length of Incubation 60

Selection Method or Assay Used ABPC (CBPC)

Electroporation Efficiency Compared to the calcium salt method, efficiency is ~100x higher

Per Cent Survival Not given

Name of Submitter Dr. Jun Naleayamer

Institution Address Okayama University, Dept. of Medicine
Virology Section
Okayama-shi 2-5-1, Okayama
JAPAN

Survey Number

040