

Gene Pulser® Electroprotocol

Cell Type	Bacterial, gram negative	Molecules Electroporated	DNA: various plasmids
Species Used	<i>Bradyrhizobium japonicum</i> ; <i>E. coli</i> , species unspecified		

Before the Pulse

Cell Growth Medium	For <i>B. japonicum</i> , see reprint; for <i>E. coli</i> , standard protocol, see Pulse Controller or <i>E. coli</i> Pulser™ Manual.	Growth Phase at Harvest	O.D. (600) =varies - usually exponentially growing
---------------------------	---	--------------------------------	--

Pre-pulse Incubation 5 min on ice

Wash Solution Sterile distilled water

The Pulse

Instruments Used Gene Pulser® apparatus

Electroporation Temperature See notes

Electroporation Medium 10 % glycerol

Cell Density 10 (9) to 10 (10) cfu / ml

Volume of Cells 40 µl

DNA Concentration Varies: 12 ng / ml to 4 µg / ml

DNA Resuspension Buffer Distilled water

Volume of DNA 2 µl

Cuvette Gap 0.2 cm

Voltage 2.5 kV

Field Strength 12.5 kV/cm

Capacitor Not given

Resistor 200 Ω (Pulse Controller)

Time Constant 5 msec

After the Pulse

Outgrowth Medium For *B. japonicum*, see reprint: YEGG

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.

Electroporation temperature: cells and cuvette are on ice - then pulsed

Publications: Gerinot, M.L., Morisseau, B.A. & T. Klapatch. 1990. Electroporation of *Bradyrhizobium japonicum* *Mol. Gen. Genet.* **221**:287

The info on this sheet is for *B. japonicum*. See reprint. For *E. coli* we just use standard conditions (see Pulse Controller or *E. coli* Pulser™ Manual).

Outgrowth Temperature	30 °C
Length of Incubation	20 hr
Selection Method or Assay Used	Various drug resistances
Electroporation Efficiency	1.8 x 10 (5) (<i>B. japonicum</i>) transformants / µg DNA
Per Cent Survival	20 to 95 %

Name of Submitter Dr. Mary Lou Guerinot

Institution Address Dartmouth College, Biological Sciences Dept.
Gilman Hall
Hanover, NH 03755

Survey Number

048