

Cell Type	Bacterial, gram positive	Molecules Electroporated	DNA: pC194 (2.9 kB), pE194 (2.9 kB), pI258 (2.9 kB), pMH109 (7.4 kB), pMH120 (7.8 kB), supercoiled.
Species Used	<i>Staphylococcus aureus</i> , RN4220		

Before the Pulse

Cell Growth Medium	B2 broth (see notes)	Growth Phase at Harvest	Mid-log; 260 Klett units, #66 red filter.
		Pre-pulse Incubation	None

Wash Solution See notes

The Pulse **Instruments Used** Gene Pulser® apparatus & Pulse Controller

Electroporation Temperature Room temperature, 20 °C

Electroporation Medium 10% glycerol

Cuvette Gap 0.1 cm

Cell Density 3 x 10 (10) cells / ml

Voltage 2.3 kV

Volume of Cells 60 µl

Field Strength 23 kV/cm

DNA Concentration 17 µg DNA / µl

Capacitor 25 µF

DNA Resuspension Buffer B2 broth

Resistor (Pulse Controller) 100 Ω

Volume of DNA 1 to 6 µl

Time Constant 2.5 msec

After the Pulse

Outgrowth Medium NYE broth: 1% casein hydrolyzate (Sigma), 0.5% yeast extract, 0.5% NaCl

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.
B2 broth = 1% casein hydrolyzate (Sigma), 2.5% yeast extract (Difco), 0.5% glucose, 0.1% KHPO₄, 0.5% NaCl pH 7.5, overnight culture diluted 1/25 in 25 ml B2 broth in 300 ml Klett flask.
Wash Solution: 3 washes deionized water; 1/5 wash 10% glycerol; 1/10 wash 10% glycerol (Wash volume refers to volume of growth medium).

Outgrowth Temperature 37 °C

Length of Incubation 2 hours

Selection Method or Assay Used NYE agar and appropriate selective agent

Electroporation Efficiency 2 x 10 (8) CFU per µg DNA

Per Cent Survival 2 %

Name of Submitter Richard A. Laddaga, Ph.D. Asst. Prof.

Institution Address Bowling Green State University
 Biological Sciences
 Bowling Green, OH 43403-0212

Survey Number

083