

Gene Pulser® Electroprotocols

* We recommend adapting this protocol to use the Gene Pulser electroporation buffer (catalog #165-2676, 165-2677), which increases cell viability and transfection efficiency in mammalian cell lines.

Cell Type	Mammalian, adherent	Molecules Electroporated	Proteins: restriction enzymes (<i>Eco</i> RI, <i>Sca</i> I, <i>Dra</i> I) and catalase
Species Used	Hamster, CHO, ovary		

Before the Pulse

Cell Growth Medium	McCoy's 5a (GIBCO/BRL, Sigma)	Growth Phase at Harvest	Exponential
		Pre-pulse Incubation	Hepes Buffered Saline (HBS)
Wash Solution	Phosphate Buffered Saline (PBS) or serum-free medium		

The Pulse

Instruments Used Gene Pulser® apparatus & Capacitance

Electroporation Temperature	4°C		
Electroporation Medium*	HEPES Buffered Saline (HBS)	Cuvette Gap	0.4 cm
Cell Density	2 x 10 ⁶ (6) cells / ml	Voltage	0.3 kV
Volume of Cells	0.8 ml	Field Strength	0.750 kV/cm
DNA Concentration	Not given	Capacitor	960 µF
DNA Resuspension Buffer	Not given	Resistor	(Pulse Controller) Ω none
Volume of DNA	Not given	Time Constant	16 to 22 msec
After the Pulse			
Outgrowth Medium	McCoy's 5A		

Outgrowth Temperature	37 °C
Length of Incubation	20 hr.
Selection Method or Assay Used	Not given
Electroporation Efficiency	Not given
Per Cent Survival	40 to 60%

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.
Cortez, F., and Ortiz, T. *Mutation Research* **246**(1):221-6
PBS: 1x = 8g NaCl, 0.2g KCl, 0.2g KH₂PO₄, 1.15g Na₂HPO₄
HBS: 10mM HEPES, pH 7.2, 150 mM NaCl, 5 mM CaCl₂

Name of Submitter Dr. Felipe Cortes

Institution Address Faculty of Biology
Cell Biology Avenida Reina Mercedes s.n.
E-41012 Sevilla
Spain

Survey Number

090