

Gene Pulser® Electroprotocols

* We recommend adapting this protocol to use the Gene Pulser electroporation buffer (catalog #165-2676, 165-2677), which increases cell viability and transfection efficiency in mammalian cell lines.

Cell Type	Mammalian, adherent, suspension	Molecules Electroporated	DNA: plasmid
Species Used	Human, HeLa, epithelial carcinoma		

Before the Pulse

Cell Growth Medium	Not given	Growth Phase at Harvest	50 to 80%
		Pre-pulse Incubation	Ice, 10 minutes
Wash Solution	Not given		

The Pulse

Instruments Used Not given

Electroporation Temperature	Room temperature		
Electroporation Medium*	Phosphate Buffered Saline	Cuvette Gap	0.4 cm
Cell Density	1 x 10 ⁷ (7) cells / ml	Voltage	0.25 kV
Volume of Cells	800 µl	Field Strength	0.625 kV/cm
DNA Concentration	1 µg / ml		
DNA Resuspension Buffer	TE Buffer (10 mM Tris, 1 mM EDTA)	Capacitor	960 µF
Volume of DNA	5 µl	Resistor	(Pulse Controller) none Ω
		Time Constant	18 to 22 msec

After the Pulse

Outgrowth Medium Not given

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.
PBS: 1x = 8g NaCl, 0.2g KCl, 0.2g KH₂PO₄, 1.15g Na₂HPO₄

Outgrowth Temperature	Not given
Length of Incubation	Not given
Selection Method or Assay Used	Not given
Electroporation Efficiency	Not given
Per Cent Survival	Not given

Name of Submitter Phang Seng Meng, Graduate Student

Institution Address Dr. Thomas Leong Laboratory
National University of Singapore
Institute Molecular Cell & Biology
Kent Ridge Crescent
Singapore 0511

Survey Number
095