

Gene Pulser® Electroprotocols

* We recommend adapting this protocol to use the Gene Pulser electroporation buffer (catalog #165-2676, 165-2677), which increases cell viability and transfection efficiency in mammalian cell lines.

Cell Type	Mammalian, suspension	Molecules Electroporated	DNA: pBLCAT2 series, 4.5 kB
Species Used	Human, HeLa, epithelial carcinoma		

Before the Pulse

Cell Growth Medium	MEM + 10% Fetal Calf Serum (GIBCO/BRL, Sigma)	Growth Phase at Harvest	Log phase
		Pre-pulse Incubation	5 minutes on ice
Wash Solution	Phosphate Buffered Saline		

The Pulse

Instruments Used Gene Pulser® apparatus & Capacitance

Electroporation Temperature	4 °C (from ice to chamber at 25 °C)		
Electroporation Medium*	Phosphate Buffered Saline	Cuvette Gap	0.4 cm
Cell Density	2 x 10 ⁶ (6) cells in cuvette	Voltage	0.250 kV
Volume of Cells	0.8 ml	Field Strength	0.625 kV/cm
DNA Concentration	2 to 20 µg		
DNA Resuspension Buffer	TE Buffer (10 mM Tris, 1 mM EDTA)	Capacitor	960 µF
Volume of DNA	2 to 20 µl	Resistor	(Pulse Controller) Ω none
After the Pulse		Time Constant	18 to 20 msec
Outgrowth Medium	MEM + 10% Fetal Calf Serum		

		Relevant Publications and/or Comments
Outgrowth Temperature	37 °C	Note: exponential values designated in parentheses.
Length of Incubation	40 hours	
Selection Method or Assay Used	CAT assay	<i>Nucleic Acid Research</i> , 18 (3): 465-470 (1990).
Electroporation Efficiency	Not determined	PBS: 1x = 8g NaCl, 0.2g KCl, 0.2g KH ₂ PO ₄ , 1.15g Na ₂ HPO ₄
Per Cent Survival	20 to 60%	

Name of Submitter Mr. Terence Chong

Institution Address National University of Singapore
Institute Molecular Cell & Biology
Kent Ridge Crecent
Singapore 0511

Survey Number
096