




Nuvia aPrime 4A Hydrophobic Anion Exchange Resin



HYDROPHOBIC INTERACTION + ANION EXCHANGE

- Wide design space
- Optimal product recovery
- Mechanical and chemical stability
- Full regulatory support

High Purity and Recovery in Downstream Purification Processes

Nuvia aPrime 4A Hydrophobic Anion Exchange Resin possesses the optimal balance of ion exchange and hydrophobic interactions to deliver simultaneous purity and yield of therapeutic proteins and monoclonal antibodies that are difficult to purify.

Benefits of Nuvia aPrime 4A Resin

- Removes multiple impurities, such as viruses, host cell proteins, aggregates, and nucleic acids in a single chromatography step
- Enables capture and polish purification without extensive feedstream conditioning, saving process time and cost
- Facilitates straightforward method development

Nuvia aPrime 4A Resin is designed with a positively charged hydrophobic ligand (Figure 1). It can operate across a wide range of salt concentrations and pH. The resin's ionic capacity and hydrophobicity are designed to facilitate selective and reversible binding of target molecules for maximal purity and recovery. The technical properties of Nuvia aPrime 4A Resin are listed in Table 1.

Nuvia aPrime 4A Resin can be used to purify both established therapeutic proteins and diverse new constructs in development, including those that lack affinity handles. It is also effective in the purification of salt- and pH-sensitive proteins with a high propensity for aggregation and/or degradation.




Fig. 1. Mixed-mode ligand for Nuvia aPrime 4A Resin.

Table 1. Properties of Nuvia aPrime 4A Resin.

Property	Description
Function	Mixed-mode, hydrophobic interaction and anion exchange
Functional group	Aromatic hydrophobic strong anion exchanger
Base matrix composition	Polycrylamide
Mean particle size	50 µm
Ionic capacity	80–120 µeq/ml
Dynamic binding capacity (DBC)*	≥50 mg/ml at 300 cm/hr
Recommended linear flow rate	50–300 cm/hr
Pressure/flow performance**	<3 bar at flow rate of 300 cm/hr
Compression factor	1.1–1.25
pH stability***	Short term: 2–14
Shipping solution	20% ethanol + 1 M NaCl
Regeneration	1 M NaOH
Sanitization	1 M NaOH
Storage conditions	20% ethanol
Storage temperature	Room temperature
Chemical stability†	1 M NaOH, 1 M HCl, 25% acetic acid, 8 M urea, 6 M guanidine HCl, 3 M NaCl, 1% Triton X-100, 20–70% ethanol, 30% isopropanol
Shelf life††	5 years

* 10% breakthrough capacity determined with 1.0 mg/ml of an acidic monoclonal antibody (mAb; pI ~6.9) in 20 mM NaPO₄, pH 7.8.

** 20 x 20 cm packed bed (1.25 compression factor).

*** No significant change in ligand density after 200 hr contact at 22°C.

† No significant change in ligand density after 1 week contact at 22°C.

†† Stored at room temperature in 20% ethanol + 1 M NaCl.

Exceptional Purity and Recovery

The Nuvia aPrime 4A bead is built to provide a good dynamic binding capacity (DBC) and excellent recovery of biomolecules under mild conditions. The 50 µm particle size offers good resolution and high productivity while keeping the column pressure relatively low.

The ionic capacity and hydrophobicity of Nuvia aPrime 4A Resin have been fine-tuned to facilitate selective and reversible binding of target molecules in order to achieve maximal purity and recovery. We tested the ability of this resin to clear high molecular weight species from an acidic mAb (pI ~6.9) feedstream. At pH 7.8, Nuvia aPrime 4A interacts with the mAb through both electrostatic and hydrophobic interactions. The mAb monomer was eluted at ~pH 6.0 with an overall recovery of ~90%. The high molecular weight impurities were retained until the resin was stripped with a pH 4.0 buffer (Figure 2).




Fig. 2. Purification of mAb S (pI ~6.9) on Nuvia aPrime 4A Resin. A solution containing ~10 mg of target mAb in 20 mM NaPO₄, pH 7.8 was injected onto a 1 ml column and eluted by a buffer at pH 6. Clear separation of the high molecular weight impurities from the monomeric fraction was observed. AU (—); pH (---). AU, absorbance unit; ND, not detected.

Straightforward Method Development

The mixed-mode nature of the Nuvia aPrime 4A ligand facilitates straightforward method development and process optimization. A simple design of experiment (DOE) protocol can also be used by method developers to optimize their process conditions. We tested the ease of method development for purification of a basic protein with a pI of ~8.45. Our goal was to find the optimal conditions for the highest monomer recovery and purity. A fractional factorial screening design with two center points was performed. Binding capacity was low under all tested conditions (Figure 3A), suggesting that flow-through purification was more appropriate for this protein. Figures 3B and 3C show the high molecular weight impurity clearance and monomer recovery in response to purification buffer composition (or purification condition).




Fig. 3. Basic protein (pI ~8.45) purification using Nuvia aPrime 4A Resin. A, effect of binding buffer pH and NaCl on the binding of the basic protein by Nuvia aPrime 4A Resin; B, effect of binding buffer pH and NaCl on monomer content in the flow-through fraction; C, effect of binding buffer pH and NaCl on monomer recovery in the flow-through fraction.

Excellent Pressure-Flow Properties

Nuvia aPrime 4A Resin is designed with a bead size optimized to achieve process-scale purification of biomolecules at high flow rates without being limited by column pressure. This provides an increase in productivity. The column pressure remains below 3 bar up to a linear velocity of 350 cm/hr (Figure 4). Nuvia aPrime 4A Resin’s fast mass transfer dynamics ensure efficient chromatography at these high flow rates, making it an ideal choice for manufacturing purifications.




Fig. 4. Pressure-flow performance of Nuvia aPrime 4A Resin. Nuvia aPrime 4A Resin slurry prepared in 1x PBS, pH 7.5 was packed into a 20 x 20 cm column by axial compression at 300 cm/hr with a compression factor of 1.25.

Consistent Binding Capacity at Fast Flow Rates

Nuvia aPrime 4A Resin is built to have superior mechanical strength and optimized pore structure to ensure that the DBC is maintained at the high flow rates required for process purification (Figure 5). This ensures that it delivers the productivity needed in downstream manufacturing.




Fig. 5. Effect of flow rate on Nuvia aPrime 4A Resin binding capacity for polyclonal human IgG (hIgG). Bio-Scale MT10 Column dimension: 1.2 x 8.8 cm; equilibration buffer: 50 mM sodium phosphate, pH 8.0. DBC, dynamic binding capacity.

Robust Performance and Recovery

Nuvia aPrime 4A Resin is built on a mechanically stable base bead. It is produced by a validated manufacturing process, which ensures batch-to-batch reproducibility. The chemical stability of Nuvia aPrime 4A allows the resin to perform consistently with minimal changes to dynamic binding capacity or recovery even after multiple exposures to NaOH (Figure 6).




Fig. 6. Stability, reusability, and recovery with Nuvia aPrime 4A Resin. Polyclonal hIgG (1.0 mg/ml) in equilibration buffer (50 mM Na₂PO₄, pH 8.0) was loaded onto a 1 ml Bio-Scale Mini Column (0.56 x 4 cm) packed with Nuvia aPrime 4A. The column was operated at 150 cm/hr. hIgG was eluted with 100 mM sodium citrate, pH 3.0. The column was cleaned in place for 30 min with 1 M NaOH. The 10% DBC and recovery of hIgG at linear velocity of 150 cm/hr was determined after every ten cycles. DBC (—◆—); recovery (---■---). DBC, dynamic binding capacity; hIgG, human IgG.

Viral Clearance

The viral clearance properties of Nuvia aPrime 4A Resin were evaluated with minute virus of mice (MVM), a nonenveloped DNA parvovirus, and xenotropic murine leukemia virus-related virus (XMuLV), an enveloped mRNA retrovirus. Nuvia aPrime 4A, a mixed-mode resin, has a wide design space and excellent performance across various salt concentrations and pH; various conditions were initially screened using spin columns to generate a design of experiment (DOE) model. Conditions selected for further screening on prepacked 5 ml Foresight™ Nuvia aPrime 4A Columns are listed in Table 2. Studies were conducted with a third party in accordance with ICH Q5A(R2)* guidelines. Purified IgG monoclonal antibody (mAb) was prepared in the respective buffer conditions and the samples were spiked with the virus stock solution as shown. Samples at ~50 mg protein/ml of resin were loaded onto 5 ml columns at 120 cm/hr, and the virus particle concentration in the flowthrough measured. The log₁₀ reduction factor under various buffer conditions is shown in Table 2.

Table 2. Viral clearance of flow-through chromatography on Nuvia aPrime 4A Resin

Virus	Buffer	Log ₁₀ Reduction Factor
MVM	10 mM NaPO ₄ , 0 mM NaCl, pH 6.0	2.23
MVM	10 mM NaPO ₄ , 150 mM NaCl, pH 8.0	>4.21
MVM	10 mM NaPO ₄ , 200 mM NaCl, pH 8.0	>4.53
MVM	10 mM Tris, 150 mM NaCl, pH 8.0	>4.70
XMuLA	10 mM NaPO ₄ , 150 mM NaCl, pH 8.0	>5.02

MVM, minute virus of mice; xMuLA, xenotropic murine leukemia virus-related virus.

* International Council for Harmonisation of Technical Requirements for Pharmaceuticals for Human Use guidance entitled Q5A(R2) Viral Safety Evaluation of Biotechnology Products Derived from Cell Lines of Human or Animal Origin.

As MVM has a pI of 6.1–6.2, even when neutrally charged as in the first experimental condition at pH 6.0, MVM was reduced over 100-fold due to hydrophobic interactions with the resin. Typically, with traditional ion exchange chromatography (IEX), minimal viral clearance is observed when operating near the pI of the virus. The results from the 5 ml column runs were greater than the value predicted by the preliminary spin column experiments (<1.50 logs predicted versus 2.23 achieved) likely due to the differences in static binding and dynamic flow conditions. Under basic pH conditions, MVM will have a net negative charge and thus be adsorbed due to the anionic mode of interaction; additionally, the higher salt concentration enhances the hydrophobic interaction. This combination of results suggests that the viral particles are bound by both electrostatic and hydrophobic interactions, further highlighting the power of utilizing more than one mode of interaction for viral clearance.

Easy Scalability from Laboratory to Bioprocess Manufacturing

Nuvia aPrime 4A is specifically designed for easy scalability to meet manufacturing demands. It is available in multiple user-friendly formats, including prepacked Foresight Columns and Plates for purification condition screening and bulk bottles for pilot-to manufacturing-scale purifications. It is backed by our regulatory support documentation and security of supply commitment.

Ordering Information

Catalog # Description

Bulk Resins

12007397	Nuvia aPrime 4A Media, 25 ml
12007396	Nuvia aPrime 4A Media, 100 ml
12007379	Nuvia aPrime 4A Media, 500 ml
12007380	Nuvia aPrime 4A Media, 5 L
12007391	Nuvia aPrime 4A Media, 10 L

Large volumes and special packaging for industrial applications are available upon request.

Prepacked Screening Tools

12009315	EconoFit Aggregate Clearance Kit , includes 1 ml each of Nuvia cPrime, Nuvia aPrime 4A, Nuvia HR-S, CHT™ XT Media
12009314	EconoFit Large Biomolecule Kit , includes 1 ml each of Nuvia Q; Nuvia HP-Q; Nuvia cPrime; Nuvia aPrime 4A; CHT Type II, 40 µm; and CHT XT Media
12009313	EconoFit mAb Kit , includes 1 ml each of Nuvia S, Nuvia Q, Nuvia cPrime, Nuvia aPrime 4A, Nuvia HR-S, and CHT XT Media
12009280	EconoFit Nuvia aPrime 4A Column , 1 x 1 ml, 7 x 25 mm
12007411	Foresight Nuvia aPrime 4A Plates , 2 x 96-well, 20 µl
12007394	Foresight Nuvia aPrime 4A RoboColumn Unit , 200 µl
12007395	Foresight Nuvia aPrime 4A RoboColumn Unit , 600 µl
12007392	Foresight Nuvia aPrime 4A Column , 1 ml
12007393	Foresight Nuvia aPrime 4A Column , 5 ml

Technical Assistance

A regulatory support file is available upon request. Bio-Rad Laboratories, Inc. is an ISO 13485 registered corporation.

Visit [bio-rad.com/NuviaaPrime4A](https://www.bio-rad.com/NuviaaPrime4A) for more information and to request samples.

Visit [bio-rad.com/ProcessResins](https://www.bio-rad.com/ProcessResins) for more information about Bio-Rad's complete line of process chromatography supports.

For additional information and technical assistance, contact your local Bio-Rad office or email our process specialists at process@bio-rad.com. In the U.S. and Canada, call 1-800-4BIORAD (1-800-424-6723).

BIO-RAD, CHT, and FORESIGHT are trademarks of Bio-Rad Laboratories, Inc. in certain jurisdictions. Nuvia aPrime 4A Media is covered by U.S. Patent Number 9,669,402 and foreign counterparts. All trademarks used herein are the property of their respective owner. © 2024 Bio-Rad Laboratories, Inc.

BIO-RAD

**Bio-Rad
Laboratories, Inc.**

Life Science
Group

Website [bio-rad.com](https://www.bio-rad.com) **USA** 1 800 424 6723 **Australia** 61 2 9914 2800 **Austria** 00 800 00 24 67 23 **Belgium** 00 800 00 24 67 23 **Brazil** 4003 0399 **Canada** 1 905 364 3435 **China** 86 21 6169 8500 **Czech Republic** 00 800 00 24 67 23 **Denmark** 00 800 00 24 67 23 **Finland** 00 800 00 24 67 23 **France** 00 800 00 24 67 23 **Germany** 00 800 00 24 67 23 **Hong Kong** 852 2789 3300 **Hungary** 00 800 00 24 67 23 **India** 91 124 4029300 **Israel** 0 3 9636050 **Italy** 00 800 00 24 67 23 **Japan** 81 3 6361 7000 **Korea** 82 080 007 7373 **Luxembourg** 00 800 00 24 67 23 **Mexico** 52 555 488 7670 **The Netherlands** 00 800 00 24 67 23 **New Zealand** 64 9 415 2280 **Norway** 00 800 00 24 67 23 **Poland** 00 800 00 24 67 23 **Portugal** 00 800 00 24 67 23 **Russian Federation** 00 800 00 24 67 23 **Singapore** 65 6415 3188 **South Africa** 00 800 00 24 67 23 **Spain** 00 800 00 24 67 23 **Sweden** 00 800 00 24 67 23 **Switzerland** 00 800 00 24 67 23 **Taiwan** 886 2 2578 7189 **Thailand** 66 2 651 8311 **United Arab Emirates** 36 1 459 6150 **United Kingdom** 00 800 00 24 67 23

