iQ-Check *E. coli* O157:H7 Kit

User Guide

Test for the real-time PCR detection of *Escherichia coli* O157:H7
in food and environmental samples

Catalog # 3578114

The Bio-Rad logo, featuring the text "BIO-RAD" in white, bold, sans-serif capital letters inside a green rounded rectangle.

Table of Contents

Section 1.	Introduction	1
Section 2.	The iQ-Check <i>E. coli</i> O157:H7 Technology.....	1
Section 3.	Kit Components	2
Section 4.	Shelf Life and Storage	2
Section 5.	Materials Required but Not Supplied	2
	Equipment.....	2
	Supplies	3
Section 6.	Safety Precautions and Recommendations for Best Results	4
Section 7.	Protocol.....	5
	A. Sample Enrichment	5
	B. Free DNA Removal Treatment.....	7
	C. DNA Extraction	7
	D. Real-Time PCR	8
	E. Data Analysis.....	9
Section 8.	Confirmation of Positive Results.....	10
Section 9.	Confirmation of Single Colonies Using iQ-Check Kit.....	11
Section 10.	Test Performance and Validations.....	11
Section 11.	References.....	12
Section 12.	Revision History	12
	Appendix — PCR Mix Calculation Guide.....	13

Section 1

Introduction

Conventional bacteriological detection methods are often long and tedious. In comparison, iQ-Check *Escherichia coli* O157:H7 is a simple and rapid qualitative test, allowing the detection of specific DNA sequences unique to *E. coli* O157:H7 found in food products and environmental samples. Using real-time polymerase chain reaction (PCR), *E. coli* O157:H7 specific DNA sequences are amplified and detected simultaneously by means of fluorescent probes. Up to 94 samples can be processed, with a minimized risk of contamination and an easy-to-use procedure. The intended users of this kit are trained laboratory personnel who are performing tests to detect *E. coli* O157:H7. The use of this test allows results to be obtained within a few hours following enrichment of a sample.

Section 2

The iQ-Check *E. coli* O157:H7 Technology


The iQ-Check *E. coli* O157:H7 Kit is a test based on gene amplification and detection by real-time PCR. Ready-to-use PCR reagents contain oligonucleotides (primers and probes) specific for *E. coli* O157:H7, as well as a DNA polymerase and nucleotides. Detection and data analysis are optimized for use with a Bio-Rad real-time PCR instrument, such as the CFX96 Touch Deep Well Real-Time PCR Detection System.

PCR is a powerful technique used to generate many copies of target DNA. During the PCR reaction, several cycles of heating and cooling allow DNA denaturation followed by primers annealing to the target region at which point DNA polymerase uses these primers and deoxynucleotide triphosphates (dNTPs) to extend the DNA, creating copies of the target DNA. These copies are called amplicons.

In real-time PCR, specific probes are used to detect DNA during the amplification step by hybridizing to the amplicons. These probes are linked to a fluorophore that fluoresces only when hybridized to the target sequence. FAM is the fluorophore linked to the probe hybridizing to the *E. coli* O157:H7-specific DNA sequence. In the absence of target DNA, no fluorescence will be detected. As the amount of amplicons increases with each round of amplification, fluorescence intensity also increases. The optical module measures this fluorescence at the annealing step during each PCR cycle while the associated software plots the fluorescence intensity versus number of cycles.

A synthetic DNA internal control is included in the reaction mix. This control is amplified with a specific probe at the same time as the *E. coli* O157:H7 target DNA sequence and detected by a second fluorophore.

This test allows the qualitative detection of *E. coli* O157:H7 in select foods and environmental samples previously enriched by culture. It includes the following main steps:

*Please refer to the user guide of the iQ-Check Free DNA Removal Solution (#10000058391) for the conditions of use.

Section 3

Kit Components

The iQ-Check *E. coli* O157:H7 Kit contains sufficient reagents for 96 tests (94 samples).

Reagent ID	Reagent	Quantity Provided, ml
A	Lysis reagent	1 bottle, 20
B	Fluorescent probes	1 tube, 0.55
C	Amplification mix	1 tube, 4.4
D	PCR negative control	1 tube, 0.5
E	PCR positive control	1 tube, 0.25
F	Lysis beads	1 bottle (17.6 g)

Section 4

Shelf Life and Storage

Once received, the kit must be stored at 2–8°C. Reagents stored at this temperature can be used until the expiration date indicated on the tubes.

Section 5

Materials Required but Not Supplied

Equipment

- Lab paddle blender for homogenizing test samples
- Incubator for sample microbiological enrichment
- Specific for extraction in sterile 1.5 ml conical screwcap tubes:
 - Benchtop centrifuge capable of 10,000–12,000 x g
 - Dry heat block at 37 ± 2°C and/or 95–100°C
- Specific for extraction in deep well plate:
 - Heating thermoshaker/agitator-incubator* capable of maintaining 37 ± 2°C and/or 95–100°C, with a mixing speed of at least 1,300 rpm
- Vortexer
- Magnetic stir plate

Section 5 Materials Required but Not Supplied

- 20, 200, and 1,000 µl micropipets
- Tips for repeat pipettors; sterile, individually packaged
- Bio-Rad real-time PCR system;* for example, the CFX96 Touch Deep Well (catalog #3600037) or CFX Opus Deepwell (catalog #17007991) Real-Time PCR Systems
- Bio-Rad iQ-Check Prep System for automated DNA extraction and PCR plate setup (catalog #3594911)

Note: We recommend using an uninterrupted power supply (UPS) with the thermal cycler and iQ-Check Prep Systems.

*Contact Bio-Rad Technical Support for information on recommended instruments.

Supplies

- Enrichment medium: buffered peptone water (for example, BPW Plus catalog #3564684, dehydrated, 500 g; 3554179, 225 ml x 6 bottles; 3555789, 2.3 L x 5 bags; 3555790, 5 L x 2 bags. BPW Standard catalog #12013259, dehydrated, 500 g; 12013258 dehydrated, 5 kg; 12013260 5L x 2 bags)
- iQ-Check Free DNA Removal Solution (catalog #3594970)
- STEC supplement (catalog #3564005)
- iQ-Check Purification Reagent (catalog #12012383)
- Specific for environmental samples
 - Environmental sponges
 - Environmental swabs
 - Neutralizing broth for sponges and swabs, such as Dey-Engley (D/E) or HiCap Neutralizing Broth, or Lethen broth
- Specific for extraction in tubes
 - 1.5 ml conical screwcap tubes, sterile (catalog #2240110XTU)
- Specific for extraction in a deep well plate
 - 96-well deep well plate (iQ-Check Deep Well Microplates, catalog #3594900)
 - Plastic sealing film (TeSeE NSP Plastic Sealing Film, catalog #3590139)
 - PCR plate sealing film (X-Pierce Films, catalog #3593977, or Pre-Pierced Plate Sealing Film, #3600040, North America only)
- Specific for iQ-Check Prep System:
 - 60 ml dilution container (catalog #12014473)
 - Filter tips (catalog # 12014486, 50 µl x 5,760; 12014483, 1,000 µl x 3,840)
 - PCR mix tubes (catalog #12016673, 5 ml x 25)
- PCR plates, tubes, sealing tape, and caps
- Sterile filter tips adaptable to 20, 200, and 1,000 µl micropipets
- Tips for Combitip pipets or equivalent repeat pipettors; sterile, individually packaged
- 2 and 5 ml sterile test tubes

Section 6 Safety Precautions and Recommendations for Best Results

- RAPID[®]*E.coli* O157:H7 agar (catalog #3564748, 100 g) supplemented with novobiocin and potassium tellurite (can be used outside the scope of AOAC validation)
- For confirmations: CT-SMAC agar, plus O157 and H7 latex tests
- Powder-free gloves
- Distilled sterile water
- Bleach, 5%
- Cleaning agent such as DNA AWAY or RNase AWAY

Section 6

Safety Precautions and Recommendations for Best Results

- This test must be performed by trained personnel
- Samples and enrichment cultures must be handled as potentially infectious material and discarded according to local rules and regulations
- All potentially infectious material should be autoclaved before disposal
 - The quality of results depends on strict compliance with Good Laboratory Practices (for example, the EN ISO 7218 standard), especially concerning PCR:
 - Never circulate laboratory equipment (pipets, tubes, etc.) from one workstation to another
 - Always use a positive control and a negative control for each series of amplification reactions
 - Do not use reagents after their expiration date
 - Vortex reagents from the kit before using them to ensure homogeneity
 - Periodically verify the accuracy and precision of pipets, as well as correct functioning of the instruments
 - Change gloves often, especially if you suspect they are contaminated
 - Clean work spaces periodically with 5% bleach and other decontaminating agents such as DNA AWAY
 - Use powder-free gloves and avoid fingerprints and writing on tube caps. Both will interfere with data acquisition
- It is strongly advised to follow the general requirements described in the standard EN ISO 22174:2005 (Microbiology of food and animal feeding stuffs — Polymerase chain reaction (PCR) for the detection of food pathogens — General requirements and definitions)
- iQ-Check *E. coli* O157:H7 Kit
 - All substances or mixtures in the test kit are classified products, according to the Globally Harmonized System (GHS). Contact with acids may cause release of toxic gases. No special precautions are necessary if used correctly. If the product is inhaled, supply fresh air and consult a doctor in case of complaints. After eye contact with the product, rinse opened eye for several minutes under running water. If the products are swallowed, induce vomiting and call for medical help

Section 7 Protocol

- iQ-Check Prep System
 - Improper use of the iQ-Check Prep System may cause personal injury or damage to the instrument. Some components may pose a risk of personal injury due to excessive heat if improperly handled. For safe use, the iQ-Check Prep System must be operated only by qualified laboratory personnel who have been appropriately trained. Servicing of instrument must be performed only by Bio-Rad field service engineers
- CFX96 Touch or CFX Opus Deep Well Real-Time PCR Detection System
 - Improper use of the CFX96 Touch or CFX Opus Deep Well Real-Time PCR Detection System may cause personal injury or damage to the instrument. Some components may pose a risk of personal injury due to excessive heat if improperly handled. For safe use, the CFX96 Touch or CFX Opus Deep Well Real-Time PCR Detection System must be operated only by qualified laboratory personnel who have been appropriately trained. Servicing of instrument must be performed only by Bio-Rad Field Service Engineers
- Enrichment
 - The user should read, understand, and follow all safety information in the instructions for the iQ-Check *E. coli* O157:H7 Kit. Retain the safety instructions for future reference. To reduce the risks associated with exposure to chemicals and biohazards, perform pathogen testing in a properly equipped laboratory under the control of trained personnel. Always follow standard laboratory safety practices, including wearing appropriate protective apparel and eye protection while handling reagents and contaminated samples. Avoid contact with the contents of the enrichment media and reagent tubes after amplification. Dispose of enriched samples according to current industry standards
 - Pathogenic *E. coli* is a Biosafety Level 2 or 3 organism. Biological samples such as enrichments have the potential to transmit infectious diseases. Follow all applicable local, state/provincial, and/or national regulations on disposal of biological waste. Wear appropriate protective equipment, which includes but is not limited to: protective eyewear, face shield, clothing/lab coat, and gloves. All work should be conducted in properly equipped facilities utilizing the appropriate safety equipment (for example, physical containment devices). Individuals should be trained in accordance with applicable regulatory and company/institution requirements before working with potentially infectious materials
 - When testing is complete, all materials and media possibly containing pathogens should be decontaminated following current industry standards for the disposal of contaminated waste (that is, autoclave for 20 min at 120°C). Consult the Safety Data Sheet for additional information and local regulations for disposal

Section 7

Protocol

A. Sample Enrichment

It is strongly recommended to read the entire protocol before starting the test.

The use of an enrichment bag with incorporated filter is highly recommended

Enrichment media must be at the appropriate incubation temperature (37°C or 41.5°C when required) before use.

Section 7 Protocol

The following table outlines the different protocols that can be used depending on the application and the scope of the validation.

The iQ-Check *E. coli* O157:H7 kit is part of the USDA MLG 5C.02 standard method and can be used following the enrichment indicated in the method.

NF VALIDATION BRD 7/15-06/08		
Scope (Matrices)	Sample Preparation	Enrichment/DNA Extraction
Raw meat (excluding poultry) ^{1,2} (25 g)	Homogenize <i>n</i> g sample in 9 x <i>n</i> ml BPW (25 g in 225 ml)	Easy II extraction - Incubate 8-16 hr at 41.5 ± 1°C - Tube/deep well format
Raw meat (excluding poultry) ^{1,2} (375 g)	Homogenize <i>n</i> g sample in 3 x <i>n</i> ml BPW (375 g in 1,125 ml)	Easy II extraction - Incubate 10–18 hr at 41.5 ± 1°C - Tube/deep well format
AOAC PTM 020801		
Scope (Matrices)	Sample Preparation	Enrichment/DNA Extraction
Raw ground beef, fresh spinach, apple cider (25 g) ¹	Homogenize <i>n</i> g sample in 9 x <i>n</i> ml prewarmed BPW (25 g in 225 ml)	Easy I extraction - Incubate 8-24 hr at 41.5 ± 1°C - Tube/deep well format
Raw ground beef (25 g) ^{1,2}	Homogenize <i>n</i> g sample in 9 x <i>n</i> ml prewarmed BPW (25 g in 225 ml)	Standard II extraction - Incubate 8-12 hr at 36 ± 1°C - Tube/deep well format
Raw ground beef and raw beef trim ^{1,3} (375 g)	Homogenize <i>n</i> g sample in 3 x <i>n</i> ml prewarmed BPW (375 g in 1,125 ml)	Easy II extraction - Incubate 8-22 hr at 41.5 ± 1°C - Tube/deep well format
Fresh spinach ^{1,3} (375 g)	Homogenize <i>n</i> g sample in 3 x <i>n</i> ml prewarmed BPW (375 g in 1,125 ml)	Easy II extraction - Incubate 10-22 hr at 41.5 ± 1°C - Tube/deep well format
Raw chicken with and without skin, mechanically separated chicken and raw ground pork (25 g) ¹	Homogenize <i>n</i> g sample in 9 x <i>n</i> ml prewarmed BPW (25 g in 225 ml)	Easy II extraction - Incubate 8-22 hr at 41.5 ± 1°C - Tube/deep well format
MicroTally swab ^{1,4}	Homogenize swab in 200 ml prewarmed BPW by stomaching or hand massaging for 30 sec	Easy II extraction - Incubate 8-22 hr at 41.5 ± 1°C - Tube/deep well format
Environmental surface samples	Homogenize swabs in 10 ml prewarmed BPW by vortexing or shaking vigorously for 10 sec	Easy II extraction - Incubate 8-24 hr at 41.5 ± 1°C - Tube/deep well format
	Homogenize sponges in 60-90 ml prewarmed BPWn by stomaching or hand massaging for 30 sec	

¹ Validation also includes the use of the "O157_H7 Fast" Application Protocol File for a reduced PCR run time. Contact your Bio-Rad representative for more information.

² Protocol for simultaneous detection of *Salmonella*.

³ Harmonized protocol with iQ-Check *Salmonella* II and iQ-Check STEC VirX Kits

⁴ Harmonized protocol with iQ-Check STEC VirX Kit

B. Free DNA Removal Treatment

The iQ-Check Free DNA Removal Solution (catalog #3594970) provides an ideal way to get rid of free DNA (outside the current scope of AOAC and NF validations). Follow Bio-Rad's recommendations in the user guide.

C. DNA Extraction

General recommendations:

1. Turn on the heat block or thermoshaker to preheat before starting the test. Set it to 95–100°C. Keep the lysis reagent in suspension while pipetting by stirring at medium speed on a magnetic stir plate.
2. In general, avoid shaking the enrichment bag and collecting large fragments of food debris. For food samples with a fatty supernatant, collect the sample just below this layer.
3. Open tubes and wells carefully to avoid any possible cross contamination.
4. Cool the deep well plate before pipetting directly through pre-pierced sealing film.
5. Use the magnetic bar to keep the lysis reagent in suspension. Pipet while it is stirring at medium speed.
6. Gently shake the lysis reagent by hand first to resuspend the resin. Then pipet while it is stirring at medium speed with the magnetic bar contained in the bottle, in order to keep it in suspension.
7. Reconstitute the final lysis reagent as follows:
 - a. Carefully pour all the contents from reagent F (lysis beads) into reagent A (lysis reagent).
 - b. Use consumables with a wide enough tip to allow pipetting of the homogenized lysis reagent.
 - c. The lysis reagent mixed with lysis beads (reagents A + F) has a shelf life of 6 months when stored at 4°C.

Standard II Protocol

1. Collect 1 ml of decanted enriched sample into a tube or in a well of a deep well plate. Seal the deep well plate with a plastic film.
2. Centrifuge tubes at 10,000–12,000 x g for 5 min. Discard the supernatant. Centrifuge deep well plates at 2,250 x g for 20 min and discard supernatant manually or using the DW40.
3. Add 200 µl of homogenized lysis reagent (reagents A + F) to the pellet and resuspend pellet by pipetting the reagent up and down. Close tubes or seal deep well plate with pre-pierced sealing film.

Note: Gently shake lysis reagent by hand first to resuspend the beads.

4. Place tube in the cell disruptor for 3 ± 1 min, or shake deep well plate in the plate agitator-incubator at 1,300–1,600 rpm at 95–100°C for 15–20 min.
5. Incubate in the appropriate heat block at 95–100°C for 10–15 min.
6. Vortex tubes at high speed and centrifuge at 10,000–12,000 x g for 5 min. Centrifuge deep well plate at 2,250 x g for 2 min.

This is the recommended stopping point for temporarily stopping the procedure.

The supernatant can be stored for up to 1 year at –20°C. Always allow it to thaw and homogenize, and then centrifuge at 10,000–12,000 x g for 5 min before reusing.

Easy I Protocol

1. Aliquot 100 µl of homogenized lysis (reagent A) to tubes or wells of a deep well plate.
2. Add 100 µl of enriched sample. Mix by pipetting up and down and close the tube with caps or seal the deep well plate with pre-pierced sealing film.
3. Incubate in the appropriate heat block at 95–100°C for 10–15 min or in the plate thermoshaker for 15–20 min at 1,300 rpm.
4. Vortex tubes at high speed.
5. If using a deep well plate, allow it to cool to room temperature (20–25°C).
6. Centrifuge tubes at 10,000–12,000 x g for at least 2 min. Centrifugation is not needed for deep well plate.

This is the recommended stopping point for temporarily stopping the procedure.

The supernatant can be stored for up to 1 year at –20°C. Always allow it to thaw and homogenize, and then centrifuge at 10,000–12,000 x g for 5 min before reusing.

Easy II Protocol

1. Aliquot 100 µl of homogenized lysis reagent (reagents A + F) to the wells of a deep well plate.

Note: Gently shake the lysis reagent by hand first to resuspend the beads.

2. Add 100 µl of enriched sample.
3. Mix the solution by pipetting up and down until homogenized.
4. Close the tubes or seal the deep well plate with pre-pierced sealing film.
5. Place tubes in the cell disruptor for 3 ± 1 min. Incubate tubes in the heat block at 95–100°C for 10–15 min.
6. If using the thermoshaker for the tube or deep well format, incubate the tubes or deep well plate in the thermoshaker under agitation at 1,300 rpm at 95–100°C for 15–20 min.
7. Vortex tubes at high speed, and then centrifuge at 10,000–12,000 x g for at least 2 min. Centrifugation is not needed for deep well plate.

This is the recommended stopping point for temporarily stopping the procedure.

The supernatant can be stored for up to 1 year at –20°C. Always allow it to thaw and homogenize, and then centrifuge at 10,000–12,000 x g for 5 min before reusing.

D. Real-Time PCR

Instrument and Software Setup

For instrument and software setup, follow the instructions in the real-time PCR system user guide for iQ-Check Kits.

PCR Mix Preparation

1. Prepare a PCR mix containing the amplification solution (reagent C) and the fluorescent probes (reagent B). The volume of PCR mix needed depends on the number of samples and controls to be analyzed. At least one positive and one negative control must be included in each PCR run. Use the pipetting table in the Appendix to find the correct volumes to use for each reagent.

Note: Use the PCR mix (reagent B + C) immediately after preparation. It is stable for 1 hr maximum at 2–8°C.

2. Pipet 45 µl of the PCR mix into each well according to your plate setup.
3. Add 5 µl of sample, reagent D (negative control), or reagent E (positive control). Do not vortex the sample before pipetting. Hermetically seal the wells of the plate or strips. It is important to avoid bubbles at the bottom of the wells by pipetting carefully. As an optional step, centrifuge the sealed PCR plate or tube strips (quick spin) to eliminate any bubbles.
4. Place the plate or strips in the thermal cycler. Be sure to place the plate with the A1 well at the upper left corner. Close the reaction module.

Run PCR

To start the PCR run, follow instructions in the real-time PCR system user guide for iQ-Check Kits.

E. Data Analysis

Data can be analyzed directly at the end of the PCR run or at a later time by opening the stored data file. Follow instructions in the corresponding CFX Manager IDE Software User Manual for opening data files and setting the data analysis parameters.

Interpreting Results

Once the data analysis parameters have been set, results are interpreted by analyzing the Cq values of each sample (the cycle at which the amplification curve crosses the threshold).

CFX Manager IDE Software allows complete automated analysis for Bio-Rad real-time PCR detection systems. Verification of the typical characteristics of the amplification curves should be performed prior to releasing results. Contact your local Bio-Rad technical support team for additional support.

Controls

Verify the positive and negative controls before interpreting sample results.

For the experiment to be valid, the controls must have the following results, as summarized in the table below. Otherwise the PCR reaction must be repeated.

	<i>E. coli</i> O157:H7 Detection (FAM channel)	Internal Control Detection (HEX channel)
Negative control	Cq = N/A*	$28 \leq Cq \leq 40$
Positive control	$26 \leq Cq \leq 36$	NA

* The software indicates a Cq value of N/A (not applicable) when the fluorescence of a sample does not rise significantly above the background noise, and hence does not cross the threshold.

If results of negative and positive controls differ from those in the Controls table (invalid control), repeat the run and analysis described in D. Real-Time PCR and E. Data Analysis in Section 7 Protocol.

Samples

A **positive** iQ-Check *E. coli* O157:H7 PCR test must show a typical amplification curve and a Cq value ≥ 10 for the FAM fluorophore.

- If the Cq value for both channels is below 10, verify the raw data curve is a regular amplification curve (with a flat baseline followed by a rapid exponential increase of fluorescence and then a flattening out). If the curve seems correct, it may be considered a positive *E. coli* O157:H7 sample

Section 8 Confirmation of Positive Results

If there is no Cq value (Cq = N/A) for FAM or if the curve is not a typical amplification curve, the internal control for that sample must then be analyzed:

- If there is no Cq value for FAM and the internal control has a Cq ≥ 28 , this sample is considered as a **negative** *E. coli* O157:H7 sample
- If the internal control also has no Cq value (Cq = N/A), this probably indicates inhibition of the PCR reaction. Dilute the sample (perform a 1:10 dilution in distilled sterile water using 10 μ l of DNA extract), use 5 μ l of the dilution for amplification, and repeat the PCR test
- If the Cq value for the internal control is <28 , it is not possible to interpret the result. Verify that the threshold was correctly placed, or that the raw data curve is a regular amplification curve. If the curve does not have a characteristic shape, it will be necessary to repeat the PCR test

Interpretation of sample results is summarized in the following table:

<i>E. coli</i> O157:H7 Detection (FAM)	Internal Control Detection (HEX channel)	Interpretation
Cq ≥ 10	NA	Positive.
Cq = N/A	Cq > 28	Negative
Cq = N/A	Cq = N/A	Inhibition*

* When both target and internal control detection give a Cq value = N/A, the sample must be diluted 1/10 and tested again.

An invalid interpretation can be given when validation criteria are not met. Check the raw data and proceed as if the sample was inhibited.

Section 8

Confirmation of Positive Results

In the context of the NF VALIDATION certified method, all positive iQ-Check *E. coli* O157:H7 results must be confirmed in the following ways:

1. Using classic tests described in the standardized methods CEN or ISO by going back to the sample.
2. Using one of the following two methods, starting from the BPW enrichment:
 - a. Isolation on CT-SMAC agar (with or without an immunoseparation step first), followed by isolation of colonies on nonselective agar, then O157 and H7 latex tests on 1–3 characteristic colonies.
 - b. Isolation on RAPID *E. coli* O157:H7 chromogenic medium, with an immunoseparation step first, followed by isolation of colonies on nonselective agar, then O157 and H7 latex tests 1–3 on characteristic colonies. Refer to the confirmation method in the RAPID *E. coli* O157:H7 chromogenic media user guide

In the case of discrepant results between the iQ-Check *E. coli* O157:H7 Kit and one of the confirmation options listed above (particularly by latex tests), follow the necessary steps to ensure valid results, for example, adding an immunoseparation step and then isolating on CT-SMAC or RAPID *E. coli* O157:H7 agar.

In the context of the AOAC PTM validated method, a positive iQ-Check *E. coli* O157:H7 result is presumed to be positive, and confirmation according to an appropriate reference method (for example, USDA MLG or FDA BAM) is recommended.

Section 9

Confirmation of Single Colonies Using iQ-Check Kit

The iQ-Check *E. coli* O157:H7 Kit may also be used for confirming single isolated *E. coli* O157:H7 colonies on agar plates.

1. Pick an isolated colony from a selective or nonselective agar plate with a toothpick, sterile loop, or other adapted consumable (for example, a pipet tip).
2. Resuspend the colony in 100 µl tryptone salt or distilled sterile water in a microcentrifuge tube. Homogenize using a vortexer.
3. Use 5 µl of the suspension with 45 µl of PCR mix (see D. Real-Time PCR in Section 7 Protocol) and follow the rest of the iQ-Check *E. coli* O157:H7 protocol for the data and results interpretation.

Note that colony confirmation can be performed in accordance with the protocol described in the USDA MLG 5C.00 method.

Section 10

Test Performance and Validations

BRD 07/15 – 06/08

ALTERNATIVE
ANALYTICAL METHODS
FOR AGRIBUSINESS

<http://nf-validation.afnor.org>

iQ-Check *E. coli* O157:H7 (Easy II Protocol) is certified NF VALIDATION as an alternative method to the reference method ISO 16654 (July 2001) and its amendment A1 (June 2017), for the detection of *E. coli* O157:H7 in raw meat (excluding poultry meats) in 25 and 375 g samples. The validation followed the protocol of the NF EN ISO 16140-2:2016 standard, and includes the use of the CFX96, CFX96 Deep Well and CFX Opus Deepwell Real-Time PCR Systems. The associated software program is CFX Manager Software IDE (V2.2 and later) and the O157_H7 Fast APF. Certificate number: BRD 07/15–06/08. Valid until: For more information about the end of validity of the NF VALIDATION certification, please refer to the certificate available on the website <http://nf-validation.afnor.org/en>

iQ-Check *E. coli* O157:H7 has been validated by the AOAC-Research Institute under the Performance Tested Method Program for detection of *E. coli* O157:H7 from raw ground beef (25 and 375 g), raw beef trim (375 g), raw boneless chicken breast without skin (25 g), raw chicken thigh with skin (25 g), mechanically separated chicken (25 g), raw ground pork (25 g), fresh spinach (25 and 375 g), apple cider (25 ml) and MicroTally swabs (4" x 4" surface area). Standard Protocol II has also been validated by AOAC-RI for the simultaneous detection of *E. coli* O157:H7 and *Salmonella* spp. in raw ground beef. A positive result with iQ-Check should be considered presumptive and it is recommended that it be confirmed by standard reference methods. Certificate number 020801.

Section 11

References

ISO 16654. Microbiology of food and animal feeding stuffs – Horizontal method for the detection of *E. coli* O157.

ISO 16140. Microbiology of the food chain – Method validation – Part 2: Protocol for the validation of alternative (proprietary) methods against a reference method.

United States Department of Agriculture, Food Safety and Inspection Service (2015). Microbiology Laboratory Guidebook. Chapter 5.09: Detection, Isolation and Identification of *Escherichia coli* O157:H7 from Meat Products and Carcass and Environmental Sponges.

United States Department of Agriculture, Food Safety and Inspection Service (2021). Microbiology Laboratory Guidebook –Chapter 5C.02. Detection, Isolation and Identification of Top Seven Shiga Toxin-Producing *Escherichia coli* (STECs) from Meat Products and Carcass and Environmental Sponges.

United States Food and Drug Administration, Center for Food Safety and Applied Nutrition (2020). Bacteriological Analytical Manual. Chapter 4a. Diarrheagenic *Escherichia coli*.

Section 12

Revision History

Release date	Document number	Change
May 2020	10000131517 Ver A	- Renewal and extension of AFNOR validation for: 25 and 375 g meat samples and “O157_H7 Fast” APF - AOAC extension - New document design and updated references, content, USDA reference method, and colony confirmation - Document number change – previous version 808473 REV D
December 2022	10000131517 Ver B	- Added new catalog numbers for consumables - Replaced protocol section with summary table - AFNOR extension: CFX Opus Deepwell Real-Time PCR System and upgrade to CFX Manager Software IDE version 3.1 - Clarification of result interpretation protocol - General content updates
March 2024	10000131517 Ver C	- AOAC extension: MicroTally swabs - Clarification of extraction protocol

Appendix — PCR Mix Calculation Guide

To find the correct volumes to use when preparing the PCR mix, add the total number of samples and controls to be analyzed and find the corresponding volumes of reagent B and reagent C in the table.

Total Number of Samples and Controls	Probes Reagent B, µl	Amplification Mix Reagent C, µl	Total Number of Samples and Controls	Probes Reagent B, µl	Amplification Mix Reagent C, µl	Total Number of Samples and Controls	Probes Reagent B µl	Amplification Mix Reagent C, µl
1	5	40	33	178	1,400	65	351	2,800
2	11	86	34	184	1,500	66	356	2,900
3	16	130	35	189	1,500	67	362	2,900
4	22	173	36	194	1,600	68	367	2,900
5	27	216	37	200	1,600	69	373	3,000
6	32	259	38	205	1,600	70	378	3,000
7	38	302	39	211	1,700	71	383	3,100
8	43	346	40	216	1,700	72	389	3,100
9	49	389	41	221	1,800	73	394	3,200
10	54	432	42	227	1,800	74	400	3,200
11	59	475	43	232	1,900	75	405	3,200
12	65	518	44	238	1,900	76	410	3,300
13	70	562	45	243	1,900	77	416	3,300
14	76	605	46	248	2,000	78	421	3,400
15	81	648	47	254	2,000	79	427	3,400
16	86	691	48	259	2,100	80	432	3,500
17	92	734	49	265	2,100	81	437	3,500
18	97	778	50	270	2,200	82	443	3,500
19	103	821	51	275	2,200	83	448	3,600
20	108	864	52	281	2,200	84	454	3,600
21	113	907	53	286	2,300	85	459	3,700
22	119	950	54	292	2,300	86	464	3,700
23	124	994	55	297	2,400	87	470	3,800
24	130	1,000	56	302	2,400	88	475	3,800
25	135	1,100	57	308	2,500	89	481	3,800
26	140	1,100	58	313	2,500	90	486	3,900
27	146	1,200	59	319	2,500	91	491	3,900
28	151	1,200	60	324	2,600	92	497	4,000
29	157	1,300	61	329	2,600	93	502	4,000
30	162	1,300	62	335	2,700	94	508	4,100
31	167	1,300	63	340	2,700	95	513	4,100
32	173	1,400	64	346	2,800	96	518	4,100

Visit bio-rad.com/iqcheck for more information.

BIO-RAD is a trademark of Bio-Rad Laboratories, Inc.

IQ-CHECK is a trademark of Bio-Rad Europe GMBH in certain jurisdictions.

All trademarks used herein are the property of their respective owner.

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