

Gene Pulser® Electroprotocols

* We recommend adapting this protocol to use the Gene Pulser electroporation buffer (catalog #165-2676, 165-2677), which increases cell viability and transfection efficiency in mammalian cell lines.

Cell Type	Mammalian	Molecules Electroporated	DNA, various linear constructs of human sequences & selectable markers.
Species Used	Mouse, A-9, derivative of mouse L cell (contains human chromosomes)		

Before the Pulse

Cell Growth Medium	DMEM (GIBCO/BRL, Sigma)	Growth Phase at Harvest	Not given
		Pre-pulse Incubation	10 min. ice

Wash Solution Phosphate Buffered Saline, without Ca++, Mg++

The Pulse

Instruments Used Not given

Electroporation Temperature Cuvette on ice just prior to pulse

Electroporation Medium* Phosphate Buffered Saline, without Ca++, Mg++ **Cuvette Gap** 0.4 cm

Cell Density 10 (7) cells / ml

Voltage 0.45 kV

Volume of Cells 0.8 ml

Field Strength 1.125 kV/cm

DNA Concentration 1 µg / 1 µl

Capacitor 500 µF and 960 µF

DNA Resuspension Buffer water

Resistor (Pulse Controller) Ω none

Volume of DNA 10 µl

Time Constant Not given

After the Pulse

Outgrowth Medium DMEM

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.

Outgrowth Temperature	37 °C
Length of Incubation	48 hours
Selection Method or Assay Used	Hygromycin, G418, or MX
Electroporation Efficiency	Low
Per Cent Survival	Not given

Name of Submitter Keith Rosenbach

Institution Address UMDNJ
Microbiology Department
185 South Orange Ave
Newark, NJ 07103

Survey Number

120