

Gene Pulser® Electroprotocols

* We recommend adapting this protocol to use the Gene Pulser electroporation buffer (catalog #165-2676, 165-2677), which increases cell viability and transfection efficiency in mammalian cell lines.

Cell Type	Mammalian, suspension	Molecules Electroporated	DNA: linearized DNA used for stable transfections.
Species Used	Mouse, X-63, myeloma [p3 X63 - Ag8.653]		

Before the Pulse

Cell Growth Medium	DMEM, 10% Fetal Calf Serum (FCS) (GIBCO/BRL, Sigma)	Growth Phase at Harvest	50 to 70% confluency
		Pre-pulse Incubation	4° C, 10 min.

Wash Solution Wash two times in electroporation buffer

The Pulse **Instruments Used** Gene Pulser® apparatus & Capacitance

Electroporation Temperature	Room temperature		
Electroporation Medium*	Phosphate Buffered Saline	Cuvette Gap	0.4 cm
Cell Density	5 x 10 (5) cells/pulse, stable transfection	Voltage	0.180 kV
Volume of Cells	0.5 ml	Field Strength	0.45 kV/cm
DNA Concentration	10 µg / pulse		
DNA Resuspension Buffer	Not given; pulse volume: 0.8 ml	Capacitor	960 µF
Volume of DNA	Not given; pulse volume: 0.8 ml	Resistor	(Pulse Controller) Ω none
After the Pulse		Time Constant	24 msec
Outgrowth Medium	DMEM, 10% Fetal Calf Serum (FCS)		

Outgrowth Temperature	37 °C	Relevant Publications and/or Comments	
Length of Incubation	48 to 72 hrs.	Note: exponential values designated in parentheses.	
Selection Method or Assay Used	G418 (stable transfections)	Note: Stable transfections generally do not use carrier DNA. Also, the level of selective agent required to kill off non-transfected cells needs to be established before transfection. The level required should kill non-transfected cells in approximately 7 days.	
Electroporation Efficiency	Not given	PBS: 1x = 8g NaCl, 0.2g KCl, 0.2g KH ₂ PO ₄ , 1.15g Na ₂ HPO ₄	
Per Cent Survival	about 50%		

Name of Submitter Jackie Beall

Institution Address University of Adelaide
Department of Biochemistry
P.O. Box 498
Adelaide 5001
Australia

Survey Number
145