

Gene Pulser® Electroprotocols

* We recommend adapting this protocol to use the Gene Pulser electroporation buffer (catalog #165-2676, 165-2677), which increases cell viability and transfection efficiency in mammalian cell lines.

Cell Type	Mammalian, suspension	Molecules Electroporated	Fluorescent dyes (SBFI and MQAE); (3)H-inositol.
Species Used	Rat, submandibular acini (secretory cells).		

Before the Pulse

Cell Growth Medium	KRH + Ca	Growth Phase at Harvest	Not applicable
		Pre-pulse Incubation	A variation of the growth medium, but high in potassium and low in sodium and calcium.
Wash Solution	Cell growth medium		

The Pulse

Instruments Used Not given

Electroporation Temperature	Room temperature		
Electroporation Medium*	Same as pre-pulse	Cuvette Gap	0.4 cm
Cell Density	10 mg protein / ml	Voltage	0.6 to 1.25 kV
Volume of Cells	10 mg protein / ml	Field Strength	1.5 to 3.125 kV/cm
DNA Concentration	1 mg protein / ml		
DNA Resuspension Buffer	Same as Pre-Pulse buffer or KRH + Ca	Capacitor	25 μ F
Volume of DNA	Not assayed	Resistor	(Pulse Controller) Ω none.
		Time Constant	0.4 msec

After the Pulse

Outgrowth Medium Not applicable

Outgrowth Temperature	Not given	Relevant Publications and/or Comments	
Length of Incubation	Not given	Note:	exponential values designated in parentheses.
Selection Method or Assay Used	Not given		**It is NOT RECOMMENDED to use high voltage with out the Pulse Controller.
Electroporation Efficiency	Not given		
Per Cent Survival	Not given		

Name of Submittor Sheryll Barker, Research Technologist

Institution Address Lovelace Medical Foundation
Research Dept.
2425 Ridgerest Dr. S.E.
Albuquerque, NM 87108

Survey Number
147