

Gene Pulser® Electroprotocols

* We recommend adapting this protocol to use the Gene Pulser electroporation buffer (catalog #165-2676, 165-2677), which increases cell viability and transfection efficiency in mammalian cell lines.

Cell Type	Mammalian, adherent, suspension	Molecules Electroporated	DNA: linearized constructs
Species Used	Human, K562, chronic myeloid leukemia; Hamster, CHO, ovary; Hybrid, rat / mouse, MEL cells.		

Before the Pulse

Cell Growth Medium	RPMI 1640, 10% Fetal Bovine Serum (GIBCO/BRL, Sigma)	Growth Phase at Harvest	3 x 10 ⁵ to 5 x 10 ⁵ cells / ml (suspension)
		Pre-pulse Incubation	10 min. at room temp in Hepes buffered saline with dextrose
Wash Solution	Phosphate Buffred Saline		

The Pulse

Instruments Used Gene Pulser® apparatus & Capacitance

Electroporation Temperature	Room temperature		
Electroporation Medium*	HEPES buffered saline with dextrose	Cuvette Gap	0.4 cm
Cell Density	2 x 10 ⁷ cells / ml	Voltage	0.2 kV
Volume of Cells	0.4 ml	Field Strength	0.5 kV/cm
DNA Concentration	Not given		
DNA Resuspension Buffer	RPMI 1640 +10% Fetal Bovine Serum	Capacitor	960 µF
Volume of DNA	Not given	Resistor	(Pulse Controller) Ω none
After the Pulse		Time Constant	25.0 msec
Outgrowth Medium	RPMI 1640+10% Fetal Bovine Serum		

Outgrowth Temperature	37 °C
Length of Incubation	Not given
Selection Method or Assay Used	G418
Electroporation Efficiency	Not given
Per Cent Survival	Not given

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.
PBS: 1x = 8g NaCl, 0.2g KCl, 0.2g KH₂PO₄, 1.15g Na₂HPO₄
HBS: 10mM HEPES,pH 7.2,150 mM NaCl, 5 mM CaCl₂

Name of Submitter No name or address given

Institution Address

Survey Number

168