

Cell Type Bacterial, gram negative

Molecules Electroporated DNA: plasmid, covalently closed, pTTQ18, pUC, pBR322, M13

Species Used *E. coli*, DH5α

Before the Pulse

Cell Growth Medium LB or 2X-YT

Growth Phase at Harvest O.D. (600) = 0.5

Pre-pulse Incubation < 30 sec

Wash Solution 10 mM HEPES, pH 7.0, water, water /10 % glycerol

The Pulse

Instruments Used Gene Pulser® apparatus

Electroporation Temperature 0 °C (ice)

Electroporation Medium water or 10% glycerol/water
See comments.

Cuvette Gap 0.2 cm

Cell Density >10 (10) cells/ml

Voltage 2.5 kV

Volume of Cells 40 µl

Field Strength 12.5 kV/cm

DNA Concentration 2.5 - 250 µg/µl cell suspension

DNA Resuspension Buffer water, 1X TE, 10 mM HEPES

Capacitor 25 µF

Volume of DNA 1 to 2 µl

Resistor 200 Ω (Pulse Controller)

After the Pulse

Time Constant 4.6 to 4.8 msec

Outgrowth Medium SOC

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.

Comments: Necessary to clean up ligations before electroporation. Ligations ethanol precipitated or Gene Clean® purified and resuspended in water or 1mM HEPES, pH 7.0. Also found high efficiencies associated with scrupulously clean apparatus and all equipment and reagents pre-cooled to 0°C. High efficiencies associated with cells grown only to O.D. (600) = 0.5.

SOC: 2% Bacto tryptone, 0.5% Bacto yeast extract, 10mM NaCl, 2.5mM KCl, 10 mM MgCl₂, 10 mM MgSO₄, 20 mM glucose.

LB: 1% Bacto tryptone, 0.5% Bacto yeast extract, 0.5% NaCl.

Outgrowth Temperature 37 °C

Length of Incubation 1 hour

Selection Method or Assay Used ampicillin

Electroporation Efficiency Not given

Per Cent Survival 10-20%, depends on cell type

Name of Submitter Dr. Glenn Lilley

Institution Address CSIRO
Div. of Biomolecular Engineering
343 Rooyal Parade
Parkville, Victoria 3052
Australia

Survey Number

007