

Gene Pulser® Electroporation

Cell Type	Bacterial, gram positive	Molecules Electroporated	DNA: plasmid, covalently closed circular form, pYT937, 5.7 kB
Species Used	<i>Mycobacterium bovis</i> , BCG		

Before the Pulse

Cell Growth Medium	Middlebrook 7H9 (Difco)+ Tween 80	Growth Phase at Harvest	O.D. (600) =0.3
---------------------------	-----------------------------------	--------------------------------	-----------------

Pre-pulse Incubation	None
-----------------------------	------

Wash Solution	10 mM Tris (pH 7.0), 20% sucrose; or 7 mM phosphate buffer, 10% sucrose
----------------------	--

The Pulse	Instruments Used	Gene Pulser® apparatus
------------------	-------------------------	------------------------

Electroporation Temperature	4 °C
------------------------------------	------

Electroporation Medium	10 mM Tris (pH 7.0), 20% sucrose
-------------------------------	----------------------------------

Cuvette Gap	0.2 cm
--------------------	--------

Cell Density	10 (9) cells / ml
---------------------	-------------------

Voltage	1.25 kV
----------------	---------

Volume of Cells	200 µl
------------------------	--------

Field Strength	6.25 kV/cm
-----------------------	------------

DNA Concentration	1 to 2 µg
--------------------------	-----------

DNA Resuspension Buffer	Middlebrook 7H9 + Tween 80
--------------------------------	----------------------------

Capacitor	25 µF
------------------	-------

Volume of DNA	10 µl
----------------------	-------

Resistor	(Pulse Controller) 600 Ω
-----------------	--------------------------

After the Pulse

Time Constant	Not given
----------------------	-----------

Outgrowth Medium	Middlebrook 7H9 + Tween 80
-------------------------	----------------------------

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.

Outgrowth Temperature	37 °C
------------------------------	-------

Length of Incubation	overnight
-----------------------------	-----------

Selection Method or Assay Used	kanamycin resistance
---------------------------------------	----------------------

Electroporation Efficiency	10 (3) to 10 (4) transformants / µg
-----------------------------------	-------------------------------------

Per Cent Survival	50%
--------------------------	-----

Name of Submitter	Dr. Hatsumi Taniguchi
--------------------------	-----------------------

Institution Address	Dept. Microbiology, Univ. Occupational & Environmental Health, 1-1 Isseigaoka, Yahatanishiku, Kitakyushu, Fukuoka 80 Japan
----------------------------	---

Survey Number	
----------------------	--

078