

Gene Pulser® Electroprotocols

* We recommend adapting this protocol to use the Gene Pulser electroporation buffer (catalog #165-2676, 165-2677), which increases cell viability and transfection efficiency in mammalian cell lines.

Cell Type	Mammalian, adherent	Molecules Electroporated	DNA: plasmids of 8.4 kb, linearized
Species Used	Hamster, CHO, ovary		

Before the Pulse

Cell Growth Medium	DMEM, 10% Fetal Bovine Serum, NEAA, Pennicillin /Streptomycin	Growth Phase at Harvest	Log
		Pre-pulse Incubation	10 min on ice

Wash Solution Berg buffer (see Ref in notes)

The Pulse

Instruments Used Gene Pulser®

Electroporation Temperature	4 °C		
Electroporation Medium*	Berg buffer (see Ref. in notes)	Cuvette Gap	0.4 cm
Cell Density	10(8) cells /0.8ml	Voltage	0 .80 kV
Volume of Cells	0.8 ml	Field Strength	2.0 kV/cm
DNA Concentration	20 µg	Capacitor	25 µF
DNA Resuspension Buffer	water	Resistor	(Pulse Controller) Ω none
Volume of DNA	20 µl	Time Constant	0 .6 to 0.8 msec
After the Pulse			
Outgrowth Medium	DMEM, 10% FBS, NEAA, Pen/Strep		

Outgrowth Temperature	37 °C
Length of Incubation	2 weeks
Selection Method or Assay Used	G418, 8-Aza adenine
Electroporation Efficiency	10 transfectants / µg DNA
Per Cent Survival	20 to 80% (varies)

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.
Ref: Chu,G., Hayakawa, H. and Berg, P. *NAR* **15**, 1131.
Berg buffer: 20 mM HEPES, pH 7.05, 137 mM NaCl, 5 mM KCl, 0.7 mM Na₂HPO₄, 6 mM dextrose.

Name of Submittor Not given

Institution Address Not given

Survey Number

089