





# Nuvia HR-S Cation Exchange Resin



## CATION EXCHANGE

- Best-in-class resolution
- High recovery of target biomolecules
- Robust polymeric matrix
- Scalable to industrial purification applications
- Full regulatory support

## Exceptional Resolution of Related Species for Bioprocessing

Nuvia HR-S Resin is a strong cation exchanger (CEX), optimized for particle size and chemistry, that provides exceptional resolution and high recovery. The resin demonstrates fast mass transfer kinetics and superlative flow characteristics, chemical stability, and scalability. Nuvia HR-S Resin is the preferred solution when process developers face challenging separations of closely related biomolecules.

### High Resolution


Improvements in upstream cell culture expression systems have produced high feed stream titers that place a larger demand on downstream purifications to separate closely related contaminants from target biomolecules. In monoclonal antibody therapeutics, high molecular weight impurities in the form of aggregates are a key challenge for process developers. Nuvia HR-S Resin delivers the resolution performance that bioprocess scientists need to address this challenge. The technical properties of Nuvia HR-S Resin are listed in Table 1.

An example of a high-resolution separation is shown in Figure 1, where a mixture of monomer and aggregate was separated on Nuvia HR-S Resin. From a starting aggregate concentration of 8.9%, the final eluate pool contained 0.11% residual aggregate with an overall monomer recovery of 82%. Furthermore, as shown in Figure 2, better resolution was obtained using the Nuvia HR-S Resin with near baseline resolution of a protein mixture compared to Resin 1, another commercially available high-resolution small-particle cation exchanger.

**Table 1. Properties of Nuvia HR-S Resin.**

| Property                                                  | Description                                                          |
|-----------------------------------------------------------|----------------------------------------------------------------------|
| Type of ion exchanger                                     | Strong cation                                                        |
| Functional group                                          | $-\text{SO}_3^-$                                                     |
| Base matrix composition                                   | Polyacrylamide                                                       |
| Mean particle size                                        | 50 $\mu\text{m}$                                                     |
| Total ionic capacity                                      | 100–180 $\mu\text{eq/ml}$                                            |
| Dynamic binding capacity*                                 | >70 mg/ml at 300 cm/hr                                               |
| Recommended linear flow rate                              | 50–200 cm/hr                                                         |
| Compression factor (settled bed volume/packed bed volume) | 1.20–1.30                                                            |
| pH stability                                              | Short term: 2–14<br>Long term: 4–13                                  |
| Shipping solution                                         | 20% ethanol + 0.1 M NaCl or 2% benzyl alcohol                        |
| Regeneration                                              | 1–2 M NaCl                                                           |
| Sanitization                                              | 0.5–1.0 M NaOH                                                       |
| Storage conditions                                        | 20% ethanol or 0.1 M NaOH                                            |
| Storage temperature                                       | Room temperature                                                     |
| Chemical stability                                        | 1.0 M NaOH (20°C), up to 5 weeks or 0.1 M NaOH (20°C), up to 5 years |
| Shelf life                                                | 5 years                                                              |


\* 10% breakthrough capacity determined with 5.0 mg/ml human IgG in 20 mM Na acetate, pH 5.0.



**B**

| Chromatography Condition | Specification                                                      |
|--------------------------|--------------------------------------------------------------------|
| Column                   | 7 x 56 mm Nuvia HR-S Resin packed bed column                       |
| Load                     | 20 ml aggregate and monomer mixture at a loading level of 46 mg/ml |
| Flow rate                | 150 cm/hr                                                          |
| Equilibration            | 40 mM Na acetate, pH 5.0                                           |
| Gradient elution         | 40 mM Na acetate, 0–0.4 M NaCl, pH 5.0                             |

**Fig. 1. Monoclonal antibody separation on Nuvia HR-S Resin.** **A**, chromatogram of the purification of monomer from aggregate; **B**, chromatography conditions. A<sub>280</sub> (—); pH (—); gradient elution (—).




**Fig. 2. Resolution performance of Nuvia HR-S Resin.** Comparison of Nuvia HR-S Resin to another commercially available small-particle cation exchanger (Resin 1). Column size, 7 x 56 mm. A protein standard mixture of myoglobin, ribonuclease, and cytochrome C was loaded onto the columns in 20 mM sodium phosphate, pH 7.0, washed, and then eluted with a gradient of 20 mM sodium phosphate, pH 7.0 + 1 M NaCl. Nuvia HR-S (—); Resin 1 (—).

**Consistent and Reproducible Performance**

Chromatography media and resins used in downstream purification are required to perform consistently, even after repeated exposure to caustic sanitization protocols. Nuvia HR-S Resin binding capacity and protein recovery remained unchanged after exposure to 840 hr of 1.0 N NaOH (equivalent to 420 cycles at 2-hr sanitization hold), thus ensuring the resin can be used repeatedly without compromising performance, as shown in Figure 3.

Produced in a validated manufacturing process under strict specifications, the cornerstone of Nuvia HR-S Resin is batch-to-batch reproducibility. With more than 50 years of experience in chromatography, Bio-Rad understands and is committed to maintaining the highest level of quality and security of supply.




**Fig. 3. Stability of Nuvia HR-S Resin.** Results of hlgG binding capacity after exposing Nuvia HR-S Resin to 1.0 N NaOH. 10% breakthrough (—); % recovery (—). BT, breakthrough; DBC, dynamic binding capacity.

**Comparative Data**

**Superior Aggregate Content and Monomer Recovery Using Nuvia HR-S**


Nuvia HR-S Resin was compared to a commercially available small particle size CEX resin (Resin 1) to determine the aggregate content and monomer recovery as fractions were progressively pooled during elution. Nuvia HR-S Resin delivered a final aggregate content of <0.3% and a recovery of >80%, while the recovery using agarose-based Resin 1 was <70% for the same final aggregate content (Figure 4).



**Fig. 4. Performance of Nuvia HR-S (■) vs. Resin 1 (◆).**

### Minimal On-Column Aggregate Generation

Six commercial resins including Nuvia HR-S were tested for the formation of on-column aggregates during mAb purification. Nuvia HR-S produced the smallest amount of on-column aggregates, as shown in Figure 5.



**Fig. 5. On-column aggregates.** 0 min hold time (---); 1,000 min hold time (—); percentage of on-column aggregate shown in plots. For further details, refer to [Bio-Rad bulletin 7000](#), Minimizing On-Column Monoclonal Antibody Aggregate Formation.

### Viral Clearance

Nuvia HR-S viral clearance was evaluated using xenotropic murine leukemia virus (XMuLV), an enveloped mRNA retrovirus. Studies were conducted with a third party in accordance with ICH Q5A(R2)\* guidelines. Purified IgG monoclonal antibody (mAb) was prepared in 20 mM Na acetate, 25 mM NaCl, pH 5.0 to a concentration of 5 mg/ml. The samples were spiked with virus stock solution and loaded at >50 mg protein/ml of resin onto a prepacked 5 ml Foresight™ Nuvia HR-S Column (Bio-Rad, catalog #7324743). The mAb was eluted using ten column volumes (CV) at 300 cm/hr with a 0–100% linear gradient of 20 mM acetate, 500 mM NaCl, pH 5.0. The log<sub>10</sub> reduction factor was 3.76 logs.

\* International Council for Harmonisation of Technical Requirements for Pharmaceuticals for Human Use guidance entitled Q5A(R2) Viral Safety Evaluation of Biotechnology Products Derived from Cell Lines of Human or Animal Origin.

### Ordering Information

| Catalog #          | Description                      |
|--------------------|----------------------------------|
| <b>Bulk Resins</b> |                                  |
| 1560511            | Nuvia HR-S Media, 25 ml          |
| 12003712           | Nuvia HR-S Media (in BA), 25 ml  |
| 1560513            | Nuvia HR-S Media, 100 ml         |
| 12003713           | Nuvia HR-S Media (in BA), 100 ml |
| 156-0515           | Nuvia HR-S Media, 500 ml         |
| 12003714           | Nuvia HR-S Media (in BA), 500 ml |
| 12018125           | Nuvia HR-S Media, 5 L            |
| 156-0517           | Nuvia HR-S Media, 10 L           |
| 12003715           | Nuvia HR-S Media (in BA), 10 L   |

BA, 2% benzyl alcohol.

### Prepacked Screening Tools

|          |                                                                                                                                         |
|----------|-----------------------------------------------------------------------------------------------------------------------------------------|
| 12009310 | <b>EconoFit Cationic IEX Kit</b> , includes 1 ml each of Macro-Prep™ High S, UNOsphere S, Nuvia S, Nuvia HR-S, and Macro-Prep CM Resins |
| 12009284 | <b>EconoFit Nuvia HR-S Column</b> , 1 x 1 ml, 7 x 25 mm                                                                                 |
| 732-4707 | <b>Foresight Nuvia HR-S Plates</b> , 2 x 96-well, 20 µl                                                                                 |
| 732-4831 | <b>Foresight Nuvia HR-S RoboColumn Unit</b> , 200 µl                                                                                    |
| 732-4832 | <b>Foresight Nuvia HR-S RoboColumn Unit</b> , 600 µl                                                                                    |
| 732-4723 | <b>Foresight Nuvia HR-S Column</b> , 1 ml                                                                                               |
| 732-4743 | <b>Foresight Nuvia HR-S Column</b> , 5 ml                                                                                               |



Technical Assistance

A regulatory support file is available upon request. Bio-Rad Laboratories, Inc. is an ISO 13485 registered corporation.

Visit [bio-rad.com/NuviaHRS](https://bio-rad.com/NuviaHRS) for more information and to request samples.

Visit [bio-rad.com/ProcessResins](https://bio-rad.com/ProcessResins) for more information about Bio-Rad's complete line of process chromatography supports.

For additional information and technical assistance, contact your local Bio-Rad office or email our process specialists at [process@bio-rad.com](mailto:process@bio-rad.com). In the U.S. and Canada, call 1-800-4BIORAD (1-800-424-6723).

BIO-RAD, FORESIGHT, and MACRO-PREP are trademarks of Bio-Rad Laboratories, Inc. in certain jurisdictions. All trademarks used herein are the property of their respective owner.  
© 2024 Bio-Rad Laboratories, Inc.



Bio-Rad  
Laboratories, Inc.

Life Science  
Group

Website [bio-rad.com](https://bio-rad.com) USA 1 800 424 6723 Australia 61 2 9914 2800 Austria 00 800 00 24 67 23 Belgium 00 800 00 24 67 23 Brazil 4003 0399  
Canada 1 905 364 3435 China 86 21 6169 8500 Czech Republic 00 800 00 24 67 23 Denmark 00 800 00 24 67 23 Finland 00 800 00 24 67 23  
France 00 800 00 24 67 23 Germany 00 800 00 24 67 23 Hong Kong 852 2789 3300 Hungary 00 800 00 24 67 23 India 91 124 4029300 Israel 0 3 9636050  
Italy 00 800 00 24 67 23 Japan 81 3 6361 7000 Korea 82 080 007 7373 Luxembourg 00 800 00 24 67 23 Mexico 52 555 488 7670  
The Netherlands 00 800 00 24 67 23 New Zealand 64 9 415 2280 Norway 00 800 00 24 67 23 Poland 00 800 00 24 67 23 Portugal 00 800 00 24 67 23  
Russian Federation 00 800 00 24 67 23 Singapore 65 6415 3188 South Africa 00 800 00 24 67 23 Spain 00 800 00 24 67 23 Sweden 00 800 00 24 67 23  
Switzerland 00 800 00 24 67 23 Taiwan 886 2 2578 7189 Thailand 66 2 651 8311 United Arab Emirates 36 1 459 6150 United Kingdom 00 800 00 24 67 23

