

## 01.09.2020 Kit components

Product code	Description	
1762510 1762500 1762600	ProteOn™ HTG & HTE Reagent Kit	
Components:		
10021442	ProteOn <sup>TM</sup> HTG Chip Activation Solution, 10mM Nickel Sulfate pH 6.0	
10021443	ProteOn <sup>TM</sup> HTG Chip Regeneration Solution, 300 mM EDTA, pH 8.5	





Printing date 01.09.2020 Revision: 11.08.2020

Hazardous according to criteria of Safe Work Australia.

#### 1 Identification

- · Product identifier
- · Trade name: ProteOn™ HTG Chip Activation Solution, 10mM Nickel Sulfate pH 6.0
- · Article number: 10021442
- · Relevant identified uses of the substance or mixture and uses advised against

No further relevant information available.

- · Application of the substance / the mixture Laboratory chemicals
- · Details of the supplier of the safety data sheet
- · Manufacturer/Supplier:

Bio-Rad Laboratories Pty. Ltd Level 5, 446 Victoria Road Gladesville NSW 2111

Australia

Tel 61-2-9914-2800 (1 800 224 354) Fax 61-2-9914-2888

- · Information department: sales.australia@bio-rad.com
- · Emergency telephone number: 61-2-9037-2994

## 2 Hazard(s) Identification

· Classification of the substance or mixture

Skin Sens. 1 H317 May cause an allergic skin reaction.
Carc. 1A H350i May cause cancer by inhalation.

STOT RE 2 H373 May cause damage to organs through prolonged or repeated exposure.

Aquatic Acute 3 H402 Harmful to aquatic life.

Aquatic Chronic 3 H412 Harmful to aquatic life with long lasting effects.

- · Label elements
- · GHS label elements The product is classified and labelled according to the Globally Harmonised System (GHS).
- · Hazard pictograms





GHS07 GHS08

- · Signal word Danger
- · Hazard-determining components of labelling:

nickel sulphate

· Hazard statements

H317 May cause an allergic skin reaction.

H350i May cause cancer by inhalation.

H373 May cause damage to organs through prolonged or repeated exposure.

H412 Harmful to aquatic life with long lasting effects.

(Contd. on page 2)





Printing date 01.09.2020 Revision: 11.08.2020

Trade name: ProteOn™ HTG Chip Activation Solution, 10mM Nickel Sulfate pH 6.0

(Contd. of page 1)

· Precautionary statements

Do not breathe dust/fume/gas/mist/vapours/spray. P260

P273 Avoid release to the environment.

Wear protective gloves/protective clothing/eye protection/face protection. P280

P308+P313 IF exposed or concerned: Get medical advice/attention.

P405 Store locked up.

P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

- · Other hazards
- · Results of PBT and vPvB assessment
- · **PBT**: Not applicable. · vPvB: Not applicable.

## 3 Composition and Information on Ingredients

- · Chemical characterisation: Mixtures
- · Description: Mixture of the substances listed below with nonhazardous additions.
- Dangerous components:

7786-81-4 nickel sulphate 0.25-<0.3 Resp. Sens. 1, H334; Muta. 2, H341; Carc. 1A, H350i; Repr. 1B, H360; STOT RE 1, H372	
Resn Sens 1 H334: Muta 2 H341: Carc 14 H350i: Renr 18 H360: STOT RF 1 H372	5-<0.3%
Aguatic Chronic 1, H410	
Acute Tox. 4, H302; Acute Tox. 4, H332; Skin Irrit. 2, H315; Skin Sens. 1, H317	

· Additional information For the wording of the listed hazard phrases refer to section 16.

### 4 First Aid Measures

- · Description of first aid measures
- · General information Immediately remove any clothing soiled by the product.
- · After inhalation

Supply fresh air and to be sure call for a doctor.

In case of unconsciousness place patient stably in side position for transportation.

· After skin contact

Immediately wash with water and soap and rinse thoroughly.

Immediately rinse with water.

- · After eye contact Rinse opened eye for several minutes under running water.
- · After swallowing Rinse mouth with water. Seek medical attention and appropriate follow-up.
- · Information for doctor
- · Most important symptoms and effects, both acute and delayed No further relevant information available.
- · Indication of any immediate medical attention and special treatment needed No further relevant information available.

### 5 Fire Fighting Measures

- Extinguishing media
- · Suitable extinguishing agents Use fire fighting measures that suit the environment.
- · Special hazards arising from the substance or mixture During heating or in case of fire poisonous gases are produced. (Contd. on page 3)



Page 3/8

## Safety Data Sheet according to WHS Regulations

Printing date 01.09.2020 Revision: 11.08.2020

Trade name: ProteOn™ HTG Chip Activation Solution, 10mM Nickel Sulfate pH 6.0

(Contd. of page 2)

- · Advice for firefighters
- · Protective equipment:

Mount respiratory protective device.

Wear self-contained respiratory protective device.

#### 6 Accidental Release Measures

· Personal precautions, protective equipment and emergency procedures

Mount respiratory protective device.

Wear protective clothing.

Environmental precautions:

Do not allow product to reach sewage system or any water course.

Inform respective authorities in case of seepage into water course or sewage system.

Dilute with plenty of water.

Do not allow to enter sewers/ surface or ground water.

· Methods and material for containment and cleaning up:

Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).

Dispose contaminated material as waste according to item 13.

Ensure adequate ventilation.

Reference to other sections

See Section 7 for information on safe handling

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

### 7 Handling and Storage

- Handling
- · Precautions for safe handling

Ensure good ventilation/exhaustion at the workplace.

Open and handle receptacle with care.

Prevent formation of aerosols.

- · Information about protection against explosions and fires: Keep respiratory protective device available.
- · Conditions for safe storage, including any incompatibilities
- Storage
- · Requirements to be met by storerooms and receptacles: According to product specification
- · Information about storage in one common storage facility: Not required.
- · Further information about storage conditions: Keep receptacle tightly sealed.
- · Specific end use(s) No further relevant information available.

### 8 Exposure controls and personal protection

· Additional information about design of technical systems: No further data; see item 7.

(Contd. on page 4)





Printing date 01.09.2020 Revision: 11.08.2020

Trade name: ProteOn™ HTG Chip Activation Solution, 10mM Nickel Sulfate pH 6.0

(Contd. of page 3)

- · Control parameters
- · Components with limit values that require monitoring at the workplace:

#### 7786-81-4 nickel sulphate

WES Long-term value: 0.1 mg/m3

Sen; Note (g); as Ni

- · Additional information: The lists that were valid during the creation were used as basis.
- · Exposure controls
- · Personal protective equipment
- · General protective and hygienic measures

Keep away from foodstuffs, beverages and feed.

Immediately remove all soiled and contaminated clothing

Wash hands before breaks and at the end of work.

Store protective clothing separately.

Avoid contact with the skin.

· Breathing equipment:

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use respiratory protective device that is independent of circulating air.

Not necessary if room is well-ventilated.

- · Protection of hands: Protective gloves.
- · Material of gloves Synthetic gloves
- · Penetration time of glove material

The exact break trough time has to be found out by the manufacturer of the protective gloves and has to be observed.

· Eye protection:

Safety glasses

Tightly sealed goggles.

0.51					
u Phi	veical	andl	Chemica	il Pro	nartias
<i> </i>		and			per aco

- · Information on basic physical and chemical properties
- · General Information
- · Appearance:

Form: Fluid
Colour: Light green
Odour: Odourless
Odour threshold: Not determined.

· pH-value at 20 °C:

· Change in condition

Melting point/freezing point: 0 °C

Initial boiling point and boiling range: undetermined

6

• Flash point: Not applicable Not determined.

· Flammability (solid, gaseous) Not applicable.

(Contd. on page 5)





Printing date 01.09.2020 Revision: 11.08.2020

Trade name: ProteOn™ HTG Chip Activation Solution, 10mM Nickel Sulfate pH 6.0

	(Contd. of page 4)
· Decomposition temperature:	Not determined.
· Self igniting:	Product is not selfigniting.
· Explosive properties:	Not determined.
· Explosion limits: Lower: Upper:	Not determined. Not determined.
· Vapour pressure at 20 °C:	23 hPa
<ul> <li>Density at 20 °C:</li> <li>Relative density</li> <li>Vapour density</li> <li>Evaporation rate</li> </ul>	1 g/cm³ Not determined. Not determined. Not determined.
· Solubility in / Miscibility with Water:	Fully miscible
· Partition coefficient: n-octanol/water:	Not determined.
· Viscosity: dynamic: kinematic:	Not determined. Not determined.
· Solvent content: Water: · Other information	99.5 % No further relevant information available.

## 10 Stability and Reactivity

- · Reactivity No further relevant information available.
- · Chemical stability
- · Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
- · Possibility of hazardous reactions No dangerous reactions known
- · Conditions to avoid No further relevant information available.
- · Incompatible materials: No further relevant information available.
- · Hazardous decomposition products: No dangerous decomposition products known

## 11 Toxicological Information

- · Information on toxicological effects
- · Acute toxicity
- Primary irritant effect:
- · Skin corrosion/irritation No irritant effect.
- · Serious eye damage/irritation No irritating effect.
- · Respiratory or skin sensitisation Sensitisation possible through skin contact.

(Contd. on page 6)



Page 6/8

## Safety Data Sheet according to WHS Regulations

Printing date 01.09.2020 Revision: 11.08.2020

Trade name: ProteOn™ HTG Chip Activation Solution, 10mM Nickel Sulfate pH 6.0

(Contd. of page 5)

· Additional toxicological information:

The product shows the following dangers according to the calculation method of the General EU Classification Guidelines for Preparations as issued in the latest version:

Irritant

· CMR effects (carcinogenity, mutagenicity and toxicity for reproduction)
Carc. 1A

#### 12 Ecological Information

- · Toxicity
- · Aquatic toxicity: No further relevant information available.
- · Persistence and degradability No further relevant information available.
- Behaviour in environmental systems:
- · Bioaccumulative potential No further relevant information available.
- · Mobility in soil No further relevant information available.
- · Ecotoxical effects:
- · Remark: Harmful to fish
- · Additional ecological information:
- · General notes:

Water hazard class 2 (German Regulation) (Self-assessment): hazardous for water.

Do not allow product to reach ground water, water course or sewage system.

Danger to drinking water if even small quantities leak into the ground.

Harmful to aquatic organisms

- · Results of PBT and vPvB assessment
- · PBT: Not applicable.
- · vPvB: Not applicable.
- · Other adverse effects No further relevant information available.

#### 13 Disposal considerations

- · Waste treatment methods
- · Recommendation

Must not be disposed of together with household garbage. Do not allow product to reach sewage system. Dispose of waste in accordance to applicable national, regional, or local regulations.

- · Uncleaned packagings:
- · Recommendation: Disposal must be made according to official regulations.
- · Recommended cleansing agent: Water, if necessary with cleansing agents.

#### 14 Transport information

- · UN-Number
- · ADG, ADN, IMDG, IATA

Not Requated

(Contd. on page 7)





Printing date 01.09.2020 Revision: 11.08.2020

Trade name: ProteOn™ HTG Chip Activation Solution, 10mM Nickel Sulfate pH 6.0

	(Contd. of pag	ge 6
· UN proper shipping name · ADG	Not Reguated Not Reguated	
· ADN, IMDG, IATA	Not Reguated	
· Transport hazard class(es)		
· ADG, ADN, IMDG, IATA · Class	Not Reguated	
· Packing group · ADG, IMDG, IATA	Not Reguated	
· Environmental hazards: · Marine pollutant:	No	
· Special precautions for user	Not applicable.	
· Transport in bulk according to Annex II of IBC Code	Marpol and the Not applicable.	
· UN "Model Regulation":	Not Reguated	

### 15 Regulatory information

- · Safety, health and environmental regulations/legislation specific for the substance or mixture
- · Australian Inventory of Chemical Substances

All ingredients are listed.

· Standard for the Uniform Scheduling of Medicines and Poisons

7786-81-4 nickel sulphate

S6

· Australia: Priority Existing Chemicals

None of the ingredients is listed.

- · Directive 2012/18/EU
- · Named dangerous substances ANNEX I None of the ingredients is listed.
- · National regulations
- · Information about limitation of use:

Workers are not allowed to be exposed to the hazardous carcinogenic materials contained in this preparation. Exceptions can be made by the authorities in certain cases.

· Technical instructions (air):

Class	Share in %
11	0.1-1

- · Water hazard class: Water hazard class 2 (Self-assessment): hazardous for water.
- · Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

- AU





Printing date 01.09.2020 Revision: 11.08.2020

Trade name: ProteOn™ HTG Chip Activation Solution, 10mM Nickel Sulfate pH 6.0

(Contd. of page 7)

### 16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

#### · Relevant phrases

H302 Harmful if swallowed.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H332 Harmful if inhaled.

H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.

H341 Suspected of causing genetic defects.

H350i May cause cancer by inhalation.

H360 May damage fertility or the unborn child.

H372 Causes damage to organs through prolonged or repeated exposure.

H410 Very toxic to aquatic life with long lasting effects.

#### Department issuing SDS: Environmental Health and Safety.

#### · Contact:

Life Science Group, Environmental Health and Safety, 2000 Alfred Nobel Drive, Hercules, California, 94547: 1(510) 741-1000

Diagnostic Group, Environmental Health and Safety, 4000 Alfred Nobel Drive, Hercules, California, 94547: 1(510) 724-7000

#### · Abbreviations and acronyms:

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

PBT: Persistent, Bioaccumulative and Toxic

vPvB. very Persistent and very Bioaccumulative

Acute Tox. 4: Acute toxicity - oral - Category 4

Skin Irrit. 2: Skin corrosion/irritation - Category 2

Resp. Sens. 1: Respiratory sensitisation - Category 1

Skin Sens. 1: Skin sensitisation - Category 1

Muta. 2: Germ cell mutagenicity – Category 2

Carc. 1A: Carcinogenicity – Category 1Ai

Repr. 1B: Reproductive toxicity - Category 1B

STOT RE 1: Specific target organ toxicity (repeated exposure) - Category 1

STOT RE 2: Specific target organ toxicity (repeated exposure) – Category 2

Aquatic Acute 3: Hazardous to the aquatic environment - acute aquatic hazard – Category 3
Aquatic Chronic 1: Hazardous to the aquatic environment - long-term aquatic hazard – Category 1

Aquatic Chronic 3: Hazardous to the aquatic environment - long-term aquatic hazard – Category 3

· \* Data compared to the previous version altered.

AU





Printing date 01.09.2020 Revision: 15.08.2020

Hazardous according to criteria of Safe Work Australia.

#### 1 Identification

- · Product identifier
- · Trade name: ProteOn™ HTG Chip Regeneration Solution, 300 mM EDTA, pH 8.5
- · Article number: 10021443
- · Relevant identified uses of the substance or mixture and uses advised against

No further relevant information available.

- · Application of the substance / the mixture Laboratory chemicals
- · Details of the supplier of the safety data sheet
- · Manufacturer/Supplier:

Bio-Rad Laboratories Pty. Ltd Level 5, 446 Victoria Road Gladesville NSW 2111

Australia

Tel 61-2-9914-2800 (1 800 224 354) Fax 61-2-9914-2888

- · Information department: sales.australia @bio-rad.com
- · Emergency telephone number: 61-2-9037-2994

## 2 Hazard(s) Identification

· Classification of the substance or mixture

Eye Irrit. 2A H319 Causes serious eye irritation.

- · Label elements
- · GHS label elements The product is classified and labelled according to the Globally Harmonised System (GHS).
- · Hazard pictograms



- · Signal word Warning
- · Hazard statements

H319 Causes serious eye irritation.

· Precautionary statements

P264 Wash thoroughly after handling.
P280 Wear eye protection / face protection.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy

to do. Continue rinsing.

P337+P313 If eye irritation persists: Get medical advice/attention.

- · Other hazards
- · Results of PBT and vPvB assessment
- · PBT: Not applicable.

(Contd. on page 2)



Page 2/7

## Safety Data Sheet according to WHS Regulations

Printing date 01.09.2020 Revision: 15.08.2020

Trade name: ProteOn™ HTG Chip Regeneration Solution, 300 mM EDTA, pH 8.5

(Contd. of page 1)

· **vPvB:** Not applicable.

## 3 Composition and Information on Ingredients

- · Chemical characterisation: Mixtures
- · Description: Mixture of the substances listed below with nonhazardous additions.
- · Dangerous components:

60-00-4 ethylenediaminetetraacetic acid

Eye Irrit. 2A, H319 Acute Tox. 5, H303 10-20%

· Additional information For the wording of the listed hazard phrases refer to section 16.

### 4 First Aid Measures

- · Description of first aid measures
- · General information No special measures required.
- · After inhalation Supply fresh air; consult doctor in case of complaints.
- · After skin contact Generally the product does not irritate the skin.
- · After eye contact Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.
- · After swallowing Rinse mouth with water. Seek medical attention and appropriate follow-up.
- · Information for doctor
- · Most important symptoms and effects, both acute and delayed No further relevant information available.
- · Indication of any immediate medical attention and special treatment needed No further relevant information available.

### 5 Fire Fighting Measures

- · Extinguishing media
- · Suitable extinguishing agents Use fire fighting measures that suit the environment.
- · Special hazards arising from the substance or mixture No further relevant information available.
- · Advice for firefighters
- · Protective equipment: Wear self-contained respiratory protective device.

#### 6 Accidental Release Measures

- · Personal precautions, protective equipment and emergency procedures Wear protective clothing.
- · Environmental precautions:

Dilute with plenty of water.

Do not allow to enter sewers/ surface or ground water.

Methods and material for containment and cleaning up:

Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).

Reference to other sections

No dangerous substances are released.

See Section 7 for information on safe handling

See Section 8 for information on personal protection equipment.

(Contd. on page 3)



Page 3/7

## Safety Data Sheet according to WHS Regulations

Printing date 01.09.2020 Revision: 15.08.2020

Trade name: ProteOn™ HTG Chip Regeneration Solution, 300 mM EDTA, pH 8.5

(Contd. of page 2)

See Section 13 for disposal information.

### 7 Handling and Storage

- · Handling
- · Precautions for safe handling No special measures required.
- · Information about protection against explosions and fires: No special measures required.
- · Conditions for safe storage, including any incompatibilities
- · Storage
- · Requirements to be met by storerooms and receptacles: According to product specification
- Information about storage in one common storage facility: Not required.
- · Further information about storage conditions: Keep receptacle tightly sealed.
- · Specific end use(s) No further relevant information available.

## 8 Exposure controls and personal protection

- · Additional information about design of technical systems: No further data; see item 7.
- · Control parameters
- · Components with limit values that require monitoring at the workplace:

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

- · Additional information: The lists that were valid during the creation were used as basis.
- Exposure controls
- · Personal protective equipment
- · General protective and hygienic measures

Keep away from foodstuffs, beverages and feed.

Immediately remove all soiled and contaminated clothing

Wash hands before breaks and at the end of work.

Avoid contact with the eyes.

Avoid contact with the eyes and skin.

- · Breathing equipment: Not necessary if room is well-ventilated.
- · Protection of hands: Protective gloves.
- · Material of gloves Synthetic gloves
- Penetration time of glove material

The exact break trough time has to be found out by the manufacturer of the protective gloves and has to be observed.

· Eye protection:

Safety glasses

Tightly sealed goggles.

(Contd. on page 4)





Printing date 01.09.2020 Revision: 15.08.2020

Trade name: ProteOn™ HTG Chip Regeneration Solution, 300 mM EDTA, pH 8.5

(Contd. of page 3)

	(Contd. or pag
Physical and Chamical Proper	tion
Physical and Chemical Propert	iles
· Information on basic physical and ch	nemical properties
· General Information	
· Appearance:	
Form:	Liquid
Colour:	Colourless
· Odour:	Odourless
· Odour threshold:	Not determined.
· pH-value at 20 °C:	8.5
· Change in condition	
Melting point/freezing point:	undetermined
Initial boiling point and boiling rang	ge: undetermined
· Flash point:	Not applicable
	Not determined.
· Flammability (solid, gaseous)	Not applicable.
Decomposition temperature:	Not determined.
· Self igniting:	Product is not selfigniting.
· Explosive properties:	Not determined.
· Explosion limits:	
Lower:	Not determined.
Upper:	Not determined.
· Vapour pressure at 20 °C:	23 hPa
Density:	Not determined
· Relative density	Not determined.
· Vapour density	Not determined.
· Evaporation rate	Not determined.
· Solubility in / Miscibility with	
Water:	Fully miscible
Partition coefficient: n-octanol/water	: Not determined.
· Viscosity:	
dynamic:	Not determined.
kinematic:	Not determined.
· Solvent content:	
Water:	88.8 %
Other information	No further relevant information available.

## 10 Stability and Reactivity

· Reactivity No further relevant information available.

(Contd. on page 5)

Page 5/7

## Safety Data Sheet according to WHS Regulations

Printing date 01.09.2020 Revision: 15.08.2020

Trade name: ProteOn™ HTG Chip Regeneration Solution, 300 mM EDTA, pH 8.5

(Contd. of page 4)

- · Chemical stability
- · Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
- · Possibility of hazardous reactions No dangerous reactions known
- · Conditions to avoid No further relevant information available.
- · Incompatible materials: No further relevant information available.
- · Hazardous decomposition products: No dangerous decomposition products known

### 11 Toxicological Information

- · Information on toxicological effects
- · Acute toxicity
- · LD/LC50 values that are relevant for classification:

60-00-4 ethylenediaminetetraacetic acid

Oral LD50 4,500 mg/kg (rat)

- · Primary irritant effect:
- · Skin corrosion/irritation Irritant to skin and mucous membranes.
- · Serious eye damage/irritation Irritating effect.
- · Respiratory or skin sensitisation No sensitising effects known.
- Additional toxicological information:

The product shows the following dangers according to the calculation method of the General EU Classification Guidelines for Preparations as issued in the latest version:

Irritant

### 12 Ecological Information

- Toxicity
- · Aquatic toxicity: No further relevant information available.
- · Persistence and degradability No further relevant information available.
- · Behaviour in environmental systems:
- · Bioaccumulative potential No further relevant information available.
- · Mobility in soil No further relevant information available.
- · Additional ecological information:
- · General notes:

Water hazard class 2 (German Regulation) (Self-assessment): hazardous for water.

Do not allow product to reach ground water, water course or sewage system.

Danger to drinking water if even small quantities leak into the ground.

- · Results of PBT and vPvB assessment
- · PBT: Not applicable.
- · vPvB: Not applicable.
- · Other adverse effects No further relevant information available.

— AU

(Contd. on page 6)





Printing date 01.09.2020 Revision: 15.08.2020

Trade name: ProteOn™ HTG Chip Regeneration Solution, 300 mM EDTA, pH 8.5

(Contd. of page 5)

## 13 Disposal considerations

- · Waste treatment methods
- · Recommendation

Must not be disposed of together with household garbage. Do not allow product to reach sewage system. Dispose of waste in accordance to applicable national, regional, or local regulations.

- Uncleaned packagings:
- · Recommendation: Disposal must be made according to official regulations.
- · Recommended cleansing agent: Water, if necessary with cleansing agents.

UN-Number		
ADG, ADN, IMDG, IATA	Not Reguated	
UN proper shipping name		
ADG	Not Reguated	
	Not Reguated	
ADN, IMDG, IATA	Not Reguated	
Transport hazard class(es)		
ADG, ADN, IMDG, IATA		
Class	Not Reguated	
Packing group		
ADG, IMDG, İATA	Not Reguated	
Environmental hazards:		
Marine pollutant:	No	
Special precautions for user	Not applicable.	
Transport in bulk according to Annex II of N	Marpol and the	
BC Code	Not applicable.	

## 15 Regulatory information

- · Safety, health and environmental regulations/legislation specific for the substance or mixture
- · Australian Inventory of Chemical Substances

All ingredients are listed.

· Standard for the Uniform Scheduling of Medicines and Poisons

60-00-4 ethylenediaminetetraacetic acid

S4

· Australia: Priority Existing Chemicals

None of the ingredients is listed.

(Contd. on page 7)





Printing date 01.09.2020 Revision: 15.08.2020

Trade name: ProteOn™ HTG Chip Regeneration Solution, 300 mM EDTA, pH 8.5

(Contd. of page 6)

- · Directive 2012/18/EU
- · Named dangerous substances ANNEX I None of the ingredients is listed.
- · National regulations
- · Water hazard class: Water hazard class 2 (Self-assessment): hazardous for water.
- · Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

### 16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· Relevant phrases

H303 May be harmful if swallowed.

H319 Causes serious eye irritation.

- · Department issuing SDS: Environmental Health and Safety.
- · Contact:

Life Science Group, Environmental Health and Safety, 2000 Alfred Nobel Drive, Hercules, California, 94547: 1(510) 741-1000

Diagnostic Group, Environmental Health and Safety, 4000 Alfred Nobel Drive, Hercules, California, 94547: 1(510) 724-7000

Abbreviations and acronyms:

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic vPvB: very Persistent and very Bioaccumulative Acute Tox. 5: Acute toxicity - oral – Category 5

Eye Irrit. 2A: Serious eye damage/eye irritation - Category 2A

· \* Data compared to the previous version altered.

AU