

# **SAFETY DATA SHEET (SDS)**

### SECTION 1: IDENTIFICATION OF PRODUCT (MIXTURE) AND SUPPLIER

Product Name: GS HIV Combo Ag/Ab EIA

**Product** 12009912 (192 tests)

**Number:** Catalog number(s) for replacement, separately purchased components that can be obtained for use with this

kit, and which are covered by this SDS include: 25260, 25261, 26181, 26182, 26220, 26221, 26222,

26223, 26224, 26225, 12009798, 12009799 and 26228 (refer to Section 2).

**Intended Use:** The GS HIV Combo Ag/Ab EIA is an enzyme immunoassay kit for the simultaneous qualitative detection

of Human Immunodeficiency Virus (HIV) p24 antigen and antibodies to HIV Type 1 (HIV-1 Groups M and O) and HIV Type 2 (HIV-2) in human serum or plasma. This kit is intended as an aid in the diagnosis of HIV-1 and/or HIV-2 infection, including acute or primary HIV-1 infection. The assay may also be used as an aid in the diagnosis of HIV-1 and/or HIV-2 infection in pediatric subjects (i.e, children as young as 2 years of age). The GS HIV Combo Ag/Ab EIA is intended for manual use and for use with the Bio-Rad

**EVOLIS** Automated Microplate System.

Results from the GS HIV Combo Ag/Ab EIA cannot be used to distinguish between the presence of HIV-1

p24 antigen, HIV-1 antibody, or HIV-2 antibody in a sample.

The GS HIV Combo Ag/Ab EIA is not intended for routine use in screening blood or plasma donors, as the effectiveness of this test for use in the screening of these donors has not been established. However, in urgent situations where traditional licensed blood donor screening tests are unavailable or their use is

impractical, this assay can be used as a blood donor screening assay.

Manufactured by: Bio-Rad Laboratories, Inc.

**Address:** 6565 185th Avenue NE

Redmond, WA 98052-5039, USA

Website: www.bio-rad.com

**Phone Number:** 1-800-2-BIORAD (1-800-224-6723); or (425) 881-8300 (daytime PT)

SDS e-mail <u>ro-sds@bio-rad.com</u>

contact:

**Technical** Bio-Rad provides a toll free line for technical assistance, available 24 hours a day, 7 days a week. In the

**Information** United States of America and Puerto Rico, call toll free 1-800-2-BIORAD (1-800-224-6723). Outside the **Contacts:** U.S.A., please contact your regional Bio-Rad office for assistance. *Refer to section 16 for non-US local* 

Bio-Rad agent contact information.

**Authorized** FRANCE: Bio-Rad

**Representative** 3 boulevard Raymond Poincaré in the European 92430 Marnes-la-Coquette

**Community:** Phone: +33 (0) 1 47 95 60 00 / Fax: +33 (0) 1 47 41 91 33

[fds-msds.fr@bio-rad.com]

**Emergency Phone** 

Number:

This SDS is listed with CHEMTREC 1-800-424-9300 or 1-703-527-3887 (US/CA) / +1-703-741-5970

(international - can be called collect). Use only in the event of a CHEMICAL EMERGENCY involving a

SPILL, LEAK, FIRE, EXPLOSION or ACCIDENT with this product.

### SECTION 2: HAZARDS IDENTIFICATION -- HAZARDOUS COMPONENTS

This test kit should be handled only by qualified personnel trained in laboratory procedures and familiar with their potential hazards. Specific warnings are given in the instructions for use. The absence of a specific warning should not be interpreted as an indication of safety. The following information is furnished for those product hazardous constituents that require regulatory control or disclosure at the concentration found in the product. Refer to Section 16 for the full text of any solely abbreviated or coded hazard statements provided below and for the key / legend to abbreviations and acronyms



Component *	Content
R1. Microwell Strip Plate (2 plates)	<ul> <li>Microwell plate coated with monoclonal antibodies to HIV p24 (mouse) and HIV-1 and HIV-2 antigens.</li> <li>Potential residue of ProClin 300 used as a production preservative (aspirated prior to drying strips).</li> <li>Tabs are labeled "JJ"</li> <li>Contains sealed pelletized desiccant packet: There are no health hazards associated with intact desiccant container; however, health hazards could result from dusts generated if the packet is cut, split or otherwise compromised and is crushed.</li> </ul>
R2. Wash Solution Concentrate (30X) 1 bottle (120 mL) Catalog No. 25261	- Sodium chloride (NaCl) [CAS# 7647-14-5, EC No 231-598-3] aqueous solution with < 2% Tween 20 (C <sub>58</sub> H <sub>114</sub> O <sub>26</sub> ) [CAS# 9005-64-5, EC No 585-580-06-X] (clear liquid). Not subject to GHS, US HCS, EC CLP and analogous global GHS-based regulatory requirements in this product mixture and concentration.
C0. Negative Control 1 vial (1.5 mL) Catalog No. 26220	<ul> <li>Normal human serum/plasma that is non-reactive for HBsAg, HIV-1 Ag and antibodies to HIV and HCV (amber / yellow liquid).</li> <li>Preserved with 0.16% ProClin 950 containing 0.016% active ingredient: 9.5-9.9% 2-methyl-4-isothiazolin-3-one (C<sub>4</sub>H<sub>5</sub>NOS); CAS# 2682-20-4, EC No 220-239-6. Not subject to GHS, US HCS, EC CLP and analogous global GHS-based regulatory requirements in this product mixture and concentration.</li> <li>Preserved with 0.005% gentamicin sulfate, CAS# 1405-41-0, EC No 215-778-9. Not subject to GHS, US HCS, EC CLP and analogous global GHS-based regulatory requirements in this product mixture and concentration.</li> </ul>
C1. HIV-1 Ab Positive Control 1 vial (1.5 mL) Catalog No. 26221	<ul> <li>Heat-treated human serum/plasma containing HIV-1 antibody, non-reactive for HBsAg and antibody to HCV (yellow or green aqueous solution).</li> <li>Preserved with 0.16% ProClin 950 containing 0.016% active ingredient: 9.5-9.9% 2-methyl-4-isothiazolin-3-one (C<sub>4</sub>H<sub>5</sub>NOS); CAS# 2682-20-4, EC No 220-239-6. Not subject to GHS, US HCS, EC CLP and analogous global GHS-based regulatory requirements in this product mixture and concentration.</li> <li>Preserved with 0.005% gentamicin sulfate, CAS# 1405-41-0, EC No 215-778-9. Not subject to GHS, US HCS, EC CLP and analogous global GHS-based regulatory requirements in this product mixture and concentration.</li> </ul>
C2. HIV-2/O Ab Positive Control 1 vial (1.5 mL) Catalog No. 26222	<ul> <li>Murine monoclonal antibody to HIV-2 and rabbit antibody to HIV-1 Group O diluted in normal human serum/plasma; non-reactive for HBsAg and antibodies to HCV (yellow or green aqueous solution).</li> <li>Preserved with 0.16% ProClin 950 containing 0.016% active ingredient: 9.5-9.9% 2-methyl-4-isothiazolin-3-one (C<sub>4</sub>H<sub>5</sub>NOS); CAS# 2682-20-4, EC No 220-239-6. Not subject to GHS, US HCS, EC CLP and analogous global GHS-based regulatory requirements in this product mixture and concentration.</li> <li>Preserved with 0.005% gentamicin sulfate, CAS# 1405-41-0, EC No 215-778-9. Not subject to GHS, US HCS, EC CLP and analogous global GHS-based regulatory requirements in this product mixture and concentration.</li> </ul>
C3. HIV Ag Positive Control 1 vial (1.5 mL) Catalog No. 26223  WARNING	<ul> <li>- Purified HIV-1 viral lysate antigen inactivated with heat and a chaotropic agent, in synthetic diluent (bovine) (blue-green aqueous liquid).</li> <li>- Preserved with 0.5% ProClin 300 (0.015% active ingredient), EC Index No 613-167-00-5 with CAS# 55965-84-9. GHS\US HCS\EC CLP Classification: WARNING; GHS07; H317, H412; P273, P280; P302 + P352, P333 + P313; P501.</li> <li>- Contains &lt; 0.1% Tartrazine [FD&amp;C Yellow #5, Acid Yellow 23; C₁6H9N₄O₂S₂•3Na]; CAS# 1934-21-0; EC 217-699-5. Not subject to GHS, US HCS, EC CLP and analogous global GHS-based regulatory requirements in this product mixture and concentration.</li> <li>- Preserved with 0.005% gentamicin sulfate, CAS# 1405-41-0, EC No 215-778-9. Not subject to GHS, US HCS, EC CLP and analogous global GHS-based regulatory requirements in this product mixture and concentration.</li> </ul>
C4. Cutoff Calibrator 3 vials (1.7 mL) Catalog No. 26224	<ul> <li>Normal human serum/plasma that is non-reactive for HBsAg, HIV-1 Ag and antibodies to HIV and HCV (amber / yellow liquid).</li> <li>Preserved with 0.16% ProClin 950 containing 0.016% active ingredient: 9.5-9.9% 2-methyl-4-isothiazolin-3-one (C<sub>4</sub>H<sub>5</sub>NOS); CAS# 2682-20-4, EC No 220-239-6. Not subject to GHS, US HCS, EC CLP and analogous global GHS-based regulatory requirements in this product mixture and concentration.</li> <li>Preserved with 0.005% gentamicin sulfate, CAS# 1405-41-0, EC No 215-778-9. Not subject to GHS, US HCS, EC CLP and analogous global GHS-based regulatory requirements in this product mixture and concentration.</li> </ul>



Component *	Content
R6. Conjugate 1 1 bottle (10 mL)	- Biotinylated polyclonal antibodies (sheep) to HIV p24 Ag in a buffer with protein stabilizers (bovine and ovine) and yellow-green dye.
Catalog No. 26225	- Preserved with <b>0.5% ProClin 300</b> (0.015% active ingredient), EC Index No 613-167-00-5 with CAS# 55965-84-9. GHS \ US HCS \ EC CLP Classification: WARNING; GHS07; H317, H412; P273, P280; P302 + P352, P333 + P313; P501.
WARNING	- Diluted < 0.1% Bromocresol Purple [C <sub>21</sub> H <sub>16</sub> Br <sub>2</sub> O <sub>5</sub> S], CAS# 115-40-2, EC No 204-087-8. Not subject to GHS, US HCS, EC CLP and analogous global GHS-based regulatory requirements in this product mixture and concentration.
	- Preserved with <b>0.005% gentamicin sulfate</b> , CAS# 1405-41-0, EC No 215-778-9. Not subject to GHS, US HCS, EC CLP and analogous global GHS-based regulatory requirements in this product mixture and concentration.
R7a. Conjugate 2 Concentrate (11X)	- Peroxidase-labeled Streptavidin and peroxidase-labeled HIV-1 and HIV-2 antigens in a buffer with protein stabilizers (bovine and caprine) and green dye.
1 vial (12 mL)  Catalog No. 12009798	- Preserved with <b>0.5% ProClin 300</b> (0.015% active ingredient), EC Index No 613-167-00-5 with CAS# 55965-84-9. GHS \ US HCS \ EC CLP Classification: WARNING; GHS07; H317, H412; P273, P280; P302 + P352, P333 + P313; P501.
<b>①</b>	- Preserved with <b>0.005% gentamicin sulfate</b> , CAS# 1405-41-0, EC No 215-778-9. Not subject to GHS, US HCS, EC CLP and analogous global GHS-based regulatory requirements in this product mixture and concentration.
WARNING	
R7b. Conjugate 2	- Buffer (pH neutral) with protein stabilizers (bovine and ovine) and red dye (red aqueous liquid).
<b>Diluent,</b> 1 bottle (120 mL)  Catalog No. 12009799	- Preserved with <b>0.5% ProClin 300</b> (0.015% active ingredient), EC Index No 613-167-00-5 with CAS# 55965-84-9. GHS \ US HCS \ EC CLP Classification: WARNING; GHS07; H317, H412; P273, P280; P302 + P352, P333 + P313; P501.
<b>(</b>	- Diluted <b>0.1% Amaranth [FD&amp;C Red #2;</b> C <sub>20</sub> H <sub>11</sub> N <sub>2</sub> O <sub>10</sub> S <sub>3</sub> •3Na], CAS# 915-67-3, EC No 213-022-2. [IARC Group 3]. Not subject to GHS, US HCS, EC CLP and analogous global GHS-based regulatory requirements in this product mixture and concentration.
WARNING	- Preserved with <b>0.005% gentamicin sulfate</b> , CAS# 1405-41-0, EC No 215-778-9. Not subject to GHS, US HCS, EC CLP and analogous global GHS-based regulatory requirements in this product mixture and concentration.
<b>R8.</b> Substrate Buffer, 1 bottle (120 mL)	- Dilute citric acid (C <sub>6</sub> H <sub>8</sub> O <sub>7</sub> , CAS# 79-92-9, EC No 201-069-1) / sodium acetate buffer, (pH ~ 4.0, clear liquid) < 5% dimethylsulfoxide [DMSO - C <sub>2</sub> H <sub>6</sub> OS], CAS# 67-68-5, EC No 200-644-3.
Catalog No. 26181	- < 0.1% hydrogen peroxide [H <sub>2</sub> O <sub>2</sub> ], CAS# 7722-84-1, EC No 231-765-0.
	Not subject to GHS, US HCS, EC CLP and analogous global GHS-based regulatory requirements in this product mixture and concentration.
<b>R9.</b> Chromogen (11X) 1 bottle (12 mL)	- ≤ <b>0.25% 3,3',5,5' tetramethylbenzidine dihydrochloride</b> [TMB- C <sub>16</sub> H <sub>20</sub> N <sub>2</sub> •2HCl], CAS# 207738-08-7, EC No 264-769-6.
Catalog No. 26182	$- \le$ <b>0.04 N hydrochloric acid</b> [ $\sim$ 0.3% HCl, CAS# 7647-01-0, EC No 231-595-7] solution (pH $\sim$ 1.5, clear liquid).
	Not subject to US HCS, EC CLP and analogous global GHS-based regulatory requirements in this product mixture and concentration.
R10. Stopping Solution 1 bottle (120 mL) Catalog No. 25260	- 1N H <sub>2</sub> SO <sub>4</sub> (4.4% w/w Sulfuric acid), CAS# 7664-93-9, EC No 231-639-5 [pH ≤ 2, clear liquid]. Severely irritating to skin, corrosive to eyes. GHS \ US HCS \ EC CLP: DANGER! GHS05; H290, H314; P280; P301 + P330 + P331, P305 + P351 + P338; P501.
DANGER!	

<sup>\*</sup> Replacement, optional and separately purchased component Catalog numbers are provided in this column where available.

Markings according to the *United Nations* (UN) Globally Harmonized System (GHS), *United States* Hazard Communication Standard (US HCS), *European Commission* (EC) 2008/1272/EC (EC CLP) guidelines and analogous GHS-based global regulations:

This product has been conservatively classified and labeled in accordance with *United Nations (UN)* GHS, *United States* Hazard Communication Standard (US HCS), related *European Commission (EC)* 2008/1272/EC (EC CLP) guidelines and applicable analogous GHS-based global regulations. The following regulated hazardous chemical concentrations are found in product component(s):

Component R10 - 1N H<sub>2</sub>SO<sub>4</sub> [4.4% w/w Sulfuric acid], CAS# 7664-93-9, EC No 231-639-5 (pH ≤ 2); severely irritating to skin, corrosive to eyes. [This STOP solution has been evaluated with the CORROSITEX test method to determine its corrosive potential and classification. The results of this testing classified this STOP solution as Class: 8, Packing

group **II** (**UN2796**)]

Comprehensive GHS Based Classification: Skin Corrosive Category 1B Serious eye damage (Category 1)



Label(s):

Signal Word: DANGER!

**Label Hazard Statements:** 

H290 May be corrosive to metals.

H314 Causes severe skin burns and eye damage.

Precautionary Statements (statements for product intended use and as codified on the product label):

P280 Wear protective gloves / protective clothing / eye protection / face protection.

P301 + P330 + P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

P303 + P361 + P353 IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with

water/shower.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to

do. Continue rinsing.

P501 This material and its container must be disposed of as hazardous waste.

Supplemental Precautionary Statements (additional precautions to consider relative to specific customer use):

P308 + P313 If exposed or concerned: Get medical advice/ attention.

P405 Store locked up.

[Source: Raw Material vendor SDS, CCOHS databases and regulatory research]

Component C3, R6, R7a and R7b - 0.5% ProClin 300 [≤ 0.015%. active ingredients – reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one ( $C_4H_4$ ClNOS; CAS# 26172-55-4, EC No 247-500-7) and 2-methyl-2H -isothiazol-3-one ( $C_4H_5$ NOS; CAS# 2682-20-4, EC No 220-239-6) (3:1)], EC Index No 613-167-00-5 with CAS# 55965-84-9.

Comprehensive GHS Based Classification: Skin Sensitizer Category 1



Label(s):

Signal Word: WARNING

<u>Label Hazard Statements:</u>

H317 May cause an allergic skin reaction.

H412 Harmful to aquatic life with long lasting effects.

Precautionary Statements (statements for product intended use and as codified on the product label):

P273 Avoid release to the environment.

P280 Wear protective gloves / protective clothing / eye protection / face protection.

P302 + P352 IF ON SKIN: Wash with plenty of soap and water.

P333 + P313 If skin irritation or rash occurs: Get medical advice/ attention.

P501 Dispose of contents and container in accordance to local, regional, national and international regulations.

Supplemental Precautionary Statements (additional precautions to consider relative to specific customer use):

P261 Avoid breathing mist / vapors/vapours / spray.

P272 Contaminated work clothing should not be allowed out of the workplace.

#### SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS

The following information is furnished for those product hazardous constituents that require regulatory control or disclosure regardless of the concentration found in the product. Note that the information here is often based on data from the chemical raw material safety data sheet and literature (LD<sub>50</sub>, exposure limits, etc.). Chemical constituents that do not require regulatory disclosure are not generally included here. This product contains a significantly diluted concentration in an aqueous solution, thus the assessment below has not considered the dilution reduction effect on the hazard. That hazard communication information is provided in Section 2 above. Some components were tested at the concentration found in the kit. In that case, the assessment is provided for the chemical dilution tested and the tested concentration will be provided at the beginning of the *Chemical Ingredient Data/Information* box. Refer to section 16 for the full text of any *Comprehensive GHS-based Classification* statements coded below, for the list of sources utilized in the assessment and for the key / legend to abbreviations and acronyms.

#### **Chemical Ingredient Data / Information**

Chemical Ingredient: Sulfuric acid

Chemical concentrations found in this product: 1 N (< 5% H<sub>2</sub>SO<sub>4</sub> in water) in R10

Data for Concentrated / 100% chemical used in the product mixture (concentration tested):

CAS#: 7664-93-9 (Conc. sulfuric acid 100%)

LD<sub>50</sub> (oral-rat): 2,140 mg/kg

EC No: 231-639-5 (100%)

LC<sub>50</sub> (inhalation-rat): 510 mg/m<sup>3</sup>/2H

Index No: 016-020-00-8 (100%) LD<sub>50</sub> (skin-rabbit): NE

Registration No: 01-2119458838-20-XXXX LC<sub>50</sub> (96 hr-fish): Gambusia affinis (Mosquito fish) – 42 mg/l

RTECS#: WS5600000 (100%) pH value: 1.2 at 5 g/L

Chemical Formula: H<sub>2</sub>SO<sub>4</sub> (100%) Molecular weight: 98.08 g/mol (100%)

Synonyms/Trade Names: Acide sulfurique ;Acido solforico; BOV; Battery acid; Dihydrogen sulfate; Dipping acid; Electrolyte acid; Mattling acid; Oil of vitriol; Schwefelsaeureloesungen; Strong inorganic acid mists containing sulfuric acid; Sulfuric

Skin corrosion/irritation: Skin - rabbit - Extremely corrosive and destructive to tissue.

Serious eye damage/eye irritation: Eyes - rabbit - Severe eye irritation

#### Raw Material GHS / US HCS / EC CLP Classification (100%):

#### DANGER!

Skin Corr. Cat. 1A, Eye Damage Cat. 1, Aquatic Acute Cat. 3, Aquatic Chron. Cat. 3

H290, H314, H412

P264, P273, P280, P301 + P330 + P331, P303 + P361 + P353, P304 + P340,

P305 + P351 + P338, P310, P363, P405, P501





#### **Chemical Ingredient Data / Information**

#### Chemical Ingredient: <u>Dimethyl-sulfoxide [DMSO]</u>

CChemical concentrations found in this product:  $\leq 5\%$  in in R8, an aqueous solution

#### Data for Concentrated / 100% chemical used in the product mixture (concentration tested):

CAS#: 67-68-5 (100%) LD<sub>50</sub> (oral-rat): 14500 mg/kg

EC No: 200-644-3 (100%)

RTECS#: PV6210000 (100%)

LC<sub>50</sub> (inhalation-rat): 1600 mg/m<sup>3</sup> (4 hr)

LD<sub>50</sub> (skin-rabbit): >5000 mg/m<sup>3</sup>

Index No: NA (100%) LC<sub>50</sub> (96 hr-fish): NE

Chemical Formula: C<sub>2</sub>H<sub>6</sub>OS (100%) Molecular weight: 78.13 g/mol (100%)

Flash point: 188-192°F / 86.7-88.9°C (100%)

Flammable limits: LEL/LFL is 3.5%; UEL/UFL is 42% vv in air

Synonyms/Trade Names: Dimethyl sulfoxide; Dimethyl sulphoxide; Dimexide; Dipirartril-tropico; DMS-70; DMS-90; DMSO; Dolicur; Domoso; Dromisol; Durasorb; A 10846; Deltan; Demeso; Demasorb; Demavet; Demsodrox; Dermasorb; Gamasol 90; Hyadur; Infiltrina; M 176, Methane, sulfinylbis-; Methylsulfinylmethane; Somipront; SQ 9453, Topsym; NSC-763; Rimso-50; Sulfinylbis(methane); Syntexan

#### Raw Material GHS / US HCS / EC CLP Classification (100%):

No Pictogram required

#### WARNING

Fla. Liq. Cat. 4

H227

P210, P280, P370 + P378, P403 + P235, P501

[Source: Raw Material vendor SDS, CCOHS databases and regulatory research]

### Chemical Ingredient: Citric acid

Chemical concentrations found in this product: < 1.5% w/v in R8

### Data for Concentrated / 100% chemical used in the product mixture (concentration tested):

CAS#: 77-92-9 (100%) LD<sub>50</sub> (oral-rat): 5400 mg/kg EC No: 201-069-1 (100%) LC<sub>50</sub> (inhalation-rat): NE RTECS#: GE7350000 (100%) LD<sub>50</sub> (skin-rabbit): >2000 mg/kg

Index No: NA (100%)

Toxicity to fish mortality  $LC_{50}$  - Leuciscus idus melanotus - 440 mg/l - 48 h - Method: OECD Test Guideline 203 Chemical Formula:  $C_6H_8O_7$  (100%) pH value: 1.8 at ca.50 g/l at 25 °C (77 °F)

Molecular weight: 192.12 g/mol (100%)

Synonyms/Trade Names: Aciletten; Anhydrous citric acid; Citretten; Citro; 2-Hydroxy-1,2,3-propanetricarboxylic acid; beta-Hydroxytricarballylic acid; Kyselina citronova

Skin corrosion/irritation: Skin - rabbit - Mild skin irritation - OECD Test Guideline 404 Serious eye damage/eye irritation: Eyes - rabbit - Irritating to eyes. - OECD Test Guideline 405

Respiratory or skin sensitization: Prolonged or repeated exposure may cause allergic reactions in certain sensitive individuals.

#### Raw Material GHS / US HCS / EC CLP Classification (100%):

#### WARNING

Acute Tox. - skn Cat. 5, Skin Irrit. Cat. 1, Eye Irrit. Cat. 2A

H313, H316, H319

P264, P280, P305 + P351 + P338, P337 + P313



#### **Chemical Ingredient Data / Information**

#### **Chemical Ingredient:** <u>Tartrazine Solution</u>

Chemical concentrations found in this product: < 0.1% in C3, an aqueous solution

Data for Concentrated / 100% chemical used in the product mixture (concentration tested):

CAS#: 1934-21-0 (100%) LD $_{50}$  (oral-mouse): 12750 mg/kg EC No: 217-699-5 (100%) LC $_{50}$  (inhalation-rat): NE RTECS#: UQ6400000 (100%) LD $_{50}$  (skin-rabbit): NE Chemical Formula:  $C_{16}H_9N_4O_9S_2\bullet 3Na$  LC $_{50}$  (96 hr-fish): NE

Molecular weight: 534.36 g/mol (100%)

Synonyms/Trade Names: Trisodium 5-hydroxy-1-(4-sulphophenyl)-4-(4-sulphophenylazo)pyrazole-3-carboxylate; Acid Yellow 23; Yellow

no. 4; FD&C Yellow #5

#### Raw Material GHS / US HCS / EC CLP Classification (100%):

#### DANGER!

Resp. Sens. Cat. 1, Skin Sens. Cat. 1

H317, H334

P261, P272, P280, P285, P302 + P352, P304 + P341, P333 + P313, P342 + P311, P363, P501

[Source: Raw Material vendor SDS, CCOHS databases and regulatory research]



#### Chemical Ingredient: ProClin 300

Chemical concentrations found in this product: <u>0.5% (0.015% active ingredient) in C3, R6, R7a and R7b</u>

#### [Potential residue dried on plates in R1]

### Data for Concentrated / 100% chemical used in the product mixture (concentration tested):

Hazardous ingredient concentration in raw material:

60-100% Glycols;

1-5% Mixture (3:1) of 5-Chloro-2-methyl-4-isothiazolin-3-one (C<sub>4</sub>H<sub>5</sub>NOS; CAS# 2682-20-4, EC# 220-239-6)

and **2-Methyl-2H -isothiazol-3-one** (C<sub>4</sub>H<sub>4</sub>ClNOS; CAS# 26172-55-4, EC# 247-500-7)

CAS#: 55965-84-9, Index No: 613-167-00-5.

RTECS#: NE

LD<sub>50</sub> (oral-rat): 862 mg/kg LD<sub>50</sub> (skin-rabbit): 2,800 mg/kg LC<sub>50</sub> (inhalation-rat): NE LD<sub>50</sub> (skin-rabbit): NE pH value: 4.1 at 100 g/L

Flash Point: 244° F / 118° C (concentrated solution)

Synonyms/Trade Names: 5-Chloro-2-methyl-4-isothiazolin-3-one solution; Kathon 300; Isothiazolinone chloride solution

Skin corrosion/irritation - rabbit - Corrosive (concentrated solution)

Serious eye damage/eye irritation - rabbit - Corrosive to eyes (concentrated solution)

Respiratory or skin sensitization - May cause allergic skin reaction (concentrated solution)

OECD TG 429 (LLNA (2000a)):  $SI \ge 3$  for all concentrations; EC3 value of  $\le 2$  (0.003%) – Positive, Skin Sens. Cat. 1A

OECD TG 429 (LLNA (2000b)): SI  $\geq$  3 from 70ppm; EC3 value of  $\leq$  2 (0.007%) – inconclusive

GPMT (2000a): OECD TG 416, GLP – Negative GPMT (2000b): OECD TG 416, GLP – Positive

Buehler (1982): GLP 9/15 animals responded to an induction concentration of 0.01% - Positive, Skin Sens. Cat. 1A

.....

Open Epicutaneous Test (2001) - Positive, Skin Sens. Cat. 1

#### Raw Material GHS / US HCS / EC CLP Classification (100%):

#### DANGER!

Acute Tox. – oral Cat. 4, Skin Corr. Cat. 1B, Eye Damage Cat. 1, Skin. Sens. Cat.1, Aquatic Acute Cat. 1, Aquatic Chron. Cat. 1 H302, H314, H317, H410

P261, P264, P270, P272, P273, P280, P301 + P312 + P330, P301 + P330 + P331,

P303 + P361 + P353, P305 + P351 + P338, P310, P333 + P313, P363, P391, P405, P501





[Catalog # 12009912]

#### Chemical Ingredient Data / Information

#### Chemical Ingredient: Hydrochloric acid

Chemical concentrations found in this product: ≤ 0.04N (< 0.4% v/v HCl) in R9

#### Data for Concentrated / 100% chemical used in the product mixture (concentration tested):

CAS#: 7647-01-0 (concentrate solution) LD<sub>50</sub> (oral, rat): 700 mg/kg (unconfirmed) EC No: 231-595-7 (concentrate solution) LC<sub>50</sub> (inhalation-rat): 3124 ppm/1H Index No: 017-002-01-X (concentrate solution) LD<sub>50</sub> (skin-rabbit): > 5010 mg/kg (unconfirmed)

RTECS#: MW4025000 (concentrate solution) Fish LC<sub>50</sub> – Bluegill/Sunfish – 282 mg/l – 48 h

Chemical Formula: HCl (concentrate solution)

Molecular weight: 36.46 g/mol (concentrate solution) pH value: highly acidic (30-50% concentrated solution)

Synonyms/Trade Names: Acide chlorhydrique; Acido cloridrico; Anhydrous hydrochloric acid; Chlorowaterstof; Chlorohydric acid; Chlorowodor; Chlorwasserstoff; Hydrochloride; Hydrogen chloride; Hydrogen chloride; Muriatic acid; Spirits of salt

LC<sub>50</sub> (male rat): 1405 ppm (4-hour exposure; head-only); cited as 2810 ppm (1-hour exposure; head-only) (30-50% solution) LC<sub>50</sub> (male rat): 1562 ppm (4-hour exposure; whole-body); cited as 3124 ppm (1-hour exposure; whole-body) (30-50% solution)

Skin corrosion/irritation: Skin - rabbit - Causes burns. (30-50% solution)

Serious eye damage/eye irritation: Eyes - rabbit - Corrosive to eyes (30-50% solution) IARC: Group 3: Not classifiable as to its carcinogenicity to humans (30-50% solution)

### Raw Material GHS / US HCS / EC CLP Classification (100%):

#### DANGER!

Skin Corr. Cat. 1B, Eye Damage Cat. 1, STOT SE Cat. 3, Met. Corr. Cat. 1 H290, H314, H335

P234, P261, P264, P271, P280, P301 + P330 + P331, P303 + P361 + P353, P304 + P340,

P305 + P351 + P338, P310, P363, P390, P403 + P233, P405, P406, P501.

[Source: Raw Material vendor SDS, CCOHS databases and regulatory research]



#### Chemical Ingredient: 3,3',5,5'-Tetramethylbenzidine, Dihydrochloride

Chemical concentrations found in this product: < 0.3% w/v in R9, an aqueous solution

### Data for Concentrated / 100% chemical used in the product mixture (concentration tested):

CAS#: 207738-08-7 (54827-17-7 Free base) (100%) LD<sub>50</sub> (ipr-mouse): 135 mg/kg

EC No: 264-769-6 (100%) LD<sub>50</sub> (oral-rat): NE RTECS#: DV2300000 (100%) LC<sub>50</sub> (inhalation-rat): NE LD<sub>50</sub> (skin-rabbit): NE Chemical Formula: C<sub>16</sub>H<sub>20</sub>N<sub>2</sub>•2HCl (100%) Molecular weight: 313.27 g/mol (100%) LC<sub>50</sub> (96 hr-fish): NE

Synonyms/Trade Names: TMB

Raw Material GHS / US HCS / EC CLP Classification (100%): Not a dangerous substance according to GHS, US HCS, EC CLP and analogous global GHS-based regulatory requirements.



#### **Chemical Ingredient Data / Information**

#### Chemical Ingredient: ProClin 950

Chemical concentrations found in this product: ≤ 0.16% (< 0.015% active ingredient) in C0, C1, C2, and C4

#### Data for Concentrated / 100% chemical used in the product mixture (concentration tested):

Hazardous ingredient concentration in raw material: the concentrated preservative contains:

5-10% of **2-methyl-4-isothiazolin-3-one** (active ingredient).

CAS#: 2682-20-4 (active ingredient) LD50 (oral-rat): No data available EC No: 220-239-6 (active ingredient) LC50 (inhalation-rat): No data available RTECS#: NE LD50 (skin-rabbit): No data available

Chemical Formula: C<sub>4</sub>H<sub>5</sub>NOS (active ingredient) pH value: 3.0-6.0

### Raw Material GHS / US HCS / EC CLP Classification (100%):

DANGER!

Acute Tox. – inhl. Cat. 3, Skin Corr. Cat. 1B, Eye Damage Cat. 1, Skin Sens. Cat. 1, Aquatic Acute Cat. 1, Aquatic Chron. Cat. 1

H314, H317, H331, H410,

P261, P264, P271, P272, P273, P280, P301 + P330 + P331, P303 + P361 + P353,

P305 + P351 + P338, P310, P403 + P233, P405, P501

[Source: Raw Material vendor SDS, CCOHS databases and regulatory research]







#### Chemical Ingredient: Gentamicin sulfate

Chemical concentrations found in this product: < 0.01% from a 50 mg/ml Solution in C0, C1, C2, C4, R6 R7a, and

R7b (an aqueous solution)

Data for Concentrated / 100% chemical used in the product mixture (concentration tested):

CAS#: 1405-41-0 (100%) LD<sub>50</sub> (oral-rat): > 5000 mg/kg EC No: 215-778-9 (100%) LC<sub>50</sub> (inhalation-rat): NE RTECS#: LY2625000 (100%) LD<sub>50</sub> (skin-rabbit): NE

Synonyms/Trade Names: Gentamicin sulfate salt; Garamycin; Gentiomycin C

#### Raw Material GHS / US HCS / EC CLP Classification (10%):

#### DANGER!

Resp. Sens. Cat. 1A, Skin Sens. Cat. 1, Rep. Tox. Cat. 2

H317, H334, H361

P261, P272, P280, P285, P302 + P352, P304 + P341, P308 + P313, P333 + P313, P342 + P311, P363, P501



<b>Biological Ingredient</b>	Data / Information
Animal proteins	This material is of animal origin (fetal bovine / and bovine / and goat/caprine / and rabbit / and monkey / and porcine / and avian / and murine / and equine / and burro / and sheep/ovine / and mouse / and pig / and rat / and guinea pig, etc.) and may be a potential contact irritant. Hazard Unknown. Handle as potentially infectious. The chemical, physical and toxicological properties have not been thoroughly investigated. Handle appropriately with the requisite Good Laboratory Practices, <i>Standard</i> and <i>Universal Precautions</i> . Dispose of this material in accordance with local, regional, national and international regulation.



<b>Biological Ingredient</b>	Data / Information
Human blood derivatives	Human source material used in the preparation of the <b>Negative Control</b> (C0) and <b>Cutoff Calibrator</b> (C4), and as a diluent for the <b>HIV-1 Ab Positive Control</b> (C1) and <b>HIV-2/O Ab Positive Control</b> (C2), has been tested and found nonreactive for Hepatitis B surface antigen (HBsAg), and antibodies to Hepatitis C virus (HCV Ab) and human immunodeficiency virus (HIV-1 and HIV-2).
	Human source material, containing HIV-1 and HIV-2 human antibody used in the preparation of the <b>HIV Ab Positive Controls</b> (C1 and C2) has been heat-treated. It has been tested and found nonreactive for Hepatitis B surface antigen (HBsAg) and antibodies to Hepatitis C virus (HCV Ab).
	The HIV-1 viral lysate used in the <b>HIV Ag Positive Control</b> (C3) has been inactivated using a chaotropic agent and heat.
	No known test method can offer complete assurance that HIV, hepatitis B or C virus or other infectious agents are absent. Moreover, patient blood samples tested with this kit represent an unknown, heightened hazard. Employ Standard and Universal Precautions when handling these reagents and all human blood or specimens. Handle as if capable of transmitting infectious disease, in a Biosafety Level 2 lab, applying the guidelines from the current CDC/NIH Biosafety in Microbiological and Biomedical Laboratories or WHO Laboratory Biosafety Manual or equivalent. Avoid splashing, spills and the generation of aerosols. Secure in secondary containment with proper biohazard labeling. Do not inhale mists or aerosols; avoid contact with skin, eyes, mucous membranes and clothing. In case of contact with eyes, immediately rinse with copious water and seek medical attention. Employ decontamination procedures with appropriate decon agent or disinfectant (typically a 1:10 dilution of household bleach, 70-80% ethanol or isopropanol, an iodophor like 0.5% Wescodyne Plus (EPA Reg. #4959-16), an o-phenylphenol/amyphenol such as 0.8% Vesphene (EPA Reg. #1043-87), or equiv.) before discarding any materials utilized or returning equipment used to general use. Dispose of this material in accordance with local, regional, national and international regulations. Handle appropriately with the requisite Good Laboratory Practices, Standard and Universal Precautions. Persons handling blood samples should have the option of receiving hepatitis B vaccination.

NA: Not Applicable.

NE: Not Established or Unknown (unable to locate data); typically for concentrate form unless otherwise specified.

#### **Related product information:**

- ♦ R Refer to section 16 for the full text of any *Comprehensive GHS-based Classification* statements, for the list of sources utilized in the assessment and for the key / legend to abbreviations and acronyms.
- No significant adverse health effects are expected by any route for the miscellaneous Tween 20 [C<sub>58</sub>H<sub>114</sub>O<sub>26</sub>, CAS #9005-64-5], hydrogen peroxide [H<sub>2</sub>O<sub>2</sub>, CAS# 7722-84-1,-≤ 0.1% v/v in R8], salts, sugars, buffers, water, animal sera, and other chemicals found in the HRP conjugate, buffers with protein stabilizers, dyes, and sodium acetate solutions, in the kit volumes and/or concentrations present [chemical or dilution is not subject to US HCS, EC CLP or other GHS-based hazard labeling].
- ♦ Component R1 contains < 0.1% of Cobalt (II) Chloride [CAS# 7646-79-9, EC No. 231-589-4], which is classified as an IARC Group 2B (possible human carcinogen) and EU Category 2 carcinogen, and silica quartz [CAS# 14808-60-7, EC No. 238-87-4], which in dust form is classified as an ACGIH Class A2 (suspected human carcinogen) and IARC Group 1 (carcinogenic to humans). These materials are pelletized and sealed in a desiccant packet within the plate pouch, which is unlikely to generate significant dust under normal conditions of use and is thus not typically considered a health hazard. However, health hazards could result from dusts generated if the packet is cut, split or otherwise compromised and a significant number of pellets were crushed to a powder form. Keep the desiccant packet intact as received in the microwell plate component package.
- ◆ According to the concept of *Universal Precautions* (29 CFR 1910.1030), all human blood and certain human body fluids must be treated as if known to be infectious for HIV, HBV and other bloodborne pathogens. No known test method can offer complete assurance that products derived from human blood will not transmit infection; thus, they should be handled as though they contain infectious agents. Furthermore, individual patient samples being tested represent a heightened, unknown hazard. Aerosolization/inhalation, contact and mucous membrane exposure should be avoided during sample and kit handling. Consider equipment that potentially comes in contact with human source material as contaminated until appropriately decontaminated.
- Do not eat, drink or smoke when using this product.
- ♦ Wear protective gloves / protective clothing / eye protection / face protection. Take off contaminated clothing and wash before reuse.



	SECTION 4: EMERGENCY FIRST AID MEASURES		
Health Effects:	Symptoms of over exposure may include headache, dizziness and congestion. Causes severe skin burns and eye damage. Severely irritating or corrosive to eyes; greater exposures can cause eye damage, including permanent impairment of vision. Risk of serious damage to eyes. May cause ingestion corrosive effects, including burning throat, mouth and stomach. Skin contact may result in dermatitis and may cause allergic skin reaction upon repeated exposure. May be toxic to developing fetus, generally at concentrations and volumes that greatly exceed that of this kit.		
Eye Contact:	Flush eyes with copious water for at least 15 minutes. Ensure adequate flushing by separating the eyelids with fingers while flushing with water. OBTAIN MEDICAL ATTENTION.		
Skin Contact:	Remove contaminated clothing. Flush skin with copious water and wash affected area with soap and water. If blood-to-blood contact occurs, or if more severe symptoms develop, consult a physician.		
Inhalation:	Remove person from exposure area to fresh air. If breathing becomes difficult, immediately call for emergency medical assistance. Treat symptomatically and supportively. Generally, this aqueous product is not a significant inhalation hazard in the kit volumes and concentrations present.		
If Swallowed:	If ingested, rinse out mouth thoroughly with water, provided the person is conscious, and OBTAIN MEDICAL ATTENTION. Rinse mouth. Do NOT induce vomiting. Call a physician or the local poison control center. Treat symptomatically and supportively. If vomiting occurs, keep head lower than hips to prevent aspiration.		
Notes to Physician	According to the OSHA Bloodborne Pathogens Standard (29 CFR 1910.1030), Universal Precautions apply. Persons handling human blood source samples should be offered hepatitis B vaccination prior to working with human source material.		

SECTION 5: FIREFIGHTING MEASURES			
Flammable Properties:	This product is aqueous based so is not considered flammable or combustible.		
Extinguishing Media:	Use extinguishing media appropriate for the surrounding fire.		
Hazardous Combustion Products:	May release toxic oxides of carbon, nitrogen and sulfur or toxic hydrogen chloride gas.		
Special Firefighting Procedures:	Conventional firefighting full protective equipment (with NIOSH-approved self-contained breathing apparatus) and procedures appropriate for the surrounding fire should be sufficient.		

### **SECTION 6: ACCIDENTAL RELEASE MEASURES**

- Avoid direct contact with skin, eyes, mucous membranes and clothing by wearing appropriate lab Personal Protective Equipment (PPE), including gloves, lab coat and eye/face protection.
- In the event of a hazardous material spill, contain the spill if it is safe to do so and immediately move to a safe area, free from potential aerosols, to decontaminate and/or safely remove any contaminated clothing, as necessary. IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. Isolate the hazard area and ventilate if appropriate. Ensure that appropriate spill cleanup materials and PPE are available and used.
- Prevent material from entering sewers, waterways or confined spaces.
- Follow established laboratory policy and applicable WHO/CDC/NIH biosafety and/or WHO/OSHA hazardous material and/or equivalent guidelines for appropriate hazardous chemical and/or biological material spill response and cleanup. Avoid release to the environment.
- Wear appropriate PPE. Immediately, and on-site if possible:
  - Decontaminate Biohazard/Human Source Material spills, which should always be treated as potentially infectious, including the area, spill materials and any contaminated surfaces or equipment. Utilize an appropriate chemical decon agent or disinfectant that is effective for the known or potential pathogens relative to the samples involved (commonly a 1:10 dilution of bleach, 70-80% Ethanol or Isopropanol, an iodophor (such as Wescodyne Plus), or a phenolic, etc.).
  - o Neutralize corrosive acidic spills with the appropriate Acid neutralization / adsorbent product.



- Clean the spill area with water and wipe dry. Spills can also be absorbed with an appropriate inert material (e.g. spill pillows, absorbent pads, etc.), which are secured in an appropriate, labeled, sealed container. Material used to absorb the spill may require hazardous material waste disposal. Infectious, Chemical and Laboratory wastes must be handled and discarded in accordance with all local, regional and national regulations.
- Refer to Sections 8 and 11 for more specifics.

	SECTION 7: HANDLING AND STORAGE INFORMATION
Handling:	This test kit should be handled only by qualified personnel trained in laboratory procedures and familiar with their potential hazards. Follow proper good laboratory practices and safety guidelines for handling chemical, biological and laboratory hazards.
	Do not smoke, eat, or drink in areas where patient samples and kit reagents are handled. Wash your hands after use. Wear appropriate personal protective equipment (PPE), including gloves, lab coat or equivalent and eye/face protection.
	Keep containers tightly closed; avoid splashing, spills and the generation of aerosols.
	Handle all human source material, specimens and equipment used to perform the operations as though they were capable of transmitting infectious disease, as per <i>Standard</i> and <i>Universal Precautions</i> .
	All personal protective equipment should be removed before leaving the work area. Refer to Section 8 for more specifics.
	Avoid release to the environment. Do not allow undiluted product hazardous chemical ingredient or large quantities of it to reach ground water or water course.
	Consult with your Environmental Health & Safety Office for assistance.
Storage:	Store according to product and label instructions (generally at 2-8 °C).
Caution, cor use.	isult accompanying documents. Read and follow all the precautions and warnings in the kit product instructions for
For in vitro	diagnostic use.

### SECTION 8: EXPOSURE CONTROL / PERSONAL PROTECTION MEASURES

Control Parameters – Component chemicals with with limit values that require monitoring at the workplace: The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

Chemical	CAS-No.	Value	Control parameter	Update	Basis	
Sulfuric acid	7664-93-9	TWA – TLV	0.2 mg/m <sup>3</sup> (thoracic fraction)	2004-01-01	USA. ACGIH Threshold Limit Values (TLV)	
		TWA – PEL	1 mg/m <sup>3</sup> *	1993-06-30	USA. Occupational Exposure Limits (OSHA) - Table Z-1 Limits for Air Contaminants	
		REL IDLH	1 mg/m <sup>3</sup> 15 mg/m <sup>3</sup>	2005-149 [SEP-2007]	USA. National Institute for Occupational Safety and Health (NIOSH)	
	laboratory a conflicting carcinogen	The value in mg/m³ is approximate  *Remarks: TLV CARCINOGENICITY DESIGNATION A2 – Suspected Human Carcinogen: Substance is carcinogenic in laboratory animals under conditions that are considered relevant to worker exposure. Available human studies are conflicting or insufficient to confirm an increased risk of cancer in exposed humans. Worker exposure to an A2 carcinogen should be controlled to levels as low as reasonably achievable below the TLV.  The A2 Carcinogenicity Designation refers to sulfuric acid contained in strong inorganic acid mists.				
Hydrochloric	7647-01-0	TLV – C	2 ppm	2007-01-01	USA. ACGIH Threshold Limit Values (TLV)	
acid		PEL – C	7 mg/m <sup>3</sup> * 5 ppm	2006-02-28	USA. Occupational Exposure Limits (OSHA) - Table Z-1 Limits for Air Contaminants	
		REL – C IDLH	7 mg/m <sup>3</sup> 5 ppm 50 ppm	2005-149 [SEP-2007]	USA. National Institute for Occupational Safety and Health (NIOSH)	
	* The value in mg/m³ is approximate. Ceiling limit is to be determined from breathing-zone air samples.					



Chemical	CAS-No.	Value	Control parameter	Update	Basis	
	<b>Remarks</b> : TLV CARCINOGENICITY DESIGNATION A4 – Not Classifiable as a Human Carcinogen: Inadequate data on which to classify the substance as a human and/or animal carcinogen.					
Hydrogen	7722-84-1	TWA – TLV	1 ppm	2007-01-01	USA. ACGIH Threshold Limit Values (TLV)	
peroxide		TWA – PEL	1.4 mg/m <sup>3</sup> * 1 ppm	1997-08-04	USA. Occupational Exposure Limits (OSHA) - Table Z-1 Limits for Air Contaminants	
		REL IDLH	1.4 mg/m <sup>3</sup> 1 ppm 75 ppm	2005-149 [SEP-2007]	USA. National Institute for Occupational Safety and Health (NIOSH)	
	* The value in mg/m <sup>3</sup> is approximate					
	<b>Remarks</b> : TLV CARCINOGENICITY DESIGNATION A3 – Animal Carcinogen: Substance is carcinogenic in labora animals under conditions that are not considered relevant to worker exposure. Available human studies and evidence suggest that the substance is not likely to cause cancer in humans except under unusual or unlikely routes or levels of exposure. Worker exposure to an A3 carcinogen should be controlled to levels as low as reasonably achievable below TLV.				xposure. Available human studies and evidence ept under unusual or unlikely routes or levels of	
Dimethyl	67-68-5	TWA-WEEL	250 ppm	2014	USA: Workplace Environmental Exposure Levels	
sulfoxide		MAK	50 ppm (160 mg/m <sup>3</sup> )	2011	GERMANY:	
	Source: CCOHS CHEMINFO 2013, RTECS September 2013 Update and Raw Material Vendor Safety Data Sheet					

Additional information: The lists that were valid during the creation were used as basis.

The following personal protective equipment (PPE) is recommended to prevent blood or other potentially infectious or hazardous materials from reaching the user's work or street clothes, skin, mouth, mucous membranes and eyes, or hazardous inhalation, under normal conditions of use and for the time during which the protective equipment is utilized:

Ventilation:	Adequate lab ventilation is required. It is recommended that users handle potentially infectious human source material / patient samples in a biological safety cabinet (BSC), expressly if aerosols might be generated.
Eye / Face Protection:	Wear ANSI approved safety glasses, goggles or face shield with safety glasses or goggles. Contact lenses should not be worn when handling lab hazards.
Protective Gloves:	Suitable gloves must be worn at all times when handling kit reagents or patient samples to provide skin protection from splash and intermittent contact. Synthetic gloves, such as Nitrile, Neoprene and Vinyl, are recommended because they are sturdy, effective and contain no natural latex ingredients associated with latex glove allergic reactions. Disposable (single use) gloves should be changed often and never be reused. Wash hands thoroughly after removing gloves.
	Guidelines for <i>Sulfuric Acid</i> , less than 30%:  RECOMMENDED (resistance to breakthrough longer than 8 hours): Butyl rubber, natural rubber, neoprene, polyethylene, polyvinyl chloride, Viton <sup>®</sup> , Viton <sup>®</sup> /Butyl rubber, Barrier (PE/PA/PE), Silver Shield/4H <sup>TM</sup> (polyethylene/ethylene vinyl alcohol), Trellchem <sup>®</sup> HPS, Trellchem <sup>®</sup> VPS, Tychem <sup>®</sup> SL (Saranex <sup>®</sup> ), Tychem <sup>®</sup> CPF 3, Tychem <sup>®</sup> F, Tychem <sup>®</sup> BR/LV, Tychem <sup>®</sup> Responder <sup>TM</sup> , Tychem <sup>®</sup> TK.  RECOMMENDED (resistance to breakthrough longer than 4 hours): Nitrile rubber.  NOT RECOMMENDED for use (resistance to breakthrough less than 1 hour): Polyvinyl alcohol.  Source: CHEMINFO 12-2013
Protective Clothing:	Wear a lab coat, clinic jacket, gown, apron and/or smock. Disposable clothing is strongly recommended when handling biohazardous material. If reusable clothing is used, procedures for handling potentially infectious laundry under the OSHA Bloodborne Pathogens Standard (29 CFR 1910.1030) are required.
Respiratory Protection:	Do not breathe mist / vapors/vapours / spray.
Other:	All personal protective equipment should be removed before leaving the work area and placed in an appropriately designated area or container for storage, processing, decontamination or disposal. Protective coverings such as plastic wrap, aluminum foil, or imperviously-backed absorbent pads used to cover equipment and/or surfaces must be removed and replaced if they become overtly contaminated.



SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES						
Appearance:	Variable, refer to section 2; generally aqueous liquids. Exceptions are the solid microtiter plate and related materials.					
Odor/odour:	No applicable information was found. Odor/odour Threshold: Not Established.					
рН:	Most of the liquid chemical components are between pH 6 and 8, Exceptions are the following acidic solutions: Substrate Buffer at pH $\sim$ 4, Chromogen at pH $\sim$ 1.5 and Stopping Solution at pH $\leq$ 2.					
Boiling point:	Undetermined.	Melting point:	Undetermined.			
Flash point:	Not Applicable. Flammable limits: LEL/LFL is Not applicable; UEL/UFL is Not applicable.					
Evaporation rate:	No applicable information was found.					
Fire hazard:	Although the components have not been tested for fire hazard and explosion data, being water-based, they are not expected to be fire hazards, but some of the kit packaging materials may burn under fire conditions.					
Vapor/vapour pressure:	No applicable information was found.					
Vapor/vapour density:	No applicable information was found.					
Relative Density:	Variable; approximately 1.					
Solubility:	The liquid chemical components are soluble in water. The acidic solutions may release heat.					
Partition coefficient (n-octanol/water):	No applicable information was found.					
Auto igniting:	Product is not known to be self-igniting.					
Decomposition temperature:	No applicable information was found.					
Viscosity:	No applicable information was found.					
Danger of explosion:	Product is not known to present an explosion hazard.					
Molecular mass:	Mixture.					
No Other Standard Characteristics applicable to the identification or hazards of the product are known.						
Source: Raw Material vendor SDS, CCOHS databases and/or regulatory research						

# SECTION 10: STABILITY AND REACTIVITY INFORMATION

NOTE: Chemical reactions that could result in a hazardous situation (e.g. generation of flammable or toxic chemicals, fire or detonation) are listed here. Although not intended to be complete, an overview of important reactions involving common chemicals is provided to assist in the development of safe work practices.

Chemical Stability / Reactivity:	Components are stable with no known inherent significant reactivity, except the acidic solutions, which may have an exothermic reaction with certain chemicals, particularly strong bases and reducing agents.
Conditions and/or Materials to Avoid:	Do not allow the acidic <i>Stop Solution</i> to come in contact with strong bases or reducing agents (may lead to a violent exothermic reaction). <i>Sulfuric Acid</i> - Bases, Halides, Organic materials, Carbides, fulminates, Nitrates, picrates, Cyanides, Chlorates, alkali halides, Zinc salts, permanganates, e.g. potassium permanganate, Hydrogen peroxide, Azides, Perchlorates., Nitromethane, phosphorous, Reacts violently with: cyclopentadiene, cyclopentanone oxime, nitroaryl amines, hexalithium disilicide, phosphorous(III) oxide, Powdered metals



Incompatible materials:	Sulfuric acid:
	Although concentrated <i>sulfuric acid</i> is referred to as an oxidizing agent in some sources, it is not a very strong oxidizing agent. The 98% acid has some oxidizing ability when hot. <i>Sulfuric acid</i> does not polymerize and does not form peroxides.
	Sulfuric acid is a very reactive substance. The concentrated acid dehydrates, or sulfonates most organic compounds. Sulfuric acid reacts vigorously, violently or explosively with many organic and inorganic chemicals including water, acrylonitrile, alkali solutions, carbides, chlorates, fulminates, nitrates, perchlorates, permanganates, picrates, powdered metals, metal acetylides or carbides, epichlorohydrin, aniline, ethylenediamine, alcohols with strong hydrogen peroxide, chlorosulfonic acid, cyclopentadiene, hydrofluoric acid, nitromethane, 4-nitrotoluene, phosphorus (III) oxide, potassium, sodium, ethylene glycol, isoprene, styrene.
	Hazardous gases, such as hydrogen, hydrogen cyanide, hydrogen sulfide and acetylene, are evolved on contact with chemicals such as metals, cyanides, sulfides and mercaptans and carbides respectively.
Hazardous Decomposition Products:	May release toxic oxides of carbon, nitrogen and sulfur or toxic hydrogen chloride gas.
Hazardous Polymerization:	Has not been reported to occur.

### SECTION 11: TOXICOLOGICAL INFORMATION -- GENERAL COMPOSITE

Refer to Sections 2 and 3 for the kit component concentrations. The composite toxicological information for this product is:

### **Acute Health Effects**

Acute Toxicity:	May be detrimental in contact with skin, if swallowed, and to eyes upon contact; in case of contact with eyes, immediately rinse with copious water and seek medical attention.	
Primary Irritant Effect:	Irritating to skin and severely irritating or corrosive to eyes, and with greater exposures can cause eye damage, including permanent impairment of vision or blindness.	
Skin Corrosivity / Metal Corrosion:	Causes severe skin burns and eye damage. The <i>Stopping Solution</i> (R10) is Corrosive, able to cause severe burns of the mucous membranes, skin and eyes. Destructive to tissue of the skin, respiratory tract, mucous membranes, and eyes; may cause permanent injury. May cause ingestion corrosive effects, including burning throat, mouth and stomach.  Note: The <i>EIA Chromogen</i> ( <i>11X</i> ) solution in this product has been evaluated with the CORROSITEX test method to determine its corrosive potential. The results of this testing classified this solution as non-corrosive.	
Serious Eye Damage / Irritation:	Causes severe eye damage. The <i>Stopping Solution</i> (R10) is severely corrosive to eyes; contact can cause eye damage, including permanent impairment of vision or blindness. The <i>Stopping Solution</i> poses a risk of serious damage to eyes. Harmful to eyes upon contact; in case of contact with eyes, immediately rinse with copious water and seek medical attention.	
STOT-Single Exposure:	No applicable information was found.	
Aspiration Hazard:	No applicable information was found.	
Other Acute Health Effects:	No significant other acute health effect known.	

### Biohazard Potential:

Inactivated Human Immunodeficiency Virus, type 1 (HIV-1), though verified to be non-infectious [inactivation evidenced by a  $\geq 99\%$  reduction in RT activity (example:  $10^{-5}$  viable  $\geq 10^{-1}$  disrupted) by standard RT assay], should be handled with *Standard* and *Universal Precautions*, as if capable of transmitting infectious disease. The HIV-1 Ab Positive Control (2) was heat-treated to inactivate the HIV. Human sera in reagents were tested and found non-reactive for Hepatitis B surface antigen and antibodies to HCV (the Negative Control (C0), HIV-1 Ag Positive Control (C3) and Cutoff Calibrator (C4) are also



non-reactive for antibodies to HIV). No known test method can offer complete assurance that HIV, hepatitis B or C virus or other infectious agents are absent. Moreover, patient blood samples tested with this kit represent an unknown, heightened hazard. Employ *Standard* and *Universal Precautions*; handle these reagents, all human blood and specimens as if capable of transmitting infectious disease, in a Biosafety Level 2 laboratory, applying the guidelines from the current CDC/NIH *Biosafety in Microbiological and Biomedical Laboratories*, WHO *Laboratory Biosafety Manual* or equivalent. Persons handling blood samples should have the option of receiving hepatitis B vaccination.

### **Chronic Toxicity**

Respiratory or Skin Sensitization:	Contains a small volume of very dilute, potentially skin-contact sensitizing preservatives, <i>ProClin</i> and <i>Gentamicin sulfate</i> (an antimicrobial biocide that is also a photosensitizer); prolonged or repeated exposure may cause allergic reaction in certain sensitive individuals. Though the potential for an allergic response is greatly reduced by the dilution, sensitization threshold is unknown; thus handle accordingly.
Carcinogenicity:	Component R1 contains < 0.1% Cobalt (II) chloride (CAS# 7646-79-9, IARC Group 2B and EU Category 2 carcinogen) and silica quartz (CAS# 14808-60-7, ACGIH class A2 and IARC Group 1 carcinogen). Keep the desiccant packet intact as received in the component package. Component R10 contains Sulfuric Acid, CAS# 7664-93-9: .IARC Group 1, The agent is Carcinogenic to Humans, NTP listed as Known to be a Human Carcinogen and ACGIH-TLV Group A2, Suspected Human Carcinogen  Note: The IARC Group and ACGIH A2 1 classifications refers specifically to sulfuric acid contained in strong inorganic acid mists are and does not apply to sulfuric acid or sulfuric acid solutions.
Germ Cell Mutagenicity:	No applicable information was found.
Reproductive hazard:	Reasonably anticipated to be a reproductive toxin. May cause harm to unborn child. <b>Gentamicin sulfate</b> is known to the State of California to cause developmental toxicity (teratogen), classified under the generic class of aminoglycosides. (Designation is for concentrated gentamicin sulfate, which is diluted to 0.005% in kit components.)
STOT-Repeated Exposure:	No applicable information was found.

Additional Toxicological Information: To the best of our knowledge, the chemical, physical and toxicological properties have NOT been thoroughly investigated for some of the component chemicals and/or mixtures.

### **SECTION 12: ECOLOGICAL INFORMATION**

This product was not tested. The following assessment is based on information for the ingredients.

Ecotoxicity:	Concentrated Sulfuric acid [CAS# 7664-93-9] *:
	Fish LC <sub>50</sub> - Gambusia affinis (Mosquito fish) – 42 mg/l - 96 h
	Concentrated Hydrochloric acid [CAS# 7647-01-0] *:
	Fish LC <sub>50</sub> - Bluegill/Sunfish – 282 mg/l - 48 h
	Concentrated 2-methyl-4-isothiazolin [CAS# 2682-20-4] **:
	Fish $LC_{50}$ – Lepomis macrochirus (Bluegill) – 0.30 mg/l [min. 0.24 mg/l, max. 0.32 mg/l] – 96 h Fish $LC_{50}$ – Oncorhynchus mykiss (rainbow trout) – 0.19 mg/l [min. 0.13 mg/l, max. 0.31 mg/l] – 96 h Daphnia $EC_{50}$ - Daphnia pulex (Water flea) – min. 0.71 mg/l, max. 1.8 mg/l – 48 h
	Concentrated <i>ProClin 300</i> (Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3- one [EC no. 220-239-6] (3:1); Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC no. 220-239-6] (3:1)) [CAS# 55965-84-9]**:
	Fish LC <sub>50</sub> – Lepomis macrochirus (US EPA FIFRA 72-1, flow-through mm) – 0.28 mg a.s. /L – 96 h
	Fish LC <sub>50</sub> – Oncorhynchus mykiss (OECD TG 203, static) – 0.22 mg a.s. /L – 96 h
	Fish LC <sub>50</sub> – Oncorhynchus mykiss (OECD TG 204, flow-through mm) – 0.09 mg a.s. /L – 14 d
	Fish 36-d NOEC – Pimephales promelas (US EPA FIFRA 72-4, flow-through mm) – 0.02-0.12 mg a.s. /L
	Fish 28-d NOEC - Oncorhynchus mykiss (OECD TG 215, semi-static nom) - 0.098 mg a.s. /L
	Daphnia EC <sub>50</sub> – Daphnia magna (US EPA 72-2, flow-through mm) – 0.16 mg a.s. /L – 48 h
	Daphnia $EC_{50}$ – Daphnia magna (US EPA 72-4, semi-static mm) – 0.10 mg a.s. /L – 21 d

[Catalog # 12009912]

	Concentrated Citric acid [CAS#: 77-92-9]*:  Toxicity to fish mortality LC <sub>50</sub> – Leuciscus idus melanotus – 440 mg/l – 48 h  Method: OECD Test Guideline 203  Toxicity to daphnia and other aquatic invertebrates: static test – Daphnia magna (Water flea) – 1,535 mg/l  - 24 h  * Source: Raw Material Vendor Safety Data Sheets, RTECS and/or CCOHS Cheminfo  **Source: PAN Pesticides Database – Chemical Studies on Aquatic Organisims [obtained 3/7/2012]
Persistence and degradability:	No information found.
Bioaccumulation potential:	No information found.
Mobility in soil:	No information found.
PBT and vPvB assessment:	No information found.
Other adverse effects:	The acidic corrosive Components <b>R8</b> ( <b>pH 4</b> ), <b>R9</b> ( <b>pH 1.5</b> ) and <b>R10</b> ( <b>pH ≤2</b> ) are hazardous for drinking water and toxic to aquatic organisms by pH modification if not neutralized.  An environmental hazard cannot be excluded in the event of unprofessional handling or disposal.

Avoid release to the environment.

### **SECTION 13: DISPOSAL CONSIDERATIONS**

Disposal of hazardous and/or laboratory wastes, product or packaging must be conducted in accordance with all applicable local, regional, national and international regulations. This section specifies the general and United States RCRA requirements. Processing, use or contamination of the kit components may change waste management requirements and options. Contact your Environmental Health & Safety Office for your specific disposal procedures.

#### **Recommended Product Disposal:**

- All *human source* and *other potentially infectious material* must be appropriately decontaminated or disposed of as infectious material; check your international, national, regional and local ordinances accordingly.
- Acidic Stopping Solution (sulfuric acid, pH ≤ 2), Chromogen (pH ~1.5), and Substrate Buffer (pH ~4.0) wastes should be neutralized to pH 6-8 for safe sewer disposal; check your local, regional, national and international, ordinances accordingly. If the final pH measures ≤ 2, it requires disposal as a corrosive material in a RCRA approved waste facility (or equivalent); the US RCRA Waste disposal Code for this waste, if not neutralized, is D002, check your international, national and regional ordinances accordingly.

Do not allow undiluted product or large quantities of it to reach ground water or water course.

**Recommended Unclean Packaging Disposal:** Dispose in accordance with all applicable local, regional, national and international regulations.

#### **SECTION 14: TRANSPORT INFORMATION**

Shipping of product, packaging and waste must be conducted in accordance with all applicable local, regional, national and international regulations. Processing, use or contamination of the kit components may change shipping requirements and options. Contact your Environmental Health & Safety Office for your specific shipping procedures.

**Recommended Unused Product Multi-Modal Transportation:** According to US DOT, IATA and UN "Model Regulations", the *STOPPING SOLUTION* in the kit must be transported as follows:

Acidic Component *Stopping Solution* in this product contains 1 N Sulfuric acid (< 5% H<sub>2</sub>SO<sub>4</sub>).

Thus any unneutralized discarded product component or waste generated from its use resulting in a corrosive liquid (pH  $\leq$  2 or a pH  $\geq$  12.5 per Method 9040 (USEPA Publication SW-846) or which corrodes Steel (NACE Standard TM-01-69)) must be transported as follows:

Proper Shipping Name and Description: Sulphuric acid [with not more than 51% acid]

Hazard Class or Division: **8** UN ID Number: **UN 2796** 

Packing group II





[Catalog # 12009912]

The EIA Chromogen (11X) solution in this product has been evaluated with the CORROSITEX test method to determine its corrosive potential and any packing group classification. The results of this testing classified this STOP solution as non-corrosive for shipping purposes.

Recommended Used Product Hazardous Waste Disposal Transportation: Air and land transportation information for discarded kit components and waste from this product when used as intended is:

**Acidic** *Chromogen* is at pH  $\sim$ 1.5 and 1N sulfuric acid *Stopping Solution* is at pH  $\leq$  2.

Thus any un-neutralized discarded product component or waste generated from its use resulting in a corrosive liquid (pH  $\leq$  2 or an pH  $\geq$  12.5 per Method 9040 (USEPA Publication SW-846) or Corrodes Steel (NACE Standard TM-01-69)) must be transported as follows:

Proper Shipping Name and Description: Waste Corrosive Liquid n.o.s.

Hazard Class or Division: 8 UN ID Number: UN 1760 Packing group III



Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code: Not applicable.

SECTION 15: REGULATORY INFORMATION				
Composite HMIS Rating:	Health: 2	Flammability: 0	Reactivity: 1	

### **Carcinogenicity Categories:**

**Component R1** contains < 0.1% **Cobalt (II) chloride** (CAS# 7646-79-9, IARC Group 2B and EU Category 2 carcinogen) and **silica quartz** (CAS# 14808-60-7, in dust form is an ACGIH class A2 and IARC Group 1 carcinogen) in a pelletized desiccant sealed packet. Keep the desiccant packet intact as received in the component package.

**Component R10** contains **1N Sulfuric Acid,** CAS# 7664-93-9: .IARC Group 1, The agent is Carcinogenic to Humans, NTP listed as Known to be a Human Carcinogen and ACGIH-TLV Group A2, Suspected Human Carcinogen

Note: The IARC Group 1 and ACGIH A2classifications refers specifically to sulfuric acid contained in strong inorganic acid mists are and does not apply to sulfuric acid or sulfuric acid solutions.

### National Regulations - Other Domestic / Foreign Laws:

*Hazard communication compliance* – This SDS contains the required information for preparation in accordance with the following GHS-based global regulations:

- United States Occupational Safety Health Administration Hazard Communication Standard 29 CFR 1910.1200 (US HCS)
- 2. Taiwan -OSHA Published National Standard CNS 15030 Classification and Labelling of Chemicals
- 3. Singapore SS 586 2:2014
- 4. Russia GOST 31340-2013, GOST 32419-2013, GOST 32423-2013, GOST 32424-2013, GOST 32425-2013, R 50.1.102-2014, R 50.1.101-2014
- 5. People's Republic of China National Standard GB/T 17519-2013, GB 30000-2013
- 6. New Zealand Hazardous Substances and New Organisms Act (HSNO)

Composite HSNO Hazard Class: Subclass 6.5 Category B (contact sensitizers)

Subclass 8.2 Category B (skin corrosive, GHS 1B)

Subclass 8.3 Category A (eye corrosive)

- 7. Mexico Standard NOM-018-STPS-2015, NMX-R-019-SCFI-2011
- 8. Korea MoEL-Public Notice 2016-19, 2013-37 Standard for classification and labeling of chemical substances and MSDS.
- 9. **Japan** Industrial Safety and Health Law (ISHL) National Standard **JIS Z7252**, **JIS Z7253**
- 10. European Commission (EC) applicable CLP related regulations (2010/453/EC, 2008/1272/EC, 2006/1907/EC etc.)
- 11. Canada Hazardous Products Regulations (HPR) / Standard Workplace Hazardous Materials Information System (WHMIS-GHS) Canadian Standard for the hazard classification criteria for this product.

Composite WHMIS Hazards: Skin Corrosion

Serious Eye Damage Skin Sensitization

- 12. Brazil Regulation ABNT NRB 14725
- 13. Australia Code of Practice Preparation of Safety Data Sheets for Hazardous Chemicals under Section 274 of the Work Health and Safety (WHS) Act.



#### 14. Analogous GHS-based global regulations

### Inventory status

Country(s) or region Inventory name	In Compliance (yes/no)*
Australia – Australian Inventory of Chemical Substances (AICS)	Yes
Canada – Domestic Substances List (DSL)	Yes
Canada – Non-Domestic Substances List (NDSL)	Yes
China – Inventory of Existing Chemical Substances in China (IECSC)	Yes
European Union – European Inventory of Existing Commercial Chemical Substances (E	EINECS)
or European List of Notified Chemical Substances (ELINCS)	Yes
Japan – Inventory of Existing and New Chemical Substances (ENCS)	Yes
Korea – Existing Chemicals List (ECL)	Yes
New Zealand – New Zealand Inventory	Yes
Philippines – Philippine Inventory of Chemicals and Chemical Substances (PICCS)	Yes
Taiwan – inventory CSNN	Yes
United States & Puerto Rico - Toxic Substances Control Act (TSCA) Inventory	Yes

<sup>\*</sup> A "Yes" indicates that all components of this product comply with the inventory requirements administered by the governing country(s)

### Regulation (EC) No. 1907/2006 (REACH):

Chemicals included in the Candidate List of Substances of Very High Concern (SVHC): None

REACH No.: A registration number is not available for this substance as the substance or its uses are exempted from registration, the annual tonnage does not require a registration or the registration is envisaged for a later registration deadline.

#### **United States SARA:**

SARA 302 (extremely hazardous substance) components: The following components are subject to reporting levels established by SARA Title III, Section 302: Sulfuric Acid, CAS# 7664-93-9; Revision Date: 2007-07-01 Hydrogen peroxide, CAS# 7722-84-1; Revision Date: 1993-04-24

SARA 313 components: The following components are subject to reporting levels established by SARA Title III, Section 313: Sulfuric Acid, CAS# 7664-93-9; Revision Date: 2007-07-01

California Proposition 65 (California Safe Drinking Water and Toxic Enforcement Act of 1986): WARNING: This Product Contains a Chemical(s) Known to the State of California to Cause Reproductive Toxicity.

Chemicals known to cause reproductive Toxicity: Gentamicin Sulfate CAS# 1405-41-0; classified under the generic class of Aminoglycosides. (Listed October 1, 1992)

### **SECTION 16: OTHER INFORMATION**

#### **Hazard statement abbreviation(s):**

Acute Tox. - oral. Acute toxicity - ingested (swallowed) Acute Tox. - skn. Acute toxicity – skin contact Acute Tox. - inhl. Acute toxicity - inhaled Resp. Sens. Respiratory sensitization Skin Sens. Skin sensitisation Skin Corr. Skin corrosion Skin Irrit. Skin irritation Eye Damage. Serious eye damage Eye Irrit. Eye irritation Rep. tox. Reproduction toxicant STOT SE

Specific target organ toxicity - single exposure

Aquatic Acute Acute aquatic toxicity Aquatic Chronic Chronic aquatic toxicity Met. Corr Corrosive to Metals Fla. Liq. Flammable liquid Cat. Category

H227 Combustible liquid. H290 May be corrosive to metals. H302 Harmful if swallowed

H313 May be harmful in contact with skin.



H314	Causes severe skin burns and eye damage.
H316	Causes mild skin irritation.
H317	May cause an allergic skin reaction.
H319	Causes serious eye irritation.
H331	Toxic if inhaled.
H334	May cause allergy or asthma symptoms or breathing difficulties if inhaled.
H335	May cause respiratory irritation.
H361	Suspected of damaging fertility or the unborn child.
H410	Very toxic to aquatic life with long lasting effects.
H412	Harmful to aquatic life with long lasting effects.
11412	Harmitul to aquatic fire with long fasting effects.
P234	Keep only in original container.
P210	Keep away from heat/sparks/open flames/hot surfaces. — No smoking
P261	Avoid breathing mist / vapors/vapours / spray.
P264	Wash skin thoroughly after handling.
P270	Do not eat, drink or smoke when using this product.
P271	Use only outdoors or in a well-ventilated area.
P272	Contaminated work clothing should not be allowed out of the workplace.
P273	Avoid release to the environment.
P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
P285	In case of inadequate ventilation wear respiratory protection.
P302 + P352	IF ON SKIN: Wash with plenty of soap and water.
P301 + P312 + P330	IF SWALLOWED: Call a POISON CENTER or doctor/ physician if you feel unwell. Rinse mouth.
P301 + P330 + P331	IF SWALLOWED: rinse mouth. Do NOT induce vomiting.
P303 + P361 + P353	IF ON SKIN (or hair): Remove/ Take off immediately all contaminated clothing. Rinse skin with water/ shower.
P304 + P340	IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
P304 + P341	IF INHALED: If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing.
	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P305 + P351 + P338 P308 + P313	IF the Execution of the
P310	Immediately call a POISON CENTER or doctor/ physician.
P333 + P313	If skin irritation or rash occurs: Get medical advice/ attention.
P337 + P313	If eye irritation persists: Get medical advice/ attention.
P342 + P311	If experiencing respiratory symptoms: Call a POISON CENTER or doctor/ physician.
P363	Wash contaminated clothing before reuse.
P370 + P378	In case of fire: Use dry sand, dry chemical or alcohol-resistant foam for extinction.
P390	Absorb spillage to prevent material damage.
P391	Collect spillage.
P403 + P233	Store in a well-ventilated place. Keep container tightly closed.
P403 + P235	Store in a well-ventilated place. Keep cool.
P405	Store locked up.
P406	Store in corrosive resistant stainless steel container with a resistant inner liner.
P501	This material and its container must be disposed of as hazardous waste.
P501	Dispose of contents / container to an approved waste disposal plant.
P501	Dispose of this material in accordance with local, regional, national and international regulation.
Caution	Contains human source material and inactivated pathogen. Handle as if capable of transmitting potentially infectious agents

This test kit should be handled only by qualified personnel trained in laboratory procedures and familiar with their potential hazards. Specific warnings are given in the instructions for use. The absence of a specific warning should not be interpreted as an indication of safety.

For in vitro diagnostic use.

Chemical safety assessment: Mixtures covered in this SDS were classified using the US HCS, EC CLP and/or UN Globally Harmonized System of Classification and Labeling of Chemicals (GHS) Fifth edition unless otherwise specified.

Sources of key data used to compile the Safety Data Sheet:

Raw Material Vendor Safety Data Sheets

United Nations (UN) Globally Harmonized System (GHS)

United States OSHA Hazard Communication Standard (US HCS) 1910.1200

Canadian Workplace Hazardous Materials Information System (WHMIS)

Mexican Standard (NOM-018-STPS-2015, NMX-R-019-SCFI-2011) [regulatory translation and summaries]

European Commission (EC) Regulations 2008/1272/EC, 2010/453/EC, 2006/1907/EC (EC CLP)

Australian Code of Practice - Preparation of Safety Data Sheets for Hazardous Chemicals (Section 274 of the Work Health and Safety Act)

New Zealand – Hazardous Substances and New Organisms Act (HSNO)

The People's Republic of China National Standard GB/T 17519-2013, GB 30000-2013 [regulatory translation if available and summaries]

Taiwan OSHA Published National Standard CNS 15030 [regulatory translation if available / summaries]

Korean MoEL-Public Notice 2016-19, 2013-37 [regulatory translation if available and summaries]



[Catalog # 12009912]

Japanese Industrial Standard JIS Z7252, JIS Z7253 [regulatory translation if available and summaries]

Registry of Toxic Effects of Chemical Substances (RTECS)

Canadian Centre for Occupational Health and Safety (CCOHS) CHEMINFO databases, etc.

International Agency for Research on Cancer (IARC)

American Conference of Governmental Industrial Hygienists (ACGIH)

Occupational Safety and Health Administration, U.S. Department of Labor (OSHA)

National Toxicity Program (NTP)

National Institute for Occupational Safety and Health (NIOSH)

World Health Organization. Laboratory Biosafety Manual

CDC/NIH Biosafety in Microbiological and Biomedical Laboratories

PAN Pesticides Database - Chemical Studies on Aquatic Organisms

Australian Inventory of Chemical Substances (AICS) Listing

California Proposition 65

#### Key / legend to abbreviations and acronyms used in the safety data sheet:

ACGIH - American Conference of Governmental Industrial Hygienists

AICS - Australian Inventory of Chemical Substances

ANSI - American National Standards Institute

CAS - Chemical Abstracts Service

CCOHS - Canadian Centre for Occupational Health and Safety

CDC - Centers for Disease Control, USA

CNS - Central Nervous System

DGSMA - Dangerous Goods Safety Management Act

DOT - Department of Transportation, USA

EC<sub>50</sub> – half maximal effective concentration

EC CLP - European Commission regulation for the Classification, Labeling and Packaging of chemical substances and mixtures

EU - European Union

GHS - Globally Harmonized System

HNOC - Hazard Not Otherwise Classified

HSNO - Hazardous Substances and New Organisms Act 1996 (New Zealand)

IARC - International Agency for Research on Cancer

IATA – International Air Transport Association

ICAO - International Civil Aviation Organization

IDLH – Immediately Dangerous to Life or Health

IMDG – International Maritime Dangerous Goods

IPCS - International Programme on Chemical Safety

ISHA – Industrial Safety and Health Act

 $LC_{50}$  – median lethal concentration, 50%

 $LD_{50}-median\ lethal\ dose,\ 50\%$ 

MSDS – Material Safety Data Sheet NIH – National Institute of Health

NIOSH – National Institute for Occupational Safety and Health

NTP - National Toxicity Program

OEL - Occupational Exposure Limit

PEL – Permissible Exposure Limit

ppm – parts per million

RTECS – Registry of Toxic Effects of Chemical Substances

SDS - Safety Data Sheet

STEL - Short Term Exposure Limit

STOT - Specific Target Organ Toxicity

TCCA - Toxic Chemical Control Act

TLV/TWA - Threshold Limit Value / Time-Weighted Average

UN - United Nations

US EPA - United States Environmental Protection Agency, USA

US HCS - Hazard Communication Standard, USA

US OSHA - Occupational Safety and Health Administration, U.S. Department of Labor

WHMIS -Workplace Hazardous Materials Information System, Canada

WHO – World Health Organization (United Nations)

Additional information: The lists that were valid during the creation were used as basis.

\_\_\_\_\_\_

This Revision: New Safety Data Sheet (SDS).

Preparation date: Refer to the date in the footer.

### **Bio-Rad Laboratories:**

Department issuing SDS: Environmental Health and Safety.

Contact for general SDS information: Seattle Operations, Environmental Health & Safety, 6565 185th Ave. NE, Redmond, WA 98052, USA, Phone: 425-881-8300 (8 am to 5 pm PT), ro-sds@bio-rad.com



Customer support contact: Clinical Diagnostics Group, 4000 Alfred Nobel Drive, Hercules, CA 94547, USA Phone: 1-800-224-6723, <a href="https://www.bio-rad.com/diagnostics">www.bio-rad.com/diagnostics</a>

*Emergency Contact* (24/7/365) – Chemtrec: 1-800-424-9300 or 1-703-527-3887 (USA/CAN) / +1-703-741-5970 (international – can be called collect). Refer below for in-country Chemtrec emergency contacts.

```
Contact Local Bio-Rad Agents for general information:
Australia, Bio-Rad Laboratories Pty. Ltd., Level 5, 446 Victoria Road, Gladesville NSW 2111 • Phone 61-2-9914-2800 • Fax 61-2-9914-2888 •
  Poison Information Centre: 13 11 26 (24 hours a day, anywhere in Australia); +61 13 11 26 • 24h/365d: 1800 039 008 • Chemtrec (24h/365d): +(61)-290372994
Austria, Bio-Rad Laboratories Ges.m.b.H., Hummelgasse 88/3-6, A-1130 Vienna • Phone 43-1-877-8901 • Fax 43-1-876-5629 •
  Poison Information Centre: +43 1 406 43 43 • Chemtrec (24h/365d): +(43)-13649237
Belgium, Bio-Rad S.A.-N.V. Winninglaan 3, BE-9140 Temse • Phone +32 (3)710-53-00 • Fax +32 (3)710-53-01 • Poison Information Centre / Belgisch Antigifcentrum:
  Brussles: +32 70 245 245; Luxembourg: 070 245 245 / 8002 5500 (every day, 24 to 24 hours) • Chemtrec (24h/365d): +(32)-28083237
Brazil, Bio-Rad Laboratórios Brasil Ltda, Rua Alfredo Albano da Costa, 100 / Distrito Industrial, Lagoa Santa - MG, CEP: 33400-000 • Phone +55 (31)3689-6600 •
  Fax +55 (31)3689-6611 • Contatos com a área de Suporte Técnico ao Cliente (CTS - Customer Technical Support): Telefones: 4003-0399 (Capital e Regiões
  Metropolitanas) e 0800-200-8900 (Outras Localidades) • Centro de Informações Toxicológicas: 0800 643 5252 •
  Telefone para emergências / Chemtrec (24h/365d): +(55)-2139581449
Canada, Bio-Rad Laboratories, Ltd., 2403 Guénette Street, Montréal, Québec H4R 2E9 • Phone 1-514-334-4372 • Fax 1-514-334-4415 •
   Poison Information Centre: Alberta: 1 800 332 1414; British Columbia: 1 800 567 8911; Manitoba: 1 855 776 4766; New Brunswick: 911;
  Newfoundland: 709 722 110); Northwest Territories: 1 800 332 1414; Nunavut: 867 979 7350 (Iqaluit), 867 983 2531 (Cambridge Bay),
  867 645 2816 (Rankin Inlet); Ontario: 1 800 268 9017; Prince Edward Island: 1 800 565 8161; Québec: 1 800 463 5060; Saskatchewan: 1 866 454 1212;
   Yukon: 867 393 8700 • Chemtrec (24h/365d): 1-800-424-9300 / 1-703-527-3887
China, Bio-Rad Laboratories Shanghai Ltd. 3rd Floor, #18 Dong Fang Road, Bldg E, Poly Plaza, Pudong, Shanghai, PRC 200120 • Phone 86-21-61698500 •
  Fax 86-21-61698599 • National Poison Control Center: +86 10 831 323 45 • Chemtrec (24h/365d): 4001-204937
Czech Republic, Bio-Rad spol. s r.o., Nad ostrovem 1119/7, 147 00 Prague 4 • Phone 420-241-430-532 • Fax 420-241-431-642 •
   Toxikologické informační středisko: +420 224 919 293, +420 224 915 402 • Chemtrec (24h/365d): +(420)-228880039
Denmark, Bio-Rad Laboratories, Symbion Science Park, Fruebjergvej 3, DK-2100 Copenhagen East • Phone +45-4452-1000 • Fax +45-4452-1001 •
  Poison Information Center / Giftlinjen Ring: +45 82 12 12 12 • Chemtrec (24h/365d): +(45)-69918573
Finland, Bio-Rad Laboratories, Linnanherrankuja 16, FIN-00950 Helsinki • Phone 358-9-804-22-00 • Fax 358-9-7597-5010 • Poison Information Centre:
  Myrkytystietokeskus: +358 9 471977 Avoinna 24 t / vrk • Chemtrec (24h/365d): +(358)-942419014
France, Bio-Rad, 3 boulevard Raymond Poincaré, 92430 Marnes-la-Coquette • Phone 33-1-47-95-60-00 • Fax 33-1-47-41-91-33 • Numéro ORFILA + 33 (0)1 45 42 59 59
  (24 heures sur 24 - 7 jours sur 7) • Chemtrec (24h/365d): +(33)-975181407
Germany, Bio-Rad Laboratories GmbH, Heidemannstrasse 164, D-80939 Munich • Phone +49-(0)89-318-840 • Fax +49-(0)89-318-84100 •
   Giftinformationszentren: Berlin: 030/19240 (Notruf); Bonn: 0228/19240 (Notruf); Erfurt: 0361/730 730; Freiburg: 0761/19240 (Notruf); Göttingen: 0551/19 240
  (Notruf); Homburg: 06841/19240 (Notruf) 06841/1628436 (Sekretariat); Mainz: 06131/19240 (Notruf); 06131-23 24 66 (Infoline);
  München: 089/19240 (Notruf) • Chemtrec (24h/365d): 0800-181-7059 oder +(49)-69643508409
Greece, Bio-Rad Laboratories M E.P.E, 2-4 Mesogeion Street, Fourth Floor 115 27 Athens • Phone 30-210-7774396 • Fax 30-210-7774376 •
   Poison Information Centre: (0030) 2107793777 (24 h/365d) • Chemtrec (24h/365d): +(30)-2111768478
Hong Kong, Bio-Rad Pacific Ltd., Unit 1101, 11/F DCH Commercial Centre, 25 Westlands Road, Quarry Bay • Phone 852-2789-3300 • Fax 852-2789-1290 •
  Poison Information Centre: +852 2772 2211 • Chemtrec (24h/365d): 800-968-793
Hungary, Bio-Rad Hungary Ltd., H-1082 Budapest, Futo Street 47-53, Hungary • Phone +36-1-459-6100 • Fax +36-1-459-6101 • Poison Information Centre:
  Információszolgáltatás akut mérgezés esetén (+36 80) 201-199 (0-24 h, díjmentesen hívható) • Chemtrec (24h/365d): +(36)-18088425
India, Bio-Rad Laboratories (India) Pvt. Ltd., Bio-Rad House, 86-87, Udyog Vihar, Phase IV, Gurgaon, Haryana 122 015 • Phone 1-800-180-1224 •
  Fax 91-124-2398115 • National Poison Information Centre: Toll Free No. – 1 800 116 117 (24h/365d) • Chemtrec (24h/365d): 000-800-100-7141
Israel, Bio-Rad Laboratories Ltd., 14 Homa Street, New Industrial Area, Rishon Le Zion 75655 • Phone 972-3-9636050 • Fax 972-3-9514129 •
  Poison Information Center: 04-7771900 (24/7) or 4 854 1900 • Chemtrec (24h/365d): +(972)-37630639
Italy, Bio-Rad Laboratories S.r.l., Via Cellini 18/A, 20090 Segrate, Milan • Phone +39-02-216091 • Fax +39-02-21609553 • Poison Information Centro / Centro
   Antiveleni: Bergamo: +39 800 883 300; Firenze: +39 55 794 7819; Foggia: +39 881 732 326; Milan: +39 02 6610 1029; Pavia: +39 38 224 444;
   Rome: +39 06 305 43 43; Turin: +39 011 663 7637 • Chemtrec (24h/365d): 800-789-767 o +(39)-0245557031
Japan, Bio-Rad Laboratories K.K., Tennoz Central Tower 20F, 2-2-24 Higashi-Shinagawa, Shinagawa-ku, Tokyo 140 0002 • Phone 81-3-6361-7070 •
  Fax 81-3-5463-8481 • Poison Information Centre: Tsukuba-City: +81 72 727 2499 +81 29 852 9999; Minoh City: +81 (0)72 727 2499 •
  Chemtrec (24h/365d): +(81)-345209637
Korea, Bio-Rad Korea Ltd., 10th Floor, Hyunjuk Building, 832-41, Gangnam-gu, Seoul 135-080 • Phone 82-2-3473-4460 • Fax 82-2-3472-7003 •
  Poison Control Information: Yongsan Army Base: 7917-5545/6001; Seoul Poison Center: 129 • Chemtree (24h/365d): +(82) 070-7686-0086 또는 00-308-13-2549
México, Bio-Rad, S.A., Avenida Eugenia 197, Piso 10-A, Colonia Narvarte, Delegación Benito Juárez, C.P. 03020 División de Diagnóstico Clínico "CDG" •
   Tel: 54887690 ó 54887670 Ext. 1031 • Correo: cts cdg mexico@bio-rad.com • Línea gratuita desde el interior de la república: 01 800 712 0190•
  Centro de Información Toxicologica: 1 800 009 2800 • Teléfono de emergencia: Chemtrec (24h/365d): 01-800-681-9531
The Netherlands, Bio-Rad Laboratories B.V., Postbus 222, 3900 AE Veenendaal • Phone +31-318-540666 • Fax +31-318-542216 • Poison Information Centre: NVIC
  (030 274 88 88), Only for the purpose of informing medical personnel in cases of acute intoxications. Uitsluitend bestemd om professionele hulpverleners te informeren
  bij acute vergiftigingen • Chemtrec (24h/365d): +(31)-858880596
New Zealand, Bio-Rad New Zealand, 189 Bush Road Unit B, Albany, Auckland • Phone 64-9-415-2280 • Fax 64-9-415-2284 •
  National Poisons Centre: +0800 764 766 • 24h/365d: 800 2436 2255 • Chemtrec (24h/365d): +(64)-98010034
Norway, Bio-Rad Laboratories, Nydalsveien 33, 0484 Oslo • Phone +47-23-38-41-30 • Fax +46(0)8-5551-2780 / +47-23-38-41-39 • Giftinformasjonen er det nasjonale
  rådgivnings- og kompetanseorganet vedrørende akutte forgiftninger og forgiftningsfare: +47 22 591300 (24h/365d) • Chemtrec (24h/365d): +(47)-21930678
Poland, Bio-Rad Polska Sp. z o.o., ul. Przyokopowa 33, 01-208 Warsaw • Phone +48-22-3319999 • Fax +48-22-3319988 • 24h/365d:+48 (22) 331 99 85 •
  [biorad@biorad.com.pl] • Poisons Information Centre: Gdansk: +48 58 682 04 04; Krakow: +48 12 411 99 99; Lodz: +48 42 63 14 724;
  Sosnowiec: +48 32 266 11 45; Warszawa: +48 22 619 66 54; Wroclaw: +48 71 343 30 08 • Chemtrec (24h/365d): +(48)-223988029
Portugal, Bio-Rad Laboratories, Lda., Edificio Prime, Ave. Quinta Grande, 53 - Fracção 3B Alfragide 26114-521 Amadora • Phone 351-21-472-7700 •
```

Fax 351-21-472-7777 • CIAV - Centro de Informação Antivenenos (Portuguese Poison Centre): 808 250 143 • Chemtrec (24h/365d): +(351)-308801773 Russian Federation, Bio-Rad Laboratorii, Russian Federation, Moscow, Varshavskoe sh., 9, Bldg., 1B • Phone: +7-495-721-1404 • Fax +7-495-721-1412 •

Singapore, Bio-Rad Laboratories (Singapore) Pte. Ltd., 27 International Business Park, #01-02 iQuest @IBP, Singapore 609924 • Phone 65-6415-3170 •

Fax 65-6415-3189 • 24/7/365: 65-6415-3188 • Poison Information Centre: 6423 9119 (24h/365d); +61 13 11 26 (Australian hotline)

Poison Information Center: Ekaterinburg: +7 343 229 98 57; Moscow: +7 495 628 1687; Saint-Petersburg: +7 921 757 3228 • Chemtrec (24h/365d): 8-800-100-6346

Chemtrec (**24h/365d**): 800-101-2201 或+(65)-31581349



[Catalog # 12009912]

South Africa, Bio-Rad Laboratories (Pty) Ltd., 34 Bolton Road, Parkwood, Johannesburg 2193 • Phone 27-11-442-85-08 • Fax 27-11-442-85-25 • Poison Information Centre: +27 824 910 160; Rondebosch: +27 (0) 861 555 777; Tygerberg: +27 (0) 861 555 777 • Chemtrec (24h/365d): 0-800-983-611

Spain, Bio-Rad Laboratories, S.A., C/ Caléndula, 95, Edificio M. Minipare II, El Soto de la Moraleja, 28109 Madrid • Phone 34-91-590-5200 • Fax 34-91-590-5211 • Servicio de Información Toxicológica: +34 91 562 0420 • Chemtrec (24h/365d): +(34)-931768545

Sweden, Bio-Rad Laboratories A.B., Box 1097, Solna Strandväg 3, SE-171 54, Solna • Phone +46-8-555-127-00 • Fax +46-8-555-127-80 • Poison Information Centre / Giftinformationscentralen: När det är akut (urgent):112 Begär giftinformation; I mindre akuta fall (less urgent): +46 (0) 10 456 6700 Direktnummer [Vi svarar på frågor om akuta förgiftningar - dygnet runt, året runt] • Chemtrec (24h/365d): +(46)-852503403

Switzerland, Bio-Rad Laboratories AG, Pra Rond 23, CH-1785 Cressier • Phone +41 (0)26-674-55-05/06 • Fax +41 (0)26-674-52-19 • Email: swiss@bio-rad.com • Toxicological Information Centre: 145 (24-h-Notfallnummer) • Chemtrec (24h/365d): +(41)-435082011

Taiwan, Bio-Rad Laboratories Taiwan Ltd., 14F-B, No. 126 Nan-King East Road, Sec. 4, Taipei, Taiwan 10546 R.O.C. • Phone 886-2-2578-7189 • Fax 886-2-2578-6890 • Chemtrec (24h/365d): 00801-14-8954

Thailand, Bio-Rad Laboratories Ltd., 1st & 2nd Floor, Lumpini I Bldg., 239/2 Rajdamri Rd., Lumpini, Pathumwan, Bangkok 10330 • Phone 662-651-8311 • Fax 662-651-8312 • Poisons Control Center: 1367, Line ID poisrequest (national callers); +66 2 419 7007 • Chemtree (24h/365d): 001-800-13-203-9987

United Kingdom, Bio-Rad Laboratories Ltd., 3rd &4th Floor, The Junction, Station Road, Watford, Hertfordshire, WD17 1ET • Phone +44 01 923 47 1301 • Fax +44 (0) 208 328 2550 • Freephone: 00 800 00 BIORAD (00 800 00 24 67 23) • Emergency number: +44 111 (NHS Emergency support) 9am-5pm Mon-Fri: • Email: techsupport.uk@bio-rad.com • Poisons Information Centre: NHS 111 / 844 892 0111 • Chemtrec (24h/365d): +(44)-870-8200418

United States: American Association of Poison Control Centers call (800) 222-1222 (24h/365d) • Chemtrec (24h/365d): 1-800-424-9300 / 1-703-527-3887

This document was developed from information obtained from reputable sources, but does not purport to be all-inclusive. The data contained herein, which is based on our present knowledge and is intended for information purposes only, shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship. Regulatory requirements are subject to change and vary from one location to another thus it is the buyer's responsibility to ensure that its activities comply with international, national, regional and local laws and regulations. Bio-Rad Laboratories makes no warranty, expressed or implied, regarding the accuracy or completeness of these data or the results to be obtained from the use thereof. Since the use of this information and the conditions of use of the product are not within the control of Bio-Rad Laboratories, it is the user's obligation to determine the suitability of the information for the intended application and use appropriate safety procedures.