

09/01/2020 Kit Components

Product code	Description
1762510 1762500 1762600	ProteOn™ HTG & HTE Reagent Kit
Components:	
10021442	ProteOn TM HTG Chip Activation Solution, 10mM Nickel Sulfate pH 6.0
10021443	ProteOn TM HTG Chip Regeneration Solution, 300 mM EDTA, pH 8.5





Printing date 09/01/2020 Reviewed on 08/11/2020

1 Identification

- · Product identifier
- · Trade name: ProteOn™ HTG Chip Activation Solution, 10mM Nickel Sulfate pH 6.0
- · Catalog or product number: 10021442
- · Application of the substance / the mixture Laboratory chemicals
- · Details of the supplier of the safety data sheet
- · Manufacturer/Supplier:

Bio-Rad Laboratories (Canada) Ltd.

1329 Meyerside Drive Mississauga, Ontario

Mississauga, Oni L5T 1C9

Canada phone: 1-800-268-0213 Fax: 1-888-913-9779

- · Information department: sales canada@bio-rad.com
- · Emergency telephone number: 1-800-424-6723 option 2

2 Hazard identification

· Classification of the substance or mixture

Respiratory Sensitizer - Category 1

Skin Sensitizer - Category 1 H317 May cause an allergic skin reaction.

Carcinogenicity - Category 1A H350 May cause cancer.

Reproductive Toxicity - Category 1B H360 May damage fertility or the unborn child.

Specific Target Organ Toxicity - Repeated Exposure - H373 May cause damage to organs through prolonged or

H334 May cause allergy or asthma symptoms or breathing

difficulties if inhaled.

Category 2 repeated exposure.

- · Label elements
- · GHS label elements The product is classified and labeled according to the Globally Harmonized System (GHS).
- · Hazard pictograms



GHS08

- · Signal word Danger
- · Hazard-determining components of labeling: nickel sulphate

· Hazard statements

H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.

H317 May cause an allergic skin reaction.

H350 May cause cancer.

H360 May damage fertility or the unborn child.

H373 May cause damage to organs through prolonged or repeated exposure.

Precautionary statements

P260 Do not breathe dust/fume/gas/mist/vapours/spray.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P284 [In case of inadequate ventilation] wear respiratory protection.

(Contd. on page 2)





Printing date 09/01/2020 Reviewed on 08/11/2020

Trade name: ProteOn™ HTG Chip Activation Solution, 10mM Nickel Sulfate pH 6.0

(Contd. of page 1)

P342+P311 If experiencing respiratory symptoms: Call a poison center/doctor.

P405 Store locked up.

P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

3 Composition/Information on ingredients

- · Chemical characterization: Mixtures
- · Description: Mixture of the substances listed below with non-hazardous additions.
- · Listing of dangerous and non-hazardous components:

·	
7786-81-4 nickel sulphate	0.25-<1% w/w *
Respiratory Sensitizer - Category 1, H334; Germ Cell Mutagenicity - Category 2, H341;	
Carcinogenicity - Category 1A, H350; Reproductive Toxicity - Category 1B, H360; S	pecific
Target Organ Toxicity - Repeated Exposure - Category 1, H372	
Acute Toxicity (Oral) - Category 4, H302; Acute Toxicity (Inhalation) - Category 4, H33.	2; Skin
Irritation - Category 2, H315; Skin Sensitizer - Category 1, H317	

^{*} Actual concentration ranges are withheld as a trade secret.

4 First-aid measures

- · Description of first aid measures
- · General information Immediately remove any clothing soiled by the product.
- · After inhalation

Supply fresh air and to be sure call for a doctor.

In case of unconsciousness place patient stably in side position for transportation.

· After skin contact

Immediately wash with water and soap and rinse thoroughly.

Immediately rinse with water.

- · After eye contact Rinse opened eye for several minutes under running water.
- · After swallowing Rinse mouth with water. Seek medical attention and appropriate follow-up.
- · Information for doctor
- · Most important symptoms and effects, both acute and delayed No further relevant information available.
- Indication of any immediate medical attention and special treatment needed No further relevant information available.

5 Fire-fighting measures

- · Extinguishing media
- · Suitable extinguishing agents Use fire fighting measures that suit the environment.
- · Special hazards arising from the substance or mixture During heating or in case of fire poisonous gases are produced.
- · Advice for firefighters
- · Protective equipment:

Mount respiratory protective device.

Wear self-contained respiratory protective device.

CA

[·] Additional information For the wording of the listed hazard phrases refer to section 16.





Printing date 09/01/2020 Reviewed on 08/11/2020

Trade name: ProteOn™ HTG Chip Activation Solution, 10mM Nickel Sulfate pH 6.0

(Contd. of page 2)

Page 3/8

6 Accidental release measures

· Personal precautions, protective equipment and emergency procedures

Mount respiratory protective device.

Wear protective clothing.

· Environmental precautions:

Do not allow product to reach sewage system or any water course.

Inform respective authorities in case of seepage into water course or sewage system.

Dilute with plenty of water.

Do not allow to enter sewers/ surface or ground water.

· Methods and material for containment and cleaning up:

Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).

Dispose contaminated material as waste according to item 13.

Ensure adequate ventilation.

· Reference to other sections

See Section 7 for information on safe handling

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

7 Handling and storage

- · Handling
- · Precautions for safe handling

Ensure good ventilation/exhaust at the workplace.

Open and handle receptacle with care.

Prevent formation of aerosols.

- · Information about protection against explosions and fires: Keep respiratory protective device available.
- · Conditions for safe storage, including any incompatibilities
- Storage
- · Requirements to be met by storerooms and receptacles: According to product specification
- Information about storage in one common storage facility: Not required.
- · Further information about storage conditions: Keep receptacle tightly sealed.
- · Specific end use(s) No further relevant information available.

8 Exposure controls/ Personal protection

- · Additional information about design of technical systems: No further data; see item 7.
- · Control parameters
- · Components with limit values that require monitoring at the workplace:

7786-81-4 nickel sulphate

EL | TWA: 0.05 mg/m³

as Ni; ACGIH A1, IARC 1

EV TWA: 0.1 mg/m³

Inhalable fraction, as Ni

(Contd. on page 4)





Printing date 09/01/2020 Reviewed on 08/11/2020

Trade name: ProteOn™ HTG Chip Activation Solution, 10mM Nickel Sulfate pH 6.0

(Contd. of page 3)

- · Additional information: The lists that were valid during the creation were used as basis.
- · Exposure controls
- · Personal protective equipment
- · General protective and hygienic measures

Keep away from foodstuffs, beverages and feed.

Immediately remove all soiled and contaminated clothing

Wash hands before breaks and at the end of work.

Store protective clothing separately.

Avoid contact with the skin.

· Breathing equipment:

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use respiratory protective device that is independent of circulating air.

Not necessary if room is well-ventilated. Adequate ventilation required.

- · Protection of hands: Protective gloves.
- · Material of gloves Synthetic gloves
- Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

· Eye protection:

Safety glasses

Lower:

Tightly sealed goggles.

9 Physical and chemical properties

· Information on basic physical and	chemical properties
General Information	
· Appearance:	
Form:	Liquid
Color:	Light green
· Odor:	Odorless
· Odor threshold:	Not determined.
· pH-value at 20 °C:	6
· Change in condition	
Melting point/Melting range:	0 ℃
Boiling point/Boiling range:	undetermined
· Flash point:	Not applicable
·	Not determined.
· Flammability (solid, gaseous)	Not applicable.
· Decomposition temperature:	Not determined.
· Auto igniting:	Product is not selfigniting.
· Danger of explosion:	Not determined.
· Explosion limits:	

Not determined.

(Contd. on page 5)





Printing date 09/01/2020 Reviewed on 08/11/2020

Trade name: ProteOn™ HTG Chip Activation Solution, 10mM Nickel Sulfate pH 6.0

		(Contd. of page 4)
Upper:	Not determined.	
· Vapor pressure at 20 °C:	23 hPa	
· Density at 20 °C:	1 g/cm³	
· Relative density	Not determined.	
· Vapor density	Not determined.	
· Evaporation rate	Not determined.	
· Solubility in / Miscibility with		
Water:	Fully miscible	
· Partition coefficient (n-octanol/w	ater): Not determined.	
· Viscosity:		
dynamic:	Not determined.	
kinematic:	Not determined.	
· Solvent content:		
Water:	99.5 %	
· Other information	No further relevant information available.	

10 Stability and reactivity

- · Reactivity No further relevant information available.
- · Chemical stability
- · Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
- · Possibility of hazardous reactions No dangerous reactions known
- · Conditions to avoid No further relevant information available.
- · Incompatible materials: No further relevant information available.
- · Hazardous decomposition products: No dangerous decomposition products known

11 Toxicological information

- · Information on toxicological effects
- · Acute toxicity:
- Primary irritant effect:
- · on the skin: No irritant effect.
- · on the eye: No irritant effect.
- · Sensitization:

Sensitization possible through inhalation. Sensitization possible through skin contact.

Additional toxicological information:

The product shows the following dangers according to the calculation method of the General EU Classification Guidelines for Preparations as issued in the latest version:

Harmful Irritant

(Contd. on page 6)





Printing date 09/01/2020 Reviewed on 08/11/2020

Trade name: ProteOn™ HTG Chip Activation Solution, 10mM Nickel Sulfate pH 6.0

(Contd. of page 5)

- · Carcinogenic categories
- · IARC (International Agency for Research on Cancer)

7786-81-4 nickel sulphate

1

· NTP (National Toxicology Program)

7786-81-4 nickel sulphate

K

12 Ecological information

- · Toxicity
- Aquatic toxicity: No further relevant information available.
- · Persistence and degradability No further relevant information available.
- · Behavior in environmental systems:
- · Bioaccumulative potential No further relevant information available.
- · Mobility in soil No further relevant information available.
- · Additional ecological information:
- · General notes:

Water hazard class 2 (Self-assessment): hazardous for water.

Do not allow product to reach ground water, water course or sewage system.

Danger to drinking water if even small quantities leak into the ground.

- Results of PBT and vPvB assessment
- · PBT: Not applicable.
- · vPvB: Not applicable.
- · Other adverse effects No further relevant information available.

13 Disposal considerations

- · Waste treatment methods
- · Recommendation

Must not be disposed of together with household garbage. Do not allow product to reach sewage system. Dispose of waste in accordance to applicable national, regional, or local regulations.

- · Uncleaned packagings:
- · Recommendation: Disposal must be made according to official regulations.
- · Recommended cleansing agent: Water, if necessary with cleansing agents.

14 Transport information

- · UN-Number
- · DOT/TDG, ADR, ADN, IMDG, IATA Not Regulated
- · UN proper shipping name
- · DOT/TDG, ADN, IMDG, IATA

Not Regulated · ADR Not Regulated

Not Regulated

(Contd. on page 7)





Printing date 09/01/2020 Reviewed on 08/11/2020

Trade name: ProteOn™ HTG Chip Activation Solution, 10mM Nickel Sulfate pH 6.0

	(Contd. of page
· Transport hazard class(es)	
· DOT, ADR, ADN, IMDG, IATA · Class	Not Regulated
· Packing group · DOT/TDG, ADR, IMDG, IATA	Not Regulated
· Environmental hazards: · Marine pollutant:	No
· Special precautions for user	Not applicable.
 Transport in bulk according to Annex II of MARPOL7 and the IBC Code 	73/78 Not applicable.
· UN "Model Regulation":	Not Regulated

15 Regulatory information

- · Safety, health and environmental regulations/legislation specific for the substance or mixture
- · SARA (Superfund Amendents and Reauthorization Act of 1986 USA)
- Section 302/304 (40CFR355.30 / 40CFR355.40):

None of the ingredients is listed.

Section 313 (40CFR372.65):

None of the ingredients is listed.

· TSCA (Toxic Substances Control Act):

All components have the value ACTIVE.

· MAK (German Maximum Workplace Concentration)

7786-81-4 nickel sulphate

· National regulations

This product has been classified in accordance with the hazard criteria of the Controlled Products Regulations and the MSDS contains all of the information required by those regulations

· Information about limitation of use:

Workers are not allowed to be exposed to the hazardous carcinogenic materials contained in this preparation. Exceptions can be made by the authorities in certain cases.

· Technical instructions (air):

Class	Share in %
II.	0.1-1

- · Water hazard class: Water hazard class 2 (Self-assessment): hazardous for water.
- · Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

__ C/

1



Page 8/8

Safety Data Sheet according to HPR, Schedule 1

Printing date 09/01/2020 Reviewed on 08/11/2020

Trade name: ProteOn™ HTG Chip Activation Solution, 10mM Nickel Sulfate pH 6.0

(Contd. of page 7)

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- · Department issuing SDS: Environmental Health and Safety.
- · Contact:

Life Science Group, Environmental Health and Safety, 2000 Alfred Nobel Drive, Hercules, California, 94547: 1(510) 741-1000

Diagnostic Group, Environmental Health and Safety, 4000 Alfred Nobel Drive, Hercules, California, 94547: 1(510) 724-7000

- Date of the latest revision of the safety data sheet 09/01/2020 / -
- · Abbreviations and acronyms:

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation IATA: International Air Transport Association

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

PBT: Persistent, Bioaccumulative and Toxic vPvB: very Persistent and very Bioaccumulative

C۸





Printing date 09/01/2020 Reviewed on 08/15/2020

1 Identification

- · Product identifier
- · Trade name: ProteOn™ HTG Chip Regeneration Solution, 300 mM EDTA, pH 8.5
- · Catalog or product number: 10021443
- · Application of the substance / the mixture Laboratory chemicals
- · Details of the supplier of the safety data sheet
- · Manufacturer/Supplier:

Bio-Rad Laboratories (Canada) Ltd.

1329 Meyerside Drive Mississauga, Ontario

Mississauga, Ontario L5T 1C9

Canada phone: 1-800-268-0213 Fax: 1-888-913-9779

- · Information department: sales_canada@bio-rad.com
- · Emergency telephone number: 1-800-424-6723 option 2

2 Hazard identification

· Classification of the substance or mixture

Eye Irritation - Category 2A H319 Causes serious eye irritation.

- · Label elements
- · GHS label elements The product is classified and labeled according to the Globally Harmonized System (GHS).
- · Hazard pictograms



GHS07

- · Signal word Warning
- · Hazard statements

H319 Causes serious eye irritation.

· Precautionary statements

P264 Wash thoroughly after handling.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy

to do. Continue rinsing.

P337+P313 If eye irritation persists: Get medical advice/attention.

3 Composition/Information on ingredients

- · Chemical characterization: Mixtures
- · Description: Mixture of the substances listed below with non-hazardous additions.
- · Listing of dangerous and non-hazardous components:

60-00-4 ethylenediaminetetraacetic acid

Eye Irritation - Category 2A, H319

10-13% w/w *

(Contd. on page 2)

^{*} Actual concentration ranges are withheld as a trade secret.





Printing date 09/01/2020 Reviewed on 08/15/2020

Trade name: ProteOn™ HTG Chip Regeneration Solution, 300 mM EDTA, pH 8.5

(Contd. of page 1)

· Additional information For the wording of the listed hazard phrases refer to section 16.

4 First-aid measures

- · Description of first aid measures
- · General information No special measures required.
- · After inhalation Supply fresh air; consult doctor in case of complaints.
- · After skin contact Generally the product does not irritate the skin.
- · After eye contact Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.
- · After swallowing Rinse mouth with water. Seek medical attention and appropriate follow-up.
- · Information for doctor
- · Most important symptoms and effects, both acute and delayed No further relevant information available.
- Indication of any immediate medical attention and special treatment needed No further relevant information available.

5 Fire-fighting measures

- · Extinguishing media
- · Suitable extinguishing agents Use fire fighting measures that suit the environment.
- · Special hazards arising from the substance or mixture No further relevant information available.
- Advice for firefighters
- · Protective equipment: Wear self-contained respiratory protective device.

6 Accidental release measures

- · Personal precautions, protective equipment and emergency procedures Wear protective clothing.
- · Environmental precautions:

Dilute with plenty of water.

Do not allow to enter sewers/ surface or ground water.

· Methods and material for containment and cleaning up:

Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).

· Reference to other sections

No dangerous substances are released.

See Section 7 for information on safe handling

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

7 Handling and storage

- · Handling
- · Precautions for safe handling No special measures required.
- · Information about protection against explosions and fires: No special measures required.
- · Conditions for safe storage, including any incompatibilities
- Storage
- · Requirements to be met by storerooms and receptacles: According to product specification
- · Information about storage in one common storage facility: Not required.

(Contd. on page 3)





Printing date 09/01/2020 Reviewed on 08/15/2020

Trade name: ProteOn™ HTG Chip Regeneration Solution, 300 mM EDTA, pH 8.5

(Contd. of page 2)

- · Further information about storage conditions: Keep receptacle tightly sealed.
- · Specific end use(s) No further relevant information available.

8 Exposure controls/ Personal protection

- · Additional information about design of technical systems: No further data; see item 7.
- · Control parameters
- · Components with limit values that require monitoring at the workplace:

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

- · Additional information: The lists that were valid during the creation were used as basis.
- · Exposure controls
- · Personal protective equipment
- · General protective and hygienic measures

Keep away from foodstuffs, beverages and feed.

Immediately remove all soiled and contaminated clothing

Wash hands before breaks and at the end of work.

Avoid contact with the eyes.

Avoid contact with the eyes and skin.

- · Breathing equipment: Not necessary if room is well-ventilated. Adequate ventilation required.
- · Protection of hands: Protective gloves.
- · Material of gloves Synthetic gloves
- Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

· Eye protection:

Safety glasses

Tightly sealed goggles.

9 Physical and chemical properties

- · Information on basic physical and chemical properties
- · General Information
- · Appearance:

Form: Liquid
Color: Colorless
Odor: Odorless
Odor threshold: Not determined.

• pH-value at 20 °C: 8.5

· Change in condition

Melting point/Melting range: undetermined Boiling point/Boiling range: undetermined

• Flash point: Not applicable Not determined.

(Contd. on page 4)





Printing date 09/01/2020 Reviewed on 08/15/2020

Trade name: ProteOn™ HTG Chip Regeneration Solution, 300 mM EDTA, pH 8.5

	(Contd. of pa
Flammability (solid, gaseous)	Not applicable.
Decomposition temperature:	Not determined.
· Auto igniting:	Product is not selfigniting.
Danger of explosion:	Not determined.
Explosion limits:	
Lower:	Not determined.
Upper:	Not determined.
Vapor pressure at 20 °C:	23 hPa
Density:	Not determined
Relative density	Not determined.
Vapor density	Not determined.
Evaporation rate	Not determined.
Solubility in / Miscibility with	
Water:	Fully miscible
Partition coefficient (n-octanol/wa	ter): Not determined.
Viscosity:	
dynamic:	Not determined.
kinematic:	Not determined.
Solvent content:	
Water:	88.8 %
· Other information	No further relevant information available.

10 Stability and reactivity

- · Reactivity No further relevant information available.
- · Chemical stability
- Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
- · Possibility of hazardous reactions No dangerous reactions known
- · Conditions to avoid No further relevant information available.
- · Incompatible materials: No further relevant information available.
- · Hazardous decomposition products: No dangerous decomposition products known

11 Toxicological information

- · Information on toxicological effects
- · Acute toxicity:
- · LD/LC50 values for hazardous components per OSHA criteria:

60-00-4 ethylenediaminetetraacetic acid

Oral LD50 4,500 mg/kg (rat)

(Contd. on page 5)



Page 5/7

Safety Data Sheet according to HPR, Schedule 1

Printing date 09/01/2020 Reviewed on 08/15/2020

Trade name: ProteOn™ HTG Chip Regeneration Solution, 300 mM EDTA, pH 8.5

(Contd. of page 4)

- · Primary irritant effect:
- · on the skin: Irritant to skin and mucous membranes.
- · on the eye: Irritant effect.
- · Sensitization: No sensitizing effects known.
- · Additional toxicological information:

The product shows the following dangers according to the calculation method of the General EU Classification Guidelines for Preparations as issued in the latest version:

Irritant

- · Carcinogenic categories
- · IARC (International Agency for Research on Cancer)

None of the ingredients is listed.

· NTP (National Toxicology Program)

None of the ingredients is listed.

12 Ecological information

- · Toxicity
- · Aquatic toxicity: No further relevant information available.
- · Persistence and degradability No further relevant information available.
- · Behavior in environmental systems:
- · Bioaccumulative potential No further relevant information available.
- · Mobility in soil No further relevant information available.
- · Additional ecological information:
- · General notes:

Water hazard class 2 (Self-assessment): hazardous for water.

Do not allow product to reach ground water, water course or sewage system.

Danger to drinking water if even small quantities leak into the ground.

- · Results of PBT and vPvB assessment
- · PBT: Not applicable.
- · vPvB: Not applicable.
- · Other adverse effects No further relevant information available.

13 Disposal considerations

- · Waste treatment methods
- · Recommendation

Must not be disposed of together with household garbage. Do not allow product to reach sewage system. Dispose of waste in accordance to applicable national, regional, or local regulations.

- · Uncleaned packagings:
- · Recommendation: Disposal must be made according to official regulations.
- · Recommended cleansing agent: Water, if necessary with cleansing agents.

(Contd. on page 6)





Printing date 09/01/2020 Reviewed on 08/15/2020

Trade name: ProteOn™ HTG Chip Regeneration Solution, 300 mM EDTA, pH 8.5

(Contd. of page 5)

UN-Number		
DOT/TDG, ADR, ADN, IMDG, IATA	Not Regulated	
UN proper shipping name		
DOT/TDG, ADN, IMDG, IATA	Not Regulated	
ADR	Not Regulated	
	Not Regulated	
Transport hazard class(es)		
DOT, ADR, ADN, IMDG, IATA		
Class	Not Regulated	
Packing group		
DOT/TDG, ADR, IMDG, IATA	Not Regulated	
Environmental hazards:		
Marine pollutant:	No	
Special precautions for user	Not applicable.	
Transport in bulk according to Annex II of MA	ARPOL73/78	_
and the IBC Code	Not applicable.	
UN "Model Regulation":	Not Regulated	

15 Regulatory information

- · Safety, health and environmental regulations/legislation specific for the substance or mixture
- · SARA (Superfund Amendents and Reauthorization Act of 1986 USA)
- · Section 302/304 (40CFR355.30 / 40CFR355.40):

None of the ingredients is listed.

· Section 313 (40CFR372.65):

None of the ingredients is listed.

· TSCA (Toxic Substances Control Act):

All components have the value ACTIVE.

· MAK (German Maximum Workplace Concentration)

None of the ingredients is listed.

· National regulations

This product has been classified in accordance with the hazard criteria of the Controlled Products Regulations and the MSDS contains all of the information required by those regulations

- · Water hazard class: Water hazard class 2 (Self-assessment): hazardous for water.
- · Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

CA



Page 7/7

Safety Data Sheet according to HPR, Schedule 1

Printing date 09/01/2020 Reviewed on 08/15/2020

Trade name: ProteOn™ HTG Chip Regeneration Solution, 300 mM EDTA, pH 8.5

(Contd. of page 6)

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- · Department issuing SDS: Environmental Health and Safety.
- · Contact:

Life Science Group, Environmental Health and Safety, 2000 Alfred Nobel Drive, Hercules, California, 94547: 1(510) 741-1000

Diagnostic Group, Environmental Health and Safety, 4000 Alfred Nobel Drive, Hercules, California, 94547: 1(510) 724-7000

- · Date of the latest revision of the safety data sheet 09/01/2020 / -
- · Abbreviations and acronyms:

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation IATA: International Air Transport Association

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

LC50: Lethal concentration, 50 percent LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic vPvB: very Persistent and very Bioaccumulative

^^