

KIT COMPONENTS

Product name: VINEO™ Extract DNA Kit

Product code: **3548100**

Description	Quantity Symbols
VINEO™ Extract DNA Kit - R1 (100 ml)	2
VINEO Extract DNA Kit - R2, 1 x 25 ml	1
VINEO Extract DNA Kit - W1 conc., 1 x 10 ml	1
VINEO Extract DNA Kit - W2, 1 x 100 ml	1 GHS07
	VINEO™ Extract DNA Kit - R1 (100 ml) VINEO Extract DNA Kit - R2, 1 x 25 ml VINEO Extract DNA Kit - W1 conc., 1 x 10 ml



Page 1/7

Safety Data Sheet

according to HPR, Schedule 1

Printing date 11/01/2020

Reviewed on 07/08/2015

1 Identification

- · Product identifier
 - Trade name: VINEO™ Extract DNA Kit R1 (100 ml)
 - · Article number: 4810A
 - · Application of the substance / the mixture Laboratory reagent or component.
- · Details of the supplier of the safety data sheet
 - Manufacturer/Supplier:

Bio-Rad

3 boulevard Raymond-Poincaré 92430 Marnes-la-Coquette Phone: +33 (0)1 47 95 60 00 Fax: +33 (0)1 47 41 91 33

e-mail: fds-msds.fr@bio-rad.com

Bio-Rad laboratories (Canada) Ltd. 2403 Guénette Montreal, Montréal (Québec) H4R 2E9 Phone: (800) 361-1808

Phone: (800) 361-1808 Fax: (514) 334-0872

· Information department:

Service technique clients

Customer support.

· Emergency telephone number: See item 16

2 Hazard identification

· Classification of the substance or mixture

The product is not classified, according to the Globally Harmonized System (GHS).

- · Label elements
 - · GHS label elements Void
 - · Hazard pictograms Void
 - · Signal word Void
 - · Hazard statements Void
- · Classification system:
 - · NFPA ratings (scale 0 4)

Health = 0

Fire = 0

Reactivity = 0

· HMIS-ratings (scale 0 - 4)

Health = 0Fire = 0Reactivity = 0

- CDN



Page 2/7

Safety Data Sheet

according to HPR, Schedule 1

Printing date 11/01/2020 Reviewed on 07/08/2015

Trade name: VINEO™ Extract DNA Kit - R1 (100 ml)

(Contd. of page 1)

3 Composition/Information on ingredients

· Chemical characterization: Mixtures

• **Description:** Mixture containing components at a non hazardous concentration.

· Dangerous components: Void

4 First-aid measures

- · Description of first aid measures
 - · General information: No special measures required.
 - After inhalation: Supply fresh air; consult doctor in case of complaints.
 - · After skin contact: Generally the product does not irritate the skin.
 - · After eye contact: Rinse opened eye for several minutes under running water.
 - After swallowing: Consult a doctor, if necessary.
 - · Information for doctor:
 - · Most important symptoms and effects, both acute and delayed No further relevant information available.
 - · Indication of any immediate medical attention and special treatment needed No further relevant information available.

5 Fire-fighting measures

- · Extinguishing media
 - Suitable extinguishing agents:

CO2, extinguishing powder or water spray. Fight larger fires with water spray.

- · Special hazards arising from the substance or mixture No further relevant information available.
- · Advice for firefighters
 - · Protective equipment: No special measures required.

6 Accidental release measures

- · Personal precautions, protective equipment and emergency procedures Wear disposable gloves.
- · Environmental precautions:

Dilute with plenty of water.

Do not allow to enter sewers/ surface or ground water.

- Methods and material for containment and cleaning up:
- Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- · Reference to other sections No dangerous substances are released.

CDN

(Contd. on page 3)





according to HPR, Schedule 1

Printing date 11/01/2020 Reviewed on 07/08/2015

Trade name: VINEO™ Extract DNA Kit - R1 (100 ml)

(Contd. of page 2)

7 Handling and storage

- · Handling:
 - Precautions for safe handling

Keep receptacles tightly sealed.

Avoid all contact with reagents, even those considered to be non-hazardous.

Follow good laboratory practices.

- Information about protection against explosions and fires: No special measures required.
- · Conditions for safe storage, including any incompatibilities
 - · Storage:
 - · Requirements to be met by storerooms and receptacles: No special requirements.
 - · Information about storage in one common storage facility: Not required.
 - Further information about storage conditions: None.
- · Specific end use(s) No further relevant information available.

8 Exposure controls/ Personal protection

- Additional information about design of technical systems: No further data; see item 7.
- · Control parameters
 - · Components with limit values that require monitoring at the workplace:

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

- · Additional information: The lists that were valid during the creation were used as basis.
- · Personal protective equipment:
 - · General protective and hygienic measures:

The usual precautionary measures for handling chemicals should be followed.

- Breathing equipment: Suitable respiratory protective device recommended.
- · Protection of hands:

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

· Eye protection: Goggles recommended during refilling.

9 Physical and chemical properties

- · Information on basic physical and chemical properties
 - · General Information
 - · Appearance:

Form: FluidColor: ColorlessOdor: Odorless

- · Change in condition
 - Melting point/Melting range: Undetermined.
 Boiling point/Boiling range: Undetermined.

(Contd. on page 4)





according to HPR, Schedule 1

Printing date 11/01/2020 Reviewed on 07/08/2015

Trade name: VINEO™ Extract DNA Kit - R1 (100 ml)

		(Contd. of page 3
· Flash point:	Not applicable.	
· Auto igniting:	Product is not selfigniting.	
· Danger of explosion:	Product does not present an explosion hazard.	
· Density:	Not determined.	
· Solubility in / Miscibility	/ with	
· Water:	Fully miscible.	
Other information	No further relevant information available.	

10 Stability and reactivity

- · Reactivity No further relevant information available.
 - · Chemical stability
 - Thermal decomposition / conditions to be avoided:

No decomposition if used according to specifications.

- · Possibility of hazardous reactions No dangerous reactions known.
- · Conditions to avoid No further relevant information available.
- · Incompatible materials: No further relevant information available.
- · Hazardous decomposition products: No dangerous decomposition products known.

11 Toxicological information

- · Information on toxicological effects
 - · Acute toxicity:
 - · Primary irritant effect:
 - on the skin: No irritant effect.
 - · on the eye: No irritating effect.
 - · Sensitization: No sensitizing effects known.
 - · Additional toxicological information:

The product is not subject to classification according to internally approved calculation methods for preparations:

When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us.

· Carcinogenic categories

IARC (International Agency for Research on Cancer)

None of the ingredients is listed.

· NTP (National Toxicology Program)

None of the ingredients is listed.

CDN





according to HPR, Schedule 1

Printing date 11/01/2020 Reviewed on 07/08/2015

Trade name: VINEO™ Extract DNA Kit - R1 (100 ml)

(Contd. of page 4)

12 Ecological information

- · Toxicity
 - · Aquatic toxicity: No further relevant information available.
- · Persistence and degradability No further relevant information available.
- Behavior in environmental systems:
 - · Bioaccumulative potential No further relevant information available.
 - · Mobility in soil No further relevant information available.
- Additional ecological information:
 - · General notes:

Water hazard class 1 (Self-assessment): slightly hazardous for water

Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

- Results of PBT and vPvB assessment
 - · **PBT**: Not applicable.
 - · vPvB: Not applicable.
- · Other adverse effects No further relevant information available.

13 Disposal considerations

- · Waste treatment methods
 - · Recommendation:

Must not be disposed of together with household garbage. Do not allow product to reach sewage system.

- · Uncleaned packagings:
 - · Recommendation: Disposal must be made according to official regulations.
 - · Recommended cleansing agent: Water, if necessary with cleansing agents.

UN-Number		
DOT/TDG, ADR, ADN, IMDG, IATA	Void	
UN proper shipping name		
DOT/TDG, ADR, ADN, IMDG, IATA	Void	
Transport hazard class(es)		
· DOT, ADR, ADN, IMDG, IATA		
Class	Void	
Packing group		
DOT/TDG, ADR, IMDG, IATA	Void	
Environmental hazards:		
· Marine pollutant:	No	

(Contd. on page 6)





according to HPR, Schedule 1

Printing date 11/01/2020 Reviewed on 07/08/2015

Trade name: VINEO™ Extract DNA Kit - R1 (100 ml)

	(Contd. of page
· Special precautions for user	Not applicable.
Transport in bulk according to Annex MARPOL73/78 and the IBC Code	II of Not applicable.
· Transport/Additional information:	Not dangerous according to the above specifications.
· UN "Model Regulation":	Void

15 Regulatory information

Safety, health and environmental regulations/legislation specific for the substance or mixture Sara

· Section 355 (extremely hazardous substances):

None of the ingredient is listed.

· Section 313 (Specific toxic chemical listings):

None of the ingredients is listed.

TSCA (Toxic Substances Control Act):

7732-18-5 water, distilled, conductivity or of similar purity

· Canadian substance listings:

· Canadian Domestic Substances List (DSL)

7732-18-5 water, distilled, conductivity or of similar purity

Canadian Ingredient Disclosure list (limit 0.1%)

None of the ingredients is listed.

· Canadian Ingredient Disclosure list (limit 1%)

None of the ingredients is listed.

· Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· Department issuing SDS:

MSDS / FDS Administr.

Bio-Rad

Marnes-la-Coquette, France

· Contact:

SOFRALAB SAS

79 av. A.A. Thévenet, BP1031 Magenta,

51319 EPERNAY Cedex Mob.: +33 6 25 10 00 26

e-mail: nsieczkowski@sofralab.com

(Contd. on page 7)



Page 7/7

Safety Data Sheet according to HPR, Schedule 1

Printing date 11/01/2020 Reviewed on 07/08/2015

Trade name: VINEO™ Extract DNA Kit - R1 (100 ml)

(Contd. of page 6)

Customer Support

Phone: +33 (0)1 47 95 62 77

Customer Support

Phone: +33 (0)1 47 95 61 85

Abbreviations and acronyms:

ICAO: International Civil Aviation Organisation

IMDG: International Maritime Code for Dangerous Goods DOT: US Department of Transportation

IATA: International Air Transport Association

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)
NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

PBT: Persistent, Bioaccumulative and Toxic vPvB: very Persistent and very Bioaccumulative



Page 1/7

Safety Data Sheet

according to HPR, Schedule 1

Printing date 11/01/2020

Reviewed on 07/08/2015

1 Identification

- · Product identifier
 - · Trade name: VINEO Extract DNA Kit R2, 1 x 25 ml
 - · Article number: 4810B
 - · Application of the substance / the mixture

Laboratory reagent or component.

In vitro diagnostic medical device or component.

- · Details of the supplier of the safety data sheet
 - Manufacturer/Supplier:

Bio-Rad

3 boulevard Raymond-Poincaré 92430 Marnes-la-Coquette Phone: +33 (0)1 47 95 60 00 Fax: +33 (0)1 47 41 91 33 e-mail: fds-msds.fr@bio-rad.com

Bio-Rad laboratories (Canada) Ltd. 2403 Guénette Montreal, Montréal (Québec) H4R 2E9

Phone: (800) 361-1808 Fax: (514) 334-0872

· Information department:

Service technique clients Customer support.

· Emergency telephone number: See item 16

2 Hazard identification

· Classification of the substance or mixture

The product is not classified, according to the Globally Harmonized System (GHS).

- · Label elements
 - · GHS label elements Void
 - · Hazard pictograms Void
 - · Signal word Void
 - · Hazard statements Void
- · Classification system:
 - · NFPA ratings (scale 0 4)

Health = 0

Fire = 0

Reactivity = 0

HMIS-ratings (scale 0 - 4)

Health = 0Fire = 0Reactivity = 0

- CDN



Page 2/7

Safety Data Sheet

according to HPR, Schedule 1

Printing date 11/01/2020 Reviewed on 07/08/2015

Trade name: VINEO Extract DNA Kit - R2, 1 x 25 ml

(Contd. of page 1)

3 Composition/Information on ingredients

· Chemical characterization: Mixtures

Description: Mixture containing components at a non hazardous concentration.

· Dangerous components: Void

4 First-aid measures

- · Description of first aid measures
 - General information: No special measures required.
 - After inhalation: Supply fresh air; consult doctor in case of complaints.
 - · After skin contact: Generally the product does not irritate the skin.
 - · After eye contact: Rinse opened eye for several minutes under running water.
 - After swallowing: Consult a doctor, if necessary.
 - · Information for doctor:
 - · Most important symptoms and effects, both acute and delayed No further relevant information available.
 - · Indication of any immediate medical attention and special treatment needed No further relevant information available.

5 Fire-fighting measures

- · Extinguishing media
 - Suitable extinguishing agents:

CO2, extinguishing powder or water spray. Fight larger fires with water spray.

- · Special hazards arising from the substance or mixture No further relevant information available.
- · Advice for firefighters
 - · Protective equipment: No special measures required.

6 Accidental release measures

- · Personal precautions, protective equipment and emergency procedures Wear disposable gloves.
- · Environmental precautions:

Dilute with plenty of water.

Do not allow to enter sewers/ surface or ground water.

- Methods and material for containment and cleaning up:
- Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- · Reference to other sections No dangerous substances are released.

CDN

(Contd. on page 3)





according to HPR, Schedule 1

Printing date 11/01/2020 Reviewed on 07/08/2015

Trade name: VINEO Extract DNA Kit - R2, 1 x 25 ml

(Contd. of page 2)

7 Handling and storage

- · Handling:
 - · Precautions for safe handling

Keep receptacles tightly sealed.

Avoid all contact with reagents, even those considered to be non-hazardous.

Follow good laboratory practices.

- Information about protection against explosions and fires: No special measures required.
- · Conditions for safe storage, including any incompatibilities
 - · Storage:
 - · Requirements to be met by storerooms and receptacles: No special requirements.
 - · Information about storage in one common storage facility: Not required.
 - Further information about storage conditions: None.
- · Specific end use(s) No further relevant information available.

8 Exposure controls/ Personal protection

- Additional information about design of technical systems: No further data; see item 7.
- · Control parameters
 - Components with limit values that require monitoring at the workplace:

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

- · Additional information: The lists that were valid during the creation were used as basis.
- · Personal protective equipment:
 - · General protective and hygienic measures:

The usual precautionary measures for handling chemicals should be followed.

- · Breathing equipment: Suitable respiratory protective device recommended.
- · Protection of hands:

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

• Eye protection: Goggles recommended during refilling.

9 Physical and chemical properties

- · Information on basic physical and chemical properties
 - · General Information
 - · Appearance:

Form: FluidColor: ColorlessOdor: Odorless

- · Change in condition
 - Melting point/Melting range: Undetermined.
 Boiling point/Boiling range: Undetermined.

(Contd. on page 4)





according to HPR, Schedule 1

Printing date 11/01/2020 Reviewed on 07/08/2015

Trade name: VINEO Extract DNA Kit - R2, 1 x 25 ml

		(Contd. of page 3)
· Flash point:	Not applicable.	
· Auto igniting:	Product is not selfigniting.	
· Danger of explosion:	Product does not present an explosion hazard.	
· Density:	Not determined.	
· Solubility in / Miscibility	with	
· Water:	Fully miscible.	
· Other information	No further relevant information available.	

10 Stability and reactivity

- · Reactivity No further relevant information available.
 - · Chemical stability
 - Thermal decomposition / conditions to be avoided:

No decomposition if used according to specifications.

- · Possibility of hazardous reactions No dangerous reactions known.
- · Conditions to avoid No further relevant information available.
- · Incompatible materials: No further relevant information available.
- · Hazardous decomposition products: No dangerous decomposition products known.

11 Toxicological information

- · Information on toxicological effects
 - · Acute toxicity:
 - · Primary irritant effect:
 - on the skin: No irritant effect.
 - · on the eye: No irritating effect.
 - · Sensitization: No sensitizing effects known.
 - · Additional toxicological information:

The product is not subject to classification according to internally approved calculation methods for preparations:

When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us.

Carcinogenic categories

· IARC (International Agency for Research on Cancer)

None of the ingredients is listed.

· NTP (National Toxicology Program)

None of the ingredients is listed.

CDN -





according to HPR, Schedule 1

Printing date 11/01/2020 Reviewed on 07/08/2015

Trade name: VINEO Extract DNA Kit - R2, 1 x 25 ml

(Contd. of page 4)

12 Ecological information

- · Toxicity
 - · Aquatic toxicity: No further relevant information available.
- · Persistence and degradability No further relevant information available.
- Behavior in environmental systems:
 - · Bioaccumulative potential No further relevant information available.
 - · Mobility in soil No further relevant information available.
- Additional ecological information:
 - · General notes:

Water hazard class 1 (Self-assessment): slightly hazardous for water

Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

- Results of PBT and vPvB assessment
 - · **PBT**: Not applicable.
 - · vPvB: Not applicable.
- · Other adverse effects No further relevant information available.

13 Disposal considerations

- · Waste treatment methods
 - · Recommendation:

Must not be disposed of together with household garbage. Do not allow product to reach sewage system.

- · Uncleaned packagings:
 - · Recommendation: Disposal must be made according to official regulations.
 - · Recommended cleansing agent: Water, if necessary with cleansing agents.

UN-Number		
	Void	
· DOT/TDG, ADR, ADN, IMDG, IATA	Volu	
UN proper shipping name		
DOT/TDG, ADR, ADN, IMDG, IATA	Void	
Transport hazard class(es)		
· DOT, ADR, ADN, IMDG, IATA		
Class	Void	
Packing group		
DOT/TDG, ADR, IMDG, IATA	Void	
Environmental hazards:		
· Marine pollutant:	No	

(Contd. on page 6)





according to HPR, Schedule 1

Printing date 11/01/2020 Reviewed on 07/08/2015

Trade name: VINEO Extract DNA Kit - R2, 1 x 25 ml

	(Contd. of page
· Special precautions for user	Not applicable.
Transport in bulk according to Annex MARPOL73/78 and the IBC Code	II of Not applicable.
· Transport/Additional information:	Not dangerous according to the above specifications.
· UN "Model Regulation":	Void

15 Regulatory information

Safety, health and environmental regulations/legislation specific for the substance or mixture Sara

· Section 355 (extremely hazardous substances):

None of the ingredient is listed.

· Section 313 (Specific toxic chemical listings):

None of the ingredients is listed.

TSCA (Toxic Substances Control Act):

7732-18-5 water, distilled, conductivity or of similar purity

· Canadian substance listings:

· Canadian Domestic Substances List (DSL)

7732-18-5 water, distilled, conductivity or of similar purity

Canadian Ingredient Disclosure list (limit 0.1%)

None of the ingredients is listed.

· Canadian Ingredient Disclosure list (limit 1%)

None of the ingredients is listed.

· Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· Department issuing SDS:

MSDS / FDS Administr.

Bio-Rad

Marnes-la-Coquette, France

· Contact:

SOFRALAB SAS

79 av. A.A. Thévenet, BP1031 Magenta,

51319 EPERNAY Cedex Mob.: +33 6 25 10 00 26

e-mail: nsieczkowski@sofralab.com

(Contd. on page 7)



Page 7/7

Safety Data Sheet according to HPR, Schedule 1

Printing date 11/01/2020 Reviewed on 07/08/2015

Trade name: VINEO Extract DNA Kit - R2, 1 x 25 ml

(Contd. of page 6)

Customer Support

Phone: +33 (0)1 47 95 62 77

Customer Support

Phone: +33 (0)1 47 95 61 85

Abbreviations and acronyms:

ICAO: International Civil Aviation Organisation

IMDG: International Maritime Code for Dangerous Goods DOT: US Department of Transportation

IATA: International Air Transport Association

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)
NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

PBT: Persistent, Bioaccumulative and Toxic vPvB: very Persistent and very Bioaccumulative



Page 1/7

Safety Data Sheet according to HPR, Schedule 1

Printing date 11/01/2020

Reviewed on 07/08/2015

1 Identification

- · Product identifier
 - · Trade name: VINEO Extract DNA Kit W1 conc., 1 x 10 ml
 - · Article number: 4810C
 - · Application of the substance / the mixture

Laboratory reagent or component.

In vitro diagnostic medical device or component.

- · Details of the supplier of the safety data sheet
 - Manufacturer/Supplier:

Bio-Rad

3 boulevard Raymond-Poincaré 92430 Marnes-la-Coquette Phone: +33 (0)1 47 95 60 00 Fax: +33 (0)1 47 41 91 33

e-mail: fds-msds.fr@bio-rad.com

Bio-Rad laboratories (Canada) Ltd. 2403 Guénette Montreal,

Montréal (Québec) H4R 2E9 Phone: (800) 361-1808

Phone: (800) 361-1808 Fax: (514) 334-0872

· Information department:

Service technique clients

Customer support.

· Emergency telephone number: See item 16

2 Hazard identification

· Classification of the substance or mixture

The product is not classified, according to the Globally Harmonized System (GHS).

- · Label elements
 - · GHS label elements Void
 - · Hazard pictograms Void
 - · Signal word Void
 - · Hazard statements Void
- · Classification system:
 - · NFPA ratings (scale 0 4)

Health = 0

Fire = 0

Reactivity = 0

· HMIS-ratings (scale 0 - 4)

Health = 0Fire = 0Reactivity = 0

- CDN



according to HPR, Schedule 1

Page 2/7
Safety Data Sheet

Printing date 11/01/2020

Reviewed on 07/08/2015

Trade name: VINEO Extract DNA Kit - W1 conc., 1 x 10 ml

(Contd. of page 1)

3 Composition/Information on ingredients

- · Chemical characterization: Mixtures
 - · Description: Mixture without a chemical hazard.
 - · Dangerous components:

Tris(hydroxymethyl)-aminomethane

3-7% w/w *

Skin Irritation - Category 2, H315; Eye Irritation - Category 2A, H319; Specific Target Organ Toxicity - Single Exposure - Category 3, H335

* Actual concentration ranges are withheld as a trade secret.

4 First-aid measures

- · Description of first aid measures
 - · General information: No special measures required.
 - · After inhalation: Supply fresh air; consult doctor in case of complaints.
 - · After skin contact: Generally the product does not irritate the skin.
 - · After eye contact: Rinse opened eye for several minutes under running water.
 - After swallowing: Consult a doctor, if necessary.
 - · Information for doctor:
 - · Most important symptoms and effects, both acute and delayed No further relevant information available.
 - · Indication of any immediate medical attention and special treatment needed No further relevant information available.

5 Fire-fighting measures

- · Extinguishing media
 - Suitable extinguishing agents:

CO2, extinguishing powder or water spray. Fight larger fires with water spray.

- · Special hazards arising from the substance or mixture No further relevant information available.
- Advice for firefighters
 - · Protective equipment: No special measures required.

6 Accidental release measures

- · Personal precautions, protective equipment and emergency procedures Wear disposable gloves.
- · Environmental precautions:

Dilute with plenty of water.

Do not allow to enter sewers/ surface or ground water.

· Methods and material for containment and cleaning up:

Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).

· Reference to other sections No dangerous substances are released.

CDN





according to HPR, Schedule 1

Printing date 11/01/2020 Reviewed on 07/08/2015

Trade name: VINEO Extract DNA Kit - W1 conc., 1 x 10 ml

(Contd. of page 2)

7 Handling and storage

- · Handling:
 - · Precautions for safe handling

Keep receptacles tightly sealed.

Avoid all contact with reagents, even those considered to be non-hazardous.

Follow good laboratory practices.

- Information about protection against explosions and fires: No special measures required.
- · Conditions for safe storage, including any incompatibilities
 - · Storage:
 - · Requirements to be met by storerooms and receptacles: No special requirements.
 - · Information about storage in one common storage facility: Not required.
 - Further information about storage conditions: None.
- · Specific end use(s) No further relevant information available.

8 Exposure controls/ Personal protection

- Additional information about design of technical systems: No further data; see item 7.
- · Control parameters
 - · Components with limit values that require monitoring at the workplace:

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

- · Additional information: The lists that were valid during the creation were used as basis.
- · Personal protective equipment:
 - · General protective and hygienic measures:

The usual precautionary measures for handling chemicals should be followed.

- · Breathing equipment: Suitable respiratory protective device recommended.
- · Protection of hands:

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

· Eye protection: Goggles recommended during refilling.

9 Physical and chemical properties

- · Information on basic physical and chemical properties
 - · General Information
 - · Appearance:

Form: Fluid
 Color: Colorless
 Odor: Odorless

- · Change in condition
 - Melting point/Melting range: Undetermined.
 Boiling point/Boiling range: Undetermined.

(Contd. on page 4)





according to HPR, Schedule 1

Printing date 11/01/2020 Reviewed on 07/08/2015

Trade name: VINEO Extract DNA Kit - W1 conc., 1 x 10 ml

		(Contd. of page 3
· Flash point:	Not applicable.	
· Auto igniting:	Product is not selfigniting.	
· Danger of explosion:	Product does not present an explosion hazard.	
· Density:	Not determined.	
· Solubility in / Miscibilit	y with	
· Water:	Fully miscible.	
· Other information	No further relevant information available.	

10 Stability and reactivity

- · Reactivity No further relevant information available.
 - · Chemical stability
 - Thermal decomposition / conditions to be avoided:

No decomposition if used according to specifications.

- · Possibility of hazardous reactions No dangerous reactions known.
- · Conditions to avoid No further relevant information available.
- · Incompatible materials: No further relevant information available.
- · Hazardous decomposition products: No dangerous decomposition products known.

11 Toxicological information

- · Information on toxicological effects
 - · Acute toxicity:
 - · Primary irritant effect:
 - on the skin: No irritant effect.
 - · on the eye: No irritating effect.
 - · Sensitization: No sensitizing effects known.
 - · Additional toxicological information:

The product is not subject to classification according to internally approved calculation methods for preparations:

When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us.

· Carcinogenic categories

· IARC (International Agency for Research on Cancer)

None of the ingredients is listed.

· NTP (National Toxicology Program)

None of the ingredients is listed.

- CDN





according to HPR, Schedule 1

Printing date 11/01/2020 Reviewed on 07/08/2015

Trade name: VINEO Extract DNA Kit - W1 conc., 1 x 10 ml

(Contd. of page 4)

12 Ecological information

- · Toxicity
 - · Aquatic toxicity: No further relevant information available.
- · Persistence and degradability No further relevant information available.
- Behavior in environmental systems:
 - · Bioaccumulative potential No further relevant information available.
 - · Mobility in soil No further relevant information available.
- Additional ecological information:
 - · General notes:

Water hazard class 1 (Self-assessment): slightly hazardous for water

Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

- Results of PBT and vPvB assessment
 - · **PBT**: Not applicable.
 - · vPvB: Not applicable.
- · Other adverse effects No further relevant information available.

13 Disposal considerations

- · Waste treatment methods
 - · Recommendation:

Must not be disposed of together with household garbage. Do not allow product to reach sewage system.

- · Uncleaned packagings:
 - · Recommendation: Disposal must be made according to official regulations.
 - · Recommended cleansing agent: Water, if necessary with cleansing agents.

UN-Number		
· DOT/TDG, ADR, ADN, IMDG, IATA	Void	
UN proper shipping name		
· DOT/TDG, ADR, ADN, IMDG, IATA	Void	
Transport hazard class(es)		
· DOT, ADR, ADN, IMDG, IATA		
Class	Void	
Packing group		
DOT/TDG, ADR, IMDG, IATA	Void	
Environmental hazards:		
· Marine pollutant:	No	

(Contd. on page 6)





according to HPR, Schedule 1

Printing date 11/01/2020 Reviewed on 07/08/2015

Trade name: VINEO Extract DNA Kit - W1 conc., 1 x 10 ml

	(Contd. of page
Special precautions for user	Not applicable.
Transport in bulk according to Annex MARPOL73/78 and the IBC Code	II of Not applicable.
Transport/Additional information:	Not dangerous according to the above specifications.
UN "Model Regulation":	Void

15 Regulatory information

· Safety, health and environmental regulations/legislation specific for the substance or mixture · Sara

Section 355 (extremely hazardous substances):

None of the ingredient is listed.

· Section 313 (Specific toxic chemical listings):

None of the ingredients is listed.

· TSCA (Toxic Substances Control Act):

Tris(hydroxymethyl)-aminomethane

1185-53-1 2-amino-2-(hydroxymethyl)propane-1,3-diolhydrochloride

7732-18-5 water, distilled, conductivity or of similar purity

· Canadian substance listings:

· Canadian Domestic Substances List (DSL)

All ingredients are listed.

· Canadian Ingredient Disclosure list (limit 0.1%)

None of the ingredients is listed.

Canadian Ingredient Disclosure list (limit 1%)

None of the ingredients is listed.

· Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

Department issuing SDS:

MSDS / FDS Administr.

Bio-Rad

Marnes-la-Coquette, France

· Contact:

SOFRALAB SAS

79 av. A.A. Thévenet, BP1031 Magenta,

51319 EPERNAY Cedex

Mob.: +33 6 25 10 00 26

(Contd. on page 7)



Page 7/7

Safety Data Sheet

according to HPR, Schedule 1

Printing date 11/01/2020 Reviewed on 07/08/2015

Trade name: VINEO Extract DNA Kit - W1 conc., 1 x 10 ml

(Contd. of page 6)

e-mail: nsieczkowski@sofralab.com

Customer Support

Phone: +33 (0)1 47 95 62 77

Customer Support

Phone: +33 (0)1 47 95 61 85

· Abbreviations and acronyms:

ICAO: International Civil Aviation Organisation

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation IATA: International Air Transport Association

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA) HMIS: Hazardous Materials Identification System (USA)

PBT: Persistent, Bioaccumulative and Toxic vPvB: very Persistent and very Bioaccumulative

CDN-



Safety Data Sheet according to HPR, Schedule 1

Reviewed on 07/08/2015 Printing date 11/01/2020

1 Identification

- · Product identifier
 - · Trade name: VINEO Extract DNA Kit W2, 1 x 100 ml
 - · Article number: 4810D
 - · Application of the substance / the mixture Laboratory reagent or component.
- · Details of the supplier of the safety data sheet
 - Manufacturer/Supplier:

Bio-Rad

3 boulevard Raymond-Poincaré 92430 Marnes-la-Coquette Phone: +33 (0)1 47 95 60 00 Fax: +33 (0)1 47 41 91 33 e-mail: fds-msds.fr@bio-rad.com

Bio-Rad laboratories (Canada) Ltd. 2403 Guénette Montreal. Montréal (Québec) H4R 2E9 Phone: (800) 361-1808

Fax: (514) 334-0872

Information department:

Service technique clients

Customer support.

· Emergency telephone number: See item 16

2 Hazard identification

· Classification of the substance or mixture



Skin Irritation - Category 2 H315 Causes skin irritation.

Eye Irritation - Category 2A H319 Causes serious eye irritation.

- · Label elements
 - · GHS label elements

The product is classified and labeled according to the Globally Harmonized System (GHS).

- Hazard pictograms GHS07
- Signal word Warning
- · Hazard statements

Causes skin irritation.

Causes serious eye irritation.

· Precautionary statements

Wash thoroughly after handling.

Wear protective gloves/protective clothing/eye protection/face protection.

If on skin: Wash with plenty of soap and water.

(Contd. on page 2)

Page 1/8





according to HPR, Schedule 1

Printing date 11/01/2020 Reviewed on 07/08/2015

Trade name: VINEO Extract DNA Kit - W2, 1 x 100 ml

(Contd. of page 1)

Specific treatment (see on this label).

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Take off contaminated clothing and wash it before reuse.

If skin irritation occurs: Get medical advice/attention.

If eye irritation persists: Get medical advice/attention.

· Classification system:

NFPA ratings (scale 0 - 4)

Health = 2

Fire = 0

Reactivity = 0

HMIS-ratings (scale 0 - 4)

Health = 2

Fire = 0

Reactivity = 0

3 Composition/Information on ingredients

- · Chemical characterization: Mixtures
 - · **Description:** Mixture containing hazardous components.

· Dangero	us components:	
64-19-7	acetic acid	7-13% w/w *
	 ♦ Flammable Liquids - Category 3, H226 ♦ Skin Corrosion - Category 1A, H314 ♦ Acute Toxicity (Dermal) – Category 4, H312 	
55965-84-9	reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1) Acute Toxicity (Oral) - Category 3, H301; Acute Toxicity (Dermal) – Category 2, H310; Acute Toxicity (Inhalation) - Category 2, H330 Skin Corrosion - Category 1C, H314; Serious Eye Damage - Category 1, H318 Aquatic Acute 1, H400; Aquatic Chronic 1, H410	≤0.1% w/w
* Actual o	Skin Sensitizer - Category 1A, H317 concentration ranges are withheld as a trade secret.	

4 First-aid measures

- · Description of first aid measures
 - · After inhalation: Supply fresh air and to be sure call for a doctor.
 - · After skin contact: Immediately wash with water and soap and rinse thoroughly.
 - · After eye contact:

Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.

· After swallowing: Immediately call a doctor.

(Contd. on page 3)



Page 3/8

Safety Data Sheet according to HPR, Schedule 1

Printing date 11/01/2020 Reviewed on 07/08/2015

Trade name: VINEO Extract DNA Kit - W2, 1 x 100 ml

(Contd. of page 2)

- · Information for doctor:
 - · Most important symptoms and effects, both acute and delayed No further relevant information available.
 - · Indication of any immediate medical attention and special treatment needed No further relevant information available.

5 Fire-fighting measures

- · Extinguishing media
 - · Suitable extinguishing agents:

Use fire fighting measures that suit the environment.

CO2, extinguishing powder or water spray. Fight larger fires with water spray.

- · Special hazards arising from the substance or mixture No further relevant information available.
- · Advice for firefighters
 - · Protective equipment: No special measures required.

6 Accidental release measures

 \cdot Personal precautions, protective equipment and emergency procedures

Wear protective clothing.

Wear disposable gloves.

- · Environmental precautions: Do not allow product to reach sewage system or any water course.
- Methods and material for containment and cleaning up:

Absorb liquid components with liquid-binding material.

Ensure adequate ventilation.

· Reference to other sections

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

7 Handling and storage

- · Handling:
 - · Precautions for safe handling

Ensure good ventilation/exhaustion at the workplace.

Prevent formation of aerosols.

Follow good laboratory practices.

- · Information about protection against explosions and fires: No special measures required.
- · Conditions for safe storage, including any incompatibilities
 - · Storage:
 - · Requirements to be met by storerooms and receptacles: No special requirements.
 - · Information about storage in one common storage facility: Not required.
 - Further information about storage conditions: Keep receptacle tightly sealed.

(Contd. on page 4)





according to HPR, Schedule 1

Printing date 11/01/2020 Reviewed on 07/08/2015

Trade name: VINEO Extract DNA Kit - W2, 1 x 100 ml

· Specific end use(s) No further relevant information available.

(Contd. of page 3)

8 Exposure controls/ Personal protection

- · Additional information about design of technical systems: No further data; see item 7.
- · Control parameters
 - Components with limit values that require monitoring at the workplace:

64-19-7 acetic acid

EL Short-term value: 15 ppm Long-term value: 10 ppm

- · Additional information: The lists that were valid during the creation were used as basis.
- · Personal protective equipment:
 - General protective and hygienic measures:

Immediately remove all soiled and contaminated clothing.

Wash hands before breaks and at the end of work.

Avoid contact with the eyes and skin.

- Breathing equipment: Suitable respiratory protective device recommended.
- · Protection of hands: Protective gloves
- · Eye protection: Safety glasses

9 Physical and chemical properties

- · Information on basic physical and chemical properties
 - · General Information
 - · Appearance:

· Form: Fluid
· Color: Colorless
· Odor: Characteristic

· Change in condition

Melting point/Melting range: Undetermined.
 Boiling point/Boiling range: Undetermined.

· Flash point: Not applicable.

• **Danger of explosion:** Product does not present an explosion hazard.

Density: Not determined.

· Solubility in / Miscibility with

· Water: Fully miscible.

• Other information No further relevant information available.

- CDN

(Contd. on page 5)



Page 5/8

Safety Data Sheet

according to HPR, Schedule 1

Printing date 11/01/2020 Reviewed on 07/08/2015

Trade name: VINEO Extract DNA Kit - W2, 1 x 100 ml

(Contd. of page 4)

10 Stability and reactivity

- · Reactivity No further relevant information available.
 - · Chemical stability
 - · Thermal decomposition / conditions to be avoided:

No decomposition if used according to specifications.

- · Possibility of hazardous reactions No dangerous reactions known.
- · Conditions to avoid No further relevant information available.
- · Incompatible materials: No further relevant information available.
- · Hazardous decomposition products: No dangerous decomposition products known.

11 Toxicological information

- · Information on toxicological effects
 - · Acute toxicity:
 - · Primary irritant effect:
 - · on the skin: Irritant to skin and mucous membranes.
 - · on the eye: Irritating effect.
 - Additional toxicological information:

The product shows the following dangers according to internally approved calculation methods for preparations:

Irritant

- · Carcinogenic categories
 - · IARC (International Agency for Research on Cancer)

None of the ingredients is listed.

· NTP (National Toxicology Program)

None of the ingredients is listed.

12 Ecological information

- · Toxicity
 - · Aquatic toxicity: No further relevant information available.
- $\cdot \mbox{\bf Persistence and degradability} \ \mbox{No further relevant information available}.$
- Behavior in environmental systems:
 - · Bioaccumulative potential No further relevant information available.
 - · Mobility in soil No further relevant information available.
- Additional ecological information:
 - · General notes: Avoid transfer into the environment.
- · Results of PBT and vPvB assessment
 - · PBT: Not applicable.
 - · vPvB: Not applicable.

(Contd. on page 6)





according to HPR, Schedule 1

Printing date 11/01/2020 Reviewed on 07/08/2015

Trade name: VINEO Extract DNA Kit - W2, 1 x 100 ml

· Other adverse effects No further relevant information available.

(Contd. of page 5)

13 Disposal considerations

- · Waste treatment methods
 - · **Recommendation:** Must be specially treated adhering to official regulations.
- · Uncleaned packagings:
 - · **Recommendation:** Disposal must be made according to official regulations.

	14 T	ransi	port	infori	mation
--	------	-------	------	--------	--------

· UN-Number

DOT/TDG
 DOT/TDG
 ADR, IMDG, IATA
 UN2790
 2790

· UN proper shipping name

· ADR 2790 ACETIC ACID SOLUTION ACETIC ACID SOLUTION

- · Transport hazard class(es)
 - · DOT/TDG (Transport dangerous goods):



· Class 8 Corrosive substances

· Label

· ADR



· Class 8 (C3) Corrosive substances

· Label

· IMDG, IATA



· Class 8 Corrosive substances

· Label

· Packing group

· DOT/TDG, ADR, IMDG, IATA

(Contd. on page 7)





according to HPR, Schedule 1

Printing date 11/01/2020 Reviewed on 07/08/2015

Trade name: VINEO Extract DNA Kit - W2, 1 x 100 ml

(Contd. of page 6)

· Environmental hazards:

· Marine pollutant:

No

· Special precautions for user

Warning: Corrosive substances

Hazard identification number (Kemler code): 80

· EMS Number: F-A,S-B

· Transport in bulk according to Annex II of

MARPOL73/78 and the IBC Code

Not applicable.

UN "Model Regulation":

UN2790, ACETIC ACID SOLUTION, 8, III

15 Regulatory information

· Safety, health and environmental regulations/legislation specific for the substance or mixture · Sara

· Section 355 (extremely hazardous substances):

None of the ingredient is listed.

· Section 313 (Specific toxic chemical listings):

None of the ingredients is listed.

· TSCA (Toxic Substances Control Act):

127-08-2 potassium acetate

64-19-7 acetic acid

7732-18-5 water, distilled, conductivity or of similar purity

· Canadian substance listings:

· Canadian Domestic Substances List (DSL)

127-08-2 potassium acetate

64-19-7 acetic acid

7732-18-5 water, distilled, conductivity or of similar purity

· Canadian Ingredient Disclosure list (limit 0.1%)

None of the ingredients is listed.

· Canadian Ingredient Disclosure list (limit 1%)

64-19-7 acetic acid

National regulations:

Class	Share in %
II	12.0

· Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

CDN



Page 8/8

Safety Data Sheet

according to HPR, Schedule 1

Printing date 11/01/2020 Reviewed on 07/08/2015

Trade name: VINEO Extract DNA Kit - W2, 1 x 100 ml

(Contd. of page 7)

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

Department issuing SDS:

MSDS / FDS Administr.

Bio-Rad

Marnes-la-Coquette, France

· Contact:

SOFRALAB SAS

79 av. A.A. Thévenet, BP1031 Magenta,

51319 EPERNAY Cedex Mob.: +33 6 25 10 00 26

e-mail: nsieczkowski@sofralab.com

Customer Support

Phone: +33 (0)1 47 95 62 77

Customer Support

Phone: +33 (0)1 47 95 61 85

Abbreviations and acronyms: ICAO: International Civil Aviation Organisation

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA) HMIS: Hazardous Materials Identification System (USA)

PBT: Persistent, Bioaccumulative and Toxic vPvB: very Persistent and very Bioaccumulative

CDN