

Gene Pulser® Electroporation

Cell Type	Bacterial, gram positive	Molecules Electroporated	DNA: covalently closed plasmids
Species Used	<i>Brevibacterium lactofermentum</i> , ATCC 13859, R31,1035 (thr-).		

Before the Pulse

Cell Growth Medium	Tryptic soy broth (TSB)	Growth Phase at Harvest	O.D. (600) =0.4 - 0.8
		Pre-pulse Incubation	Not given
Wash Solution	Not given		

The Pulse

Instruments Used Gene Pulser® apparatus

Electroporation Temperature	Room temperature		
Electroporation Medium	10% glycerol, distilled water	Cuvette Gap	0.2 cm
Cell Density	Approx. 10 (10) cells / ml	Voltage	2.5 kV
Volume of Cells	1 to 5 µl	Field Strength	12.5 kV/cm
DNA Concentration	0.1 to 1 µg		
DNA Resuspension Buffer	TE buffer (10 mM Tris, 1 mM EDTA, pH8.0) or water	Capacitor	25 µF
Volume of DNA	1 to 5 µl	Resistor	200 Ω (Pulse Controller)
After the Pulse		Time Constant	2.8 to 4.5 msec
Outgrowth Medium	TSB		

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.

Outgrowth Temperature	30 °C
Length of Incubation	1 hour
Selection Method or Assay Used	Minimal media plus antibiotics
Electroporation Efficiency	Maximum 10 (5) to10 (6) transformants / µg DNA
Per Cent Survival	Not done

Name of Submitter Carmen Guerrero Arroyo, Becario PFPI

Institution Address Universidad de Leo,
Dept. of Microbiology
Campus de Vegazana S/N CP 24071
Leon, Spain

Survey Number
057