

Cell Type Bacterial, gram positive

Molecules Electroporated DNA: linear; plasmids: 5 to 15 kb, supercoiled or relaxed.

Species Used *Enterococcus hirae*

Before the Pulse

Cell Growth Medium M-17 (very rich broth)

Growth Phase at Harvest O.D. (600) = late log

Pre-pulse Incubation none

Wash Solution water

The Pulse

Instruments Used Gene Pulser® apparatus

Electroporation Temperature 25 °C

Electroporation Medium water

Cuvette Gap 0.1 cm

Cell Density 1- 2 x 10 (10) cells / ml

Voltage 2.5 kV

Volume of Cells 120 µl

Field Strength 25 kV/cm

DNA Concentration 1 pg to 1 µg

Capacitor 25 µF

DNA Resuspension Buffer water

Resistor (Pulse Controller) See Comments

Volume of DNA 1 µl

Time Constant 65 msec

After the Pulse

Outgrowth Medium M-17

Relevant Publications and/or Comments

Note: exponential values designated in parentheses.

Resistor Setting: 3000 Ω serial resistor
References: *Biochemie* **72**: 279-83 (1990),
TIBS 15: 175-77 (1990).

Outgrowth Temperature 37 °C

Length of Incubation 30 min

Selection Method or Assay Used erythromycin resistance

Electroporation Efficiency 5 x 10 (6) transformants / µg DNA

Per Cent Survival 70%

Name of Submitter Marc Solioz, Dr.

Institution Address Clinical Pharmacology, University of Berne
Murteust. 35
3007 Berne, Switzerland

Survey Number
061