



Droplet Digital PCR Assay Site

ddPCR™ Assays are the only assays expertly designed specifically for Droplet Digital PCR (ddPCR). Bio-Rad offers a comprehensive assay portfolio, including validated ddPCR Assay Kits, a library of validated ddPCR Assays, and ddPCR Assays designed instantly on our Digital Assay Site.

Take advantage of 24/7 service to find or design and buy an assay in three simple steps:

1. Choose your application or search by Unique Assay ID.
2. Enter target or sequence information.
3. Pick a size and add to your Bio-Rad cart.

Visit [bio-rad.com/web/ddPCRdigital-assays](https://www.bio-rad.com/web/ddPCRdigital-assays) to find your ddPCR Assay.





ddPCR Mutation Detection Assays

ddPCR Mutation Detection Assays offer superior sensitivity for detection of rare single nucleotide polymorphisms (SNPs) and allelic variants.

- Wet-lab validated assays are available for more than 200 targets
- Wet-lab validated assays are guaranteed to work when used as instructed
- Validation data are available for review and download
- Assays are designed instantly for targets/sequences for which validated assays are not available
- Suitable for use with circulating tumor DNA, liquid biopsy, and FFPE and other tissue samples
- Orderable in convenient single-tube format (mutant and wild type assays in a single tube)
- Positive control sequences provided



ddPCR Copy Number Determination Assays

ddPCR Copy Number Variation Assays offer high precision and resolution to discriminate between small changes in copy number.

- Wet-lab validated assays are available for more than 700 targets
- Wet-lab validated centromeric reference assays are available for ~60 targets
- Validation data are available for review and download
- Assays are designed instantly for targets/sequences for which validated assays are not available
- Orderable in FAM or HEX
- Reference assays available



ddPCR Genome Edit Detection Assays

ddPCR Assays for detection of homology directed repair (HDR) and non-homologous end joining (NHEJ) events in user specified sequences allow rapid assessment of HDR and NHEJ edits generated by CRISPR-Cas9 or other genome editing tools.

- Detect editing events present at frequencies of $\leq 0.5\%$
- Provide absolute quantification of genome editing events from as little as 5 ng of total gDNA (cells from single well of a 96-well plate)
- Distinguish between homozygous and heterozygous edits in clonal populations
- Results are easily analyzable using QuantaSoft™ Analysis Pro Software



ddPCR Expert Design Assays

ddPCR Expert Design Assays encompass a wide variety of assays, including assays for translocations, mutation detection, and multiplex mutation screening.

- Assays are available for *EML4-ALK* mRNA fusions, *RET* mRNA fusions, *ROS1* mRNA fusions, *TERT C228T* promoter mutation, *TERT C250T* promoter mutation, *ESR1* multiplex mutation detection, and more
- Orderable using a Unique Assay ID corresponding to a specific formulation of primers and probes
- Validation is recommended prior to use

New assays are added regularly. Contact your Bio-Rad representative for updates.



Bio-Rad
Laboratories, Inc.

Life Science
Group

Web site bio-rad.com **USA** 1 800 424 6723 **Australia** 61 2 9914 2800 **Austria** 43 1 877 89 01 177 **Belgium** 32 (0)3 710 53 00 **Brazil** 55 11 3065 7550
Canada 1 905 364 3435 **China** 86 21 6169 8500 **Czech Republic** 420 241 430 532 **Denmark** 45 44 52 10 00 **Finland** 358 09 804 22 00
France 33 01 47 95 69 65 **Germany** 49 89 31 884 0 **Hong Kong** 852 2789 3300 **Hungary** 36 1 459 6100 **India** 91 124 4029300
Israel 972 03 963 6050 **Italy** 39 02 216091 **Japan** 81 3 6361 7000 **Korea** 82 2 3473 4460 **Mexico** 52 555 488 7670 **The Netherlands** 31 (0)318 540 666
New Zealand 64 9 415 2280 **Norway** 47 23 38 41 30 **Poland** 48 22 331 99 99 **Portugal** 351 21 472 7700 **Russia** 7 495 721 14 04
Singapore 65 6415 3188 **South Africa** 27 (0) 861 246 723 **Spain** 34 91 590 5200 **Sweden** 46 08 555 12700 **Switzerland** 41 026674 55 05
Taiwan 886 2 2578 7189 **Thailand** 66 2 651 8311 **United Arab Emirates** 971 4 8187300 **United Kingdom** 44 020 8328 2000

