





Superior Hepatitis Testing from Bio-Rad

**Color-coded Testing Procedures Provide
Safety and Ease-of-Use**

BIO-RAD



Comprehensive Product Line to Meet Your Hepatitis Testing Needs

Sensitive and Specific Assays to Detect Every Stage of Hepatitis Infection

- MONOLISA™ Anti-HAV EIA*
- MONOLISA™ Anti-HAV IgM EIA*
- MONOLISA™ Anti-HBs EIA* and MONOLISA™ Anti-HBs Calibrator Kit* for Quantification of Antibodies
- MONOLISA™ Anti-HBc EIA*
- MONOLISA™ Anti-HBc IgM EIA*
- GS HBsAg EIA 3.0 – Outstanding sensitivity with less than 0.1 ng/mL of ad/ay HBsAg**

* This assay has not been FDA cleared or approved for the screening of blood or plasma donors. Performance of the assay has not been established with cord blood, neonatal specimens, cadaver specimens, heat-inactivated specimens, or body fluids other than serum or plasma, such as saliva, urine, amniotic fluid, or pleural fluid.

** Calibrated versus PHA 805 commercial panel.




Fast incubation time:
2.5 hours



Standardized procedure
for all Bio-Rad assays



Common reagents for
all Bio-Rad kits



Manual and automated use



Microplate strip identification



Color-coded reagent
and kit labels





Bar-coded reagents




Color changes during
protocol steps



Controls with color-
coded caps



Hepatitis A Serological Profile:
Example of Acute Hepatitis




Hepatitis B Serological Profile:
Example of Acute Hepatitis

All Bio-Rad Hepatitis EIA Assays Have Standardized Procedure and Color Changes

Example: MONOLISA™ Anti-HBc EIA

Test Format



HBc Ag Recombinant Proteins
 Anti-HBc Ab
 Anti-IgG/IgM Conjugate

* Reactive sample after 30 minutes incubation.

BIO-RAD

Bio-Rad Hepatitis Family of Products

Ordering Information

Catalog No.	Description	
72496	MONOLISA™ Anti-HAV EIA	192 tests
72495	MONOLISA™ Anti-HAV IgM EIA	192 tests
25220	MONOLISA™ Anti-HBs EIA	192 tests
25219	MONOLISA™ Anti-HBs Calibrator Kit	20 calibrators
26186	MONOLISA™ Anti-HBc EIA	192 tests
26174	MONOLISA™ Anti-HBc IgM EIA	192 tests
32591	GS HBsAg EIA 3.0	480 tests
32592	GS HBsAg EIA 3.0	960 tests
32593	GS HBsAg EIA 3.0	4800 tests
930740	Ortho® HCV EIA V 3.0	480 tests
930600	Ortho® RIBA™ HCV 3.0 SIA	30 tests
For Confirmatory Neutralization of GS HBsAg EIA 3.0 Reactive Specimens:		
32594	GS HBsAg Confirmatory Assay 3.0	25 tests



MONOLISA
Anti-HBc EIA



Bio-Rad
Laboratories

For further information, please contact the Bio-Rad office nearest you or visit our website at www.bio-rad.com/diagnostics.

Clinical
Diagnostics Group

Website www.bio-rad.com/diagnostics U.S. 1-800-2BIO-RAD Australia 61-2-9914-2800 Austria 43-1-877-8901 Belgium 32-9-385-5511 Brazil 5521-3237-9400
Canada 1-514-334-4372 China 86-21-64260808 Czech Republic 420-241-430-532 Denmark +45-4452-1000 Finland 358-9-804-22-00 France 33-1-47-95-60-00
Germany +49-(0)89-318-840 Greece 30-210-7774396 Hong Kong 852-2789-3300 Hungary 36-1-455-8800 India 91-124-4029300 Israel 972-3-9636050 Italy +39-02-216091
Japan 81-3-6361-6902 Korea 82-2-3473-4460 Mexico +52-(555)-488-7670 The Netherlands +31-318-540666 New Zealand 64-9-415-2280 Norway 47-23-38-41-30
Poland 48-22-3319999 Portugal 351-21-472-7700 Russia 7-495-721-14-04 Singapore 65-6415-3188 South Africa 27-11-442-85-08 Spain 34-91-590-5200
Sweden 46-8-555-127-00 Switzerland 41-61-717-95-55 Thailand 662-651-8311 United Kingdom +44-(0)20-8328-2000