




Tumor Marker Controls

Reliable, independent QC for Tumor Marker analysis





Precision monitoring for Tumor Marker analysis

According to the World Health Organization (WHO), cancer is one of the leading causes of death globally, surpassed only by cardiovascular disease.¹ Cancer may affect almost any part of the body producing a variety of biomarkers. As efficacy of treatment becomes increasingly dependent on early detection, cancer screening to detect these biomarkers is becoming common practice around the globe.

The practice of serum tumor marker analysis is routinely performed to screen, aid in diagnosis, monitor disease progression, response to therapy or recurrence.

Bio-Rad's family of Tumor Marker Controls offer a comprehensive collection of tumor marker analytes designed to monitor a variety of cancer related procedures. These independent human based controls are similar to patient samples for reliable monitoring on most major testing platforms used in cancer screening. Bio-Rad quality

controls are combined with expert service and support to help your lab streamline QC workflows and increase confidence in patient results.

Bio-Rad Laboratories offers a vision and a commitment to quality control that is unsurpassed in the clinical laboratory industry. Using our comprehensive product line of independent controls is an important step toward helping you improve the reliability of your test results.

¹"Non Communicable Diseases" Fact Sheet, World Health Organization, accessed April 2017, <http://www.who.int/mediacentre/factsheets/fs355/en/>.

Liquichek Tumor Marker Control



4 NEW
Analytes

Standard vials (2 mL fill size)

- An assayed, multi-analyte control
- Includes low PSA levels and high levels for cancer antigens and Ferritin
- Now with HE4, HER-2/neu, ProGRP* and SCC*
- Human based
- 2 year shelf life when stored frozen at -20°C to -70°C
- 30 day open-vial stability at 2-8°C for most analytes



Compact vials (2 mL fill size)

- 2D barcode for quick and easy product insert and vial labeling scanning
- Compatible with applicable automated platforms such as Siemens Dimension Vista and Abbott Alinity
- Human based
- 10 day open-vial (punctured) stability for all analytes when stored at 2-8°C
- 2 year shelf life when stored at -20°C to -50°C

InteliQ Tumor Marker Control



Load-and-go tubes (2 mL fill size)

- Ready-to-use, barcoded tube helps automate your daily QC
- Compatible with major automated chemistry and immunoassay platforms
- Liquid, human serum based
- 2 year shelf life at -20°C to -70°C
- Open-vial stability
 - Refrigerator storage: 10 days at 2-8°C for most analytes (all instruments, off-board)**
 - Instrument storage: 20 days at 2-8°C for most analytes (Siemens Atellica, on-board)

Lyphocheck Tumor Marker Plus Control



Standard vials (2 mL fill size)

- A comprehensive, assayed, multi-analyte lyophilized control
- Ideal for monitoring both routine and esoteric tumor marker testing
- Human serum based
- 3 year shelf life at 2-8°C
- 14 day reconstituted stability at 2-8°C for most analytes


Analytes

	Liquid (Standard vial, 2 mL)	Liquid (Compact vial, 2 mL)	InteliQ	Lyophilized
β-2-Microglobulin	●	●	●	●
ACTH				●
Aldosterone				●
Alphafetoprotein (AFP)	●	●	●	●
CA 15-3	●	●	●	●
CA 19-9	●	●	●	●
CA 27.29	●		●	●
CA 50*				●
CA 72-4*	●		●	●
CA 125	●	●	●	●
Calcitonin				●
Carcinoembryonic Antigen (CEA)	●	●	●	●
CASA* ¹				●
Cyfra 21-1*	●		●	●
Ferritin	●	●	●	●
hCG	●	●	●	●
hCG-β Subunit				●
HE4	●		●	
HER2/neu	●		●	
Insulin Growth Factor I (IGF-I)	●		●	
Neuron Specific Enolase (NSE)*				●
PAP	●		●	●
ProGRP*	●		●	
Prolactin	●	●	●	●
PSA (Free)	●	●	●	●
PSA (Total)	●	●	●	●
S-100*				●
SCC*	●		●	
Thyroglobulin (Tg)	●		●	●

*International use only. Not available for diagnostic use in the United States.

¹No claim is made for expected values or stability.

**Please refer to the package insert for platform-specific stability and performance claims. Refer to myInserTs.com or the package insert of currently available lots for specific analyte and stability claims.



Unity QC Data Management Solutions

Bio-Rad Quality Controls are part of the Unity Program with its industry-leading features for benchmarking laboratory performance.

- Benefit from peer data generated from more than 40,000 instruments worldwide
- Utilize real-time access to peer group comparisons to help with troubleshooting
- Consolidate your entire laboratory into one QC system
- Facilitate run-time decisions with optimized Westgard rules
- Facilitate compliance with simple tools to review rule violations and document corrective actions
- Define your own quality requirements and track key quality performance indicators

You can consolidate your entire laboratory into one QC system through our Unity Program. This sophisticated software includes best-in-class tools to support releasing results, benchmarking laboratory performance, and designing your QC system.

Ordering Information

Cat # Description

Liquichek Tumor Marker Control

Standard Vials

547	Level 1	6 x 2 mL
548	Level 2	6 x 2 mL
549	Level 3	6 x 2 mL
548X	Trilevel MiniPak (1 of each level)	3 x 2 mL

Compact Vials

27114	Level 1	6 x 2 mL
27115	Level 2	6 x 2 mL
27116	Level 3	6 x 2 mL
27115X	Trilevel MiniPak (1 of each level)	3 x 2 mL

Cat # Description

InteliQ Tumor Marker Control

12008289	Level 1	6 x 2 mL
12008290	Level 2	6 x 2 mL
12008291	Level 3	6 x 2 mL
12008292	Trilevel MiniPak (1 of each level)	3 x 2 mL

Lyphocheck Tumor Marker Plus Control

367	Level 1	6 x 2 mL
368	Level 2	6 x 2 mL
369	Level 3	6 x 2 mL
368X	Trilevel MiniPak (1 of each level)	3 x 2 mL

BIO-RAD

**Bio-Rad
Laboratories, Inc.**

**For further information, please contact the Bio-Rad office nearest you
or visit our website at www.qcnet.com/quality-controls/tumor-markers**

**Clinical
Diagnostics Group**

Website www.bio-rad.com/diagnostics **Australia** +61 (2) 9914 2800 **Austria** +43 (0) 1 877 89 01 9 **Belgium** +32 (0) 3 710 53 00 **Brazil** +55 11 3065 7550 **Canada** +1 514 334 4372 **China** +86 21 6169 8500 **Czech Republic** +420 241 431 660 **Denmark** +45 44 52 10 00 **Finland** +358 9 804 22 00 **France** +33 (0) 1 47 95 60 00 **Germany** +49 (0) 89 31884 393 **Greece** +30 210 7774396 **Hong Kong** +85 2 2789 3300 **Hungary** +36 1 459 6190 **India** +91 124 4029300 **Israel** +972 03 963 6025 **Italy** +39 024 94 86 600 **Japan** +81 3 6361 7070 **Korea** +82 080 007 7373 **Mexico** +52 (55) 5488 7670 **The Netherlands** +31 (0)318 540 666 **New Zealand** +64 (9)415 2280 **Norway** +47 23 38 41 30 **Poland** +48 22 331 99 99 **Portugal** +351 21 47 27 700 **Russia** +7 495 721 1404 **Singapore** +65 6415 3170 **South Africa** +27 11 442 8508 **Spain** +34 91 490 6580 **Sweden** +46 844 98053 **Switzerland** +41 (0) 61 717 9555 **Taiwan** +886 (2) 2578-7189 **Thailand** (662) 651 8311 **United Kingdom** +44 (0)1923 471301

BIO-RAD, INTELIQ, LIQUICHEK and LYPHOCHEK are trademarks of Bio-Rad Laboratories, Inc. in certain jurisdictions.

© 2021 Bio-Rad Laboratories, Inc. Printed in USA QSD21-0372 Q-1556 Rev. 08/2021

Printed on recycled paper
with soy-based inks