

Liquichek Opiate Control



Liquichek Opiate Control



Targeted for analytes tested in pain management or opioid treatment and compliance programs.

- Contains eight drugs or drug metabolites, at concentrations above and below assay cutoffs
- This control has been tested by GC/MS, LC/MS/MS and other methods to confirm the recovery of drug concentration
- Representative values for screening and confirmatory methods provided
- Human urine
- 3 year shelf-life at -20°C to -70°C
- 30 day open-vial stability at 2-8°C

Gravimetric Reference Values

Analyte (ng/mL)	Level 1	Level 2
6-Monoacetylmorphine (6-MAM, Heroin Metabolite)	7.5	12.5
Buprenorphine	5	15
EDDP (Methadone Metabolite)	50	150
Fentanyl	0.25	0.75
Meperidine*	100	300
Methadone	150	450
Oxycodone	50	150
Tramadol	100	300

*International use only.

EQAS

An independent, external assessment of performance in comparison to your peers.

[QCNet.com/eqas](https://www.qcnet.com/eqas)

Independent QC

Ongoing, proactive, unbiased daily QC that helps identify errors as they occur or begin to trend.

[QCNet.com/independentqc](https://www.qcnet.com/independentqc)

Unity

QC Data Management tools that help you create a strategy to reduce risk and streamline QC workflow.

[QCNet.com/datamanagement](https://www.qcnet.com/datamanagement)

Ordering Information

Catalog No. Description

12005786	Level 1	6 x 5 mL
12005787	Level 2	6 x 5 mL
12002594	Level 1 MiniPak.	1 x 5 mL
12002595	Level 2 MiniPak.	1 x 5 mL

Refer to myinserts.com or the package insert of currently available lots for specific analyte and stability claims.

BIO-RAD

**Bio-Rad
Laboratories, Inc.**

For further information, please contact the Bio-Rad office nearest you or visit our website at www.bio-rad.com/qualitycontrol

Clinical
Diagnostics Group

Website www.bio-rad.com/diagnostics **U.S.** 1 800 224 6723 **Australia** +61 (2) 9914 2800 **Austria** +43 (0) 1 877 89 01 9 **Belgium** +32 (0) 3 710 53 00 **Brazil** +55 11 3065 7550 **Canada** +1 514 334 4372 **China** +86 21 6169 8500 **Czech Republic** +420 241 431 660 **Denmark** +45 44 52 10 00 **Finland** +358 9 804 22 00 **France** +33 (0) 1 47 95 60 00 **Germany** +49 (0) 89 31884 393 **Greece** +30 210 7774396 **Hong Kong** +85 2 2789 3300 **Hungary** +36 1 459 6190 **India** +91 124 4029300 **Israel** +972 03 963 6025 **Italy** +39 024 94 86 600 **Japan** +81 3 6361 7070 **Mexico** +52 (55) 5488 7670 **The Netherlands** +31 (0) 318 540 666 **New Zealand** +64 (9)415 2280 **Norway** +47 23 38 41 30 **Poland** +48 22 331 99 99 **Portugal** +351 21 47 27 700 **Republic of Korea** +82 080 007 7373 **Russia** +7 (495) 721-14-04 **Singapore** +65 6415 3170 **South Africa** +27 11 442 8508 **Spain** +34 91 490 6580 **Sweden** +46 844 98053 **Switzerland** +41 (0) 61 717 9555 **Taiwan** +886 (2) 2578-7189 **Thailand** (662) 651 8311 **United Kingdom** +44 (0) 1923 471301

LIQUICHEK is a trademark of Bio-Rad Laboratories, Inc. in certain jurisdictions.

© 2024 Bio-Rad Laboratories, Inc. Printed in USA CDG24-0275 Q-1627 Rev. 11/2024

Printed on recycled paper with soy-based inks