

QC Data Management Solutions

Introduce your laboratory to a whole new world of Unity data management solutions





Commitment to Quality Control

Using our comprehensive product line of independent controls is an important step toward helping you improve the reliability of your laboratory results. The ability to successfully manage and interpret your quality control results is also essential to producing solid, dependable laboratory outcomes.

Bio-Rad Laboratories offers a vision and a commitment to quality control that is unsurpassed in the clinical laboratory industry. Today, your laboratory can be connected to the Unity Interlaboratory Program, the world's largest community of quality control users. Improve patient test results and cost effectiveness with quality control products, quality control data management solutions, industry-leading interlaboratory comparisons and personalized service from Bio-Rad.



Regulatory Requirements

"The laboratory shall monitor its performance of examination methods, by comparison with results of other laboratories."

– ISO 15189:2022(E), Subclause 7.3.7.3



The Advantages of Participating in the Unity Interlaboratory Program

Detecting Errors in Your Test System

Without a high quality interlaboratory program, your laboratory may not be aware of gradual or sudden changes in your test system that may be caused by events such as reagent or calibrator reformulations, standardization changes, or instrument software changes. An interlaboratory program may offer early awareness of shifts and trends to help avoid costly test repeats and unnecessary troubleshooting.

Your Proficiency Program is Not Enough

Proficiency surveys provide a snapshot in time and should not be used as a substitute for daily quality control with an interlaboratory program.

- Acceptable performance on the day of the proficiency testing does not guarantee testing reliability every day
- Participation in an interlaboratory program is the recommended way to improve confidence in test system performance every day, especially with an interlaboratory program that offers an option to receive instant interlaboratory feedback on-demand

Some advantages for participating in an interlaboratory QC program are:

- Verifying that a laboratory is producing QC results that are consistent with other laboratories using the same measurement procedure, and thus demonstrating that the laboratory is using the measurement procedure correctly

Detecting and identifying bias in a measurement procedure

- Bias can be caused by events such as reagent or calibrator lot changes or reformulations, changes in calibration traceability to reference systems, or instrument software changes. Comparison of an individual laboratory's QC result to a peer group mean value can identify a trend or shift, or ascertain if other laboratories are experiencing the same changes."

Supplementing PT/EQA programs

- PT/EQA programs verify performance at a point in time. Acceptable performance on the day of PT/EQA testing does not guarantee testing reliability every day because errors in a measurement procedure can occur at any time. In addition, the interlaboratory QC data can be used to investigate a failed or questionable PT/EQA result.
 - CLSI, C24, 4th ed. Statistical Quality Control for Quantitative Measurement Procedures: Principles and Definitions
 - 7.2 Using Interlaboratory Quality Control to Assess a Quality Control Program.

Unity Solutions

Bio-Rad Laboratories offers several options for participation in the Unity Interlaboratory Program. Refer to the product comparison grid below to help determine which solution is the most appropriate for your laboratory.

Internet-based solutions such as UnityWeb, Unity Next Peer QC, and Unity Real Time online eliminate the need to install and update software locally in your laboratory, and reduce the amount of support necessary from on-site IT staff.

Desktop software solutions, like Unity Real Time, are appropriate if your internet connection is not sufficient, or if you prefer maintaining local software.

	Basic User Options			Advanced User Options	
	Unity Next Peer QC	UnityWeb	Unity Real Time LT	Unity Real Time Online	Unity Real Time
Platform Design					
Web Service	X	X		X	
Desktop Software			X		X
Bio-Rad Hosted Database	X	X		X	
Single Platform for Data Entry & Peer Reports	X				
Digital Sign & Save for Peer Reports	X				
Unity Interlaboratory Reports (Peer Reports)					
Monthly Evaluation	On-demand	Monthly or Instant QC	Monthly or Instant QC	Monthly or Instant QC	Monthly or Instant QC
Lab Comparison	On-demand	Monthly or Instant QC	Monthly or Instant QC	Monthly or Instant QC	Monthly or Instant QC
Lab Histogram	On-demand	Monthly	Monthly	Monthly	Monthly
Bias and Imprecision		Monthly	Monthly	Monthly	Monthly
Statistical Profile		Monthly	Monthly	Monthly	Monthly
Qualitative Report		Monthly	Monthly	Monthly	Monthly
Affiliated Report	Q3-2025	Monthly	Monthly	Monthly	Monthly
Worldwide report	On-demand	Monthly	Monthly	Monthly	Monthly
Manufacturer Report	On-demand	Monthly	Monthly	Monthly	Monthly
QC Design and Charts					
Westgard Rules	X	X	X	X	X
Levey-Jennings Chart	X	X	X	X	X
Multi-Levey-Jennings Chart		X		X	X
Additional Features					
Westgard Advisor (Subscription sold separately)		X		X	X
Bench Review		X	X	X	X
Supervisor Review				X	X
Analytical Goals				X	X
Dynamic Data Set Comparisons				X	X
Connectivity (Sold separately)					
WebConnect		X		X	
UnityConnect		X	X	X	X
Unity Next File Upload	X				



Basic



Advanced



Desktop based




Web based




Mobile

Want to know more? Discover the power of Unity Solutions at www.bio-rad.com/qc-datamanagement




Basic User Options



Unity Real Time LT

Entry-Level Desktop QC Data Management Solution


- Basic QC rules, charts and reports
- Configurable Levey-Jennings charts
- Upgrade easily to Unity Real Time for more advanced tools and features
- Upload QC data points from an LIS, middleware or instrument (optional)



UnityWeb

Entry-Level Online QC Data Management Solution

- No software to install and maintain
- Use basic QC rules, charts and reports
- Upgrade easily to Unity Real Time online for more advanced tools and features
- Upload QC data points from an LIS, middleware or instrument (optional)




Unity Next Peer QC

Entry-Level Online QC Data Management Solution with Instant Access to Peer Comparison

- Latest generation of Unity with user-friendly interface
- No software to install and maintain
- Use basic QC rules, Levey-Jennings charts and reports
- Instantly access peer reports for peer comparison and QC troubleshooting
- Upload QC data points from LIS, middleware, or instrument
- Submit QC data and download reports in one location



Connectivity



Basic

WebConnect



File uploads to UnityWeb & Unity Real Time Online Software

A web-based connectivity solution that allows laboratories to easily upload QC data from LIS systems, middleware and/or instruments directly into Unity Real Time online or UnityWeb

- Eliminate manual keying of data
- Standard LIS QC reports can be used
- No software to install

Advanced

UnityConnect



A connectivity solution that allows QC data from LIS systems, middleware and/or instruments to be quickly and easily imported into Unity software, enabling interlaboratory program comparison

- Eliminate manual keying of QC data
- Standard LIS QC reports can be used
- Automatic data import and real-time connection
- Capture and analyze QC data economically from laboratory instruments not connected to the LIS
- Easy-to-install optional interface hardware
- Choose between serial and ethernet connection

Basic


Unity Next File Upload



File uploads to Unity Next Software

- A data upload solution allowing file uploads of QC data, specifically for Unity Next software
- Upload QC data from LIS, middleware, and/or instruments to Unity Next
- No software to install


Advanced User Options



Unity Real Time

Expert QC Data Management Solution for Desktop Users


- Facilitate regulatory compliance under CLIA and ISO 15189
- Improve real-time bench and supervisor QC data review
- Implement best QC rules when used with Westgard Advisor
- Run validation with comprehensive audit trails
- Advanced charts and reports for data analysis
- Reduce non-essential retests with Analytical Goal options
- Upload QC data points from an LIS, middleware or instrument (optional)
- RiLiBÄK Advisor module available to comply with German regulations



Unity Real Time online

Expert QC Data Management Solution for Online Users


- Facilitate regulatory compliance under CLIA and ISO 15189
- Improve real-time bench and supervisor QC data review
- Implement best QC rules when used with Westgard Advisor
- Run validation with comprehensive audit trails
- Advanced charts and reports for data analysis
- Reduce non-essential retests with Analytical Goal options
- Upload QC data points from an LIS, middleware or instrument (optional)



Westgard Advisor

Automatic QC Rules Selection Engine

- Recommend and automatically apply best QC rules with patented technology
- Easy step-by-step automatic rule selection capabilities
- Reduce false rejections and desensitization to false error flags
- Save time and money by reducing unnecessary repeats and troubleshooting
- Improve laboratory test quality with optimally selected QC rules
- Available as an optional module with Unity QC data management solutions



Unity Alert

Continuously Monitor QC Status

- Notifications of missing QC runs can help to ensure proper monitoring of instrument performance
- Notifications of new violations can be useful for troubleshooting and monitoring problem analytes
- Alert rules are easy to set up and customize for specific instruments or tests
- Receive emails remotely to stay informed of QC issues in your laboratory
- Color-coded displays help draw attention to the most urgent QC issues
- Notifications are provided even when Unity Real Time is not in use
- Available as an optional module with Unity Real Time



Bio-Rad Mission: Control



Take the guesswork out of your QC strategy

Employing a robust quality control system is an essential step for any laboratory hoping to provide reliable results for appropriate patient care. Unfortunately for clinical laboratories, selecting the right QC rules and frequency can be challenging, and questions often arise.

- How often should controls be run for specific assays— is there an optimal frequency?
- How can QC be managed to minimize the risk of releasing incorrect patient results?

Bio-Rad Mission: Control is the first objective risk management solution for quality control in laboratories. The intuitive interface allows users to focus on patient safety by assessing the risk of reporting incorrect patient results in their laboratory.

Bio-Rad Mission: Control allows you to:

- Determine your laboratory's Risk Management Index (RMI)
- Evaluate your test method performance based on your QC results
- Quantify the likelihood and severity of reporting incorrect patient results
- Assess how much you can lower your risk, based on instrument performance
- Create a QC Plan by confidently selecting rules and frequency that fit your risk tolerance

Minimize your risk of reporting unreliable patient results

A patient risk approach to Quality Control is a progressive step toward improving quality assurance and patient care. Bio-Rad can help you understand what this means for your laboratory and how to take steps to implement an effective approach for Quality Control within your analytical test phase.

**For more information,
visit www.qcnet.com/missioncontrol**




Regulatory Requirements

"IQC shall be performed at a frequency that is based on the stability and robustness of the examination method and the risk of harm to the patient from an erroneous result."

– ISO 15189:2022(E), Subclause 7.3.7.2

QCNet is the Leading QC Portal on the Internet

Learn more about the many QC tools and resources from Bio-Rad at www.QCNet.com



Quick Links to the most-used features

Unity Reports


- Get your Unity Reports in just a few clicks.

My elnserts

- Access product inserts

EQAS

- View Reports
- Access EQAS Online
- Training/Tutorials
- View Product Information




Access country-specific information

At the bottom of the QCNet home page, click on the country name shown to access the pop-up menu, then select the desired country*.

*Not all countries are available.

Explore the drop-down menus to find more online QC tools on QCNet



Ordering Information

Cat # Description

Unity Real Time LT

805-1 Unity Real Time LT Annual subscription
825i Unity Remote Installation & Training. Single use

UnityWeb

870-1 UnityWeb Annual subscription

WebConnect

870-W1 WebConnect 2.0 Annual subscription
825i Unity Remote Installation & Training. Single use

UnityConnect

11000009 UnityConnect Software. Annual subscription
858i Additional Software Training. Single use

Unity Real Time

804-1 Unity Real Time Annual subscription
825i Unity Remote Installation & Training. Single use

Unity Real Time online

804-W1 Unity Real Time online. Annual subscription

Cat # Description

Westgard Advisor

811-1 Westgard Advisor. Annual subscription

Unity Alert

806-1 Unity Alert. Annual subscription
825i Unity Remote Installation & Training. Single use

Bio-Rad Mission: Control

1100002 Mission: Control. Annual Subscription
1100005 Mission: Control. . . Subscription and Service Bundle

Unity Next Peer QC

11000517 Unity Next Peer QC (5 users). . . Annual Subscription
11000506 Unity Next Onboarding Single Use

Unity Next File Upload

11001897 Unity Next File Upload Annual Subscription

EQAS

An independent, external assessment of laboratory performance in comparison to your peers.

[QCNet.com/eqas](https://qcnets.com/eqas)

Independent QC

Ongoing, proactive, unbiased daily QC that helps identify errors as they occur or begin to trend.

[QCNet.com/independentqc](https://qcnets.com/independentqc)

Unity

QC Data Management tools that help you create a strategy to reduce risk and streamline QC workflow.

[QCNet.com/datamanagement](https://qcnets.com/datamanagement)

Bio-Rad, Mission: Control, Unity, Unity Real Time and UnityWeb are trademarks of Bio-Rad Laboratories, Inc. in certain jurisdictions. All trademarks used herein are the property of their respective owner.



**Bio-Rad
Laboratories, Inc.**

Clinical Diagnostics Group

Website www.bio-rad.com/diagnostics U.S. 1 800 224 6723 Australia 61 2 9914 2800 Austria 43 0 1 877 89 01 9 Belgium 32 0 3 710 53 00
Brazil 55 11 3065 7550 Canada 1 514 334 4372 China 86 21 61698500 Czech Republic 420 241 431 660 Denmark 45 44 52 10 00
Finland 358 9 804 22 00 France 33 0 1 47 95 60 00 Germany 49 0 89 31884 393 Greece 30 210 7774396 Hong Kong 852 27893300
Hungary 36 1 459 6190 India 91 124 4029300 Israel 972 03 963 6025 Italy 39 024 94 86 600 Japan 81 3 6361 7070 Mexico 52 55 5488 7670
The Netherlands 31 0 318 540 666 New Zealand 64 9 415 2280 Norway 47 23 38 41 30 Poland 48 22 331 99 99 Portugal 351 214 727 700
Republic of Korea 82 080 007 7373 Russia 7 495 721-14-04 Singapore 65 6415 3170 South Africa 27 11 442 8508 Spain 34 91 49 06 580
Sweden 46 844 98053 Switzerland 41 0 61 717 9555 Taiwan 886 2 25787189 Thailand 66 2 6518311 United Kingdom 44 0 1923 471301

