

Stocking Up?

Reagents and Supplies Quick Reference Guide

discover.bio-rad.com



BIO-RAD

Chromatography Selection Chart

Ion Exchange Chromatography

- General biomolecule purification
- Applicable to all antibody purifications
- High capacity
- Good capture step
- Separation by charge



Packed Columns/Cartridges

- UNO® Q
- UNO S
- Econo-Pac® High Q
- Econo-Pac UNOsphere™ S
- Econo-Pac UNOsphere Q
- Econo-Pac DEAE
- Econo-Pac High S
- Econo-Pac CM



Bottled Media

- Macro-Prep® 25 Q
- Macro-Prep 25 S
- Macro-Prep High Q
- Macro-Prep High S
- Macro-Prep CM
- Macro-Prep DEAE
- UNOsphere Q
- UNOsphere S

Size Exclusion Chromatography

- Molecular weight determination
- Applicable to all antibody purifications
- High resolution
- Separation by molecular size
- Desalting or buffer exchange



Packed Columns/Cartridges

- Econo-Pac P-6
- Bio-Spin® P-6
- Micro Bio-Spin™ P-6
- Bio-Sil®/Bio-Silect®



Bottled Media

- Bio-Gel® P

Hydroxyapatite Chromatography

- Unique selectivity
- Very high purity can often be obtained in one step
- Alternative to protein A for antibody purification



Packed Columns/Cartridges

- Bio-Scale™ CHT-I
- Econo-Pac CHT-I
- Econo-Pac CHT-II



Bottled Media

- CHT™ ceramic hydroxyapatite Type I, 20, 40, 80 µm
- CHT ceramic hydroxyapatite Type II, 20, 40, 80 µm
- Bio-Gel HT gel
- Bio-Gel HTP gel

Affinity Chromatography

- High selectivity
- High purity in a single step
- Excellent capture technique
- Separation by biological interaction



Packed Columns/Cartridges

- Econo-Pac protein A
- Econo-Pac Blue
- Econo-Pac DEAE Blue



Bottled Media — Ready-to-Use

- Affi-Prep® protein A
- Affi-Gel® protein A
- Affi-Gel Blue
- DEAE Affi-Gel Blue
- Affi-Gel boronate
- Affi-Prep polymyxin



Bottled Media — Activated

- Affi-Prep Hz
- Affi-Gel Hz
- Affi-Gel 10
- Affi-Prep 10
- Affi-Gel 15
- Affi-Gel boronate
- Affi-Gel 102

Hydrophobic Interaction Chrom.

- Higher-resolution alternative to ammonium sulfate precipitation
- Gentle elution conditions
- Separation by hydrophobicity of biomolecules



Packed Columns/Cartridges

- Econo-Pac t-butyl HIC
- Econo-Pac methyl HIC



Bottled Media

- Macro-Prep t-butyl
- Macro-Prep methyl

Chromatography Media and Columns

Catalog #	Description	Size/ Quantity	Storage Temp.
Ion Exchange Bottled Media			
153-0021	Macro-Prep 25 Q Support*	50 ml	Ambient
153-0031	Macro-Prep 25 S Support*	50 ml	Ambient
156-0040	Macro-Prep High Q Support**	100 ml	Ambient
156-0030	Macro-Prep High S Support**	100 ml	Ambient
156-0020	Macro-Prep DEAE Support**	100 ml	Ambient
156-0070	Macro-Prep CM Support**	100 ml	Ambient
156-0101	UNOsphere Q Support***	25 ml	Ambient
156-0111	UNOsphere S Support***	25 ml	Ambient

* Available in 50 ml, 200 ml, 1 L, and 5 L

** Available in 25 ml, 100 ml, 500 ml, 5 L, and 10 L

*** Available in 25 ml, 100 ml, 500 ml, and 10 L

Prepacked Ion Exchange Cartridges

732-0026	Macro-Prep High Q Support, 1 x 5 ml Econo-Pac cartridge	1 cartridge	Ambient
732-0027	Macro-Prep High Q Support, 5 x 5 ml Econo-Pac cartridges	5 cartridges	Ambient
732-0028	Macro-Prep High Q Support, 5 x 1 ml Econo-Pac cartridges	5 cartridges	Ambient
732-0066	Macro-Prep High S Support, 1 x 5 ml Econo-Pac cartridge	1 cartridge	Ambient
732-0067	Macro-Prep High S Support, 5 x 5 ml Econo-Pac cartridges	5 cartridges	Ambient
732-0068	Macro-Prep High S Support, 5 x 1 ml Econo-Pac cartridges	5 cartridges	Ambient
732-0001	Macro-Prep CM Support, 1 x 5 ml Econo-Pac cartridge	1 cartridge	Ambient
732-0003	Macro-Prep CM Support, 5 x 1 ml Econo-Pac cartridges	5 cartridges	Ambient
732-0005	Macro-Prep CM Support, 5 x 5 ml Econo-Pac cartridges	5 cartridges	Ambient
732-0200	UNOsphere Q Support, 1 x 5 ml Econo-Pac cartridge	1 cartridge	Ambient
732-0201	UNOsphere Q Support, 5 x 5 ml Econo-Pac cartridge	5 cartridges	Ambient
732-0202	UNOsphere Q Support, 5 x 1 ml Econo-Pac cartridge	5 cartridges	Ambient
732-0210	UNOsphere S Support, 1 x 5 ml Econo-Pac cartridge	1 cartridge	Ambient
732-0211	UNOsphere S Support, 5 x 5 ml Econo-Pac cartridge	5 cartridges	Ambient
732-0212	UNOsphere S Support, 5 x 1 ml Econo-Pac cartridge	5 cartridges	Ambient
732-0007	Macro-Prep DEAE Support, 1 x 5 ml Econo-Pac cartridge	1 cartridge	Ambient
732-0008	Macro-Prep DEAE Support, 5 x 5 ml Econo-Pac cartridge	5 cartridges	Ambient
732-0009	Macro-Prep DEAE Support, 5 x 1 ml Econo-Pac cartridge	5 cartridges	Ambient

Prepacked Ion Exchange Columns

720-0001	UNO Q1 Column, 7 x 35 mm	1 column	Ambient
720-0003	UNO Q6 Column, 12 x 53 mm	1 column	Ambient
720-0005	UNO Q12 Column, 15 x 68 mm	1 column	Ambient
720-0009	UNO Q Polishing Column, 4.6 x 10 mm	1 column	Ambient
720-0011	UNO Q1R Column, 7 x 35 mm	1 column	Ambient
720-0013	UNO Q6R Column, 12 x 53 mm	1 column	Ambient
720-0015	UNO Q12R Column, 15 x 68 mm	1 column	Ambient
720-0021	UNO S1 Column, 7 x 35 mm	1 column	Ambient
720-0023	UNO S6 Column, 12 x 53 mm	1 column	Ambient
720-0025	UNO S12 Column, 15 x 68 mm	1 column	Ambient
720-0029	UNO S Polishing Column, 4.6 x 10 mm	1 column	Ambient

Catalog #	Description	Size/ Quantity	Storage Temp.
720-0031	UNO S1R Column, 7 x 35 mm	1 column	Ambient
720-0033	UNO S6R Column, 12 x 53 mm	1 column	Ambient
720-0035	UNO S12R Column, 15 x 68 mm	1 column	Ambient

Size Exclusion Bottled Media

150-4114	Bio-Gel P-2 Gel, fine	100 g	Ambient
150-4115	Bio-Gel P-2 Gel, fine	500 g	Ambient
150-4118	Bio-Gel P-2 Gel, extra fine	100 g	Ambient
150-4120	Bio-Gel P-4 Gel, medium	100 g	Ambient
150-4124	Bio-Gel P-4 Gel, fine	100 g	Ambient
150-4128	Bio-Gel P-4 Gel, extra fine	100 g	Ambient
150-4130	Bio-Gel P-6 Gel, medium	100 g	Ambient
150-4134	Bio-Gel P-6 Gel, fine	100 g	Ambient
150-4138	Bio-Gel P-6 Gel, extra fine	100 g	Ambient
150-0738	Bio-Gel P-6 DG Gel	100 g	Ambient
150-4140	Bio-Gel P-10 Gel, medium	100 g	Ambient
150-4144	Bio-Gel P-10 Gel, fine	100 g	Ambient
150-4150	Bio-Gel P-30 Gel, medium	100 g	Ambient
150-4154	Bio-Gel P-30 Gel, fine	100 g	Ambient
150-4160	Bio-Gel P-60 Gel, medium	100 g	Ambient
150-4164	Bio-Gel P-60 Gel, fine	100 g	Ambient
150-4170	Bio-Gel P-100 Gel, medium	100 g	Ambient
150-4174	Bio-Gel P-100 Gel, fine	100 g	Ambient

Prepacked Size Exclusion (SEC) Columns and Cartridges

732-2010	Econo-Pac 10DG Desalting Columns	30 columns	Ambient
732-0011	Bio-Gel P-6 Support, 1 x 5 ml Econo-Pac cartridge	1 cartridge	Ambient
732-0015	Bio-Gel P-6 Support, 5 x 5 ml Econo-Pac cartridges	5 cartridges	Ambient

Prepacked Size Exclusion (SEC) HPLC Stainless Steel Columns

125-0060	Bio-Sil SEC 125-5 Column, 300 x 7.8 mm	1 column	Ambient
125-0062	Bio-Sil SEC 250-5 Column, 300 x 7.8 mm	1 column	Ambient
125-0064	Bio-Sil SEC 400-5 Column, 300 x 7.8 mm	1 column	Ambient
125-0072	Bio-Sil SEC 125 Guard, 80 x 7.8 mm	1 guard	Ambient
125-0073	Bio-Sil SEC 250 Guard, 80 x 7.8 mm	1 guard	Ambient
125-0074	Bio-Sil SEC 400 Guard, 80 x 7.8 mm	1 guard	Ambient

Prepacked Size Exclusion (SEC) HPLC PEEK Biocompatible Columns

125-0475	Bio-Silect SEC 125-5 Column, 300 x 7.8 mm	1 column	Ambient
125-0476	Bio-Silect SEC 250-5 Column, 300 x 7.8 mm	1 column	Ambient
125-0477	Bio-Silect SEC 400-5 Column, 300 x 7.8 mm	1 column	Ambient
125-0478	Bio-Silect SEC 125 Guard, 50 x 7.8 mm	1 guard	Ambient
125-0479	Bio-Silect SEC 250 Guard, 50 x 7.8 mm	1 guard	Ambient
125-0480	Bio-Silect SEC 400 Guard, 50 x 7.8 mm	1 guard	Ambient

(continues)

Chromatography Media and Columns (cont.)

Catalog #	Description	Size/ Quantity	Storage Temp.
Hydroxyapatite Bottled Media*			
157-0020	CHT Ceramic Hydroxyapatite, Type I, 20 µm	100 g	Ambient
157-0040	CHT Ceramic Hydroxyapatite, Type I, 40 µm	100 g	Ambient
157-0080	CHT Ceramic Hydroxyapatite, Type I, 80 µm	100 g	Ambient
157-2000	CHT Ceramic Hydroxyapatite, Type II, 20 µm	100 g	Ambient
157-4000	CHT Ceramic Hydroxyapatite, Type II, 40 µm	100 g	Ambient
157-8000	CHT Ceramic Hydroxyapatite, Type II, 80 µm	100 g	Ambient
Prepacked Hydroxyapatite Cartridges			
732-0082	CHT Type I Support, 1 x 5 ml Econo-Pac cartridge	1 cartridge	Ambient
732-0084	CHT Type I Support, 5 x 1 ml Econo-Pac cartridges	5 cartridges	Ambient
732-0086	CHT Type I Support, 5 x 5 ml Econo-Pac cartridges	5 cartridges	Ambient
732-0081	CHT Type II Support, 1 x 5 ml Econo-Pac cartridge	1 cartridge	Ambient
732-0083	CHT Type II Support, 5 x 1 ml Econo-Pac cartridges	5 cartridges	Ambient
732-0085	CHT Type II Support, 5 x 5 ml Econo-Pac cartridges	5 cartridges	Ambient
Prepacked Hydroxyapatite Columns			
751-0021	Bio-Scale CHT2 Type I Column, 7 x 52 mm	1 column	Ambient
751-0023	Bio-Scale CHT5 Type I Column, 10 x 64 mm	1 column	Ambient
751-0025	Bio-Scale CHT10 Type I Column, 12 x 88 mm	1 column	Ambient
751-0027	Bio-Scale CHT20 Type I Column, 15 x 113 mm	1 column	Ambient
751-0029	Top-Off Resin Kit CHT Type I, 1 ml	1 kit	Ambient
Hydrophobic Interaction Bottled Media			
156-0080	Macro-Prep Methyl HIC Support	100 ml	Ambient
156-0090	Macro-Prep t-Butyl HIC Support	100 ml	Ambient
Prepacked Hydrophobic Interaction Cartridges			
732-0051	Macro-Prep Methyl HIC Support, 1 x 5 ml Econo-Pac cartridge	1 cartridge	Ambient
732-0053	Macro-Prep Methyl HIC Support, 5 x 1 ml Econo-Pac cartridges	5 cartridges	Ambient
732-0055	Macro-Prep Methyl HIC Support, 5 x 5 ml Econo-Pac cartridges	5 cartridges	Ambient
732-0056	Macro-Prep t-Butyl HIC Support, 1 x 5 ml Econo-Pac cartridge	1 cartridge	Ambient
732-0057	Macro-Prep t-Butyl HIC Support, 5 x 5 ml Econo-Pac cartridges	5 cartridges	Ambient
732-0058	Macro-Prep t-Butyl HIC Support, 5 x 1 ml Econo-Pac cartridges	5 cartridges	Ambient
Activated Affinity Bottled Media			
153-6099	Affi-Gel 10 Gel	25 ml	-20°C
156-0002	Affi-Prep 10 Gel	25 ml	-20°C
153-6051	Affi-Gel 15 Gel	25 ml	-20°C
153-6098	Affi-Gel 10/Affi-Gel 15 Combination	2 x 25 ml each	-20°C

* Also available in 10 g, 1 kg, and 5 kg

Catalog #	Description	Size/ Quantity	Storage Temp.
Ready-to-Use Affinity Media			
153-6047	Affi-Gel Hz Hydrazide Gel	25 ml	4°C
153-6060	Affi-Gel Hz Immunoaffinity Kit	1 kit	4°C
156-0016	Affi-Prep Hz Hydrazide Gel	25 ml	4°C
153-6154	Affi-Gel Protein A Agarose	50 ml	4°C
153-6159	Affi-Gel Protein A MAPS® II Kit	1 kit	4°C
156-0005	Affi-Prep Protein A Support	25 ml	4°C
153-6160	Affi-Prep MAPS II Kit	1 kit	4°C
153-6173	Affi-Gel Heparin Gel	40 ml	4°C
153-7301	Affi-Gel Blue Gel, 50–100 mesh	100 ml	4°C
153-7302	Affi-Gel Blue Gel, 100–200 mesh	100 ml	4°C
153-7304	CM Affi-Gel Blue Gel	100 ml	4°C
153-7307	DEAE Affi-Gel Blue Gel	100 ml	4°C
156-0010	Affi-Prep Polymyxin Support	25 ml	4°C
Prepacked Affinity Columns and Cartridges			
732-0031	DEAE Affi-Gel Blue Support, 1 x 5 ml Econo-Pac cartridge	1 cartridge	Ambient
732-0035	DEAE Affi-Gel Blue Support, 5 x 5 ml Econo-Pac cartridges	5 cartridges	Ambient
732-0091	Affi-Prep Protein A Support, 1 x 5 ml Econo-Pac cartridge	1 cartridge	Ambient
732-0093	Affi-Prep Protein A Support, 5 x 5 ml Econo-Pac cartridges	5 cartridges	Ambient
732-0101	Affi-Gel Blue Support, 1 x 5 ml Econo-Pac cartridge	1 cartridge	Ambient
732-0105	Affi-Gel Blue Support, 5 x 5 ml Econo-Pac cartridges	5 cartridges	Ambient
732-2022	Econo-Pac Protein A Columns	5 columns	Ambient
Chromatography Standards			
125-0561	Protein Standard for Anion Exchange	6 vials	4°C
125-0562	Protein Standard for Cation Exchange	6 vials	4°C
125-1901	Protein Standard for Size Exclusion Standard	6 vials	4°C
151-1905	Protein Standard for Hydrophobic Interaction	6 vials	4°C
Empty Chromatography Columns			
731-1550	Poly-Prep® Chromatography Columns, empty	50 columns	Ambient
732-1010	Econo-Pac Chromatography Columns, empty	50 columns	Ambient
751-0081	Bio-Scale MT2 Column, 7 x 52 mm	1 column	Ambient
751-0083	Bio-Scale MT5 Column, 10 x 64 mm	1 column	Ambient
751-0085	Bio-Scale MT10 Column, 12 x 88 mm	1 column	Ambient
751-0087	Bio-Scale MT20 Column, 15 x 113 mm	1 column	Ambient

Chromatography Media and Columns (cont.)

Catalog #	Description	Size/ Quantity	Storage Temp.
Econo-Column® Chromatography Columns*			
731-0003	Econo-Column Funnels, 250 ml	5 funnels	Ambient
737-0506	Econo-Column Chromatography Column, 0.5 x 5 cm	5 columns	Ambient
737-0511	Econo-Column Chromatography Column, 0.5 x 10 cm	5 columns	Ambient
737-0706	Econo-Column Chromatography Column, 0.7 x 5 cm	5 columns	Ambient
737-0711	Econo-Column Chromatography Column, 0.7 x 10 cm	5 columns	Ambient
737-0716	Econo-Column Chromatography Column, 0.7 x 15 cm	5 columns	Ambient
737-0721	Econo-Column Chromatography Column, 0.7 x 20 cm	5 columns	Ambient
737-1006	Econo-Column Chromatography Column, 1 x 5 cm	5 columns	Ambient
737-1011	Econo-Column Chromatography Column, 1 x 10 cm	5 columns	Ambient
737-1012	Econo-Column Chromatography Column, 1 x 10 cm	2 columns	Ambient
737-1021	Econo-Column Chromatography Column, 1 x 20 cm	5 columns	Ambient
737-1022	Econo-Column Chromatography Column, 1 x 20 cm	2 columns	Ambient
737-1031	Econo-Column Chromatography Column, 1 x 30 cm	5 columns	Ambient
737-1506	Econo-Column Chromatography Column, 1.5 x 5 cm	5 columns	Ambient
737-1511	Econo-Column Chromatography Column, 1.5 x 10 cm	5 columns	Ambient
737-1512	Econo-Column Chromatography Column, 1.5 x 10 cm	2 columns	Ambient
737-1516	Econo-Column Chromatography Column, 1.5 x 15 cm	5 columns	Ambient
737-1517	Econo-Column Chromatography Column, 1.5 x 15 cm	2 columns	Ambient
737-1521	Econo-Column Chromatography Column, 1.5 x 20 cm	5 columns	Ambient
737-1522	Econo-Column Chromatography Column, 1.5 x 20 cm	2 columns	Ambient
737-1531	Econo-Column Chromatography Column, 1.5 x 30 cm	5 columns	Ambient
737-1532	Econo-Column Chromatography Column, 1.5 x 30 cm	2 columns	Ambient
737-1552	Econo-Column Chromatography Column, 1.5 x 50 cm	2 columns	Ambient
737-1576	Econo-Column Chromatography Column, 1.5 x 75 cm	2 columns	Ambient
737-2511	Econo-Column Chromatography Column, 2.5 x 10 cm	5 columns	Ambient
737-2512	Econo-Column Chromatography Column, 2.5 x 10 cm	2 columns	Ambient
737-2521	Econo-Column Chromatography Column, 2.5 x 20 cm	5 columns	Ambient
737-2522	Econo-Column Chromatography Column, 2.5 x 20 cm	2 columns	Ambient
737-2532	Econo-Column Chromatography Column, 2.5 x 30 cm	2 columns	Ambient
737-2551	Econo-Column Chromatography Column, 2.5 x 50 cm	2 columns	Ambient
737-5011	Econo-Column Chromatography Column, 5 x 10 cm	1 column	Ambient
737-5021	Econo-Column Chromatography Column, 5 x 20 cm	1 column	Ambient
737-5031	Econo-Column Chromatography Column, 5 x 30 cm	1 column	Ambient

Catalog #	Description	Size/ Quantity	Storage Temp.
Jacketed Econo-Column Columns			
737-6108	Standard Jacketed Column, 0.7 x 15 cm	1 column	Ambient
737-6116	Standard Jacketed Column, 1 x 15 cm	1 column	Ambient
737-6131	Standard Jacketed Column, 1 x 30 cm	1 column	Ambient
737-6151	Standard Jacketed Column, 1.5 x 50 cm	1 column	Ambient
737-6201	Open-Ended Jacketed Column, with 2 flow adaptors, 1 x 30 cm	1 column	Ambient
Flow Adaptors			
738-0014	Flow Adaptor, 1.0 cm column ID, 1–7 cm functional length	1 adaptor	Ambient
738-0015	Flow Adaptor, 1.0 cm column ID, 1–14 cm functional length	1 adaptor	Ambient
738-0016	Flow Adaptor, 1.5 cm column ID, 1–14 cm functional length	1 adaptor	Ambient
738-0017	Flow Adaptor, 2.5 cm column ID, 1–14 cm functional length	1 adaptor	Ambient
738-0018	Flow Adaptor, 5.0 cm column ID, 1–14 cm functional length	1 adaptor	Ambient
738-0019	Econo-Pac Flow Adaptor, 1.5 cm column ID	1 adaptor	Ambient

* For longer column lengths, refer to the Life Science Research Products catalog.

ReadyStrip™ IPG Strip Selection Chart

	pH Gradients								Relative Focusing Power*					ReadyStrip IEF Buffer					
	3	4	5	6	7	8	9	10	7 cm	11 cm	17 cm	18 cm	24 cm	3–10	7–10	3.9–5.1	4.7–5.9	5.5–6.7	6.3–8.3
Broad Range																			
3–10									1.0x	1.6x	2.4x	2.6x	3.4x	●					
3–10 nonlinear (NL)									1.0x	1.6x	2.4x	2.6x	3.4x	●					
Narrow Range																			
3–6									2.3x	3.7x	5.7x	6.0x	8.0x	●					
5–8									2.3x	3.7x	5.7x	6.0x	8.0x	●					
7–10									2.3x	3.7x	5.7x	6.0x	—		●				
4–7									2.3x	3.7x	5.7x	6.0x	8.0x	●					
Micro Range																			
3.9–5.1									5.8x	9.2x	14.2x	15.0x	—			●			
4.7–5.9									5.8x	9.2x	14.2x	15.0x	—				●		
5.5–6.7									5.8x	9.2x	14.2x	15.0x	—					●	
6.3–8.3									3.5x	5.5x	8.5x	9.0x	—						●

* Relative Focusing Power

Relative focusing power expresses the enhanced resolution expected in the first dimension when using IPG strips of different lengths or pH ranges. The 7 cm pH 3–10 IPG strip is arbitrarily assigned a baseline focusing power of 1.0 in order to calculate the relative focusing powers of the other strips, as shown in the table above. The calculations are done as follows:

- **When two strips have the same pH range but differ in length** — the ratio of the strip lengths gives the relative focusing power of the longer strip. Thus, compared to a 7 cm pH 3–10 strip, the 11 cm pH 3–10 strip has a relative focusing power of $11 \text{ cm} / 7 \text{ cm} = 1.6$
- **When two strips have the same length but differ in pH range** — the ratio of the pH ranges (in number of pH units) gives the relative focusing power of the strip with the smaller range. Thus, compared to a 7 cm pH 3–10 strip, the 7 cm pH 5–8 strip has a relative focusing power of $(7 \text{ pH units}) / (3 \text{ pH units}) = 2.3$

ReadyStrip IPG Strips and Related Buffers

Catalog #	Description	Size/ Quantity	Storage Temp.
7 cm IPG Strips			
163-2000	ReadyStrip IPG Strips, pH 3–10, 7 cm	12 strips	–20°C
163-2001	ReadyStrip IPG Strips, pH 4–7, 7 cm	12 strips	–20°C
163-2002	ReadyStrip IPG Strips, pH 3–10 nonlinear, 7 cm	12 strips	–20°C
163-2003	ReadyStrip IPG Strips, pH 3–6, 7 cm	12 strips	–20°C
163-2004	ReadyStrip IPG Strips, pH 5–8, 7 cm	12 strips	–20°C
163-2005	ReadyStrip IPG Strips, pH 7–10, 7 cm	12 strips	–20°C
163-2028	ReadyStrip IPG Strips, pH 3.9–5.1, 7 cm	12 strips	–20°C
163-2029	ReadyStrip IPG Strips, pH 4.7–5.9, 7 cm	12 strips	–20°C
163-2030	ReadyStrip IPG Strips, pH 5.5–6.7, 7 cm	12 strips	–20°C
163-2031	ReadyStrip IPG Strips, pH 6.3–8.3, 7 cm	12 strips	–20°C
11 cm IPG Strips			
163-2014	ReadyStrip IPG Strips, pH 3–10, 11 cm	12 strips	–20°C
163-2015	ReadyStrip IPG Strips, pH 4–7, 11 cm	12 strips	–20°C
163-2016	ReadyStrip IPG Strips, pH 3–10 nonlinear, 11 cm	12 strips	–20°C
163-2017	ReadyStrip IPG Strips, pH 3–6, 11 cm	12 strips	–20°C
163-2018	ReadyStrip IPG Strips, pH 5–8, 11 cm	12 strips	–20°C
163-2019	ReadyStrip IPG Strips, pH 7–10, 11 cm	12 strips	–20°C
163-2024	ReadyStrip IPG Strips, pH 3.9–5.1, 11 cm	12 strips	–20°C
163-2025	ReadyStrip IPG Strips, pH 4.7–5.9, 11 cm	12 strips	–20°C
163-2026	ReadyStrip IPG Strips, pH 5.5–6.7, 11 cm	12 strips	–20°C
163-2027	ReadyStrip IPG Strips, pH 6.3–8.3, 11 cm	12 strips	–20°C
17 cm IPG Strips			
163-2007	ReadyStrip IPG Strips, pH 3–10, 17 cm	12 strips	–20°C
163-2008	ReadyStrip IPG Strips, pH 4–7, 17 cm	12 strips	–20°C
163-2009	ReadyStrip IPG Strips, pH 3–10 nonlinear, 17 cm	12 strips	–20°C
163-2010	ReadyStrip IPG Strips, pH 3–6, 17 cm	12 strips	–20°C
163-2011	ReadyStrip IPG Strips, pH 5–8, 17 cm	12 strips	–20°C
163-2012	ReadyStrip IPG Strips, pH 7–10, 17 cm	12 strips	–20°C
163-2020	ReadyStrip IPG Strips, pH 3.9–5.1, 17 cm	12 strips	–20°C
163-2021	ReadyStrip IPG Strips, pH 4.7–5.9, 17 cm	12 strips	–20°C
163-2022	ReadyStrip IPG Strips, pH 5.5–6.7, 17 cm	12 strips	–20°C
163-2023	ReadyStrip IPG Strips, pH 6.3–8.3, 17 cm	12 strips	–20°C

Catalog #	Description	Size/ Quantity	Storage Temp.
18 cm IPG Strips			
163-2032	ReadyStrip IPG Strips, pH 3–10, 18 cm	12 strips	–20°C
163-2033	ReadyStrip IPG Strips, pH 3–10 nonlinear, 18 cm	12 strips	–20°C
163-2034	ReadyStrip IPG Strips, pH 4–7, 18 cm	12 strips	–20°C
163-2035	ReadyStrip IPG Strips, pH 3–6, 18 cm	12 strips	–20°C
163-2036	ReadyStrip IPG Strips, pH 5–8, 18 cm	12 strips	–20°C
163-2037	ReadyStrip IPG Strips, pH 7–10, 18 cm	12 strips	–20°C
163-2038	ReadyStrip IPG Strips, pH 3.9–5.1, 18 cm	12 strips	–20°C
163-2039	ReadyStrip IPG Strips, pH 4.7–5.9, 18 cm	12 strips	–20°C
163-2040	ReadyStrip IPG Strips, pH 5.5–6.7, 18 cm	12 strips	–20°C
163-2041	ReadyStrip IPG Strips, pH 6.3–8.3, 18 cm	12 strips	–20°C
24 cm IPG Strips			
163-2042	ReadyStrip IPG Strips, pH 3–10, 24 cm	12 strips	–20°C
163-2043	ReadyStrip IPG Strips, pH 3–10 nonlinear, 24 cm	12 strips	–20°C
163-2044	ReadyStrip IPG Strips, pH 4–7, 24 cm	12 strips	–20°C
163-2045	ReadyStrip IPG Strips, pH 3–6, 24 cm	12 strips	–20°C
163-2046	ReadyStrip IPG Strips, pH 5–8, 24 cm	12 strips	–20°C
Bio-Lyte® IEF Buffers			
163-2094	100x Bio-Lyte 3–10 Ampholyte	1 ml	4°C
163-2093	100x ReadyStrip 7–10 Buffer	1 ml	4°C
163-2098	100x ReadyStrip 3.9–5.1 Buffer	1 ml	4°C
163-2097	100x ReadyStrip 4.7–5.9 Buffer	1 ml	4°C
163-2096	100x ReadyStrip 5.5–6.7 Buffer	1 ml	4°C
163-2095	100x ReadyStrip 6.3–8.3 Buffer	1 ml	4°C
Rehydration/Sample Buffers			
163-2106	ReadyPrep 2-D Starter Kit Rehydration/Sample Buffer 1	10 ml	4°C
163-2083	ReadyPrep Rehydration/Sample Buffer 1	10 ml	Ambient
Equilibration Buffers			
163-2107	ReadyPrep 2-D Starter Kit Rehydration/Sample Buffer I, with DTT	10 ml	4°C
163-2108	ReadyPrep 2-D Starter Kit Rehydration/Sample Buffer II, without DTT	20 ml	4°C

Criterion™ XT Precast Gels

Catalog #	Description	Size/ Quantity	Storage Temp.	Shelf Life
Criterion XT Bis-Tris Gels				
345-0111	Criterion XT Bis-Tris Gel, 10%, 12+2 well, 45 µl	1 gel	Ambient	12 mo
345-0112	Criterion XT Bis-Tris Gel, 10%, 18-well, 30 µl	1 gel	Ambient	12 mo
345-0113	Criterion XT Bis-Tris Gel, 10%, 26-well, 15 µl	1 gel	Ambient	12 mo
345-0114	Criterion XT Bis-Tris Gel, 10%, prep+2 well, 11 cm IPG strip	1 gel	Ambient	12 mo
345-0115	Criterion XT Bis-Tris Gel, 10%, prep+2 well, 11 cm IPG strip	1 gel	Ambient	12 mo
345-0116	Criterion XT Bis-Tris Gel, 10%, IPG well, 11 cm IPG strip	1 gel	Ambient	12 mo
345-0117	Criterion XT Bis-Tris Gel, 12%, 12+2 well, 45 µl	1 gel	Ambient	12 mo
345-0118	Criterion XT Bis-Tris Gel, 12%, 18-well, 30 µl	1 gel	Ambient	12 mo
345-0119	Criterion XT Bis-Tris Gel, 12%, 26-well, 15 µl	1 gel	Ambient	12 mo
345-0120	Criterion XT Bis-Tris Gel, 12%, prep+2 well, 800 µl	1 gel	Ambient	12 mo
345-0121	Criterion XT Bis-Tris Gel, 12%, IPG+1 well, 11 cm IPG strip	1 gel	Ambient	12 mo
345-0122	Criterion XT Bis-Tris Gel, 12%, IPG well, 11 cm IPG strip	1 gel	Ambient	12 mo
345-0123	Criterion XT Bis-Tris Gel, 4–12%, 12+2 well, 45 µl	1 gel	Ambient	12 mo
345-0124	Criterion XT Bis-Tris Gel, 4–12%, 18-well, 30 µl	1 gel	Ambient	12 mo
345-0125	Criterion XT Bis-Tris Gel, 4–12%, 26-well, 15 µl	1 gel	Ambient	12 mo
345-0126	Criterion XT Bis-Tris Gel, 4–12%, prep+2 well, 800 µl	1 gel	Ambient	12 mo
345-0127	Criterion XT Bis-Tris Gel, 4–12%, IPG+1 well, 11 cm IPG strip	1 gel	Ambient	12 mo
345-0128	Criterion XT Bis-Tris Gel, 4–12%, IPG well, 11 cm IPG strip	1 gel	Ambient	12 mo
Criterion XT Tris-Acetate Gels				
345-0129	Criterion XT Tris-Acetate Gel, 3–8%, 12+2 well, 45 µl	1 gel	4°C	8 mo
345-0130	Criterion XT Tris-Acetate Gel, 3–8%, 18-well, 30 µl	1 gel	4°C	8 mo
345-0131	Criterion XT Tris-Acetate Gel, 3–8%, 26-well, 15 µl	1 gel	4°C	8 mo
345-0132	Criterion XT Tris-Acetate Gel, 3–8%, prep+2 well, 800 µl	1 gel	4°C	8 mo
345-0133	Criterion XT Tris-Acetate Gel, 3–8%, IPG+1 well, 11 cm IPG strip	1 gel	4°C	8 mo
345-0134	Criterion XT Tris-Acetate Gel, 3–8%, IPG well, 11 cm IPG strip	1 gel	4°C	8 mo
345-0135	Criterion XT Tris-Acetate Gel, 7%, 12+2 well, 45 µl	1 gel	4°C	8 mo
345-0136	Criterion XT Tris-Acetate Gel, 7%, 18-well, 30 µl	1 gel	4°C	8 mo
345-0137	Criterion XT Tris-Acetate Gel, 7%, 26-well, 15 µl	1 gel	4°C	8 mo
345-0138	Criterion XT Tris-Acetate Gel, 7%, prep+2 well, 800 µl	1 gel	4°C	8 mo
345-0139	Criterion XT Tris-Acetate Gel, 7%, IPG+1 well, 11 cm IPG strip	1 gel	4°C	8 mo
345-0140	Criterion XT Tris-Acetate Gel, 7%, IPG well, 11 cm IPG strip	1 gel	4°C	8 mo
Criterion XT Buffers and Reagents				
161-0788	XT MOPS Running Buffer, 20x	500 ml	Ambient	
161-0789	XT MES Running Buffer, 20x	500 ml	Ambient	
161-0790	XT Tricine Running Buffer, 20x	500 ml	Ambient	
161-0791	XT Sample Buffer, 4x	10 ml	Ambient	
161-0792	XT Reducing Agent, 20x	1 ml	Ambient	
161-0793	XT MOPS Buffer Kit	1 kit	Ambient	
161-0796	XT MES Buffer Kit	1 kit	Ambient	
161-0797	XT Tricine Buffer Kit	1 kit	Ambient	

Criterion Precast Gels

Catalog #	Description	Size/ Quantity	Storage Temp.	Shelf Life
Criterion Tris-HCl Gels				
345-0001	Criterion Tris-HCl Gel, 5%, 12+2 well, 45 µl	1 gel	4°C	12 wk
345-0002	Criterion Tris-HCl Gel, 5%, 18-well, 30 µl	1 gel	4°C	12 wk
345-0003	Criterion Tris-HCl Gel, 5%, 26-well, 15 µl	1 gel	4°C	12 wk
345-0004	Criterion Tris-HCl Gel, 5%, prep+2 well, 800 µl	1 gel	4°C	12 wk
345-0005	Criterion Tris-HCl Gel, 7.5%, 12+2 well, 45 µl	1 gel	4°C	12 wk
345-0006	Criterion Tris-HCl Gel, 7.5%, 18-well, 30 µl	1 gel	4°C	12 wk
345-0007	Criterion Tris-HCl Gel, 7.5%, 26-well, 15 µl	1 gel	4°C	12 wk
345-0008	Criterion Tris-HCl Gel, 7.5%, prep+2 well, 800 µl	1 gel	4°C	12 wk
345-0009	Criterion Tris-HCl Gel, 10%, 12+2 well, 45 µl	1 gel	4°C	12 wk
345-0010	Criterion Tris-HCl Gel, 10%, 18-well, 30 µl	1 gel	4°C	12 wk
345-0011	Criterion Tris-HCl Gel, 10%, 26-well, 15 µl	1 gel	4°C	12 wk
345-0012	Criterion Tris-HCl Gel, 10%, prep+2 well, 800 µl	1 gel	4°C	12 wk
345-0013	Criterion Tris-HCl Gel, 10%, 11 cm IPG well	1 gel	4°C	12 wk
345-0014	Criterion Tris-HCl Gel, 12.5%, 12+2 well, 45 µl	1 gel	4°C	12 wk
345-0015	Criterion Tris-HCl Gel, 12.5%, 18-well, 30 µl	1 gel	4°C	12 wk
345-0016	Criterion Tris-HCl Gel, 12.5%, 26-well, 15 µl	1 gel	4°C	12 wk
345-0017	Criterion Tris-HCl Gel, 12.5%, prep+2 well, 800 µl	1 gel	4°C	12 wk
345-0018	Criterion Tris-HCl Gel, 12.5%, 11 cm IPG well	1 gel	4°C	12 wk
345-0019	Criterion Tris-HCl Gel, 15%, 12+2 well, 45 µl	1 gel	4°C	12 wk
345-0020	Criterion Tris-HCl Gel, 15%, 18-well, 30 µl	1 gel	4°C	12 wk
345-0021	Criterion Tris-HCl Gel, 15%, 26-well, 15 µl	1 gel	4°C	12 wk
345-0022	Criterion Tris-HCl Gel, 15%, prep+2 well, 800 µl	1 gel	4°C	12 wk
345-0023	Criterion Tris-HCl Gel, 18%, 12+2 well, 45 µl	1 gel	4°C	12 wk
345-0024	Criterion Tris-HCl Gel, 18%, 18-well, 30 µl	1 gel	4°C	12 wk
345-0025	Criterion Tris-HCl Gel, 18%, 26-well, 15 µl	1 gel	4°C	12 wk
345-0026	Criterion Tris-HCl Gel, 18%, prep+2 well, 800 µl	1 gel	4°C	12 wk
345-0027	Criterion Tris-HCl Gel, 4–15%, 12+2 well, 45 µl	1 gel	4°C	12 wk
345-0028	Criterion Tris-HCl Gel, 4–15%, 18-well, 30 µl	1 gel	4°C	12 wk
345-0029	Criterion Tris-HCl Gel, 4–15%, 26-well, 15 µl	1 gel	4°C	12 wk
345-0030	Criterion Tris-HCl Gel, 4–15%, prep+2 well, 800 µl	1 gel	4°C	12 wk
345-0031	Criterion Tris-HCl Gel, 4–15%, 11 cm IPG well	1 gel	4°C	12 wk
345-0032	Criterion Tris-HCl Gel, 4–20%, 12+2 well, 45 µl	1 gel	4°C	12 wk
345-0033	Criterion Tris-HCl Gel, 4–20%, 18-well, 30 µl	1 gel	4°C	12 wk
345-0034	Criterion Tris-HCl Gel, 4–20%, 26-well, 15 µl	1 gel	4°C	12 wk
345-0035	Criterion Tris-HCl Gel, 4–20%, prep+2 well, 800 µl	1 gel	4°C	12 wk
345-0036	Criterion Tris-HCl Gel, 4–20%, 11 cm IPG well	1 gel	4°C	12 wk

(continues)

Criterion Precast Gels (cont.)

Catalog #	Description	Size/ Quantity	Storage Temp.	Shelf Life
Criterion Tris-HCl Gels (cont.)				
345-0037	Criterion Tris-HCl Gel, 8–16%, 12+2 well, 45 µl	1 gel	4°C	12 wk
345-0038	Criterion Tris-HCl Gel, 8–16%, 18-well, 30 µl	1 gel	4°C	12 wk
345-0039	Criterion Tris-HCl Gel, 8–16%, 26-well, 15 µl	1 gel	4°C	12 wk
345-0040	Criterion Tris-HCl Gel, 8–16%, prep+2 well, 800 µl	1 gel	4°C	12 wk
345-0041	Criterion Tris-HCl Gel, 8–16%, 11 cm IPG well	1 gel	4°C	12 wk
345-0042	Criterion Tris-HCl Gel, 10–20%, 12+2 well, 45 µl	1 gel	4°C	12 wk
345-0043	Criterion Tris-HCl Gel, 10–20%, 18-well, 30 µl	1 gel	4°C	12 wk
345-0044	Criterion Tris-HCl Gel, 10–20%, 26-well, 15 µl	1 gel	4°C	12 wk
345-0045	Criterion Tris-HCl Gel, 10–20%, prep+2 well, 800 µl	1 gel	4°C	12 wk
345-0046	Criterion Tris-HCl Gel, 10–20%, 11 cm IPG well	1 gel	4°C	12 wk
345-0101	Criterion Tris-HCl Gel, 10%, 11 cm IPG+1 well	1 gel	4°C	12 wk
345-0102	Criterion Tris-HCl Gel, 12.5%, 11 cm IPG+1 well	1 gel	4°C	12 wk
345-0103	Criterion Tris-HCl Gel, 4–15%, 11 cm IPG+1 well	1 gel	4°C	12 wk
345-0104	Criterion Tris-HCl Gel, 4–20%, 11 cm IPG+1 well	1 gel	4°C	12 wk
345-0105	Criterion Tris-HCl Gel, 8–16%, 11 cm IPG+1 well	1 gel	4°C	12 wk
345-0106	Criterion Tris-HCl Gel, 10.5–14%, 11 cm IPG+1 well	1 gel	4°C	12 wk
345-0107	Criterion Tris-HCl Gel, 10–20%, 11 cm IPG+1 well	1 gel	4°C	12 wk
Criterion TBE Gels				
345-0047	Criterion TBE Gel, 5%, 12+2 well, 45 µl	1 gel	4°C	12 wk
345-0048	Criterion TBE Gel, 5%, 18-well, 30 µl	1 gel	4°C	12 wk
345-0049	Criterion TBE Gel, 5%, 26-well, 15 µl	1 gel	4°C	12 wk
345-0050	Criterion TBE Gel, 5%, prep+2 well, 800 µl	1 gel	4°C	12 wk
345-0051	Criterion TBE Gel, 5%, 12+2 well, 45 µl	1 gel	4°C	12 wk
345-0052	Criterion TBE Gel, 5%, 18-well, 30 µl	1 gel	4°C	12 wk
345-0053	Criterion TBE Gel, 5%, 26-well, 15 µl	1 gel	4°C	12 wk
345-0054	Criterion TBE Gel, 5%, prep+2 well, 800 µl	1 gel	4°C	12 wk
345-0055	Criterion TBE Gel, 5%, 12+2 well, 45 µl	1 gel	4°C	12 wk
345-0056	Criterion TBE Gel, 5%, 18-well, 30 µl	1 gel	4°C	12 wk
345-0057	Criterion TBE Gel, 5%, 26-well, 15 µl	1 gel	4°C	12 wk
345-0058	Criterion TBE Gel, 5%, prep+2 well, 800 µl	1 gel	4°C	12 wk
345-0059	Criterion TBE Gel, 5%, 12+2 well, 45 µl	1 gel	4°C	12 wk
345-0060	Criterion TBE Gel, 5%, 18-well, 30 µl	1 gel	4°C	12 wk
345-0061	Criterion TBE Gel, 5%, 26-well, 15 µl	1 gel	4°C	12 wk
345-0062	Criterion TBE Gel, 5%, prep+2 well, 800 µl	1 gel	4°C	12 wk

Catalog #	Description	Size/ Quantity	Storage Temp.	Shelf Life
Criterion Tris-Tricine/Peptide Gels				
345-0063	Criterion Tris-Tricine/Peptide Gel, 16.5%, 12+2 well, 45 µl	1 gel	4°C	12 wk
345-0064	Criterion Tris-Tricine/Peptide Gel, 16.5%, 18-well, 30 µl	1 gel	4°C	12 wk
345-0065	Criterion Tris-Tricine/Peptide Gel, 16.5%, 26-well, 15 µl	1 gel	4°C	12 wk
345-0066	Criterion Tris-Tricine/Peptide Gel, 16.5%, prep+2 well, 800 µl	1 gel	4°C	12 wk
345-0067	Criterion Tris-Tricine/Peptide Gel, 10–20%, 12+2 well, 45 µl	1 gel	4°C	12 wk
345-0068	Criterion Tris-Tricine/Peptide Gel, 10–20%, 18-well, 30 µl	1 gel	4°C	12 wk
345-0069	Criterion Tris-Tricine/Peptide Gel, 10–20%, 26-well, 15 µl	1 gel	4°C	12 wk
345-0070	Criterion Tris-Tricine/Peptide Gel, 10–20%, prep+2 well, 800 µl	1 gel	4°C	12 wk
Criterion IEF Gels				
345-0071	Criterion IEF Gel, pH 3–10, 12+2 well, 45 µl	1 gel	4°C	24 wk
345-0072	Criterion IEF Gel, pH 3–10, 18-well, 30 µl	1 gel	4°C	24 wk
345-0073	Criterion IEF Gel, pH 3–10, 26-well, 15 µl	1 gel	4°C	24 wk
345-0074	Criterion IEF Gel, pH 5–8, prep+2 well, 800 µl	1 gel	4°C	24 wk
345-0075	Criterion IEF Gel, pH 5–8, 12+2 well, 45 µl	1 gel	4°C	24 wk
345-0076	Criterion IEF Gel, pH 5–8, 18-well, 30 µl	1 gel	4°C	24 wk
345-0077	Criterion IEF Gel, pH 5–8, 26-well, 15 µl	1 gel	4°C	24 wk
345-0078	Criterion IEF Gel, pH 5–8, prep+2 well, 800 µl	1 gel	4°C	24 wk
Criterion Zymogram Gels				
345-0079	Criterion Zymogram Gel, 10%, gelatin, 12+2 well, 45 µl	1 gel	4°C	8 wk
345-0080	Criterion Zymogram Gel, 10%, gelatin, 18-well, 30 µl	1 gel	4°C	8 wk
345-0081	Criterion Zymogram Gel, 10%, gelatin, 26-well, 15 µl	1 gel	4°C	8 wk
345-0082	Criterion Zymogram Gel, 12.5%, casein, 12+2 well, 45 µl	1 gel	4°C	8 wk
345-0083	Criterion Zymogram Gel, 12.5%, casein, 18-well, 30 µl	1 gel	4°C	8 wk
345-0084	Criterion Zymogram Gel, 12.5%, casein, 26-well, 15 µl	1 gel	4°C	8 wk
Criterion TBE-Urea Gels				
345-0085	Criterion TBE-Urea Gel, 5%, 12+2 well, 45 µl	1 gel	4°C	8 wk
345-0086	Criterion TBE-Urea Gel, 5%, 18-well, 30 µl	1 gel	4°C	8 wk
345-0087	Criterion TBE-Urea Gel, 5%, 26-well, 15 µl	1 gel	4°C	8 wk
345-0088	Criterion TBE-Urea Gel, 10%, 12+2 well, 45 µl	1 gel	4°C	8 wk
345-0089	Criterion TBE-Urea Gel, 10%, 18-well, 30 µl	1 gel	4°C	8 wk
345-0090	Criterion TBE-Urea Gel, 10%, 26-well, 15 µl	1 gel	4°C	8 wk
345-0091	Criterion TBE-Urea Gel, 15%, 12+2 well, 45 µl	1 gel	4°C	8 wk
345-0092	Criterion TBE-Urea Gel, 15%, 18-well, 30 µl	1 gel	4°C	8 wk
345-0093	Criterion TBE-Urea Gel, 15%, 26-well, 15 µl	1 gel	4°C	8 wk

Ready Gel Precast Gels

Catalog #	Description	Size/ Quantity	Storage Temp.	Shelf Life
Ready Gel Tris-HCl Gels				
161-1100	Ready Gel Tris-HCl Gel, 7.5%, 10-well, 30 µl	1 gel	4°C	12 wk
161-1101	Ready Gel Tris-HCl Gel, 10%, 10-well, 30 µl	1 gel	4°C	12 wk
161-1102	Ready Gel Tris-HCl Gel, 12%, 10-well, 30 µl	1 gel	4°C	12 wk
161-1103	Ready Gel Tris-HCl Gel, 15%, 10-well, 30 µl	1 gel	4°C	12 wk
161-1104	Ready Gel Tris-HCl Gel, 4–15%, 10-well, 30 µl	1 gel	4°C	12 wk
161-1105	Ready Gel Tris-HCl Gel, 4–20%, 10-well, 30 µl	1 gel	4°C	12 wk
161-1106	Ready Gel Tris-HCl Gel, 10–20%, 10-well, 30 µl	1 gel	4°C	12 wk
161-1118	Ready Gel Tris-HCl Gel, 7.5%, 15-well, 15 µl	1 gel	4°C	12 wk
161-1119	Ready Gel Tris-HCl Gel, 10%, 15-well, 15 µl	1 gel	4°C	12 wk
161-1120	Ready Gel Tris-HCl Gel, 12%, 15-well, 15 µl	1 gel	4°C	12 wk
161-1121	Ready Gel Tris-HCl Gel, 15%, 15-well, 15 µl	1 gel	4°C	12 wk
161-1122	Ready Gel Tris-HCl Gel, 4–15%, 15-well, 15 µl	1 gel	4°C	12 wk
161-1123	Ready Gel Tris-HCl Gel, 4–20%, 15-well, 15 µl	1 gel	4°C	12 wk
161-1124	Ready Gel Tris-HCl Gel, 10–20%, 15-well, 15 µl	1 gel	4°C	12 wk
161-1136	Ready Gel Tris-HCl Gel, 7.5%, prep well, 450 µl	1 gel	4°C	12 wk
161-1137	Ready Gel Tris-HCl Gel, 10%, prep well, 450 µl	1 gel	4°C	12 wk
161-1138	Ready Gel Tris-HCl Gel, 12%, prep well, 450 µl	1 gel	4°C	12 wk
161-1139	Ready Gel Tris-HCl Gel, 15%, prep well, 450 µl	1 gel	4°C	12 wk
161-1140	Ready Gel Tris-HCl Gel, 4–15%, prep well, 450 µl	1 gel	4°C	12 wk
161-1141	Ready Gel Tris-HCl Gel, 4–20%, prep well, 450 µl	1 gel	4°C	12 wk
161-1142	Ready Gel Tris-HCl Gel, 10–20%, prep well, 450 µl	1 gel	4°C	12 wk
161-1154	Ready Gel Tris-HCl Gel, 7.5%, 10-well, 50 µl	1 gel	4°C	12 wk
161-1155	Ready Gel Tris-HCl Gel, 10%, 10-well, 50 µl	1 gel	4°C	12 wk
161-1156	Ready Gel Tris-HCl Gel, 12%, 10-well, 50 µl	1 gel	4°C	12 wk
161-1157	Ready Gel Tris-HCl Gel, 15%, 10-well, 50 µl	1 gel	4°C	12 wk
161-1158	Ready Gel Tris-HCl Gel, 4–15%, 10-well, 50 µl	1 gel	4°C	12 wk
161-1159	Ready Gel Tris-HCl Gel, 4–20%, 10-well, 50 µl	1 gel	4°C	12 wk
161-1160	Ready Gel Tris-HCl Gel, 10–20%, 10-well, 50 µl	1 gel	4°C	12 wk
161-1172	Ready Gel Tris-HCl Gel, 7.5%, 12-well, 20 µl	1 gel	4°C	12 wk
161-1173	Ready Gel Tris-HCl Gel, 10%, 12-well, 20 µl	1 gel	4°C	12 wk
161-1174	Ready Gel Tris-HCl Gel, 12%, 12-well, 20 µl	1 gel	4°C	12 wk
161-1175	Ready Gel Tris-HCl Gel, 15%, 12-well, 20 µl	1 gel	4°C	12 wk
161-1176	Ready Gel Tris-HCl Gel, 4–15%, 12-well, 20 µl	1 gel	4°C	12 wk
161-1177	Ready Gel Tris-HCl Gel, 4–20%, 12-well, 20 µl	1 gel	4°C	12 wk
161-1178	Ready Gel Tris-HCl Gel, 10–20%, 12-well, 20 µl	1 gel	4°C	12 wk
161-1190	Ready Gel Tris-HCl Gel, 7.5%, 9-well, 30 µl	1 gel	4°C	12 wk
161-1191	Ready Gel Tris-HCl Gel, 10%, 9-well, 30 µl	1 gel	4°C	12 wk
161-1192	Ready Gel Tris-HCl Gel, 12%, 9-well, 30 µl	1 gel	4°C	12 wk
161-1193	Ready Gel Tris-HCl Gel, 15%, 9-well, 30 µl	1 gel	4°C	12 wk
161-1194	Ready Gel Tris-HCl Gel, 4–15%, 9-well, 30 µl	1 gel	4°C	12 wk
161-1195	Ready Gel Tris-HCl Gel, 4–20%, 9-well, 30 µl	1 gel	4°C	12 wk
161-1196	Ready Gel Tris-HCl Gel, 10–20%, 9-well, 30 µl	1 gel	4°C	12 wk

Catalog #	Description	Size/ Quantity	Storage Temp.	Shelf Life
161-1210	Ready Gel Tris-HCl Gel, 5%, 10-well, 30 µl	1 gel	4°C	12 wk
161-1211	Ready Gel Tris-HCl Gel, 5%, 15-well, 15 µl	1 gel	4°C	12 wk
161-1212	Ready Gel Tris-HCl Gel, 5%, prep well, 450 µl	1 gel	4°C	12 wk
161-1213	Ready Gel Tris-HCl Gel, 5%, 10-well, 50 µl	1 gel	4°C	12 wk
161-1214	Ready Gel Tris-HCl Gel, 5%, 12-well, 20 µl	1 gel	4°C	12 wk
161-1215	Ready Gel Tris-HCl Gel, 5%, 9-well, 30 µl	1 gel	4°C	12 wk
161-1216	Ready Gel Tris-HCl Gel, 18%, 10-well, 30 µl	1 gel	4°C	12 wk
161-1217	Ready Gel Tris-HCl Gel, 18%, 15-well, 15 µl	1 gel	4°C	12 wk
161-1218	Ready Gel Tris-HCl Gel, 18%, prep well, 450 µl	1 gel	4°C	12 wk
161-1219	Ready Gel Tris-HCl Gel, 18%, 10-well, 50 µl	1 gel	4°C	12 wk
161-1220	Ready Gel Tris-HCl Gel, 18%, 12-well, 20 µl	1 gel	4°C	12 wk
161-1221	Ready Gel Tris-HCl Gel, 18%, 9-well, 30 µl	1 gel	4°C	12 wk
161-1222	Ready Gel Tris-HCl Gel, 8–16%, 10-well, 30 µl	1 gel	4°C	12 wk
161-1223	Ready Gel Tris-HCl Gel, 8–16%, 15-well, 15 µl	1 gel	4°C	12 wk
161-1224	Ready Gel Tris-HCl Gel, 8–16%, prep well, 450 µl	1 gel	4°C	12 wk
161-1225	Ready Gel Tris-HCl Gel, 8–16%, 10-well, 50 µl	1 gel	4°C	12 wk
161-1226	Ready Gel Tris-HCl Gel, 8–16%, 12-well, 20 µl	1 gel	4°C	12 wk
161-1227	Ready Gel Tris-HCl Gel, 8–16%, 9-well, 30 µl	1 gel	4°C	12 wk
161-1390	Ready Gel Tris-HCl Gel, 10%, 7 cm IPG well	1 gel	4°C	12 wk
161-1391	Ready Gel Tris-HCl Gel, 12%, 7 cm IPG well	1 gel	4°C	12 wk
161-1392	Ready Gel Tris-HCl Gel, 4–15%, 7 cm IPG well	1 gel	4°C	12 wk
161-1393	Ready Gel Tris-HCl Gel, 4–20%, 7 cm IPG well	1 gel	4°C	12 wk
161-1394	Ready Gel Tris-HCl Gel, 8–16%, 7 cm IPG well	1 gel	4°C	12 wk
161-1395	Ready Gel Tris-HCl Gel, 10–20%, 7 cm IPG well	1 gel	4°C	12 wk
Ready Gel Tris-Tricine Gels				
161-1107	Ready Gel Tris-Tricine Gel, 16.5%, 10-well, 30 µl	1 gel	4°C	12 wk
161-1108	Ready Gel Tris-Tricine Gel, 10–20%, 10-well, 30 µl	1 gel	4°C	12 wk
161-1125	Ready Gel Tris-Tricine Gel, 16.5%, 15-well, 15 µl	1 gel	4°C	12 wk
161-1126	Ready Gel Tris-Tricine Gel, 10–20%, 15-well, 15 µl	1 gel	4°C	12 wk
161-1143	Ready Gel Tris-Tricine Gel, 16.5%, prep well, 450 µl	1 gel	4°C	12 wk
161-1144	Ready Gel Tris-Tricine Gel, 10–20%, prep well, 450 µl	1 gel	4°C	12 wk
161-1161	Ready Gel Tris-Tricine Gel, 16.5%, 10-well, 50 µl	1 gel	4°C	12 wk
161-1162	Ready Gel Tris-Tricine Gel, 10–20%, 10-well, 50 µl	1 gel	4°C	12 wk
161-1179	Ready Gel Tris-Tricine Gel, 16.5%, 12-well, 20 µl	1 gel	4°C	12 wk
161-1180	Ready Gel Tris-Tricine Gel, 10–20%, 12-well, 20 µl	1 gel	4°C	12 wk
161-1197	Ready Gel Tris-Tricine Gel, 16.5%, 9-well, 30 µl	1 gel	4°C	12 wk
161-1198	Ready Gel Tris-Tricine Gel, 10–20%, 9-well, 30 µl	1 gel	4°C	12 wk
Ready Gel TBE Gels				
161-1109	Ready Gel TBE Gel, 5%, 10-well, 30 µl	1 gel	4°C	12 wk
161-1110	Ready Gel TBE Gel, 10%, 10-well, 30 µl	1 gel	4°C	12 wk
161-1127	Ready Gel TBE Gel, 5%, 15-well, 15 µl	1 gel	4°C	12 wk

(continues)

Ready Gel Precast Gels (cont.)

Catalog #	Description	Size/ Quantity	Storage Temp.	Shelf Life
Ready Gel TBE Gels (cont.)				
161-1128	Ready Gel TBE Gel, 10%, 15-well, 15 µl	1 gel	4°C	12 wk
161-1163	Ready Gel TBE Gel, 5%, 10-well, 50 µl	1 gel	4°C	12 wk
161-1164	Ready Gel TBE Gel, 10%, 10-well, 50 µl	1 gel	4°C	12 wk
161-1181	Ready Gel TBE Gel, 5%, 12-well, 20 µl	1 gel	4°C	12 wk
161-1182	Ready Gel TBE Gel, 10%, 12-well, 20 µl	1 gel	4°C	12 wk
161-1199	Ready Gel TBE Gel, 5%, 9-well, 30 µl	1 gel	4°C	12 wk
161-1200	Ready Gel TBE Gel, 10%, 9-well, 30 µl	1 gel	4°C	12 wk
161-1228	Ready Gel TBE Gel, 15%, 10-well, 30 µl	1 gel	4°C	12 wk
161-1229	Ready Gel TBE Gel, 15%, 15-well, 15 µl	1 gel	4°C	12 wk
161-1231	Ready Gel TBE Gel, 15%, 10-well, 50 µl	1 gel	4°C	12 wk
161-1232	Ready Gel TBE Gel, 15%, 12-well, 20 µl	1 gel	4°C	12 wk
161-1233	Ready Gel TBE Gel, 15%, 9-well, 30 µl	1 gel	4°C	12 wk
161-1234	Ready Gel TBE Gel, 4–20%, 10-well, 30 µl	1 gel	4°C	12 wk
161-1235	Ready Gel TBE Gel, 4–20%, 15-well, 15 µl	1 gel	4°C	12 wk
161-1237	Ready Gel TBE Gel, 4–20%, 10-well, 50 µl	1 gel	4°C	12 wk
161-1238	Ready Gel TBE Gel, 4–20%, 12-well, 20 µl	1 gel	4°C	12 wk
161-1239	Ready Gel TBE Gel, 4–20%, 9-well, 30 µl	1 gel	4°C	12 wk

Ready Gel IEF Gels

161-1111	Ready Gel IEF Gel, pH 3–10, 10-well, 30 µl	1 gel	4°C	24 wk
161-1112	Ready Gel IEF Gel, pH 5–8, 10-well, 30 µl	1 gel	4°C	24 wk
161-1129	Ready Gel IEF Gel, pH 3–10, 15-well, 15 µl	1 gel	4°C	24 wk
161-1130	Ready Gel IEF Gel, pH 5–8, 15-well, 15 µl	1 gel	4°C	24 wk
161-1147	Ready Gel IEF Gel, pH 3–10, prep well, 450 µl	1 gel	4°C	24 wk
161-1148	Ready Gel IEF Gel, pH 5–8, prep well, 450 µl	1 gel	4°C	24 wk
161-1165	Ready Gel IEF Gel, pH 3–10, 10-well, 50 µl	1 gel	4°C	24 wk
161-1166	Ready Gel IEF Gel, pH 5–8, 10-well, 50 µl	1 gel	4°C	24 wk
161-1183	Ready Gel IEF Gel, pH 3–10, 12-well, 20 µl	1 gel	4°C	24 wk
161-1184	Ready Gel IEF Gel, pH 5–8, 12-well, 20 µl	1 gel	4°C	24 wk
161-1201	Ready Gel IEF Gel, pH 3–10, 9-well, 30 µl	1 gel	4°C	24 wk
161-1202	Ready Gel IEF Gel, pH 5–8, 9-well, 30 µl	1 gel	4°C	24 wk

Ready Gel Zymogram Gels

161-1113	Ready Gel Zymogram Gel, 10%, gelatin, 10-well, 30 µl	1 gel	4°C	8 wk
161-1114	Ready Gel Zymogram Gel, 12%, casein, 10-well, 30 µl	1 gel	4°C	8 wk
161-1131	Ready Gel Zymogram Gel, 10%, gelatin, 15-well, 15 µl	1 gel	4°C	8 wk
161-1132	Ready Gel Zymogram Gel, 12%, casein, 15-well, 15 µl	1 gel	4°C	8 wk
161-1167	Ready Gel Zymogram Gel, 10%, gelatin, 10-well, 50 µl	1 gel	4°C	8 wk
161-1168	Ready Gel Zymogram Gel, 12%, casein, 10-well, 50 µl	1 gel	4°C	8 wk
161-1185	Ready Gel Zymogram Gel, 10%, gelatin, 12-well, 20 µl	1 gel	4°C	8 wk
161-1186	Ready Gel Zymogram Gel, 12%, casein, 12-well, 20 µl	1 gel	4°C	8 wk

Catalog #	Description	Size/ Quantity	Storage Temp.	Shelf Life
161-1203	Ready Gel Zymogram Gel, 10%, gelatin, 9-well, 30 µl	1 gel	4°C	8 wk
161-1204	Ready Gel Zymogram Gel, 12%, casein, 9-well, 30 µl	1 gel	4°C	8 wk
Ready Gel TBE-Urea Gels				
161-1115	Ready Gel TBE-Urea Gel, 5%, 10-well, 30 µl	1 gel	4°C	8 wk
161-1116	Ready Gel TBE-Urea Gel, 10%, 10-well, 30 µl	1 gel	4°C	8 wk
161-1117	Ready Gel TBE-Urea Gel, 15%, 10-well, 30 µl	1 gel	4°C	8 wk
161-1133	Ready Gel TBE-Urea Gel, 5%, 15-well, 15 µl	1 gel	4°C	8 wk
161-1134	Ready Gel TBE-Urea Gel, 10%, 15-well, 15 µl	1 gel	4°C	8 wk
161-1135	Ready Gel TBE-Urea Gel, 15%, 15-well, 15 µl	1 gel	4°C	8 wk
161-1151	Ready Gel TBE-Urea Gel, 5%, prep well, 450 µl	1 gel	4°C	8 wk
161-1152	Ready Gel TBE-Urea Gel, 10%, prep well, 450 µl	1 gel	4°C	8 wk
161-1153	Ready Gel TBE-Urea Gel, 15%, prep well, 450 µl	1 gel	4°C	8 wk
161-1187	Ready Gel TBE-Urea Gel, 5%, 12-well, 20 µl	1 gel	4°C	8 wk
161-1188	Ready Gel TBE-Urea Gel, 10%, 12-well, 20 µl	1 gel	4°C	8 wk
161-1189	Ready Gel TBE-Urea Gel, 15%, 12-well, 20 µl	1 gel	4°C	8 wk
161-1205	Ready Gel TBE-Urea Gel, 5%, 9-well, 30 µl	1 gel	4°C	8 wk
161-1206	Ready Gel TBE-Urea Gel, 10%, 9-well, 30 µl	1 gel	4°C	8 wk
161-1207	Ready Gel TBE-Urea Gel, 15%, 9-well, 30 µl	1 gel	4°C	8 wk

PROTEAN® Precast Gels

Catalog #	Description	Size/ Quantity	Storage Temp.	Shelf Life
PROTEAN II Ready Gel Precast Gels				
161-1450	PROTEAN II Ready Gel Tris-HCl Gel, 10%, 17 cm IPG well	1 gel	4°C	12 wk
161-1451	PROTEAN II Ready Gel Tris-HCl Gel, 12%, 17 cm IPG well	1 gel	4°C	12 wk
161-1452	PROTEAN II Ready Gel Tris-HCl Gel, 10–20%, 17 cm IPG well	1 gel	4°C	12 wk
161-1453	PROTEAN II Ready Gel Tris-HCl Gel, 8–16%, 17 cm IPG well	1 gel	4°C	12 wk
161-1454	PROTEAN II Ready Gel Tris-HCl Gel, 10%, 15-well	1 gel	4°C	12 wk
161-1455	PROTEAN II Ready Gel Tris-HCl Gel, 12%, 15-well	1 gel	4°C	12 wk
161-1456	PROTEAN II Ready Gel Tris-HCl Gel, 10–20%, 15-well	1 gel	4°C	12 wk
161-1457	PROTEAN II Ready Gel Tris-HCl Gel, 8–16%, 15-well	1 gel	4°C	12 wk
161-1458	PROTEAN II Ready Gel Tris-HCl Gel, 10%, 20-well	1 gel	4°C	12 wk
161-1459	PROTEAN II Ready Gel Tris-HCl Gel, 12%, 20-well	1 gel	4°C	12 wk
161-1460	PROTEAN II Ready Gel Tris-HCl Gel, 10–20%, 20-well	1 gel	4°C	12 wk
161-1461	PROTEAN II Ready Gel Tris-HCl Gel, 8–16%, 20-well	1 gel	4°C	12 wk
PROTEAN® Plus Gel				
161-1511	PROTEAN Plus Tris-HCl Gel, 12%, 24 cm IPG well	1 gel	4°C	12 wk


DNA Standards

Catalog #	Description	Size/ Quantity	Storage Temp.
20 bp Rulers			
170-8351	EZ Load™ 20 bp Molecular Ruler, 20–1,000 bp	500 µl	4°C
170-8201	20 bp Molecular Ruler, 20–1,000 bp	250 µl	4°C
100 bp Rulers			
170-8352	EZ Load 100 bp Molecular Ruler, 100–1,000 bp	500 µl	4°C
170-8202	100 bp Molecular Ruler, 100–1,000 bp	250 µl	4°C
170-8353	EZ Load 100 bp PCR Molecular Ruler, 100–3,000 bp	500 µl	4°C
170-8206	100 bp PCR Molecular Ruler, 100–3,000 bp	200 µl	4°C
500 bp Rulers			
170-8354	EZ Load 500 bp Molecular Ruler, 500–8,000 bp	500 µl	4°C
170-8203	500 bp Molecular Ruler, 500–8,000 bp	200 µl	4°C
1 kb Rulers			
170-8355	EZ Load 1 kb Molecular Ruler, 1–15 kb	500 µl	4°C
170-8204	1 kb Molecular Ruler, 1–15 kb	200 µl	4°C
Mass Rulers			
170-8356	EZ Load Precision Molecular Mass Ruler, 10–100 ng	500 µl	4°C
170-8207	Precision Molecular Mass Ruler, 10–100 ng	250 µl	4°C
Other Molecular Rulers			
170-8205	2.5 kb Molecular Ruler	400 µl	4°C
170-8200	AmpliSize® Molecular Ruler, 50–2,000 bp	250 µl	4°C
Pulsed Field Standards			
170-3624	CHEF DNA Size Standard, 5 kb ladder	20 µg	4°C
170-3707	CHEF DNA Size Standard, 8–48 kb ladder	20 µg	4°C
170-3635	CHEF DNA Size Standard, lambda ladder	5 plugs	4°C
170-3605	CHEF DNA Size Marker, <i>S. cerevisiae</i> chromosomes	5 plugs	4°C
170-3667	CHEF DNA Size Marker, <i>H. wingei</i> chromosomes	5 plugs	4°C
170-3633	CHEF DNA Size Marker, <i>S. pombe</i> chromosomes	5 plugs	4°C
Fluorescent Rulers			
170-8216	100 bp Fluorescein-Labeled Ruler, 100–1,000 bp	100 µl	4°C
170-8217	100 bp Texas Red-Labeled Ruler, 100–1,000 bp	100 µl	4°C
Molecular Markers			
170-8361	EZ Load™ HT Molecular Marker, 100 bp–2 kb	1.6 ml	4°C
170-8362	EZ Load HT Molecular Marker, 500 bp–10 kb	1.6 ml	4°C

Protein Standards

Catalog #	Description	Size/ Quantity	Storage Temp.
Unstained Standards			
161-0326	Polypeptide SDS-PAGE Standards, 3.4–38.6 kD	200 µl	–20°C
161-0304	SDS-PAGE Standards, low range, 14.4–97.4 kD	200 µl	–20°C
161-0303	SDS-PAGE Standards, high range, 45–200 kD	200 µl	–20°C
161-0317	SDS-PAGE Standards, broad range, 6.5–200 kD	200 µl	–20°C
161-0314	Silver Stain SDS-PAGE Standards, low range, 14.4–97.4 kD	200 µl	–20°C
161-0315	Silver Stain SDS-PAGE Standards, high range, 45–200 kD	200 µl	–20°C
161-0332	SDS-PAGE Standards, SYPRO Orange, broad range	200 µl	–20°C
161-0320	2-D SDS-PAGE Standards	500 µl	–20°C
161-0310	IEF Standards, pI range 4.45–9.6	250 µl	–20°C
161-0306	Biotinylated SDS-PAGE Standards, low range	250 µl	–20°C
161-0308	Biotinylated SDS-PAGE Standards Kit, low range, avidin-AP	250 µl	–20°C
161-0307	Biotinylated SDS-PAGE Standards Kit, low range, avidin-HRP	250 µl	–20°C
161-0311	Biotinylated SDS-PAGE Standards, high range	250 µl	–20°C
161-0312	Biotinylated SDS-PAGE Standards Kit, high range, avidin-HRP	250 µl	–20°C
161-0313	Biotinylated SDS-PAGE Standards Kit, high range, avidin-AP	250 µl	–20°C
161-0319	Biotinylated SDS-PAGE Standards, broad range	250 µl	–20°C
161-0321	Biotinylated SDS-PAGE Standards Kit, broad range, avidin-HRP	250 µl	–20°C
161-0322	Biotinylated SDS-PAGE Standards Kit, broad range, avidin-AP	250 µl	–20°C
161-0363	Precision Plus Protein™ Unstained Standards	1 ml	–20°C
161-0378	Precision Plus Protein Standard Plugs, unstained	24 plugs	4°C
Prestained Standards			
161-0305	Prestained SDS-PAGE Standards, low range, 19–107 kD	500 µl	–20°C
161-0309	Prestained SDS-PAGE Standards, high range, 47–205 kD	500 µl	–20°C
161-0318	Prestained SDS-PAGE Standards, broad range, 7–205 kD	500 µl	–20°C
161-0325	Kaleidoscope™ Polypeptide Standards, 3.4–38.6 kD	500 µl	–20°C
161-0324	Kaleidoscope Prestained Standards, 6.9–205 kD	500 µl	–20°C
161-0373	Precision Plus Protein All Blue Standards, 10–250 kD	500 µl	–20°C
161-0374	Precision Plus Protein Dual Color Standards, 10–250 kD	500 µl	–20°C
161-0375	Precision Plus Protein Kaleidoscope Standards, 10–250 kD	500 µl	–20°C

Guide to Blot Detection Reagents and Sensitivities



Blotting Detection and Reagents

Catalog #	Description	Size/ Quantity	Storage Temp.
Colorimetric Detection (Horseradish Peroxidase)			
170-8235	Opti-4CN Substrate Kit	1 kit	4°C
170-8236	Opti-4CN Goat Anti-Rabbit Detection Kit	1 kit	4°C
170-8237	Opti-4CN Goat Anti-Mouse Detection Kit	1 kit	4°C
170-8238	Amplified Opti-4CN Substrate Kit	1 kit	4°C
170-8239	Amplified Opti-4CN Goat Anti-Rabbit Detection Kit	1 kit	4°C
170-8240	Amplified Opti-4CN Goat Anti-Mouse Detection Kit	1 kit	4°C
170-6463	Immun-Blot Assay Kit, Goat Anti-Rabbit IgG (H + L) HRP Conjugate	1 kit	-20°C
170-6464	Immun-Blot Assay Kit, Goat Anti-Mouse IgG (H + L) HRP Conjugate	1 kit	-20°C
170-6465	Immun-Blot Assay Kit, Goat Anti-Human IgG (H + L) HRP Conjugate	1 kit	-20°C
170-6466	Immun-Blot Assay Kit, Protein A-HRP Conjugate and 4CN	1 kit	-20°C
170-6467	Immun-Blot Assay Kit, Protein G-HRP Conjugate and 4CN	1 kit	-20°C
170-6534	HRP Color Development Reagent, 4CN	5 g	-20°C
170-6535	HRP Color Development Reagent, DAB	5 g	-20°C
Colorimetric Detection (Alkaline Phosphatase)			
170-6460	Immun-Blot Assay Kit, Goat Anti-Rabbit IgG (H + L) AP Conjugate	1 kit	-20°C
170-6461	Immun-Blot Assay Kit, Goat Anti-Mouse IgG (H + L) AP Conjugate	1 kit	-20°C
170-6462	Immun-Blot Assay Kit, Goat Anti-Human IgG (H + L) AP Conjugate	1 kit	-20°C
170-6412	Immun-Blot Assay Kit, Goat-Anti Rabbit IgG (H + L)-Biotin	1 kit	-20°C
170-6403	Blotting Grade Biotinylated Alkaline Phosphatase	1 ml	-20°C
170-6408	Blotting Grade Streptavidin	1.0 mg	-20°C
170-6539	AP Color Development Reagent, BCIP	300 mg	-20°C
170-6532	AP Color Development Reagent, NBT	600 mg	-20°C
Chemiluminescent Detection (Horseradish Peroxidase)			
170-5040	Immun-Star HRP Substrate	500 ml	4°C
170-5041	Immun-Star HRP Substrate	100 ml	4°C
170-5042	Immun-Star Goat-Anti Rabbit (GAR)-HRP Detection Reagents	500 ml	4°C
170-5043	Immun-Star Goat-Anti Mouse (GAM)-HRP Detection Reagents	500 ml	4°C
170-5044	Immun-Star Goat-Anti Mouse (GAM)-HRP Detection Kit	1 kit	4°C
170-5045	Immun-Star Goat-Anti Rabbit (GAR)-HRP Detection Kit	1 kit	4°C
170-5046	Immun-Star Goat-Anti Rabbit (GAR)-HRP Conjugate	2 ml	4°C
170-5047	Immun-Star Goat-Anti Mouse (GAM)-HRP Conjugate	2 ml	4°C
Chemiluminescent Detection (Alkaline Phosphatase)			
170-5010	Immun-Star Goat-Anti Mouse (GAM)-AP Detection Kit	1 kit	4°C
170-5011	Immun-Star Goat-Anti Rabbit (GAR)-AP Detection Kit	1 kit	4°C
170-5012	Immun-Star AP Substrate Pack	1 pack	4°C
170-5013	Immun-Star Goat Anti-Mouse (GAM)-AP Intro Kit	1 kit	4°C
170-5014	Immun-Star Goat Anti-Rabbit (GAR)-AP Intro Kit	1 kit	4°C
170-5018	Immun-Star AP Substrate	125 ml	4°C
170-5015	Blotting Reagents Pack	1 pack	4°C

Catalog #	Description	Size/ Quantity	Storage Temp.
Antibody and Substrate Conjugates			
170-6431	Horseradish Peroxidase Conjugate Substrate Kit	1 kit	-20°C
170-6515	Goat Anti-Rabbit IgG-HRP (H + L)	2 ml	-20°C
170-6516	Goat Anti-Mouse IgG-HRP (H + L)	2 ml	-20°C
172-1050	Goat Anti-Human IgG-HRP (H + L)	2 ml	-20°C
170-6522	Protein A-HRP Conjugate	1 ml	-20°C
170-6425	Protein G-HRP Conjugate	1 ml	-20°C
170-6528	Avidin-HRP Conjugate	2 ml	-20°C
170-6432	Alkaline Phosphatase Conjugate Substrate Kit	1 kit	-20°C
170-6518	Goat Anti-Rabbit IgG-AP (H + L)	1 ml	-20°C
170-6520	Goat Anti-Mouse IgG-AP (H + L)	1 ml	-20°C
170-6521	Goat Anti-Human IgG-AP (H + L)	1 ml	-20°C
170-6401	Biotinylated Goat Anti-Rabbit IgG (H + L)	1 ml	-20°C
170-6533	Avidin-AP	1 ml	-20°C
170-3554	Streptavidin-AP	0.5 ml	-20°C
170-6529	NHS-Biotin	4 ml	-20°C
170-6403	Blotting Grade Biotinylated AP	1 ml	-20°C
161-0380	Precision Protein StrepTactin-HRP Conjugate	0.3 ml	4°C
161-0382	Precision Protein StrepTactin-AP Conjugate	0.3 ml	4°C
Blotting Reagents and Buffers			
170-6435	Blot Processing Buffer, 10x TBS	1 L	Ambient
161-0780	10x PBS	1 L	Ambient
161-0781	10% Tween 20	1 L	Ambient
170-6531	Tween 20, EIA grade	100 ml	Ambient
161-0783	1x PBS/1% Casein	1 L	4°C
161-0782	1x TBS/1% Casein	1 L	4°C
170-6404	Blotting Grade Blocker, nonfat dry milk	300 g	Ambient
170-6537	Gelatin, EIA grade	200 g	Ambient
170-3932	Mini Trans-Blot® Filter Paper, thick, 7.5 x 10 cm	50 sheets	Ambient
170-3965	Extra-Thick Blot Paper, 7.5 x 10 cm	60 sheets	Ambient
Blot Stains			
170-6517	Enhanced Colloidal Gold Total Protein Detection Kit	1 kit	4°C
170-6527	Colloidal Gold Total Protein Stain	500 ml	4°C
170-3127	SYPRO Ruby Protein Blot Stain	200 ml	Ambient
Background Remover for X-Ray Film			
170-5020	De-Expose™ Background Remover Solution, 10x	250 ml	4°C

Blotting Membranes

Catalog #	Description	Size/ Quantity	Storage Temp.
Nitrocellulose			
162-0146	Nitrocellulose Membrane, 0.2 µm, 7 x 8.4 cm	10 sheets	Ambient
162-0147	Nitrocellulose Membrane, 0.2 µm, 13.5 x 16.5 cm	10 sheets	Ambient
162-0150	Nitrocellulose Membrane, 0.2 µm, 20 x 20 cm	5 sheets	Ambient
162-0112	Nitrocellulose Membrane, 0.2 µm, 30 cm x 3.5 m	1 roll	Ambient
162-0168	Nitrocellulose Membrane, 0.2 µm, 8.5 x 13.5 cm	10 sheets	Ambient
162-0252	Nitrocellulose Membrane, 0.2 µm, 26.5 x 28 cm	10 sheets	Ambient
162-0145	Nitrocellulose Membrane, 0.45 µm, 7 x 8.4 cm	10 sheets	Ambient
162-0117	Nitrocellulose Membrane, 0.45 µm, 9 x 12 cm	10 sheets	Ambient
162-0114	Nitrocellulose Membrane, 0.45 µm, 15 x 9.2 cm	10 sheets	Ambient
162-0116	Nitrocellulose Membrane, 0.45 µm, 15 x 15 cm	10 sheets	Ambient
162-0148	Nitrocellulose Membrane, 0.45 µm, 11.5 x 16 cm	10 sheets	Ambient
162-0113	Nitrocellulose Membrane, 0.45 µm, 20 x 20 cm	5 sheets	Ambient
162-0115	Nitrocellulose Membrane, 0.45 µm, 30 cm x 3.5 m	1 roll	Ambient
162-0167	Nitrocellulose Membrane, 0.45 µm, 8.5 x 13.5 cm	10 sheets	Ambient
162-0251	Nitrocellulose Membrane, 0.45 µm, 26.5 x 28 cm	10 sheets	Ambient
Nitrocellulose Membrane Sandwiches			
162-0232	Nitrocellulose/Filter Paper Sandwich, 0.2 µm, 8.5 x 13.5 cm	20 pack	Ambient
162-0233	Nitrocellulose/Filter Paper Sandwich, 0.2 µm, 8.5 x 13.5 cm	50 pack	Ambient
162-0212	Nitrocellulose/Filter Paper Sandwich, 0.2 µm, 7 x 8.5 cm	20 pack	Ambient
162-0213	Nitrocellulose/Filter Paper Sandwich, 0.2 µm, 7 x 8.5 cm	50 pack	Ambient
162-0234	Nitrocellulose/Filter Paper Sandwich, 0.45 µm, 8.5 x 13.5 cm	20 pack	Ambient
162-0235	Nitrocellulose/Filter Paper Sandwich, 0.45 µm, 8.5 x 13.5 cm	50 pack	Ambient
162-0214	Nitrocellulose/Filter Paper Sandwich, 0.45 µm, 7 x 8.5 cm	20 pack	Ambient
162-0215	Nitrocellulose/Filter Paper Sandwich, 0.45 µm, 7 x 8.5 cm	50 pack	Ambient
Supported Nitrocellulose			
162-0095	Supported Nitrocellulose Membrane, 0.2 µm, 7 x 8.4 cm	10 sheets	Ambient
162-0097	Supported Nitrocellulose Membrane, 0.2 µm, 30 cm x 3 m	1 roll	Ambient
162-0071	Supported Nitrocellulose Membrane, 0.2 µm, 8.5 x 13.5 cm	10 sheets	Ambient
162-0253	Supported Nitrocellulose Membrane, 0.2 µm, 26.5 x 28 cm	10 sheets	Ambient
162-0090	Supported Nitrocellulose Membrane, 0.45 µm, 7 x 8.4 cm	10 sheets	Ambient
162-0094	Supported Nitrocellulose Membrane, 0.45 µm, 30 cm x 3 m	1 roll	Ambient
162-0070	Supported Nitrocellulose Membrane, 0.45 µm, 8.5 x 13.5 cm	10 sheets	Ambient
162-0254	Supported Nitrocellulose Membrane, 0.45 µm, 26.5 x 28 cm	10 sheets	Ambient
PVDF Membranes			
162-0174	Immun-Blot PVDF Membrane, 7 x 8.4 cm	10 sheets	Ambient
162-0175	Immun-Blot PVDF Membrane, 10 x 15 cm	10 sheets	Ambient
162-0176	Immun-Blot PVDF Membrane, 20 x 20 cm	10 sheets	Ambient
162-0177	Immun-Blot PVDF Membrane, 26 cm x 3.3 m	1 roll	Ambient


Catalog #	Description	Size/ Quantity	Storage Temp.
162-0186	Sequi-Blot™ PVDF Membrane, 7 x 8.4 cm	10 sheets	Ambient
162-0180	Sequi-Blot PVDF Membrane, 10 x 15 cm	10 sheets	Ambient
162-0181	Sequi-Blot PVDF Membrane, 15 x 15 cm	10 sheets	Ambient
162-0182	Sequi-Blot PVDF Membrane, 20 x 20 cm	10 sheets	Ambient
162-0185	Sequi-Blot PVDF Membrane, 20 x 20 cm	3 sheets	Ambient
162-0184	Sequi-Blot PVDF Membrane, 26 cm x 3.3 m	1 roll	Ambient
162-0255	Immun-Blot PVDF Membrane, 25 x 28 cm	10 sheets	Ambient
162-0256	Sequi-Blot PVDF Membrane, 25 x 28 cm	10 sheets	Ambient
PVDF Membrane Sandwiches			
162-0238	Immun-Blot PVDF/Filter Paper Sandwich, 8.5 x 13.5 cm	20 pack	Ambient
162-0239	Immun-Blot PVDF/Filter Paper Sandwich, 8.5 x 13.5 cm	50 pack	Ambient
162-0236	Sequi-Blot PVDF/Filter Paper Sandwich, 8.5 x 13.5 cm	20 pack	Ambient
162-0237	Sequi-Blot PVDF/Filter Paper Sandwich, 8.5 x 13.5 cm	50 pack	Ambient
162-0218	Immun-Blot PVDF/Filter Paper Sandwich, 7 x 8.5 cm	20 pack	Ambient
162-0219	Immun-Blot PVDF/Filter Paper Sandwich, 7 x 8.5 cm	50 pack	Ambient
162-0216	Sequi-Blot PVDF/Filter Paper Sandwich, 7 x 8.5 cm	20 pack	Ambient
162-0217	Sequi-Blot PVDF/Filter Paper Sandwich, 7 x 8.5 cm	50 pack	Ambient
Zeta-Probe® Membranes			
162-0153	Zeta-Probe Membrane, 9 x 12 cm	15 sheets	Ambient
162-0154	Zeta-Probe Membrane, 10 x 15 cm	15 sheets	Ambient
162-0155	Zeta-Probe Membrane, 15 x 15 cm	15 sheets	Ambient
162-0156	Zeta-Probe Membrane, 15 x 20 cm	15 sheets	Ambient
162-0157	Zeta-Probe Membrane, 20 x 20 cm	15 sheets	Ambient
162-0158	Zeta-Probe Membrane, 20 x 25 cm	3 sheets	Ambient
162-0165	Zeta-Probe Membrane, 20 cm x 3.3 m	1 roll	Ambient
162-0159	Zeta-Probe Membrane, 30 cm x 3.3 m	1 roll	Ambient
162-0206	Zeta-Probe Membrane, 7 x 10 cm	15 sheets	Ambient
Zeta-Probe® GT Membranes			
162-0190	Zeta-Probe GT Membrane, 9 x 12 cm	15 sheets	Ambient
162-0191	Zeta-Probe GT Membrane, 10 x 15 cm	15 sheets	Ambient
162-0192	Zeta-Probe GT Membrane, 15 x 15 cm	15 sheets	Ambient
162-0193	Zeta-Probe GT Membrane, 15 x 20 cm	15 sheets	Ambient
162-0194	Zeta-Probe GT Membrane, 20 x 20 cm	15 sheets	Ambient
162-0195	Zeta-Probe GT Membrane, 20 x 25 cm	3 sheets	Ambient
162-0197	Zeta-Probe GT Membrane, 20 cm x 3.3 m	1 roll	Ambient
162-0196	Zeta-Probe GT Membrane, 30 cm x 3.3 m	1 roll	Ambient
162-0208	Zeta-Probe GT Membrane, 7 x 10 cm	15 sheets	Ambient

Gel Reagents, Stains, and Buffers

Catalog #	Description	Size/ Quantity	Storage Temp.
Gel-Casting Reagents			
161-0101	Acrylamide, 99.9%	500 g	Ambient
161-0123	Acrylamide/Bis, 19:1	150 g	4°C
161-0124	Acrylamide/Bis, 29:1	150 g	Ambient
161-0125	Acrylamide/Bis, 37.5:1	150 g	Ambient
161-0201	Bis Crosslinker	50 g	Ambient
161-0144	40% Acrylamide/Bis Solution, 19:1	500 ml	4°C
161-0146	40% Acrylamide/Bis Solution, 29:1	500 ml	4°C
161-0148	40% Acrylamide/Bis Solution, 37.5:1	500 ml	4°C
161-0154	30% Acrylamide/Bis Solution, 19:1	500 ml	4°C
161-0156	30% Acrylamide/Bis Solution, 29:1	500 ml	4°C
161-0158	30% Acrylamide/Bis Solution, 37.5:1	500 ml	4°C
161-0140	40% Acrylamide Solution	500 ml	4°C
161-0142	2% Bis Solution	500 ml	4°C
161-0302	SDS	1 kg	Ambient
161-0801	TEMED	50 ml	Ambient
161-0611	Dithiothreitol	5 g	Ambient
161-0700	Ammonium Persulfate	10 g	Ambient
161-0710	2-Mercaptoethanol	25 ml	Ambient
161-0712	Tricine	100 g	Ambient
161-0718	Glycine	1 kg	Ambient
161-0719	Tris	1 kg	Ambient
161-0720	Sucrose	1 kg	Ambient
161-0729	EDTA	500 g	Ambient
161-0731	Urea	1 kg	Ambient
161-0751	Boric Acid	1 kg	Ambient
Gel Stains			
161-0786	Bio-Safe™ Coomassie Stain	1 L	Ambient
161-0435	Coomassie Brilliant Blue R-250 Staining Solutions Kit	1 kit	Ambient
161-0440	Zinc Stain and Destain Kit	1 kit	Ambient
161-0449	Silver Stain Plus™ Kit	1 kit	Ambient
161-0443	Silver Stain Kit	1 kit	4°C
170-3120	SYPRO Orange Fluorescent Protein Stain	500 µl	4°C
170-3125	SYPRO Ruby Protein Gel Stain	1 L	Ambient
170-3122	Radiant® Red Fluorescent RNA Stain	10 ml	4°C
161-0470	Copper Stain and Destain Kit	1 kit	Ambient
161-0433	Ethidium Bromide Solution	10 ml	Ambient
161-0430	Ethidium Bromide Tablets	10 tablets	Ambient
161-0480	Dodeca™ Silver Stain Kit, large	10 L	Ambient
161-0481	Dodeca Silver Stain Kit, small	7 L	Ambient

Catalog #	Description	Size/ Quantity	Storage Temp.
Premixed Buffers			
161-0798	Resolving Gel Buffer, 1.5 M Tris-HCl, pH 8.8	1 L	Ambient
161-0799	Stacking Gel Buffer, 0.5 M Tris-HCl, pH 6.8	1 L	Ambient
161-0732	10x Tris/Glycine/SDS	1 L	Ambient
161-0772	10x Tris/Glycine/SDS	5 L cube	Ambient
161-0734	10x Tris/Glycine	1 L	Ambient
161-0771	10x Tris/Glycine	5 L cube	Ambient
161-0733	10x Tris/Boric Acid/EDTA	1 L	Ambient
161-0770	10x Tris/Boric Acid/EDTA	5 L cube	Ambient
161-0741	10x Tris/Boric Acid/EDTA, extended range	1 L	Ambient
161-0744	10x Tris/Tricine/SDS	1 L	Ambient
161-0743	50x Tris/Acetic Acid/EDTA	1 L	Ambient
161-0773	50x Tris/Acetic Acid/EDTA	5 L cube	Ambient
161-0778	10x Tris/CAPS	1 L	Ambient
161-0775	20x SSC	5 L cube	Ambient
161-0774	20x SSC	1 L	Ambient
161-0737	Laemmli Sample Buffer	30 ml	Ambient
161-0738	Native Sample Buffer	30 ml	4°C
161-0739	Tricine Sample Buffer	30 ml	Ambient
161-0768	TBE/Urea Sample Buffer	30 ml	4°C
161-0767	Nucleic Acid Sample Buffer, 5x	10 ml	Ambient
161-0764	Zymogram Sample Buffer	30 ml	Ambient
161-0765	Zymogram Renaturation Buffer	125 ml	Ambient
161-0766	Zymogram Development Buffer	125 ml	Ambient
161-0416	SDS Solution, 10% (w/v)	250 ml	Ambient
161-0418	SDS Solution, 20% (w/v)	1 L	Ambient
161-0763	IEF Sample Buffer	30 ml	Ambient
161-0761	10x IEF Anode Buffer	250 ml	Ambient
161-0762	10x IEF Cathode Buffer	250 ml	Ambient
Criterion XT Buffers and Reagents			
161-0788	XT MOPS Running Buffer, 20x	500 ml	Ambient
161-0789	XT MES Running Buffer, 20x	500 ml	Ambient
161-0790	XT Tricine Running Buffer, 20x	500 ml	Ambient
161-0791	XT Sample Buffer, 4x	10 ml	Ambient
161-0792	XT Reducing Agent, 20x	1 ml	Ambient
161-0793	XT MOPS Buffer Kit	1 kit	Ambient
161-0796	XT MES Buffer Kit	1 kit	Ambient
161-0797	XT Tricine Buffer Kit	1 kit	Ambient

Sample Preparation Chart



Protein Sample Preparation

Catalog #	Description	Size/ Quantity	Storage Temp.
ReadyPrep Kits			
163-2130	ReadyPrep 2-D Cleanup Kit	1 kit	Ambient
163-2100	ReadyPrep Sequential Extraction Kit	1 kit	Ambient
163-2090	ReadyPrep Reduction-Alkylation Kit	1 kit	Ambient
163-2089	ReadyPrep Protein Extraction Kit (Cytoplasmic/Nuclear)	1 kit	Ambient
163-2088	ReadyPrep Protein Extraction Kit (Membrane I)	1 kit	Ambient
163-2087	ReadyPrep Protein Extraction Kit (Signal)	1 kit	Ambient
163-2086	ReadyPrep Protein Extraction Kit (Total Protein)	1 kit	Ambient
163-2085	ReadyPrep Protein Extraction Kit (Soluble/Insoluble)	1 kit	Ambient
163-2084	ReadyPrep Protein Extraction Kit (Membrane II)	1 kit	Ambient
ReadyPrep Components and Related Products			
163-2102	ReadyPrep Sequential Extraction Kit Reagent 1	1 vial	Ambient
163-2103	ReadyPrep Sequential Extraction Kit Reagent 2	1 vial	Ambient
163-2104	ReadyPrep Sequential Extraction Kit Reagent 3	1 vial	Ambient
163-2105	ReadyPrep 2-D Starter Kit	1 kit	4°C
163-2106	ReadyPrep 2-D Starter Kit Rehydration/Sample Buffer	10 ml	Ambient
163-2107	ReadyPrep 2-D Starter Kit Equilibration Buffer I, with DTT	20 ml	Ambient
163-2108	ReadyPrep 2-D Starter Kit Equilibration Buffer II, without DTT	20 ml	Ambient
163-2110	<i>E. coli</i> Protein Sample, lyophilized	2.7 mg	4°C
163-2111	ReadyPrep Overlay Agarose	50 ml	4°C
163-2092	PROTEAN® Plus Overlay Agarose	125 ml	4°C
163-2091	ReadyPrep Proteomic Grade Water	500 ml	Ambient
163-2083	ReadyPrep 2-D Rehydration/Sample Buffer 1	10 ml	Ambient
Column-Based Purification			
732-6713	Aurum Serum Protein Mini Kit	2 preps	4°C
732-6701	Aurum Serum Protein Mini Kit	10 preps	4°C
732-6712	Aurum Affi-Gel Blue, Mini Kit	2 pack	4°C
732-6708	Aurum Affi-Gel Blue, Mini Columns	25 pack	4°C
732-6709	Aurum Affi-Gel Blue, Mini Columns	100 pack	4°C
732-2022	Econo-Pac Protein A Columns	5 columns	4°C
732-2020	Econo-Pac Protein A Kit	1 kit	4°C
732-2027	Econo-Pac Serum IgG Purification Kit	1 kit	4°C
732-6711	Aurum CEX Mini Kit	2 pack	4°C
732-6702	Aurum CEX Mini Kit	10 pack	4°C
732-6703	Aurum CEX Mini Columns	25 pack	4°C
732-6704	Aurum CEX Mini Columns	100 pack	4°C
732-6710	Aurum AEX Mini Kit	2 pack	4°C
732-6705	Aurum AEX Mini Kit	10 pack	4°C
732-6706	Aurum AEX Mini Columns	25 pack	4°C
732-6707	Aurum AEX Mini Columns	100 pack	4°C
732-0096	Protein Purification Sampler Pack V	1 pack	4°C
732-0094	Protein Purification Sampler Pack I	1 pack	4°C

Micro Bio-Spin and Bio-Spin Columns

Refer to Prepacked Spin Columns for Sample Cleanup in the Nucleic Acid Sample Preparation section.

Protein Assays

Catalog #	Description	Size/ Quantity	Storage Temp.
Protein Assays and Cuvettes			
500-0201	Quick Start™ Bradford Protein Assay Kit 1, bovine serum albumin standard	1 kit	4°C
500-0202	Quick Start Bradford Protein Assay Kit 2, bovine serum albumin standard set	1 kit	4°C
500-0203	Quick Start Bradford Protein Assay Kit 3, bovine γ -globulin standard	1 kit	4°C
500-0204	Quick Start Bradford Protein Assay Kit 4, bovine γ -globulin standard set	1 kit	4°C
500-0205	Quick Start Bradford 1x Dye Reagent	1 L	4°C
500-0206	Quick Start Bovine Serum Albumin Standard	2 mg/ml	4°C
500-0207	Quick Start Bovine Serum Albumin Standard Set	7 standards	4°C
500-0208	Quick Start Bovine γ -Globulin Standard	2 mg/ml	4°C
500-0209	Quick Start Bovine γ -Globulin Standard Set	7 standards	4°C
500-0001	Bio-Rad Protein Assay Kit I, bovine γ -globulin standard	1 kit	4°C
500-0002	Bio-Rad Protein Assay Kit II, bovine serum albumin standard	1 kit	4°C
500-0006	Bio-Rad Protein Assay Dye Reagent Concentrate	450 ml	4°C
500-0111	DC™ Protein Assay Kit I, bovine γ -globulin standard	1 kit	4°C
500-0112	DC Protein Assay Kit II, bovine serum albumin standard	1 kit	4°C
500-0116	DC Protein Assay Reagents Package	1 package	4°C
500-0121	RC DC™ Protein Assay Kit I, bovine γ -globulin standard	1 kit	4°C
500-0122	RC DC Protein Assay Kit II, bovine serum albumin standard	1 kit	4°C
500-0005	Protein Standard I, bovine γ -globulin	1 bottle	4°C
500-0007	Protein Standard II, bovine serum albumin	1 bottle	4°C
223-9950	Standard Disposable Polystyrene Cuvettes, 3.5 ml	100 cuvettes	Ambient
223-9955	Semimicrovolume Disposable Polystyrene Cuvettes, 1.5 ml	100 cuvettes	Ambient

Nucleic Acid Sample Preparation

Catalog #	Description	Size/ Quantity	Storage Temp.
Products for Nucleic Acid Isolation			
732-6100	Quantum Prep Plasmid Miniprep Kit	1 kit	Ambient
732-6120	Quantum Prep Midiprep Kit	1 kit	Ambient
732-6130	Quantum Prep Maxiprep Kit	1 kit	Ambient
732-6165	Freeze 'N Squeeze Spin Columns	25 columns	Ambient
732-6166	Freeze 'N Squeeze Spin Columns	100 columns	Ambient
732-6030	InstaGene™ Matrix	20 ml	4°C
732-6820	Aurum Total RNA Mini Kit	50 preps	Ambient
732-6800	Aurum Total RNA 96 Kit	2 x 96-well preps	Ambient
732-6370	AquaPure RNA Isolation Kit	100 preps	Ambient
732-6371	AquaPure RNA Blood Kit	100 preps	Ambient
732-6390	Master Blaster™ RNA Extraction Reagent	100 preps	Ambient
732-6340	AquaPure Genomic DNA Isolation Kit	100 preps	Ambient
732-6343	AquaPure Genomic DNA Tissue Kit	100 preps	Ambient
732-6345	AquaPure Genomic DNA Blood Kit	100 preps	Ambient
732-6470	Aurum Vacuum Manifold	1 manifold	Ambient
732-6400	Aurum Plasmid Mini Kit	100 preps	Ambient
732-6440	Aurum Plasmid Midi Kit	20 preps	Ambient
732-6450	Aurum Plasmid Maxi Kit	10 preps	Ambient
732-6460	Aurum Plasmid 96 Kit	2 x 96-well preps	Ambient
Products for Sequencing Cleanup			
732-6500	SEQueaky Kleen H ₂ O Dye Terminator Removal Kit	2 x 96-well, 4 collection plates	4°C
732-6505	SEQueaky Kleen H ₂ O Bulk Kit	50 x 96-well, 100 collection plates	4°C
732-6510	SEQueaky Kleen H ₂ O Adaptors	Set of 2	Ambient
732-6530	SEQueaky Kleen H ₂ O Starter Kit, for half-skirt plates	1 kit	Ambient
732-6531	SEQueaky Kleen H ₂ O Starter Kit, for unskirted plates	1 kit	Ambient
732-6532	SEQueaky Kleen H ₂ O Starter Kit, for full-skirt plates	1 kit	Ambient
732-6533	SEQueaky Kleen H ₂ O Starter Kit, for segmented plates	1 kit	Ambient
732-6540	SEQueaky Kleen H ₂ O Refill Kit, for half-skirt plates	1 kit	Ambient
732-6541	SEQueaky Kleen H ₂ O Refill Kit, for unskirted plates	1 kit	Ambient
732-6542	SEQueaky Kleen H ₂ O Refill Kit, for full-skirt plates	1 kit	Ambient
732-6543	SEQueaky Kleen H ₂ O Refill Kit, for segmented plates	1 kit	Ambient
732-6520	SEQueaky Kleen H ₂ O Half-Skirt Plates	50 pack	Ambient
732-6521	SEQueaky Kleen H ₂ O Unskirted Plates	50 pack	Ambient
732-6522	SEQueaky Kleen H ₂ O Full-Skirt Plates	50 pack	Ambient
732-6523	SEQueaky Kleen H ₂ O Segmented Plates	50 pack	Ambient

Catalog #	Description	Size/ Quantity	Storage Temp.
PCR Kleen Spin Columns			
732-6300	PCR Kleen Spin Columns	25 columns	Ambient
732-6301	PCR Kleen Spin Columns	100 columns	Ambient
Prepacked Spin Columns for Sample Cleanup*			
732-6221	Micro Bio-Spin 6 Columns, in Tris buffer	25 columns	4°C
732-6222	Micro Bio-Spin 6 Columns, in Tris buffer	100 columns	4°C
732-6225	Micro Bio-Spin 6 Columns, in Tris buffer	1,000 columns	4°C
732-6200	Micro Bio-Spin 6 Columns, in SSC buffer	25 columns	4°C
732-6201	Micro Bio-Spin 6 Columns, in SSC buffer	100 columns	4°C
732-6205	Micro Bio-Spin 6 Columns, in SSC buffer	1,000 columns	4°C
732-6223	Micro Bio-Spin 30 Columns, in Tris buffer	25 columns	4°C
732-6224	Micro Bio-Spin 30 Columns, in Tris buffer	100 columns	4°C
732-6226	Micro Bio-Spin 30 Columns, in Tris buffer	1,000 columns	4°C
732-6250	Micro Bio-Spin 30 Columns, in Tris buffer, RNase-free	25 columns	4°C
732-6251	Micro Bio-Spin 30 Columns, in Tris buffer, RNase-free	100 columns	4°C
732-6202	Micro Bio-Spin 30 Columns, in SSC buffer	25 columns	4°C
732-6203	Micro Bio-Spin 30 Columns, in SSC buffer	100 columns	4°C
732-6206	Micro Bio-Spin 30 Columns, in SSC buffer	1,000 columns	4°C
732-6227	Bio-Spin 6 Columns, in Tris buffer	25 columns	4°C
732-6228	Bio-Spin 6 Columns, in Tris buffer	100 columns	4°C
732-6229	Bio-Spin 6 Columns, in Tris buffer	1,000 columns	4°C
732-6231	Bio-Spin 30 Columns, in Tris buffer	25 columns	4°C
732-6232	Bio-Spin 30 Columns, in Tris buffer	100 columns	4°C
732-6233	Bio-Spin 30 Columns, in Tris buffer	1,000 columns	4°C
732-6002	Bio-Spin 6 Columns, in SSC buffer	25 columns	4°C
732-6006	Bio-Spin 30 Columns, in SSC buffer	25 columns	4°C
732-2010	Econo-Pac 10DG Desalting Columns	30 columns	Ambient

* Bio-Spin and Micro Bio-Spin columns can also be used for protein sample preparation.

Agarose Powders and Agarose Precast Gels

Catalog #	Description	Size/ Quantity	Storage Temp.
Agarose Powder			
161-3101	Certified™ Molecular Biology Agarose	125 g	Ambient
161-3104	Certified PCR Agarose	125 g	Ambient
161-3106	Certified Low Range Ultra Agarose	25 g	Ambient
161-3109	Certified Megabase Agarose	125 g	Ambient
161-3112	Certified Low-Melt Agarose	125 g	Ambient
161-3114	Certified PCR Low-Melt Agarose	125 g	Ambient
Mini ReadyAgarose™ Gels, TBE			
161-3001	ReadyAgarose 0.8% TBE Gel, 8-well	1 gel	4–22°C
161-3002	ReadyAgarose 0.8% TBE Gel, 8-well, with ethidium bromide	1 gel	4–22°C
161-3003	ReadyAgarose 1.0% TBE Gel, 8-well	1 gel	4–22°C
161-3004	ReadyAgarose 1.0% TBE Gel, 8-well, with ethidium bromide	1 gel	4–22°C
161-3005	ReadyAgarose 3.0% TBE Gel, 8-well	1 gel	4–22°C
161-3006	ReadyAgarose 3.0% TBE Gel, 8-well, with ethidium bromide	1 gel	4–22°C
161-3007	ReadyAgarose 0.8% TBE Gel, 12-well	1 gel	4–22°C
161-3008	ReadyAgarose 0.8% TBE Gel, 12-well, with ethidium bromide	1 gel	4–22°C
161-3009	ReadyAgarose 1.0% TBE Gel, 12-well	1 gel	4–22°C
161-3010	ReadyAgarose 1.0% TBE Gel, 12-well, with ethidium bromide	1 gel	4–22°C
161-3011	ReadyAgarose 3.0% TBE Gel, 12-well	1 gel	4–22°C
161-3012	ReadyAgarose 3.0% TBE Gel, 12-well, with ethidium bromide	1 gel	4–22°C
Mini ReadyAgarose Gels, TAE			
161-3013	ReadyAgarose 0.8% TAE Gel, 8-well	1 gel	4–22°C
161-3014	ReadyAgarose 0.8% TAE Gel, 8-well, with ethidium bromide	1 gel	4–22°C
161-3015	ReadyAgarose 1.0% TAE Gel, 8-well	1 gel	4–22°C
161-3016	ReadyAgarose 1.0% TAE Gel, 8-well, with ethidium bromide	1 gel	4–22°C
161-3017	ReadyAgarose 3.0% TAE Gel, 8-well	1 gel	4–22°C
161-3018	ReadyAgarose 3.0% TAE Gel, 8-well, with ethidium bromide	1 gel	4–22°C
161-3019	ReadyAgarose 0.8% TAE Gel, 12-well	1 gel	4–22°C
161-3020	ReadyAgarose 0.8% TAE Gel, 12-well, with ethidium bromide	1 gel	4–22°C
161-3021	ReadyAgarose 1.0% TAE Gel, 12-well	1 gel	4–22°C
161-3022	ReadyAgarose 1.0% TAE Gel, 12-well, with ethidium bromide	1 gel	4–22°C
161-3023	ReadyAgarose 3.0% TAE Gel, 12-well	1 gel	4–22°C
161-3024	ReadyAgarose 3.0% TAE Gel, 12-well, with ethidium bromide	1 gel	4–22°C
161-3057	ReadyAgarose 1.0% TAE Gel, 2 x 8-well	1 gel	4–22°C
Wide Mini ReadyAgarose Gels, TBE			
161-3025	ReadyAgarose 0.8% TBE Gel, 20-well	1 gel	4–22°C
161-3026	ReadyAgarose 0.8% TBE Gel, 20-well, with ethidium bromide	1 gel	4–22°C
161-3027	ReadyAgarose 1.0% TBE Gel, 20-well	1 gel	4–22°C
161-3028	ReadyAgarose 1.0% TBE Gel, 20-well, with ethidium bromide	1 gel	4–22°C
161-3029	ReadyAgarose 3.0% TBE Gel, 20-well	1 gel	4–22°C
161-3030	ReadyAgarose 3.0% TBE Gel, 20-well, with ethidium bromide	1 gel	4–22°C

Catalog #	Description	Size/ Quantity	Storage Temp.
161-3031	ReadyAgarose 0.8% TBE Gel, 32-well	1 gel	4–22°C
161-3032	ReadyAgarose 0.8% TBE Gel, 32-well, with ethidium bromide	1 gel	4–22°C
161-3033	ReadyAgarose 1.0% TBE Gel, 32-well	1 gel	4–22°C
161-3034	ReadyAgarose 1.0% TBE Gel, 32-well, with ethidium bromide	1 gel	4–22°C
161-3035	ReadyAgarose 3.0% TBE Gel, 32-well	1 gel	4–22°C
161-3036	ReadyAgarose 3.0% TBE Gel, 32-well, with ethidium bromide	1 gel	4–22°C
161-3037	ReadyAgarose 1.0% TBE Gel, 2 x 32-well	1 gel	4–22°C
161-3038	ReadyAgarose 1.0% TBE Gel, 2 x 32-well, with ethidium bromide	1 gel	4–22°C
161-3039	ReadyAgarose 3.0% TBE Gel, 2 x 32-well	1 gel	4–22°C
161-3040	ReadyAgarose 3.0% TBE Gel, 2 x 32-well, with ethidium bromide	1 gel	4–22°C
Wide Mini ReadyAgarose Gels, TAE			
161-3041	ReadyAgarose 0.8% TAE Gel, 20-well	1 gel	4–22°C
161-3042	ReadyAgarose 0.8% TAE Gel, 20-well, with ethidium bromide	1 gel	4–22°C
161-3043	ReadyAgarose 1.0% TAE Gel, 20-well	1 gel	4–22°C
161-3044	ReadyAgarose 1.0% TAE Gel, 20-well, with ethidium bromide	1 gel	4–22°C
161-3045	ReadyAgarose 3.0% TAE Gel, 20-well	1 gel	4–22°C
161-3046	ReadyAgarose 3.0% TAE Gel, 20-well, with ethidium bromide	1 gel	4–22°C
161-3047	ReadyAgarose 0.8% TAE Gel, 32-well	1 gel	4–22°C
161-3048	ReadyAgarose 0.8% TAE Gel, 32-well, with ethidium bromide	1 gel	4–22°C
161-3049	ReadyAgarose 1.0% TAE Gel, 32-well	1 gel	4–22°C
161-3050	ReadyAgarose 1.0% TAE Gel, 32-well, with ethidium bromide	1 gel	4–22°C
161-3051	ReadyAgarose 3.0% TAE Gel, 32-well	1 gel	4–22°C
161-3052	ReadyAgarose 3.0% TAE Gel, 32-well, with ethidium bromide	1 gel	4–22°C
161-3053	ReadyAgarose 1.0% TAE Gel, 2 x 32-well	1 gel	4–22°C
161-3054	ReadyAgarose 1.0% TAE Gel, 2 x 32-well, with ethidium bromide	1 gel	4–22°C
161-3055	ReadyAgarose 3.0% TAE Gel, 2 x 32-well	1 gel	4–22°C
161-3056	ReadyAgarose 3.0% TAE Gel, 2 x 32-well, with ethidium bromide	1 gel	4–22°C
ReadyAgarose™ 96 Plus Gels, TBE			
161-3060	ReadyAgarose 96 Plus 1% TBE Gel, with ethidium bromide	1 gel	4–22°C
161-3062	ReadyAgarose 96 Plus 3% TBE Gel, with ethidium bromide	1 gel	4–22°C
ReadyAgarose 96 Plus Gels, TAE			
161-3063	ReadyAgarose 96 Plus 1% TAE Gel, with ethidium bromide	1 gel	4–22°C
161-3065	ReadyAgarose 96 Plus 3% TAE Gel, with ethidium bromide	1 gel	4–22°C

Gene Transfer

Catalog #	Description	Size/ Quantity	Storage Temp.
Electro-Competent Cells			
170-3330	EP-Max™10B Electro-Competent Cells	5 x 0.10 ml	-70°C
170-3331	EP-Max10B F' Electro-Competent Cells	5 x 0.10 ml	-70°C
170-3332	EP-Max10B T1 -Resistant Electro-Competent Cells	5 x 0.10 ml	-70°C
Chemi-Competent Cells			
170-3340	C-Max™5□ Chemi-Competent Cells	5 x 0.20 ml	-70°C
170-3342	C-Max5□ Chemi-Competent Cells	10 x 0.05 ml	-70°C
170-3341	C-Max5□ F' Chemi-Competent Cells	5 x 0.20 ml	-70°C
170-3343	C-Max5□ F' Chemi-Competent Cells	10 x 0.05 ml	-70°C
Electroporation Cuvettes			
165-2083	Gene Pulser®/MicroPulser™ Cuvettes, 0.1 cm gap	5 cuvettes	Ambient
165-2089	Gene Pulser/MicroPulser Cuvettes, 0.1 cm gap	50 cuvettes	Ambient
165-2093	Gene Pulser/MicroPulser Cuvettes, 0.1 cm gap	500 cuvettes	Ambient
165-2082	Gene Pulser/MicroPulser Cuvettes, 0.2 cm gap	5 cuvettes	Ambient
165-2086	Gene Pulser/MicroPulser Cuvettes, 0.2 cm gap	50 cuvettes	Ambient
165-2092	Gene Pulser/MicroPulser Cuvettes, 0.2 cm gap	500 cuvettes	Ambient
165-2081	Gene Pulser/MicroPulser Cuvettes, 0.4 cm gap	5 cuvettes	Ambient
165-2088	Gene Pulser/MicroPulser Cuvettes, 0.4 cm gap	50 cuvettes	Ambient
165-2091	Gene Pulser/MicroPulser Cuvettes, 0.4 cm gap	500 cuvettes	Ambient
Lipofection			
170-3350	TransFectin™ Lipid Reagent	0.5 ml	4°C
170-3351	TransFectin Lipid Reagent	1.0 ml	4°C
170-3352	TransFectin Lipid Reagent	5 x 1.0 ml	4°C
Fluorescent Assay Kits			
170-2480	DNA Quantitation Kit (Hoechst)	1 kit	4°C
170-3150	FluorAce™ □-Galactosidase Reporter Assay Kit	1 kit	4°C
170-3151	FluorAce □-Glucuronidase Reporter Assay Kit	1 kit	4°C

Amplification

Catalog #	Description	Size/ Quantity	Storage Temp.
Reagents			
170-8870	iTaq™ DNA Polymerase	250 units	-20°C
170-8875	iTaq DNA Polymerase	5,000 units	-20°C
170-8872	MgCl ₂ Solution, 50 mM	1.25 ml	Ambient
170-8874	dNTP Mix, 10 mM each dNTP	200 ml	-20°C
170-8860	iQ Supermix	100 reactions	-20°C
170-8862	iQ Supermix	500 reactions	-20°C
170-8864	iQ Supermix	1,000 reactions	-20°C
170-8880	iQ™ SYBR® Green Supermix	100 reactions	-20°C
170-8882	iQ SYBR Green Supermix	500 reactions	-20°C
170-8884	iQ SYBR Green Supermix	1,000 reactions	-20°C
170-8850	iTaq SYBR Green Supermix With ROX	200 reactions	-20°C
170-8851	iTaq SYBR Green Supermix With ROX	500 reactions	-20°C
170-8852	iTaq SYBR Green Supermix With ROX	1,000 reactions	-20°C
170-8890	iScript™ cDNA Synthesis Kit	25 reactions	-20°C
170-8891	iScript cDNA Synthesis Kit	100 reactions	-20°C
170-8892	iScript One-Step RT-PCR Kit With SYBR Green	50 reactions	-20°C
170-8893	iScript One-Step RT-PCR Kit With SYBR Green	200 reactions	-20°C
170-8894	iScript One-Step RT-PCR Kit for Probes	50 reactions	-20°C
170-8895	iScript One-Step RT-PCR Kit for Probes	200 reactions	-20°C
Micro Test Tubes and Plates			
223-9441	96-Well 0.2 ml Thin-Wall PCR Plates	25 plates	Ambient
223-9442	96-Well 0.2 ml PCR Plate Mats	5 mats	Ambient
223-9443	96-Well 0.2 ml PCR Plate Caps	300 strips	Ambient
223-9444	Optical-Quality Sealing Tape	100 sheets	Ambient
223-9466	0.6 ml Thin-Wall Tubes, natural	1,000 tubes	Ambient
223-9469	0.2 ml Thin-Wall PCR 8-Tube Strips, natural	120 strips	Ambient
223-9472	0.2 ml Thin-Wall PCR 8-Cap Strips, natural	120 strips	Ambient
223-9473	0.2 ml Thin-Wall Tubes, natural	1,000 tubes	Ambient
223-9476	0.2 ml Thin-Wall PCR 12-Tube Strips, natural	80 strips	Ambient
223-9477	0.2 ml Thin-Wall PCR 12-Cap Strips, natural	80 strips	Ambient
223-9495	384-Well PCR Plates	25 plates	Ambient
223-9496	384-Well PCR Plate Mats	25 mats	Ambient

Competent Cell Product Guide

Application	Helpful Mutation/Marker	Benefit	Host Cell
Genomic library construction	<i>mrr-hsdRMS-mcrBC</i> <i>tonA</i>	Allows cloning of methylated DNA Confers resistance to bacteriophage infection	EP-Max EP-Max T1-Resistant
cDNA library construction	<i>tonA</i>	Confers resistance to bacteriophage infection	EP-Max
Plasmid rescue/propagation	<i>recA</i> <i>endA</i>	Reduced risk of recombination events Higher-quality, higher-copy plasmid preps	EP-Max EP-Max/C-Max
ssDNA synthesis	F' episome	Allows infection by bacteriophage M13	EP-Max F', C-Max F'
Selection for transformants	Tn10 <i>lacZ</i> Δ <i>M15</i>	Tetracycline resistance (allows selection of F' hosts) Allows blue/white screening by β -complementation	EP-Max F', C-Max F' EP-Max, C-Max
Phage display	F' episome, <i>supE</i>	Required for growth of amber mutants	C-Max F'
<i>phoA</i> insertions, <i>phoA</i> fusion analysis	<i>phoA</i>	Absence of alkaline phosphatase activity	C-Max
Expression of toxic proteins	<i>lacI^q</i> gene	Permits regulated expression from <i>lac</i> , <i>tac</i> , and <i>trc</i> promoters	EP-Max F', C-Max F'

Place Magnet here

Bio-Rad Representative Business Card

Trademarks

Affi-Gel, Affi-Prep, AmpliSize, AquaPure, Aurum, Bio-Gel, Bio-Lyte, Bio-Safe, Bio-Scale, Bio-Sil, Bio-Silect, Bio-Spin, Biotin-Blot, Certified, Chelex, CHT, C-Max, Criterion, Criterion XT, *DC*, De-Expose, Dodeca, Econo-Column, Econo-Pac, EP-Max, EZ Load, EZ Load HT, FluorAce, Freeze 'N Squeeze, Gene Pulser, Immun-Blot, Immun-Star, InstaGene, iQ, iScript, iTaq, Kaleidoscope, Macro-Prep, MAPS, Master Blaster, Micro Bio-Spin, MicroPulser, Mini Trans-Blot, Opti-4CN, PCR Kleen, Poly-Prep, Precision Plus Protein, PROTEAN, PROTEAN Plus, Quantum Prep, Quick Start, Radiant, *RC DC*, ReadyAgarose, ReadyAgarose 96 Plus, Ready Gel, ReadyPrep, ReadyStrip, SEQueaky Kleen, Sequi-Blot, Silver Stain Plus, TransFectin, UNO, UNOsphere, Zeta-Probe, and Zeta-Probe GT are trademarks or registered trademarks of Bio-Rad Laboratories, Inc.

Bio-Rad Laboratories, Inc. is licensed by Molecular Probes, Inc. to sell SYPRO products for research use only, under US patent 5,616,502. Bio-Rad Laboratories, Inc. is licensed by Molecular Probes, Inc. to sell reagents containing SYBR Green I for use in real-time PCR, for research purposes only.

Coomassie is a trademark of BASF Aktiengesellschaft. PEEK is a trademark of Victrex PLC. StrepTactin is a trademark of Institut für Bioanalytik GmbH. SYBR, SYPRO, and Texas Red are trademarks of Molecular Probes, Inc. Tween is a trademark of ICI Americas, Inc.

BIO-RAD

**Bio-Rad
Laboratories, Inc.**

Life Science
Group

Web site www.bio-rad.com **USA** (800) 4BIORAD **Australia** 02 9914 2800 **Austria** (01)-877 89 01 **Belgium** 09-385 55 11 **Brazil** 55 21 2527 3454 **Canada** (905) 712-2771 **China** (86-21) 63052255 **Czech Republic** + 420 2 41 43 05 32 **Denmark** 44 52 10 00 **Finland** 09 804 22 00 **France** 01 47 95 69 65 **Germany** 089 318 84-177 **Hong Kong** 852-2789-3300 **Hungary** 36 1 455 8800 **India** (91-124)-6398112/113/114, 6450092/93 **Israel** 03 951 4127 **Italy** 39 02 216091 **Japan** 03-5811-6270 **Korea** 82-2-3473-4460 **Latin America** 305-894-5950 **Mexico** 55-52-00-05-20 **The Netherlands** 0318-540666 **New Zealand** 64 9 415 2280 **Norway** 23 38 41 30 **Poland** + 48 22 331 99 99 **Portugal** 351-21-472-7700 **Russia** 7 095 721 1404 **Singapore** 65-6415 3188 **South Africa** 00 27 11 4428508 **Spain** 34 91 590 5200 **Sweden** 08 555 12700 **Switzerland** 061 717-9555 **Taiwan** (8862) 2578-7189/2578-7241 **United Kingdom** 020 8328 2000