

ELECTROPHORESIS

NuPAGE Bis-Tris Precast Gels to Mini-PROTEAN™ TGX and Mini-PROTEAN TGX Stain-Free Precast Gels

Catalog Number Conversion Chart




NuPAGE Bis-Tris Precast Gels (Box of 10)					Mini-PROTEAN TGX Precast Gels (Box of 10)			Mini-PROTEAN TGX Stain-Free Precast Gels (Box of 10)		
Description		Well Volume	Catalog #		Catalog #	Well Volume	Description	Catalog #	Well Volume	Description
10% (MOPS buffer)	NuPAGE Novex 10% Bis-Tris, 1.0 mm, 9-well	28 µl	NP0307BOX	→	4561043	30 µl	Mini-PROTEAN TGX, 12%, 10-well	4568043	30 µl	Mini-PROTEAN TGX Stain-Free, 12%, 10-well
	NuPAGE Novex 10% Bis-Tris, 1.0 mm, 10-well	25 µl	NP0301BOX	→	4561043	30 µl	Mini-PROTEAN TGX, 12%, 10-well	4568043	30 µl	Mini-PROTEAN TGX Stain-Free, 12%, 10-well
	NuPAGE Novex 10% Bis-Tris, 1.5 mm, 10-well	37 µl	NP0315BOX	→	4561044	50 µl	Mini-PROTEAN TGX, 12%, 10-well	4568044	50 µl	Mini-PROTEAN TGX Stain-Free, 12%, 10-well
	NuPAGE Novex 10% Bis-Tris, 1.0 mm, 12-well	20 µl	NP0302BOX	→	4561045	20 µl	Mini-PROTEAN TGX, 12%, 12-well	4568045	20 µl	Mini-PROTEAN TGX Stain-Free, 12%, 12-well
	NuPAGE Novex 10% Bis-Tris, 1.0 mm, 15-well	15 µl	NP0303BOX	→	4561046	15 µl	Mini-PROTEAN TGX, 12%, 15-well	4568046	15 µl	Mini-PROTEAN TGX Stain-Free, 12%, 15-well
	NuPAGE Novex 10% Bis-Tris, 1.5 mm, 15-well	25 µl	NP0316BOX	→	4561046	15 µl	Mini-PROTEAN TGX, 12%, 15-well	4568046	15 µl	Mini-PROTEAN TGX Stain-Free, 12%, 15-well
	NuPAGE Novex 10% Bis-Tris, 1.0 mm, 1-well	700 µl	NP0304BOX	→	4561041	7 cm strip/450 µl	Mini-PROTEAN TGX, 12%, IPG/prep	4568041	7 cm strip/450 µl	Mini-PROTEAN TGX Stain-Free, 12%, IPG/prep
12% (MOPS buffer)	NuPAGE Novex 12% Bis-Tris, 1.0 mm, 10-well	25 µl	NP0341BOX	→	4569033	30 µl	Mini-PROTEAN TGX Any kD, 10-well	4568123	30 µl	Mini-PROTEAN TGX Any kD Stain-Free, 10-well
	NuPAGE Novex 12% Bis-Tris, 1.0 mm, 12-well	20 µl	NP0342BOX	→	4569035	20 µl	Mini-PROTEAN TGX Any kD, 12-well	4568125	20 µl	Mini-PROTEAN TGX Any kD Stain-Free, 12-well
	NuPAGE Novex 12% Bis-Tris, 1.0 mm, 15-well	15 µl	NP0343BOX	→	4569036	15 µl	Mini-PROTEAN TGX Any kD, 15-well	4568126	15 µl	Mini-PROTEAN TGX Any kD Stain-Free, 15-well
	NuPAGE Novex 12% Bis-Tris, 1.0 mm, 17-well	15 µl	NP0349BOX	→	4569036	15 µl	Mini-PROTEAN TGX Any kD, 15-well	4568126	15 µl	Mini-PROTEAN TGX Any kD Stain-Free, 15-well
	NuPAGE Novex 12% Bis-Tris, 1.0 mm, 1-well	700 µl	NP0344BOX	→	4569031	7 cm strip/450 µl	Mini-PROTEAN TGX Any kD, IPG/prep	4568121	7 cm strip/450 µl	Mini-PROTEAN TGX Any kD Stain-Free, IPG/prep
4–12% (MOPS buffer)	NuPAGE Novex 4–12% Bis-Tris, 1.0 mm, 9-well	28 µl	NP0327BOX	→	4561093	30 µl	Mini-PROTEAN TGX, 4–20%, 10-well	4568093	30 µl	Mini-PROTEAN TGX Stain-Free, 4–20%, 10-well
	NuPAGE Novex 4–12% Bis-Tris, 1.0 mm, 10-well	25 µl	NP0321BOX	→	4561093	30 µl	Mini-PROTEAN TGX, 4–20%, 10-well	4568093	30 µl	Mini-PROTEAN TGX Stain-Free, 4–20%, 10-well
	NuPAGE Novex 4–12% Bis-Tris, 1.5 mm, 10-well	37 µl	NP0335BOX	→	4561094	50 µl	Mini-PROTEAN TGX, 4–20%, 10-well	4568094	50 µl	Mini-PROTEAN TGX Stain-Free, 4–20%, 10-well
	NuPAGE Novex 4–12% Bis-Tris, 1.0 mm, 12-well	20 µl	NP0322BOX	→	4561095	20 µl	Mini-PROTEAN TGX, 4–20%, 12-well	4568095	20 µl	Mini-PROTEAN TGX Stain-Free, 4–20%, 12-well
	NuPAGE Novex 4–12% Bis-Tris, 1.0 mm, 15-well	15 µl	NP0323BOX	→	4561096	15 µl	Mini-PROTEAN TGX, 4–20%, 15-well	4568096	15 µl	Mini-PROTEAN TGX Stain-Free, 4–20%, 15-well
	NuPAGE Novex 4–12% Bis-Tris, 1.5 mm, 15-well	25 µl	NP0336BOX	→	4561096	15 µl	Mini-PROTEAN TGX, 4–20%, 15-well	4568096	15 µl	Mini-PROTEAN TGX Stain-Free, 4–20%, 15-well
	NuPAGE Novex 4–12% Bis-Tris, 1.0 mm, 17-well	15 µl	NP0329BOX	→	4561096	15 µl	Mini-PROTEAN TGX, 4–20%, 15-well	4568096	15 µl	Mini-PROTEAN TGX Stain-Free, 4–20%, 15-well
	NuPAGE Novex 4–12% Bis-Tris, 1.0 mm, 1-well	700 µl	NP0324BOX	→	4561091	7 cm strip/450 µl	Mini-PROTEAN TGX, 4–20%, IPG/prep	4568091	7 cm strip/450 µl	Mini-PROTEAN TGX Stain-Free, 4–20%, IPG/prep

Mini-PROTEAN TGX and TGX Stain-Free Precast Gels

The Mini-PROTEAN TGX (Tris-Glycine eXtended) and TGX Stain-Free Gels deliver accelerated run times with an extended shelf life of 12 months

- Run times as short as 15 min
- Transfer in as little as 15 min
- Formulated for use with Laemmli buffers
- Inexpensive buffer system, low running costs
- Robust system for difficult samples
- Bottom-open cassette design for ease of gel handling and blotting
- Compatible with Mini-PROTEAN Tetra and Mini-PROTEAN 3 Dodeca Systems — ability to run up to 12 mini gels at once
- Compatible with Bio-Rad's Stain-Free Western Blotting Workflow



Side-by-side comparison of migration patterns for unstained standards with a broad range of molecular weights run on NuPAGE (■) and Mini-PROTEAN TGX and TGX Stain-Free Precast Gels (■).

Related Literature

1658100	Mini-PROTEAN Precast Gels Instruction Manual and Application Guide	5911	Mini-PROTEAN TGX Gel: A Versatile and Robust Laemmli-Like Precast Gel for SDS-PAGE
5535	Mini-PROTEAN Electrophoresis System Brochure	3668	ChemiDoc™ Go Imaging System Brochure
5871	Mini-PROTEAN TGX Precast Gels Product Information Sheet	7284	GelDoc™ Go Imaging System Brochure
5910	Mini-PROTEAN TGX Precast Gel: A Gel for SDS-PAGE with Improved Stability — Comparison with Standard Laemmli Gels	2998	Protein Standards Selection Guide
		2414	The Little Book of Standards
		6040	A Guide to Polyacrylamide Gel Electrophoresis and Detection
		2895	Protein Blotting Guide

Visit [bio-rad.com/MiniPROTEAN](https://www.bio-rad.com/MiniPROTEAN) for more information.

Still unsure what gels would work best for your experiment? Try the [Free Western Blotting Selector Tool](#).

BIO-RAD, CHEMIDOC, GELDOC, and MINI-PROTEAN are trademarks of Bio-Rad Laboratories, Inc. in certain jurisdictions. All trademarks used herein are the property of their respective owner. TGX Stain-Free Precast Gels are covered by U.S. Patent Numbers 7,569,130 and 8,007,646. © 2024 Bio-Rad Laboratories, Inc.